



Laboratory Procedure Manual

Analyte: **Vitamin B₆ (pyridoxal 5'-phosphate, 4-pyridoxic acid)**

Matrix: **Serum**

Method: **HPLC with fluorescence detection**

Method No:

Revised:

as performed by:

Inorganic Toxicology and Nutrition Branch
Division of Laboratory Sciences
National Center for Environmental Health

contact:

[Dr. Eric J. Sampson, Director](#)
[Division of Laboratory Sciences](#)

Important Information for Users

CDC periodically refines these laboratory methods. It is the responsibility of the user to contact the person listed on the title page of each write-up before using the analytical method to find out whether any changes have been made and what revisions, if any, have been incorporated.

Public Release Data Set Information

This document details the Lab Protocol for testing the items listed in the following table:

File Name	Variable Name	SAS Label
VITB6_D	LBX4PA	4-pyridoxic acid (nmol/L)
	LBXPLP	Pyridoxal 5'-phosphate (nmol/L)

1. Summary of Test Principle and Clinical Relevance

a. Clinical relevance

First identified in 1934 by György as the preventative factor for rat pellagra [1], vitamin B₆ presently refers to a collective of six biologically interconvertible 3-hydroxy-2-methylpyridine compounds: pyridoxal (PL); pyridoxine (PN); pyridoxamine (PM) and their respective 5'-phosphates (PLP, PNP, PMP) [2]. Of these compounds, PLP is the primary biologically active form of vitamin B₆, serving as a coenzyme in numerous biologic processes that include tryptophan-niacin conversion (kynureninase), heme synthesis (δ -aminolevulinic acid synthetase), gluconeogenesis (glycogen phosphorylase), neurotransmitter synthesis, and amino acid metabolism (various transaminases and decarboxylases) [3]. In addition, PLP acts as a coenzyme in the catabolism of homocysteine to cystathionine (cystathionine β -synthase) and ultimately cysteine (cystathionine γ -lyase) [3]. As a result, vitamin B₆ plays a role in lowering homocysteine levels and its associated risk with atherosclerosis and coronary heart disease [4,5], although the clinical evidence of this is somewhat inconsistent [6–8]. Additionally, low PLP status may be associated with an increased risk for stroke and transient ischemic attack, irrespective of homocysteine levels [9].

Clinical vitamin B₆ status typically is assessed by measuring the level of one or more of the B₆ vitamers in vivo. Plasma or serum PLP concentration traditionally has been used as a status indicator because it is the active coenzyme form of vitamin B₆. Circulating PLP levels were also believed to correlate with indicators of long-term B₆ storage in the body, such as the amount of glycogen phosphorylase-bound PLP in skeletal muscle [10] but more recent evidence suggests that this is not the case [11,12]. As a result, the use of PLP alone as a status indicator has been contested, and the measurement of PLP along with other B₆ vitamers and indicators has been proposed [13]. An example is the measurement of the end product of vitamin B₆ catabolism, 4-pyridoxic acid (4-PA), in either urine [14] or plasma/serum as an indicator of short-term status.

b. Test principle

Serum is prepared for analysis by first combining the specimen 1:1 with a 5% solution of metaphosphoric acid to precipitate proteins out of the sample. After vortex mixing and centrifuging, the sample supernatant is combined with an approximately equivalent amount of dichloromethane to extract lipids from the sample, and vortex mixed and centrifuged a second time. The top (aqueous) layer from the sample is then filtered through a syringe and ready for high-performance liquid chromatographic (HPLC) analysis.

Vitamin B₆, in the form of PLP, and the metabolite 4-PA are measured by reversed-phase HPLC using fluorometric detection at 325 nm excitation and 425 nm emission. Post-column introduction of a sodium chlorite derivatization reagent is incorporated into the HPLC system to improve the PLP signal. Quantitation is based on analyte peak area interpolated against a five-point calibration curve obtained from aqueous standards.

2. Safety Precautions

Consider all plasma or serum specimens potentially positive for infectious agents including HIV and the hepatitis B virus. We recommend the hepatitis B vaccination series for all analysts working with whole blood and/or plasma. Observe universal precautions; wear protective gloves, laboratory coats, and safety glasses during all steps of this method. Discard any residual sample material by autoclaving after analysis is completed. Place disposable plastic, glass, and paper (pipet tips, autosampler vials, gloves, etc.) that contact plasma in a biohazard autoclave bag and keep these bags in appropriate containers until sealed and autoclaved. Wipe down all work surfaces with 10% bleach solution when work is finished.

Handle acids and bases with extreme care; they are caustic and toxic. Handle organic solvents only in a well-ventilated area or, as required, under a chemical fume hood.

Reagents and solvents used in this study include those listed in Section 6. Material safety data sheets (MSDSs) for these chemicals are readily accessible as hard copies in the lab. If needed, MSDS for other chemicals can be viewed at <http://www.ilpi.com/msds/index.html> or at <http://intranet.cdc.gov/ohs>.

3. Computerization; Data System Management

- a. Raw data collection, peak integration and analyte quantitation is performed using the HPLC operating software (Agilent ChemStation). After a run is complete and any additional corrections by the analyst are made, a Microsoft Excel results file (containing the patient data as well as the QC data) is generated from the HPLC operating software and electronically transferred to Q:\LabDeviceOutput\NCEHNUTR10\ on the NCEH/DLS Local Area Network (LAN). The analyst also prepares a hardcopy of the results from the HPLC operating software that is to be circulated through the QA/QC process. This includes: a table summarizing the metrics of the analyte signals for each calibrator, QC and patient sample; calibration curves and equations; and the chromatogram for each calibrator, QC and patient sample. The analyst then imports the Excel file containing the run data into the Microsoft Access Network Database (FrontEnds) in the method described below:

Step 1 – Analyst – Import data file into ACCESS:

Double click the ACCESS icon on desktop, password entry required

[Add Sample Results to Database] (under Batch & X-Batch)

[Import Instrument Data File] - Enter information (instrument, assay, date, time, analyst, study)

[Import] – In “select data file” window, choose A: and import file number assigned. Check that sample ID's are recognized.

[Transfer]

Step 2 – Analyst – Review run in ACCESS:

[Run Review] (under Batch & X-Batch) – Select assay

[Show runs] – Cursor to desired run, enter sample set name and comments

[QC Results] – Review QC results for transmission errors and whether they pass the 2S limits

[Print Report] [Back]

[Sample Results] – Review patient results to assure proper information transmission, enter appropriate comment codes on flagged samples

[Set Final] results that are ready to be reported

[Set Reviewed]

[Print Report] [Back]

Step 3 – Analyst – Send email and run folder to QA Officer:

An e-mail is sent to the QA Officer including the following run information: Analysis date, Instrument, Study, Groups, File name, Batch ID, Run #, and QC Status. Noteworthy comments are included in the email. All printouts including raw data are submitted in a run folder to the QA Officer who reviews the Bench QC data via the ACCESS database as described below.

Step 4 – QA Officer – Review Bench QC via ACCESS:

Double click the ACCESS icon on desktop, password entry required

[Export QC to SAS] (under Batch & X-Batch) – Select Assay, Date range and Controls

[Make QC Data Infile] – Save file to I: appropriate subfolder for archival

[Run SAS] – SAS will automatically open, [go], review each generated plot, print QC cover page and standard deviation plot, [Back]

[Run Review] (under Batch & X-Batch) – Select assay

[Show runs]

[Sample Results]

[Set Batch QC] – accept or reject

[Set Reviewed]

Forward email from Analyst to Second QA reviewer (for Blind QC review) specifying Bench QC status of the run.

Step 5 – Second QA Reviewer – Review Blind QC and other parameters in ACCESS:

Double click the ACCESS icon on desktop, password entry required

[Run Review] (under Batch & X-Batch) – Select assay, then desired run

[Blind QC Results] – Review whether Blind QC results pass the 2S limits

[Print Report] [Back]

Check other parameters if applicable (i.e., background, calibration curve, repeat values, replicates, signal intensity)

[Set RQC] – accept or reject

Verify that appropriate comment codes have been applied and that final values have been set correctly

[Set Reviewed]

Forward email from QA Officer to Supervisor specifying Blind QC status of the run and other relevant comments.

Step 6 – Supervisor – Approval and Export of Results via ACCESS:

Double click the ACCESS icon on desktop, password entry required
[Run Review] (under Batch & X-Batch) – Select assay, then desired run
Perform final review of Bench and Blind QC status, comment codes, repeat results
[Set Ready] – Final results will be set ready to be exported
[Set Reviewed]
[Export/Report Results] (under Study Functions) – Select study, select analytes/panel, use selected panel
[Generate Excel Spreadsheet] – Review file on I:\To be transmitted
[Generate Export Text File and Set Results Exported] – Review file on I:\To be transmitted
FTP file to Westat
Send Westat an email that file was transmitted
After min. of 1 day, move transmittal file from I:\To be transmitted to I:\Transmitted Data\Appropriate Year Folder.

For NHANES, data is transmitted electronically several times weekly to Westat's ISIS computer system, and transferred from there to NCHS. Abnormal values are confirmed by the analyst, and codes for missing data are entered by the analyst and are transmitted as part of the data file to the Westat ISIS computer, and are eventually forwarded to NCHS. Westat also prepares the abnormal report notifications for the NCHS Survey Physician.

- b. Backup of the daily data containing all raw data files and result files for each run are the responsibility of the analyst. These files are typically backed up to the CDC network on a daily basis and periodically backed up on CD. All sample, QC, and calibration data are stored on the CDC network and are the responsibility of the analyst. Files stored on the network or CDC mainframe are automatically backed up nightly by DLS LAN support staff and CDC Data Center Staff, respectively.
- c. Documentation for data system maintenance is contained in printed copies of data records, as well as in "system log" files on the local hard drives used for the archival of data.

4. Specimen Collection, Storage, and Handling Procedures; Criteria for Specimen Rejection

- a. For best results, a fasting sample should be obtained.
- b. Specimens for vitamin B₆ analysis may be fresh or frozen serum. In the absence of serum, heparinized plasma may be used. EDTA plasma should be avoided due to a potential interference with the PLP signal.

- c. A 350- μ L sample of plasma or serum is preferable to allow for repeat analyses; a volume of 100 μ L is required for analysis.
- d. The appropriate amount of plasma or serum is dispensed into a Nalgene 2.0-mL cryovial or other plastic screw-capped vial labeled with the participant's ID.
- e. Specimens collected in the field are frozen, and then shipped on dry ice by overnight carrier. Frozen samples are stored at -70°C. Samples are stable for at least 5 years if stored at -70°C. They can withstand at least 5 freeze/thaw cycles [15].
- f. Specimens generally arrive frozen. Refrigerated samples may be used provided they are kept cold and brought promptly (within 2 hours) from the site of collection.
- g. Specimens that have been through more than five freeze-thaw cycles, been refrigerated for more than one week, or undergone hemolysis may give inaccurate results.
- h. Specimen handling conditions are outlined in the Policies and Procedures Manual of DLS (copies are available in the Nutritional Laboratory and the electronic copy of this file is located at Q: /ITN/Nutrition Laboratory/CLIA). The protocol discusses collection and transport of specimens and the special equipment required. In general, plasma should be transported and stored at no more than -20°C. Samples thawed and refrozen less than five times are not compromised. If there is more than one analyte of interest in the specimen and it needs to be divided, the appropriate amount of blood or plasma should be transferred into a sterile Nalgene cryovial labeled with the participant's ID.

5. Procedures for Microscopic Examinations; Criteria for Rejection of Inadequately Prepared Slides

Not applicable for this procedure

6. Preparation of Reagents, Calibration (Standards), Controls, and All Other Materials; Equipment and Instrumentation

a. Reagent Preparation

Prepare all reagents using deionized water with resistance of at least 15 megaOhm-cm.

1) 5% metaphosphoric acid (MPA) solution

CAUTION! *Metaphosphoric acid is corrosive. It is a contact and inhalation hazard. Consult MSDS before use.*

Verify the actual MPA content (%) of the solid reagent you are using. This is stated on the label, and is usually around 35%. Prepare the 5% MPA solution using the following formula:

$$\text{MPA required(solid, g)} = \left(\frac{\text{volume of solution being made (mL)}}{\text{actual MPA content of reagent (\%)}} \right) \times 5 \%$$

For example, if the actual MPA content of the reagent was 35% and you wanted to make 30 mL of reagent, you would add 4.29 g of the MPA reagent to 30 mL of water.

Prepare 30 mL of 5% MPA according to the formula above by adding the appropriate amount of MPA reagent to 30 mL of water in a 50 mL polypropylene screw-top centrifuge tube. Mix vigorously using a vortex mixer until the solid reagent dissolves. Refrigerate immediately and store at 4 °C at all times. Solutions maintained at 4 °C are good for 5 days.

2) 50 mM phosphate buffer, pH 3.2, 0.2% acetonitrile (mobile phase)

Weigh 6.9 g of sodium phosphate, monobasic, monohydrate (NaH₂PO₄·H₂O) into a 1-L glass beaker or bottle. Add 1000 mL of water using a graduated cylinder, and add 2 mL of acetonitrile using a 1-mL micropipette twice. Mix the solution using a magnetic stir bar until all of the sodium phosphate dissolves. Continue stirring the solution and measure its pH using a digital pH meter set to read to two decimal places. The pH should be around 4.6. Slowly add concentrated (85%) orthophosphoric acid dropwise until the pH of the solution is 3.2 (3.18–3.22 is acceptable). Vacuum filter the solution through a 0.45 μm polyvinylidene fluoride filter into a clean 1-L amber glass bottle and cap tightly. This solution will be good for 5 days.

3) 2 g/L sodium chlorite solution (post-column derivatization solution)

CAUTION!!! Sodium chlorite is a strong oxidant. Contact with combustible material may cause fire. Do not heat sodium chlorite with organic matter or strong reducing agents. Contact with acids releases toxic fumes. Sodium chlorite is toxic in contact with the skin and if swallowed. Consult MSDS before use.

Weigh 2.0 g of sodium chlorite into a 1-L glass beaker or bottle. Add 1000 mL of water, and mix the solution with a magnetic stir bar until all of the sodium chlorite dissolves. Vacuum filter the solution through a 0.45 μm polyvinylidene fluoride filter into a clean 1-L amber glass bottle and cap tightly. This solution will be good for two weeks.

b. Preparation of standards

NOTE: Vitamin B₆ is light sensitive! Perform all sample preparation under yellow lighting and in the absence of any direct sunlight. Use amber colored glassware and plasticware whenever possible.

1) Individual stock solutions

(a) Pyridoxal 5'-phosphate (PLP) – 100 ppm:

Using an electronic balance readable to 0.1 mg, transfer ~10 mg of PLP into a weighing boat and record the mass. Quantitatively transfer the PLP to a 100-mL

volumetric flask and dilute with some water. Mix the solution thoroughly until the 4-PA dissolves. Dilute to the mark with water.

(b) 4-Pyridoxic acid (4-PA) – 100 ppm:

Using an electronic balance readable to 0.1 mg, transfer ~10 mg of 4-PA into a weighing boat and record the mass. Quantitatively transfer the 4-PA to a 100-mL volumetric flask and dilute with some water. Mix the solution thoroughly until the 4-PA dissolves. Dilute to the mark with water. *NOTE: 4-PA dissolves very slowly in water. A sonicating bath can be used to help bring the 4-PA into solution. Avoid exposing the solution to excessive heat when doing this.*

2) Formulas to calculate the concentration (μmol/L) of the stock solutions

(a) PLP:

If the free form of PLP is used:

$$PLP \text{ stock solution } (\mu\text{mol} / L) = \frac{(\text{mass of PLP used (mg)}) \times (1000)}{(\text{final volume (L)}) \times (247.14 \text{ g / mol})}$$

If the monohydrate form of PLP is used (PLP-H₂O):

$$PLP \text{ stock solution } (\mu\text{mol} / L) = \frac{(\text{mass of PLP used (mg)}) \times (1000)}{(\text{final volume (L)}) \times (265.2 \text{ g / mol})}$$

(b) 4-PA:

$$4\text{-PA stock solution } (\mu\text{mol} / L) = \frac{(\text{mass of PLP used (mg)}) \times (1000)}{(\text{final volume (L)}) \times (183.16 \text{ g / mol})}$$

3) Working solutions

(a) PLP – 100 ppb:

Prepare a 1:1000 dilution of the stock solution by transferring 1000 μL of the stock solution with a micropipette to a 1-L volumetric flask and diluting to the mark with water. Mix thoroughly. Aliquot the solution (~1.5 mL) into pre-labeled 2-mL cryogenic vials. Approximately 200 vials will be needed for one year of analyses. Store vials at –70°C. These solutions are stable for at least one year.

(b) 4-PA – 100 ppb:

Prepare a 1:1000 dilution of the stock solution by transferring 1000 μL of the stock solution with a micropipette to a 1-L volumetric flask and diluting to the mark with water. Mix thoroughly. Aliquot the solution (~1.5 mL) into pre-labeled 2-mL cryogenic vials. Approximately 200 vials will be needed for one year of analyses. Store vials at –70°C. These solutions are stable for at least one year.

4) Formulas to calculate the concentration (nmol/L) of the working solutions

(a) PLP:

$$PLP \text{ working solution (nmol/L)} = \frac{(PLP \text{ stock conc.}(\mu\text{mol/L})) \times (\text{vol. of stock used(mL)})}{\text{final volume (L)}}$$

(b) 4-PA:

$$4\text{-PA working solution (nmol/L)} = \frac{(4\text{-PA stock conc.}(\mu\text{mol/L})) \times (\text{vol. of stock used(mL)})}{\text{final volume (L)}}$$

5) Calibration standards

When a run is performed thaw one vial of the PLP working solution and one vial of the 4-PA working solution. Mix the working solutions thoroughly prior to use. Using a micropipette, combine 550 μL of the PLP working standard and 450 μL of the 4-PA working standard in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer. This will be referred to as the **mixed standard** in steps (a) through (e) that follow.

(a) Level 1 calibration standard (~4 nmol/L PLP, 4-PA):

Using a micropipette, combine 10 μL of the mixed standard and 590 μL of water in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer.

(b) Level 2 calibration standard (~11 nmol/L PLP, 4-PA):

Using a micropipette, combine 30 μL of the mixed standard and 570 μL of water in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer.

(c) Level 3 calibration standard (~37 nmol/L PLP, 4-PA):

Using a micropipette, combine 100 μL of the mixed standard and 500 μL of water in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer.

(d) Level 4 calibration standard (~75 nmol/L PLP, 4-PA):

Using a micropipette, combine 200 μL of the mixed standard and 400 μL of water in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer.

(e) Level 5 calibration standard (~225 nmol/L PLP, 4-PA):

Using a micropipette, combine 600 μL of the mixed standard and 0 μL of water in a 1.5-mL microcentrifuge tube with a snap-top closure. Close the tube and mix its contents thoroughly using a vortex mixer.

(f) Calibration blank:

Using a micropipette, transfer 600 μL of water to a 1.5-mL microcentrifuge tube with a snap-top closure and close.

6) Formula to calculate the concentrations (nmol/L) of the calibration standards

$$PLP \text{ calibration std. (nmol/L)} = \frac{(\text{vol. mixed std. used } (\mu\text{L})) \times (PLP \text{ working conc. (nmol/L)}) \times 550}{600 \times 1000}$$

$$4\text{-PA calibration std. (nmol/L)} = \frac{(\text{vol. mixed std. used } (\mu\text{L})) \times (4\text{-PA working conc. (nmol/L)}) \times 450}{600 \times 1000}$$

The exact calculated concentrations of the calibration standards should be used in the instrument method (see sections 8.c.(5) and the Appendix).

c. Preparation of Quality Control Materials

The low QC pool is prepared by selecting and pooling serum that contains low levels of PLP (~15 nmol/L) and 4-PA (~10 nmol/L). The medium QC pool is prepared by selecting and pooling serum that contains moderate levels of PLP (~40 nmol/L) and 4-PA (~20 nmol/L). The high QC pool is prepared by selecting and pooling serum that contains high levels of PLP (~100 nmol/L) and 4-PA (~80 nmol/L). If no serum with elevated PLP and/or 4-PA levels is available, spiking the serum with known amounts of PLP or 4-PA may be a useful alternative.

All serum pools are filtered through gauze before being dispensed. Serum (400 μL) is aliquoted into pre-labeled 2.0-mL Nalgene cryovials, capped, and frozen. The QC pools are stored at -70 °C and are stable for at least 2 years. Means plus range limits for all pools are established by analyzing duplicates for at least 20 consecutive runs.

d. Other Materials

- (1) Acetonitrile, HPLC-grade (Burdick & Jackson, Muskegon, MI, USA).
- (2) Methanol, HPLC-grade (Burdick & Jackson, Muskegon, MI, USA).
- (3) Deionized water, 0.45 μm filtered lab source (Aqua Solutions, Jasper, GA, USA).
- (4) Pyridoxal 5'-phosphate (Sigma, St. Louis, MO, USA).
- (5) 4-pyridoxic acid (Sigma, St. Louis, MO, USA).
- (6) Dichloromethane (J.T. Baker, Phillipsburg, NJ, USA).
- (7) Metaphosphoric acid (HPO_3) (Mallinkrodt, Paris, KY, USA).
- (8) Orthophosphoric acid (H_3PO_4) (Fisher, Fair Lawn, NJ, USA).
- (9) Sodium phosphate, monobasic, monohydrate ($\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$) (Sigma, St. Louis, MO, USA).
- (10) Sodium chlorite (NaClO_2) (Aldrich, Milwaukee, WI, USA).
- (11) Serum samples (for QC pools) (Tennessee Blood Service, Memphis, TN, USA).
- (12) Hypersil BDS C₁₈ HPLC column, 150 \times 3-mm, 5- μm particle, 130-Å pore (Thermo Hypersil-Keystone, Bellefonte, PA, USA).

- (13) SecurityGuard HPLC guard cartridge holder, 4 × 2-mm (Phenomenex, Torrance, CA, USA).
- (14) SecurityGuard HPLC C₁₈ guard cartridges, 4 × 2-mm (Phenomenex, Torrance, CA, USA).
- (15) 1.5-mL polypropylene microcentrifuge tubes (VWR, Suwanee, GA, USA).
- (16) 2.0-mL polypropylene cryovials (Nalgene, Rochester, NY, USA).
- (17) 5.75" Disposable glass pasteur pipettes (Kimble, Toledo, OH).
- (18) 1-cc luer slip polypropylene syringes (Becton Dickinson, Franklin Lakes, NJ, USA).
- (19) 0.22- μ m × 4-mm polyvinyl difluoride syringe filters (Millipore, Billerica, MA, USA).
- (20) 0.22- μ m × x-mm polyvinyl difluoride vacuum solvent filters (Millipore, Billerica, MA, USA).
- (21) 2-mL, 12 mm × 32 mm amber HPLC vials (National Scientific, Duluth, GA, USA).
- (22) 11 mm PTFE/silicone HPLC vial caps (National Scientific, Duluth, GA, USA).
- (23) 200- μ L glass HPLC inserts (Kimble, Vineland, NJ, USA).
- (24) 100-1000 μ L air-displacement pipette tips (Eppendorf) (Brinkmann, Westbury, NY, USA).
- (25) 10-100 μ L air-displacement pipette tips (Eppendorf) (Brinkmann, Westbury, NY, USA).
- (26) 5-mL repeater pipette tips (Eppendorf Combi-Tip Plus) (Brinkmann, Westbury, NY, USA).
- (27) 10–100 μ L positive-displacement pipette tips (Gilson Pipetteman) (Rainin, Woburn, MA, USA).
- (28) 50–250 μ L positive-displacement pipette tips (Gilson Pipetteman) (Rainin, Woburn, MA, USA).
- (29) Nitrile laboratory gloves (Best Manufacturing, Menlo, GA, USA).
- (30) Various glass beakers, volumetric flasks, graduated cylinders, and bottles, class A glassware.

e. Instrumentation

- (1) Agilent Technologies HPLC System (Agilent, Palo Alto, CA, USA).
 - (a) G1311A quaternary pump (×2)
 - (b) G1329/30A thermostatted autosampler
 - (c) G1316A thermostatted column compartment (×2)
 - (d) G1321 fluorescence detector with 8- μ L flow cell
 - (e) ChemStation workstation and software (version A.09.01 or higher)
- (2) Microfiltration assembly, 47 mm, with 1 L solvent bottle and fritted glass support (Kimble-Kontes, Vineland, NJ, USA).

- (3) Vortex mixer, single sample (Daigger, Vernon Hills, IL).
- (4) Vortex mixer, multiple samples (Glas-Col, Terre Haute, IN).
- (5) Magnetic stir plate (Fisher, Fair Lawn, NJ, USA).
- (6) Refrigerated microcentrifuge (Labnet, Edison, NJ, USA).
- (7) Digital pH/temperature meter (Beckman, Fullerton, CA, USA).
- (8) Sonicating bath (VWR, Suwanee, GA, USA).
- (9) Digital balance, 0.001 g accuracy (Mettler Toledo, Columbus, OH, USA).
- (10) Digital balance, 1 mg accuracy (Mettler Toledo, Columbus, OH, USA).
- (11) 100–1000 μ L adjustable air-displacement micropipette (Eppendorf) (Brinkmann Instruments Co., Westbury, NY)
- (12) 10–100 μ L adjustable air-displacement micropipette (Eppendorf) (Brinkmann Instruments Co., Westbury, NY)
- (13) 10–100 μ L adjustable positive-displacement micropipette (Gilson Pipetteman) (Rainin, Woburn, MA, USA).
- (14) 50 –250 μ L adjustable positive-displacement micropipette (Gilson Pipetteman) (Rainin, Woburn, MA, USA).
- (15) Digiflex diluter (Micromedic Division, ICN Biomedical, Costa Mesa, CA).

7. Calibration and Calibration Verification Procedures

Spiked recovery studies for this method [15] showed mean recoveries (\pm SD) of $97 \pm 4\%$, $100 \pm 3\%$ and $101 \pm 3\%$ for serum samples with PLP levels of 12 nmol/L, 20 nmol/L and 33 nmol/L, respectively. For 4-PA, mean recoveries (\pm SD) were $97 \pm 3\%$ and $101 \pm 3\%$ at concentrations of 11 nmol/L and 21 nmol/L, respectively. In all cases, the amount of analyte spiked was approximately equal to the endogenous levels of the serum.

At the beginning of each run, five calibrators and a blank are prepared as described in section 6(b)5(a)–(f) of this document. Sample preparation for the calibrators and blank is identical to the sample preparation used for patient samples and QCs.

Calculation of calibration curves for PLP and 4-PA, and interpolation of standards, QCs and patient samples against these equations is performed in ChemStation software. After determining the peak areas for PLP and 4-PA in the standards, a five-point calibration curve is generated by performing a non-weighted linear regression (not forced through zero). At the end of each run, the two calibrators are reanalyzed as unknowns. The measured concentrations of these calibrators must agree within 15% of their set values.

NIST reference materials are not available for vitamin B₆ assays. External quality assessment programs do not exist for vitamin B₆. An inter-laboratory round-robin study involving 10 internationally-recognized laboratories [16] was organized by our laboratory in 2004. Results from our laboratory were found to satisfy optimal performance requirements in terms of method bias and imprecision when evaluated against

empirically-derived performance criteria based on biological variation for vitamin B₆ [17]. Proficiency testing is conducted in-house at least twice a year.

8. Procedure Operating Instructions; Calculations; Interpretation of Results

a. Preliminaries

- (1) Allow frozen patient serum specimens, QC serum specimens, and PLP and 4-PA working solution aliquots to thaw and reach ambient temperature.
- (2) If needed, prepare the 5% MPA solution as indicated in section 6(a)(1) and refrigerate. If already available, keep the 5% MPA solution refrigerated until ready to use.
- (3) Turn on the refrigerated microcentrifuge and set the thermostat to 4 °C.
- (4) Turn on all components of the HPLC system including the computer. Start an online session of the ChemStation software (double-click the “online” ChemStation icon).
 - (a) Click on “View” → “Method and Run Control” to open the run control instrument view.
 - (b) Click on “Instrument” → “More Injector...” → “Sample Thermostat. Turn the autosampler thermostat on by clicking “on” and set its temperature to 4 °C.

b. Sample preparation

A typical run consists of one reagent blank, 5 calibrators, 3 QC samples, 60 patient samples, and 3 QC samples.

Use positive displacement pipettes wherever possible for the preparation of the calibration standards. Use positive displacement pipettes or an automated positive displacement pipettor for transferring the preliminary assay volume (100 µL) of standards, QCs and patient samples into the microcentrifuge tubes containing 100 µL of 5% MPA (steps (3) and (4) below). There are no volume-dependent steps after this point.

NOTE: Vitamin B₆ is light sensitive! Perform all sample preparation under yellow lighting and in the absence of any direct sunlight. Use amber colored glassware and plasticware whenever possible.

- (1) Prepare the five calibrators and blank as described in section 6(b)5(a)–(f).
- (2) Label two sets of 1.5-mL microcentrifuge tubes, one tube per sample in each set, and place them in racks.
- (3) Using a repeater pipette (positive displacement recommended), transfer 100 µL of the cold 5% MPA solution to each of the microcentrifuge tubes in the first set.
- (4) Using a manual pipette or automated pipettor (positive displacement recommended), transfer 100 µL of the blank, standards, QCs and patient samples into their respective microcentrifuge tubes in the first set.
- (5) Tightly cap all of the microcentrifuge tubes in the first set and mix the tubes for 5 min at 2000 rpm or greater with a multitube vortexer.
- (6) Transfer the first set of tubes to a refrigerated microcentrifuge set to 4 °C and centrifuge at 18,000 g for 15 min.
- (7) While the first set of tubes are in the centrifuge, transfer 150 µL of dichloromethane (CH₂Cl₂) to each of the microcentrifuge tubes in the second set.

CAUTION!!! Dichloromethane (CH₂Cl₂) is highly toxic and potentially carcinogenic. It is harmful or fatal if swallowed, inhaled, or absorbed through the skin. Dichloromethane affects the central nervous and cardiovascular systems, liver and blood. It is an irritant to the skin, eyes, and respiratory tract. Because risk of cancer depends on level and duration of exposure, appropriate precautions to minimize exposure must be taken. ALWAYS WORK WITH DICHLOROMETHANE IN A CHEMICAL FUMEHOOD. WEAR APPROPRIATE PPE. Consult MSDS before use

- (8) Once the centrifuge has stopped, carefully remove the samples, trying not to disturb the precipitate that has accumulated at the bottom of the tubes, and place them in racks.
- (9) For each tube in the first set:
 - (a) Draw off the supernatant (the clear liquid above the precipitate) using a micropipette or glass Pasteur pipette.
 - (b) Transfer the liquid to its corresponding tube in the second set.
- (10) Once finished, dispose of the first set of tubes.
- (11) Tightly cap all of the microcentrifuge tubes in the second set and mix the tubes for 2 min at 2000 rpm or greater with a multitube vortexer.
- (12) Transfer the second set of tubes to a refrigerated microcentrifuge set to 4 °C and centrifuge at 18,000 g for 15 min.

- (13) Once the centrifuge has stopped, carefully remove the samples, trying not to disturb the separation of the two liquid layers in the tubes, and place them in racks.
- (14) For each tube in the second set,
 - (a) Draw off the top liquid layer (aqueous layer) with a 1-cc luer slip polypropylene syringe.
 - (b) Firmly place a syringe filter on the end of the syringe.
 - (c) *SLOWLY* push the liquid through the filter into an amber HPLC vial with a 200- μ L reduced volume insert.
- (15) Cap the HPLC vials and gently tap them to ensure that there is no air trapped at the bottom of the insert.
- (16) If the HPLC is ready, transfer the HPLC vials to their appropriate positions in the autosampler tray; otherwise, refrigerate the samples until the HPLC is ready for the analysis.

c. Instrument Preparation

(1) Solvent Reservoir Filling

The HPLC system has two quaternary pumps, identified as Pump 1 and Pump 2 in the ChemStation operating software. Pump 1 controls solvent delivery for the actual HPLC separation, while Pump 2 is responsible for delivering the postcolumn derivatization reagent. Each pump has four 1-L solvent reservoirs designated A through D. The solvent reservoirs should be filled as follows:

For Pump 1:

- A: 50 mM phosphate buffer, pH 3.2, 0.2% acetonitrile
- B: Acetonitrile
- C: Methanol
- D: Water

For Pump 2:

- A: Water
- B: 2 g/L sodium chlorite solution
- C: Methanol
- D: not used

Make sure that all reservoir bottles and solvent lines are clean and free of foreign matter (particulates, algae, etc.). If dirty, clean or replace these components as necessary.

Fill all reservoirs with the appropriate solvent or reagent. For the 50 mM phosphate buffer (Pump 1, Reservoir A), use an amber colored reservoir. Once the reservoirs have been filled, replace the solvent lines, making sure that they are in the correct reservoirs and reach the bottom of each bottle. Then click on "Instrument" → "More Pump (1 or 2)..." → "Bottle Filling" and enter the amount of solvent in each reservoir. This will enable the HPLC system to track how much solvent is left in each reservoir

and take preventative measures (*i.e.*, shutdown the system) before running out of solvent.

(2) Solvent Priming

If any of the solvent lines on Pump 1 or Pump 2 need to be primed (*i.e.*, if there are signs of air present in a solvent line, or the instrument has not been used for several days), the following procedure should be performed for each line as necessary.

- (a) Click on “Instrument” → “Set Up Pump (1 or 2)”.
- (b) Set the flow rate to 4.0 mL/min.
- (c) Set the reservoir (A, B, C, or D) to 100%.
- (d) Click “OK”.
- (e) Open the purge valve (the large black knob on the front of the pump) by turning it at least one full turn counter-clockwise.
- (f) Click on “Instrument” → “More Pump (1 or 2)” → “Control...”
- (g) Turn pump “on”.
- (h) Click “OK”.
- (i) Allow the pump to run for at least 5 min. After 5 min, check the solvent line for any signs of air bubbles. Allow the pump to run until no more air bubbles appear.
- (j) Click on “Instrument” → “More Pump (1 or 2)” → “Control...”
- (k) Turn pump “off”.
- (l) Click “OK”.
- (m) Close the purge valve.

After priming all the solvent lines (or as needed), prime Pump 1 with 95% B (acetonitrile), 5% D (water) for 30 s. Do the same for Pump 2 with 95% C (methanol), 5% B (water).

NOTE: Prime each solvent line separately. When priming, avoid switching directly from solvent lines containing high amounts of dissolved salts (buffer, postcolumn reagent) to lines containing 100% organic solvents (acetonitrile, methanol) and vice versa. Doing so may cause precipitates to form in the solvent lines and the pump. To prevent this, briefly purge (30 s) a solvent line containing water between purging lines containing dissolved salts and organic solvents.

(3) Guard Cartridge Replacement

The HPLC system uses a guard column to protect the integrity of the analytical column by blocking out particulates and other potential contaminants. The guard

column uses a holder that connects directly to the inlet of the analytical column and disposable drop-in cartridges. The guard cartridge should be replaced before each new sample run. Always wear clean, powder-free gloves when working with the guard cartridge holder. Replace the guard cartridge as follows:

- (a) Disconnect the HPLC inlet line to the guard cartridge holder. Then, disconnect the cartridge from the analytical column.
- (b) Place the guard cartridge holder on a bench with the outlet tip pointing upward.
- (c) Using the wrenches provided with the holder kit, disconnect the top (female) portion of the guard cartridge holder from the bottom (male) part.
- (d) Remove the guard cartridge from the male part of the holder.
- (e) Insert a new replacement cartridge into the male part of the holder. Make sure that the gold portion of the cartridge is facing down, and that the silver portion of the cartridge is facing up.
- (f) Reconnect the female portion of the guard cartridge holder to the male portion. Tighten the parts together firmly with the wrenches provided.
- (g) Reconnect the HPLC inlet line to the guard cartridge holder, but do not connect it to the analytical column.
- (h) Run 95% acetonitrile / 5% water through the guard cartridge holder at 2 mL/min for approximately 30 s. To do this, follow the steps in section c (2) above for Pump 1, with Reservoir B set to 95%, Reservoir D set to 5%, and a flow rate of 2 mL/min. Do not open the purge valve. Collect the flow coming out of the guard cartridge holder in a small beaker and dispose of this waste appropriately.
- (i) Reconnect the guard cartridge holder to the analytical column. Do not use wrenches to do this!

(4) Analytical Column Replacement

The analytical column should be replaced when its chromatographic performance has been compromised. Problem areas include, but are not limited to, sudden changes and/or poor reproducibility in the following:

- Peak shapes
- Retention times
- Resolution
- Column backpressure

Analytical column replacement is at the sole discretion of the operator, but it is recommended that the column is replaced after 2000 sample injections, even if it is still apparently in good working order. To replace the analytical column, disconnect the guard cartridge holder and the HPLC outlet line from the old column and reconnect a new column in its place. New analytical columns are shipped in 95% acetonitrile / 5% water, so there is no need to purge the column as is done when the guard cartridge is replaced.

(5) Setting Up and Running a Sequence

Once the HPLC system has been primed and the guard and analytical columns have been replaced (if needed), a sequence is set up for running the samples. A sequence can be used to perform all steps of the analysis including system equilibration and cleanup steps. Information is entered in the form of a sequence table – a line-by-line set of analysis steps to be entered by the analyst. Rows in the sequence table are numbered 1 through n and indicate the order in which the analysis steps will be executed. Columns in the sequence table represent fields for variables that can be set for each analysis step. The fields that will be used are:

Location:

This field indicates the position of the sample vial in the autosampler tray. Valid field entries for vial positions are 1–99. Position 91 is reserved for a wash vial. If the field is left blank, that line of the sequence table will be executed without injecting a sample.

Sample Name:

This field indicates the name of the sample being injected. Standard, QC and patient sample IDs are to be entered in this field.

Method Name:

This field reveals a drop-down menu in each line from which a method used to acquire/analyze the sample for a given sequence line can be selected. The methods that will be used are:

B6-002: The analysis method for running standards, QCs and patient samples.

STARTUP: Performs all startup and equilibration steps to prepare the HPLC system for analysis.

SHUTDOWN: Performs all cleanup steps following an analysis.

STANDBY: Prepares the HPLC system to be shut down.

The instrumental settings for each of these methods appear in the Appendix.

Inj/Location:

Indicates the number of times the sample is to be injected. Default is 1 unless instructed otherwise.

Sample Type:

This field reveals a drop-down menu in which the sample type can be designated. Calibration standards run at the beginning of a sequence are designated as “Calibration”, all other samples (including when the calibration standards are re-injected at the end of the run) are designated as “Sample”.

Cal Level:

This drop-down menu field designates a point on a calibration curve. It is active on the sequence line if “Calibration” is selected in the Sample Type field. The level 1 calibration standard is designated as “1”; level 2 is designated “2”, and so on.

Update RF, Update RT:

This drop-down menu field determines if the recalibration should update the calibration response factors and retention times. Both the Update RF and Update RT fields should be set to "Replace" when the Sample Type field is "Calibration".

Multiplier:

This is a scaling factor applied to the calculated results before they are reported. All standards, QCs and samples are given a value of "1". If a sample is further diluted (e.g., if it is being reanalyzed because it exceeded the calibration range) the multiplier must reflect the degree to which the sample has been diluted (e.g., a 1:4 (5x) dilution of serum with water must have a multiplier of 5).

A run sequence is set up as follows:

- (a) Click on "Sequence" → "New Sequence". Alternatively, you can use an existing sequence as a template and edit it by clicking on "Sequence" → "Load Sequence..." and selecting the appropriate sequence file.
- (b) Click on "Sequence" → "Save Sequence as..." Name the sequence using the date in YYMMDD format.
- (c) Click on "Sequence" → "Sequence Parameters..." to open the Sequence Parameters dialog.
 - (i) Under the "Data File" heading, click on "Prefix/Counter".
 - (ii) In the "Prefix" field, enter the date using the format YYMMDD.
 - (iii) In the "Counter" field, enter "00".
 - (iv) In the "Subdirectory" field, enter the name of the appropriate subdirectory for the data files to be collected (e.g., NH04, etc.).
 - (v) Under the "Shutdown" heading, check "Post Sequence Cd/Macro".
 - (vi) Select "Standby" from the drop-down menu.
 - (vii) Click on "OK".
- (d) Click on "Sequence" → "Sequence Table". This will display the sequence table.
- (e) On sequence line 1, run the startup method:
 - (i) Location: Leave blank.
 - (ii) Sample Name: "Startup"
 - (iii) Method Name: "STARTUP"
 - (iv) Inj/Location: "1"
 - (v) Sample Type: "Sample"

- (f) On sequence lines 2, 3, and 4, run the blank, the level 3 calibration standard, and the medium QC, respectively. These three injections are system checks to verify that the instrument is functioning properly:
- (i) Location: Specify the positions of each sample in the autosampler tray.
 - (ii) Sample Name: Name these samples “Blank Test”, “Cal Test”, and “QC Test”, respectively.
 - (iii) Method Name: “B6-002”
 - (iv) Inj/Location: “1”
 - (v) Sample Type: “Sample”
 - (vi) Multiplier: “1”
- (g) On sequence lines 5–10, run the blank and the five calibration standards. These injections are used to generate the calibration curves for the analysis:
- (i) Location: Specify the positions of the blank and the calibration standards in the autosampler tray.
 - (ii) Sample Name: “blank” for the blank, “Standard 1” for the level 1 calibration standard, “Standard 2” for the level 2 calibration standard, etc. (*Note*: different names can be used for the blank and standards, as long as they match their respective names in the FrontEnds database).
 - (iii) Method Name: “B6-002”
 - (iv) Inj/Location: “1”
 - (v) Sample Type: “Sample” for the blank; “Calibration” for all 5 calibration standards.
 - (vi) Cal Level: “1” for the level 1 calibration standard, “2” for the level 2 calibration standard, etc.
 - (vii) Update RF: “Replace” for all calibration standards.
 - (viii) Update RT: “None” for all calibration standards.
 - (ix) Multiplier: “1”.
- (h) On the subsequent sequence lines, run the following samples in this order:
- First set of QCs (low, medium, high)
 - Patient samples (recommended maximum of 60)
 - Second set of QCs (low, medium, high)
 - Blank and calibration standards
- These samples will be quantified according to the calibration curve generated from the calibration standards:
- (i) Location: Specify the positions of the first set of QCs, the patient samples, the second set of QCs, and the blank and calibration standards

- (ii) Sample Name: For the QCs and patient samples, enter the sample ID exactly as it appears on the specimen vial. If the vials have barcodes, use the barcode reader to avoid data entry errors. For the blank and standards at the end of the sequence, use the same names as before: “blank” for the blank, “Standard 1” for the level 1 calibration standard, “Standard 2” for the level 2 calibration standard, etc. (*Note*: different names can be used for the blank and standards, as long as they match their respective names in the FrontEnds database).
 - (iii) Method Name: “B6-002”.
 - (iv) Inj/Location: “1”.
 - (v) Sample Type: “Sample” (including the calibration standards at the end – do not use “Calibration”).
 - (vi) Multiplier: “1” (or otherwise if sample is further diluted).
- (i) On the two sequence lines following, run the SHUTDOWN and STANDBY methods:
 - (i) Location: Leave blank.
 - (ii) Sample Name: “shutdown” and “standby”, respectively.
 - (iii) Method Name: “SHUTDOWN” and “STANDBY”, respectively.
 - (iv) Inj/Location: “1”.
 - (v) Sample Type: “Sample”.
 - (j) Click on “OK”.
 - (k) Click on “Sequence” → “Save Sequence as...” Name the sequence using the date in YYMMDD format.
 - (l) Click on “Method” → “Load Method...” and load the method “STARTUP”.
 - (m) Click on “Instrument” → “More Pump (1 or 2)...” → “Bottle Filling”. Verify that the solvent reservoirs contain sufficient amounts to perform the analysis and that their entered values are correct.
 - (n) Click on “Instrument” → “System On” to turn on all pumps, thermostats and the detector.
 - (o) Click on “Run Control” → “Run Sequence”. The sequence will now start running.

d. Data Review

Once the sequence has finished, the peak integration has to be reviewed and calibration equations updated before it can be exported for QC verification. This is done by reviewing the batch file generated when the sequence is run.

- (1) Batch Review: Preliminary Settings

- (a) Start an offline session of the ChemStation software (double-click the “offline” ChemStation icon).
 - (b) Click on “View” → “Data Analysis” to open the data analysis view.
 - (c) Click on “Batch” → “Load Batch...” Look for a batch file (.B file extension) with the same name as the sequence containing the samples you wish to process and load the file.
 - (d) A dialog window will appear listing the files (by sample name) that were run in the sequence. Check the samples to be reviewed. This will typically be all files except those created at startup and shutdown (“Startup”, “Cleanup”, “Standby”) and files created from the system checks (“Blank Test”, “Cal Test”, “QC Test”).
 - (e) Click on “OK”.
 - (f) Click on “Batch” → “Options...” Go to the “Processing” tab.
 - (i) Under the heading “During processing of a run the following steps will be executed:”, activate the “Integration” and “Identification/Quantitation” options. Deactivate the “Print individual report...” option for now. Reports will be printed in a later step.
 - (ii) Set the pause between processing of runs to 0 seconds.
 - (iii) Under the “Signal display” heading, activate both the “Freeze X-axis” and “Freeze y-axis” options. This allows a specified zoom level to be used when reviewing the integration.
 - (g) Go to the “Report Options” tab.
 - (i) Under the “Destination” heading, activate “screen” and “file” options.
 - (ii) Set the file prefix to the date from the sequence and batch (use YYMMDD format), and set the file type to “.xls”.
 - (h) Click on “OK”.
- (2) Integration:
- (a) Click on the “Select Integration Task” icon in the top left corner of the screen.
 - (i) Go to the run table at the bottom of the screen and highlight “Standard 05” (the highest calibrator) by clicking on it.
 - (ii) Click on “Integration” → “Integration Events...”
 - (iii) Verify that the settings in the integration events table match those that appear in Appendix for B6-002.
 - (iv) Examine the chromatogram for the highest calibrator and verify that the peaks for PLP (retention time ~ 3.2 min) and 4-PA (retention time ~ 5.6 min) are being integrated properly (*i.e.*, verify that the entire peak is being integrated and that no extraneous signals are being integrated in the process). If necessary, adjust the on/off times in the events table, but do

not set the integration excessively wide as this may result in additional peaks being integrated in the samples.

- (v) Close the integration events table by clicking on the “Exit and save events to method” icon in the top left corner of the table.
 - (b) Return to the run table at the bottom of the screen and highlight the first sample. Visually examine the integration. If you need to examine an area more closely, click on the “Zoom in on window” icon above the chromatogram and click and drag over the region of interest. Double-clicking the chromatogram will zoom out.
 - (i) If any peaks other than PLP and 4-PA have been integrated, these can be removed by clicking on the “Remove integrated peaks...” icon above the chromatogram and clicking on the peak to be removed. Note: leaving integrated peaks other than PLP and 4-PA in the chromatogram will have no deleterious effects on the sample processing. Removing them is done simply to clean up chromatogram appearance.
 - (ii) If the peaks for PLP or 4-PA does not appear to be integrated properly, you can remove the integration baseline as described in step (ii) above, and then draw in the baseline manually by clicking on the “Draw baseline and integrate” icon above the chromatogram. Do this only as a last resort.
 - (iii) If it is acceptable, click on the “Mark changes for current run ... for saving...” icon just above the run table. The line in the run table will now be marked with an “s” indicating that the integration has been saved.
 - (c) Click on the green “down” arrow (↓) to advance to the next sample and repeat steps (i) through (iii) above until all samples have been reviewed and saved.
- (3) Calibration:
- (a) Click on the “Select Calibration Task” icon in the top left corner of the screen.
 - (b) Verify that the retention times and concentrations for the PLP and 4-PA standards in the calibration table are correct.
 - (c) Click on “Batch” → “Update Calibration”. This will update the calibration curves and equations based on the integration performed in section (2) above.
 - (d) Check the calibration curves and equations for PLP and 4-PA. The calibration is acceptable if $r \geq 0.999$.
- (4) Report Printing and Export:
- (a) Click on “File” → “Printer Setup...” and set the paper orientation to “Landscape”.
 - (b) “Click on “Batch” → “Start” to start an automated review of the samples using the updated calibration information from section (3) above. This will take a few minutes to complete.

- (c) Click on “Batch” → “Output Batch Report”. A summary of the sample results will appear on the screen and a Microsoft Excel file of the results will be generated. Print a copy of the batch report displayed on the screen (use the “Print Batch Report” icon) and close the display.
- (d) Click on “File” → “Printer Setup...” and set the paper orientation to “Portrait”.
- (e) Click on “File” → “Print” → “Calibration Tables and Curves” to print the calibration data.
- (f) Click on “Batch” → “Options...”. Go to the “Processing” tab. Under the heading “During processing of a run the following steps will be executed:” activate the “Print individual report...” option. Click on “OK”.
- (g) Click on “Batch” → “Start”. This will start an automated review of the samples again, only this time a summary report with a chromatogram will be printed for each sample.
- (h) Assemble the printouts from the run into a single document in following order:
 - Batch report
 - Calibration tables and curves
 - Sample reports and chromatograms
- (i) Retrieve the Excel file of the batch report from the directory specified in section 6.c. (5) (c) (iv).
- (j) Proceed with the QA/QC review process as described in section 3.a.

e. System Maintenance

The system maintenance consists of following components:

- (1) Per run:
 - (a) Inspection of the solvent reservoirs, inlet lines and filters. Replenishment or replacement of solvents as needed. Cleaning or replacement of components as needed.
 - (b) Replacement of the guard column cartridge.
 - (c) Replacement of the wash vial (position 91 of autosampler tray).
 - (d) Inspection and replacement of column fittings, tubing (as needed).
- (2) Per annum:
 - (a) Preventative maintenance (performed by a service engineer):
 - (i) Pump pressure and leak tests
 - (ii) Replacement of pump seals, inlet filters, inlet valves (if needed).
 - (iii) Fluorescence detector dark current and wavelength calibration test.
 - (iv) Replacement of fluorescence bulb (if needed).
- (3) Other:

(a) Analytical column replacement (after 2000 injections or as needed depending on chromatographic performance).

i. Special Method Notes

If the HPLC is to remain idle for a week or longer, purge the solvent lines containing buffers and dissolved salts with water followed by methanol to prevent salt precipitation or algae growth from occurring in the HPLC components.

j. CDC Modifications

This method is based on the method described by Rybak et al. [15].

9. Reportable Range of Results

This method is linear for PLP and 4-PA in the range of 4–250 nmol/L. Samples with PLP or 4-PA concentrations greater than the concentration of the highest calibrator (~220 nmol/L) are diluted such that the appropriate analyte falls within its linear range. This method has a total coefficient of variation in the range of 3.6–6.7% for PLP and 3.7–5.6% for 4-PA [15].

10. Quality Control (QC) Procedures

a. Blind Quality Controls

Blind QC specimens are inserted prior to the arrival of the samples in the Inorganic Toxicology and Nutrition Branch. These specimens are prepared at two levels so as to emulate the patient samples; the labels used are identical to those used for patient samples. One blind QC specimen randomly selected for concentration is included at a randomly selected location in every 20 specimens analyzed.

b. Bench Quality Controls

Bench QC specimens are prepared from three plasma pools, which represent low, intermediate, and high levels of PLP and 4-PA in serum, based on 10th, 50th and 90th percentiles observed for these analytes in a nonrepresentative population subset [15]. These pools are prepared in the same manner as patient samples and analyzed in duplicate as part of each run.

The results from the pools are checked after each run. The system is declared “in control” if all three QC results are within 2s limits and the run is accepted. If one of the three QC results is outside the 2s limits then apply rules below and reject if any condition is met – the run is then declared “out of control”:

1_{3s} Any of the three QC results are outside the 3s limit

2_{2s} Two of the three QC results in the run are outside the 2s limit (same side of mean)

R_{4s} Sequential QC results (either within the run or across runs) are outside the 2s limit on the opposite sides of the mean

10_x Ten sequential QC results (across pools and across runs) are on the same side of the mean

A QC program written in SAS is available from the DLS Quality Assurance Officer and should be used to apply these rules to QC data and generate Shewhart QC charts. No results for a given analyte are to be reported from an analytical run that has been declared "out of control" for that analyte as assessed by internal (bench) QC.

The initial limits are established by analyzing pool material in 20 consecutive runs and then are reevaluated quarterly. When necessary, limits are updated to include more runs.

While a study is in progress, QC results are stored in the ACCESS database. For runs that are not imported into ACCESS (exception, research-type runs), QC results are stored electronically in the analyte-specific folder on Q:\ITN\Nutrition Lab\Data handling\QC results in Excel. A hardcopy of the QC results from each run is also kept by the analyst.

11. Remedial Action if Calibration or QC Systems Fail to Meet Acceptable Criteria

Check to make sure that the hardware is functioning properly, and verify that all instrumental settings are correct. Check the autosampler to make sure the injections are being made as programmed.

Look for sample preparation errors.

Check the calibrations of the pipettes.

If the steps outlined above do not result in correction of the "out of control" values for QC materials, consult the supervisor for other appropriate corrective actions. Do not report analytical results for runs not in statistical control.

12. Limitations of Method; Interfering Substances and Conditions

If serum is not available, heparinized plasma may be used as a substitute [15]. EDTA plasma should be avoided due to an interference with the signal for PLP.

Increases in the following activities may increase (↑) or decrease (↓) serum PLP and 4-PA levels as indicated [13]:

Use of vitamin B₆ supplements: ↑ PLP, ↑ 4-PA

Protein intake: ↓ PLP

Glucose: ↓ PLP

Plasma volume: ↑ PLP

Physical activity: ↑ PLP

13. Reference Ranges (Normal Values)

Based on literature data [15], means, ranges and percentiles for this method are shown in the Appendix.

14. Critical Call Results (“Panic Values”)

Though PLP values of 30 nmol/L [13] and 20 nmol/L [18] have been proposed as cutoff values for normal status, there is no consensus on a critical call value for PLP in serum. There is no critical call value for 4-PA in serum.

15. Specimen Storage and Handling During Testing

Specimens are allowed to reach room temperature during preparation. Once the samples are ready to run, the prepared samples are placed into the autosampler tray. The unused portion of the patient specimen is returned to the freezer.

16. Alternate Methods for Performing Test of Storing Specimens if Test System Fails

Because the analysis of serum vitamin B₆ is inherently complex and challenging, there are no acceptable alternative methods of analysis in the NHANES laboratory. If the analytical system fails, we recommend that the extracted specimens be stored at ≤ -20 °C until the analytical system is restored to functionality.

17. Test Result Reporting System; Protocol for Reporting Critical Calls (If Applicable)

Test results are reported to the collaborating agency at a frequency and by a method determined by the study coordinator. Generally, data from this analysis are compiled with results from other analyses and sent to the responsible person at the collaborating agency as an ASCII text file or Excel file, either through electronic mail or on a diskette.

For NHANES 1999+, all data are reported electronically periodically to the Westat ISIS computer and then are transferred to NCHS. For some smaller studies, hard copies of a data report are sent, as well as the results in electronic format.

18. Transfer or Referral of Specimens; Procedures for Specimen Accountability and Tracking

The Microsoft Access database is used to keep records and track specimens for NHANES 1999+. If serum vitamin B₆ analyses are used for smaller, non-NHANES studies, records are kept on files in Q:\ITN\Nutrition Lab on the DLS LAN.

We recommend that records, including related QA/QC data, be maintained for 10 years after completion of the NHANES study. Only numerical identifiers should be used (e.g., case ID numbers). All personal identifiers should be available only to the medical supervisor or project coordinator. Residual serum from these analyses for non-NHANES studies may be discarded at the request of the principal investigator, or may be transferred to the CDC CASPIR facility for use by other investigators. Very little residual material will be available after NHANES analyses are completed, and these vials may be routinely autoclaved.

The exact procedure used to track specimens varies with each study and is specified in the study protocol or the interagency agreement for the study. Copies of these documents are kept by the supervisor. In general, when specimens are received, the specimen ID number is entered into a database and the specimens stored in a freezer at -70°C. The specimen ID is read off of the vial by a barcode reader attached to the


computer used to prepare the electronic specimen table for the analytical system. When the analyses are completed, the DIF file containing the electronic copy of the results is loaded into the database, and the analytical results are linked to the database by ID number. The analyst is responsible for keeping a notebook containing the ID numbers of specimens prepared incorrectly, those with labeling problems, and those with abnormal results, together with information about these discrepancies.

19. Summary Statistics and QC Graphs

a. 4-Pyridoxic Acid

Summary Statistics for 4-Pyridoxic Acid by Lot


Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
ALS03420	9	3/2/2005	4/13/2005	10.494	0.477	4.5
AMS03421	9	3/2/2005	4/13/2005	21.044	0.759	3.6
AHS03422	9	3/2/2005	4/13/2005	21.083	0.530	2.5
BLS03420_V02	19	4/19/2005	7/14/2005	11.021	0.764	6.9
BHS03422_V02	19	4/19/2005	7/14/2005	21.295	0.944	4.4
BMS03421_V02	19	4/19/2005	7/14/2005	21.579	0.744	3.4
CLS03420_V03	142	7/18/2005	10/22/2007	11.089	0.748	6.7
CHS03422_V03	142	7/18/2005	10/22/2007	20.731	0.979	4.7
CMS03421_V03	142	7/18/2005	10/22/2007	20.951	0.986	4.7
3LS07420_V03	13	2/20/2008	9/12/2008	10.446	1.013	9.7
3MS07421_V03	13	2/20/2008	9/12/2008	22.081	0.956	4.3
3HS07422_V03	13	2/20/2008	9/12/2008	50.623	1.233	2.4



b. Pyridoxal 5'-Phosphate

Summary Statistics for Pyridoxal 5'-Phosphate by Lot

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
ALS03420	8	3/2/2005	4/18/2005	11.706	0.629	5.4
AMS03421	9	3/2/2005	6/13/2005	21.061	1.032	4.9
AHS03422	8	3/2/2005	4/18/2005	35.744	1.281	3.6
BLS03420_V02	22	4/19/2005	7/14/2005	11.884	0.780	6.6
BMS03421_V02	22	4/19/2005	7/14/2005	21.409	0.985	4.6
BHS03422_V02	22	4/19/2005	7/14/2005	36.014	1.441	4.0
CLS03420_V03	151	7/18/2005	10/29/2007	11.264	0.753	6.7
CMS03421_V03	151	7/18/2005	10/29/2007	20.327	0.999	4.9
CHS03422_V03	151	7/18/2005	10/29/2007	34.426	1.456	4.2
3LS07420_V03	10	2/15/2008	9/12/2008	9.945	0.829	8.3
3MS07421_V03	10	2/15/2008	9/12/2008	37.195	1.061	2.9
3HS07422_V03	10	2/15/2008	9/12/2008	63.660	1.581	2.5



References

ACKNOWLEDGMENTS

We gratefully acknowledge the contributions of:

- Michael E. Rybak, Ph.D. for developing and validating this method and preparing the manuscript for this chapter
- Donna LaVoie, M.T. (A.S.C.P.) and Janette Collins-Mitchell, M.S., R.N., M.T. (N.C.A.) for their assistance in the validation and day-to-day performance of this method
- Christine M. Pfeiffer, Ph.D. for the assistance and advice provided

Appendix (Tables)

Table 1. Method parameters: B6–002

<i>Injector setup –</i>				
Injection	Injection with needle wash			
Injection volume	10 µL			
Wash vial	Position 91			
Optimization	Prefetch sample vial 12.50 min after injection			
Stop time	As pump (13.00 min)			
Post time	Off			
Draw speed	200 µL/min			
Eject speed	200 µL/min			
Draw position	2.0 mm			
<i>Pump 1 setup –</i>				
Flow	0.7 mL/min			
Stop time	13.00 min			
Post time	Off			
Solvents:				
A: 50 mm phosphate buffer, pH 3.2, 0.2% acetonitrile	100%			
B: acetonitrile	Off			
C: methanol	0%			
D: water	Off			
Pressure limits:				
Max	400 bar			
Min	0 bar			
Timetable:				
3.50 min	0.0% B	0.0% C	0.0% D	0.700 mL/min
4.00 min	0.0% B	30.0% C	0.0% D	0.700 mL/min
8.50 min	0.0% B	30.0% C	0.0% D	0.700 mL/min
9.00 min	0.0% B	0.0% C	0.0% D	0.700 mL/min
11.00 min	0.0% B	0.0% C	0.0% D	0.700 mL/min
<i>Pump 2 setup –</i>				
Flow	0.1 mL/min			
Stop time	13.00 min			
Post time	Off			
Solvents:				
A: water	Off			
B: 2g/L sodium chlorite	100%			
C: methanol	Off			
D: not used	Off			
Pressure limits:				
Max	20 bar			
Min	0 bar			
<i>Column thermostat 1 setup –</i>				
Temperature	35.0 °C			
Stop time	No limit			
Post time	No limit			
Enable analysis	When temperature is within setpoint +/- 0.8 °C			
<i>Column thermostat 2 setup –</i>				
Temperature	75 °C			
Stop time	No limit			
Post time	No limit			
Enable analysis	When temperature is within setpoint +/- 1.2 °C			

Table 1. Method parameters: B6–002 (continued...)

Fluorescence signals –	
Excitation	325 nm
Emission	425 nm
Stop time	7.00 min
Post time	off
Peak width	> 0.2 min (4 s, standard)
PMT-gain	18
Integration events –	
Initial	Slope sensitivity = 1 Peak width = 0.04 Area reject = 1 Height reject = 1.7 Shoulders = OFF
0.000 min	Integration = ON
0.001 min	Integration = OFF
2.700 min	Integration = ON
3.700 min	Integration = OFF
5.100 min	Integration = ON
6.200	Integration = OFF
Calibration settings –	
Use Sample Data	From data file
Amount Units	nmol/L
Default RT windows:	
Reference peaks	0.25 min + 5.00%
Other Peaks	0.25 min + 5.00%
Default calibration curve:	
Type	Linear
Origin	Include
Weight	Equal
If peaks missing	Partial calibration
Calculate uncalibrated peaks	No
Calibration table* –	
PLP:	
Level 1	3.710 nmol/L
Level 2	11.130 nmol/L
Level 3	37.096 nmol/L
Level 4	74.192 nmol/L
Level 5	222.580 nmol/L
4-PA:	
Level 1	3.808 nmol/L
Level 2	11.422 nmol/L
Level 3	38.073 nmol/L
Level 4	76.146 nmol/L
Level 5	228.440 nmol/L

* Amounts shown are examples based on 404.6 µmol/L PLP and 507.8 µmol/L 4-PA working standards. Equations in section 6.b.[(2)–(6)] and actual concentrations of working standards should be used to determine these values.

Table 2. Method parameters: STARTUP*Pump 1 setup –*

Flow	0.7 mL/min			
Stop time	75.00 min			
Post time	Off			
Solvents:				
A: 50 mm phosphate buffer, pH 3.2, 0.2% acetonitrile	0%			
B: acetonitrile	95%			
C: methanol	0%			
D: water	5%			
Pressure limits:				
Max	400 bar			
Min	0 bar			
Timetable:				
25.00 min	95.0% B	0.0% C	5.0% D	0.700 mL/min
30.00 min	0.0% B	30.0% C	70.0% D	0.700 mL/min
55.00 min	0.0% B	30.0% C	70.0% D	0.700 mL/min
60.00 min	0.0% B	0.0% C	0.0% D	0.700 mL/min
75.00 min	0.0% B	0.0% C	0.0% D	0.700 mL/min

Pump 2 setup –

Flow	0.1 mL/min			
Stop time	75.00 min			
Post time	Off			
Solvents:				
A: water	0%			
B: 2g/L sodium chlorite	0%			
C: methanol	100%			
D: not used	Off			
Pressure limits:				
Max	20 bar			
Min	0 bar			
Timetable:				
25.00 min	0% B	100 %C	0% D	0.100 mL/min
30.00 min	0% B	0 %C	0% D	0.100 mL/min
55.00 min	0% B	0% C	0% D	0.100 mL/min
60.00 min	100% B	0% C	0% D	0.100 mL/min
75.00 min	100% B	0% C	0% D	0.100 mL/min

Column thermostat 1 setup –

Temperature	35.0 °C
Stop time	No limit
Post time	No limit
Enable analysis	With any temp.

Column thermostat 2 setup –

Temperature	35 °C
Stop time	No limit
Post time	No limit
Timetable	
30.00 min	35 °C
45.00 min	55 °C
60.00 min	75 °C
Enable analysis	With any temp

Table 2. Method parameters: STARTUP (continued...)*Fluorescence signals –*

Excitation	325 nm
Emission	425 nm
Stop time	0.01 min
Post time	off
Peak width	> 0.2 min (4 s, standard)
PMT-gain	18

Table 3. Method parameters: SHUTDOWN

<i>Pump 1 setup –</i>				
Flow	0.7 mL/min			
Stop time	120.00 min			
Post time	Off			
Solvents:				
A: 50 mm phosphate buffer, pH 3.2, 0.2% acetonitrile	100%			
B: acetonitrile	0%			
C: methanol	0%			
D: water	0%			
Pressure limits:				
Max	400 bar			
Min	0 bar			
Timetable:				
5.00 min	0.0% B	0.0% C	100.0% D	0.700 mL/min
25.00 min	0.0% B	0.0% C	100.0% D	0.700 mL/min
30.00 min	95.0% B	0.0% C	5.0% D	0.700 mL/min
60.00 min	95.0% B	0.0% C	5.0% D	0.700 mL/min
<i>Pump 2 setup –</i>				
Flow	0.1 mL/min			
Stop time	120.00 min			
Post time	Off			
Solvents:				
A: water	0%			
B: 2g/L sodium chlorite	100%			
C: methanol	0%			
D: not used	Off			
Pressure limits:				
Max	20 bar			
Min	0 bar			
Timetable:				
5.00 min	0.0% B	0.0% C	0.0% D	0.100 mL/min
25.00 min	0.0% B	0.0% C	0.0% D	0.100 mL/min
30.00 min	0.0% B	100.0% C	0.0% D	0.100 mL/min
60.00 min	0.0% B	100.0% C	0.0% D	0.100 mL/min
<i>Column thermostat 1 setup –</i>				
Temperature	Not controlled			
Stop time	No limit			
Post time	No limit			
<i>Column thermostat 2 setup –</i>				
Temperature	Not controlled			
Stop time	No limit			
Post time	No limit			
<i>Fluorescence signals –</i>				
Excitation	325 nm			
Emission	425 nm			
Stop time	0.01 min			
Post time	off			
Peak width	> 0.2 min (4 s, standard)			
PMT-gain	18			

Table 4. Method parameters: STANDBY

<i>Pump 1 setup –</i>	
Flow	0.7 mL/min
Stop time	1.00 min
Post time	Off
Solvents:	
A: 50 mm phosphate buffer, pH 3.2, 0.2% acetonitrile	Off
B: acetonitrile	95%
C: methanol	Off
D: water	5%
Pressure limits:	
Max	400 bar
Min	0 bar
<i>Pump 2 setup –</i>	
Flow	0.1 mL/min
Stop time	1.00 min
Post time	Off
Solvents:	
A: water	Off%
B: 2g/L sodium chlorite	Off
C: methanol	100%
D: not used	Off
Pressure limits:	
Max	20 bar
Min	0 bar
<i>Column thermostat 1 setup –</i>	
Temperature	Not controlled
Stop time	No limit
Post time	No limit
<i>Column thermostat 2 setup –</i>	
Temperature	Not controlled
Stop time	No limit
Post time	No limit
<i>Fluorescence signals –</i>	
Excitation	325 nm
Emission	425 nm
Stop time	0.01 min
Post time	off
Peak width	> 0.2 min (4 s, standard)
PMT-gain	18

Table 5. Serum PLP and 4-PA means, ranges and percentiles^{†‡}

Analyte	n	Concentration found (nmol/L)							
		Geometric mean	95% CLs	Range	Percentile				
					10th	25th	50th	75th	90th
PLP	303	42.4	38.8, 46.4	7.6–1260	16.9	25.0	40.1	64.1	103
4-PA	303	27.3	24.4, 30.6	6.2–3140	11.0	14.8	21.8	36.8	84.2

[†] Data taken from reference [15].

[‡] Values were determined using a nonrepresentative population subset and should be used for general orientation purposes only, not as representing the general U.S. population.

- [1] P. Gyorgy, Vitamin B2 and the pellagra-like dermatitis of rats, *Nature* 244 (1934) 448–449.
- [2] IUPAC-IUB Commission on Biochemical Nomenclature, Nomenclature for vitamin B6 and related compounds, *Eur. J. Biochem.* 40 (1973) 325–327.
- [3] J.E. Leklem, Vitamin B6, in: R.B. Rucker, J.W. Suttie D.B. McCormick, L.J. Machlin (Eds.), *Handbook of Vitamins* (3rd ed.), Marcel Dekker Inc., New York, 2001, pp. 339–396.
- [4] L. Chasan-Taber, J. Selhub, I.H. Rosenberg, M.R. Malinow, P. Terry, P.V. Tishler, W. Willett, C.H. Hennekens, M.J. Stampfer, A prospective study of folate and vitamin B6 and risk of myocardial infarction in U.S. physicians, *J. Am. Coll. Nutr.* 15 (1996) 136–143.
- [5] M.C. McKinley, H. McNulty, J. McPartlin, J.J. Strain, K. Pentieva, M. Ward, D.G. Weir, J.M. Scott, Low-dose vitamin B-6 effectively lowers fasting plasma homocysteine in healthy elderly persons who are folate and riboflavin replete, *Am. J. Clin. Nutr.* 73 (2001) 759–764.
- [6] L. Brattström, B. Israelsson, B. Norrving, D. Bergqvist, J. Thorne, B. Hultberg, A. Hamfelt, Impaired homocysteine metabolism in early-onset atherosclerotic and peripheral occlusive arterial disease. Effects of pyridoxine and folic acid treatment, *Atherosclerosis* 81 (1990) 51–60.
- [7] J.B. Ubbink, W.J.H. Vermaak, A. van der Merwe, P.J. Becker, R. Delport, H.C. Potgieter, Vitamin requirements for the treatment of hyperhomocysteinemia in humans, *J. Nutr.* 124 (1994) 1927–1933.
- [8] J. Dierkes, M. Kroesen, K. Pietrzik, Folic acid and vitamin B-6 supplementation and plasma homocysteine concentrations in healthy young women, *Int. J. Vitam. Nutr. Res.* 68 (1998) 98–103.
- [9] P.J. Kelly, V.E. Shih, J.P. Kistler, M. Barron, H. Lee, R. Mandell, K.L. Furie, Low vitamin B6 but not homocyst(e)ine is associated with increased risk of stroke and transient ischemic attack in the era of folic acid grain fortification, *Stroke* 34 (2003) E51–E54.

-
- [10] L. Lumeng, M.P. Ryan, T.K. Li, Validation of the diagnostic value of plasma pyridoxal-5'-phosphate measurements in vitamin B6 nutrition in the rat, *J. Nutr.* 108 (1978) 545–553.
- [11] A.L. Black, B.M. Guirard, E.E. Snell, The behavior of muscle phosphorylase as a reservoir for vitamin B6 in the rat, *J. Nutr.* 108 (1978) 670–677.
- [12] S.P. Coburn, P.J. Ziegler, D.L. Costill, J.D. Mahuren, W.J. Fink, W.E. Schaltenbrand, T.A. Pauly, D.R. Pearson, P.S. Conn, T.R. Guilarte, Response of vitamin B-6 content of muscle to changes in vitamin B-6 intake in men, *Am. J. Clin. Nutr.* 53 (1991) 1436–1442.
- [13] J.E. Leklem, Vitamin B-6: a status report, *J. Nutr.* 120 (1990) 1503S–1507S.
- [14] J. Zemleni, Pharmacokinetics of vitamin B6 supplements in humans, *J. Am. Coll. Nutr.* 14 (1995) 579–586.
- [15] M.E. Rybak, C.M. Pfeiffer, Clinical analysis of vitamin B6: Determination of pyridoxal 5'-phosphate and 4-pyridoxic acid in human serum by reversed-phase high-performance liquid chromatography with chlorite postcolumn derivatization, *Anal. Biochem.* 333 (2004) 336-344.
- [16] M.E. Rybak, R.B. Jain, C.M. Pfeiffer, Clinical vitamin B6 analysis: an interlaboratory comparison of pyridoxal 5'-phosphate measurements in serum. *Clin. Chem* 51 (2005) 1223–1231.
- [17] C.G. Fraser, P.H. Petersen, J.-C. Libeer, C. Ricos, Proposals for setting generally applicable quality goals solely based on biology, *Ann. Clin. Biochem.* 34 (1997) 8–12.
- [18] A. Lui, L. Lumeng, G.R. Aronoff, T.-K. Li, Relationship between body store of vitamin B6 and plasma pyridoxal-P clearance: Metabolic balance studies in humans, *J. Lab. Clin. Med.* 106 (1985) 491–497.