



Laboratory Procedure Manual

Analytes: **Antimony, Arsenic, Barium, Beryllium, Cadmium, Cesium, Cobalt, Lead, Manganese, Molybdenum, Platinum, Strontium, Thallium, Tin, Tungsten, and Uranium**

Matrix: **Urine**

Method: **Urine Multi-Element ICP-DRC-MS**
Renamed from "Inductively Coupled Plasma-Mass Spectrometry (ICP-DRC-MS)"

Method No: **3018.6-02 (15 element panel) and 3018A.4-02 (total arsenic)**

As performed by: Inorganic and Radiation Analytical Toxicology
Division of Laboratory Sciences
National Center for Environmental Health

Contact: Dr. Kathleen L. Caldwell
Phone: 770-488-7990
Fax: 770-488-4097
Email: KCaldwell@cdc.gov

James L. Pirkle, M.D., Ph.D.
Director, Division of Laboratory Sciences

Important Information for Users

The Centers for Disease Control and Prevention (CDC) periodically refines these laboratory methods. It is the responsibility of the user to contact the person listed on the title page of each write-up before using the analytical method to find out whether any changes have been made and what revisions, if any, have been incorporated.

Public Release Data Set Information

This document details the Lab Protocol for testing the items listed in the following table:

This method file describes measurements of UM_I, UMS_I, UTAS_I, and UTASS_I. One method was used to measure the urinary metals, urinary metals for smokers, total arsenic and total arsenic for smokers. However, these results are released as 4 separate data files.

| Data File Name | Variable Name | SAS Label |
|-------------------|---------------|------------------------------|
| UTAS_I UTASS_I | URXUAS | Urinary Total Arsenic (µg/L) |
| UM_I UMS_I | URXUBA | Barium, urine (µg/L) |
| | URXUCD | Cadmium, urine (µg/L) |
| | URXUCO | Cobalt, urine (µg/L) |
| | URXUCS | Cesium, urine (µg/L) |
| | URXUMN | Manganese, urine (µg/L) |
| | URXUMO | Molybdenum, urine (µg/L) |
| | URXUPB | Lead, urine (µg/L) |
| | URXUSN | Tin, urine (µg/L) |
| | URXUSR | Strontium, urine (µg/L) |
| | URXUSB | Antimony, urine (µg/L) |
| | URXUTL | Thallium, urine (µg/L) |
| | URXUTU | Tungsten, urine (µg/L) |
| | URXUUR | Uranium, urine (µg/L) |

1) Clinical relevance & summary of test principle

a. Clinical relevance:

These methods are used to achieve rapid and accurate quantification of elements of toxicological and nutritional interest including Antimony (Sb), Barium (Ba), Beryllium (Be), Cadmium (Cd), Cesium (Cs), Cobalt (Co), Lead (Pb), Manganese (Mn), Molybdenum (Mo), Platinum (Pt), Strontium (Sr), Thallium (Tl), Tin (Sn), Tungsten (W), and Uranium (U) or Arsenic (As). Use these methods to screen urine when people are suspected to be acutely exposed to these elements or to evaluate chronic environmental or other non-occupational exposure. [1-4].

b. Test principle:

Inductively coupled plasma mass spectrometry (ICP-MS) is a multi-element analytical technique capable of trace level elemental analysis [1-4]. When used with dynamic reaction cell technology (DRC) the technique is referred to as ICP-DRC-MS. This ICP-DRC-MS method is used to measure either arsenic, a 15 element panel (antimony, barium, beryllium, cadmium, cesium, cobalt, lead, manganese, molybdenum, platinum, strontium, thallium, tin, tungsten, and uranium), or any subgroup of these.

Liquid samples are introduced into the ICP through a nebulizer and spray chamber carried by a flowing argon stream. By coupling radio-frequency power into flowing argon, plasma is created in which the predominant species are positive argon ions and electrons and has a temperature of 6000-8000 K. The sample passes through a region of the plasma and the thermal energy atomizes the sample and then ionizes the atoms. The ions, along with the argon, enter the mass spectrometer through an interface that separates the ICP (at atmospheric pressure, ~760 torr) from the mass spectrometer (operating at a pressure of 10^{-5} torr). The ions pass through a focusing region, the dynamic reaction cell, the quadrupole mass filter, and finally are counted in rapid sequence at the detector allowing individual isotopes of an element to be determined. The dynamic reaction cell operates in one of two modes. In 'standard' mode, the cell is not pressurized and ions pass through the cell to the quadrupole mass filter unaffected. In 'DRC' mode, the cell is pressurized with a gas, which will collide or react with the incoming ions to either eliminate an interfering ion or change the ion of interest to a new mass, which is free from interference. In this method, the instrument is operated in DRC mode when analyzing for cadmium, manganese and arsenic, but in standard mode when analyzing for all of the other analytes. For arsenic, the reaction cell is pressurized with a mixture of hydrogen (10%) and argon (90%) or 100% argon which causes the breakup of the $^{40}\text{Ar}^{35}\text{Cl}^+$ ion which would otherwise interfere with detection of ^{75}As at m/z 75. When analyzing for cadmium in biomonitoring applications, the reaction cell is pressurized with oxygen. The $^{98}\text{Mo}^{16}\text{O}^+$ ions, which would normally interfere with detection of ^{114}Cd at m/z 114, react with the oxygen in the cell creating $^{98}\text{Mo}^{16}\text{O}_2^+$ and $^{98}\text{Mo}^{16}\text{O}_3^+$ at masses that no longer represent interference to low level ^{114}Cd analysis. When low level Cd analysis is not the principle purpose of the analysis (i.e. emergency response situations), Cd analysis in vented (standard) mode will reduce analytical time and still yield quantitative results that are suitable for the identification of elevated exposures as long as the results are interpreted with the caveat that the $^{98}\text{Mo}^{16}\text{O}_2^+$ interference on ^{114}Cd is not eliminated. The DRC is also pressurized with oxygen gas when analyzing for ^{55}Mn . The $^{39}\text{K}^{16}\text{O}^+$ ions, which

would normally interfere with the detection of ^{55}Mn at m/z 55, react with the oxygen in the cell and no longer represent interference to ^{55}Mn analysis. In DRC mode, the voltage applied by the axial field technology (AFT) and the additional axial push from spectator ions also stable in the DRC bandpass serve to keep the ions moving axially through the pressurized DRC chamber where they would normally slow down due to collisions which result in loss of momentum. Gold is added to the diluent to normalize the spectator ion population in the DRC cell, which could otherwise change significantly between low concentration and high concentration samples, a phenomenon called crosstalk [3, 5, 6]. Electrical signals resulting from the detection of ions are processed into digital information that is used to indicate first the intensity of the ions and then the concentration of the element. This method was originally based on the method by Mulligan et al. [7]. The DRC portions of the method are based on work published by Tanner et al. [2, 3]. Urine samples are diluted 1+ 9 with 2% (v/v) concentrated nitric acid (and 1.5% ethanol in the case of arsenic). The diluent for the 15-element panel contains iridium (Ir), rhodium (Rh) for multi-internal standardization. The diluent for arsenic contains gallium (Ga) for internal standardization. Nitric acid is used for the purpose of solubilizing and stabilizing metals in solution. Internal standards are a constant concentration in all blanks, calibrators and samples. Monitoring the instrument signal ratio of a metal to its internal standard allows correction for instrument noise and drift, and sample-to-sample matrix differences. Ethanol is used when analyzing for arsenic in biomonitoring situations for the purpose of providing a constant amount of signal enhancement (carbon effect) across all blanks, calibrators, and samples.

2) Limitations of method; interfering substances and conditions

a. Interferences addressed by this method

i. Argon chloride ($^{40}\text{Ar}^{35}\text{Cl}$) on arsenic (^{75}As): The dynamic reaction cell is used in this method to diminish the presence of the argon chloride ($^{40}\text{Ar}^{35}\text{Cl}$) interference on arsenic at m/z 75 [8] which is common to urine analysis by ICP-MS (see Section 1.b for an explanation of this process). The dynamic reaction cell gas used for this purpose is a mixture of hydrogen (10%) in argon; however, 100% argon can also be used when the mixture is unavailable.

ii. Correction & elimination of interferences (^{114}Sn , $^{98}\text{Mo}^{16}\text{O}$) on cadmium (^{114}Cd).

1. Mathematical correction for tin (^{114}Sn) interference:

The correction equation ($-0.026826 \cdot \text{Sn}118$) is used in the "Equations" tab of the method to correct the counts observed as m/z 114 to exclude counts due to ^{114}Sn .

2. Elimination of molybdenum oxide ($^{98}\text{Mo}^{16}\text{O}$) interference using DRC:

The dynamic reaction cell is used in this method to eliminate interference from molybdenum oxide ($^{98}\text{Mo}^{16}\text{O}$) onto cadmium at m/z 114 [9]. Oxygen (100%) is the gas used in the dynamic reaction cell for this purpose. When low level Cd analysis is not the principle purpose of the analysis (i.e. emergency response situations), Cd analysis in vented (standard) mode will reduce analytical time and still yield quantitative results that are suitable for the identification of elevated exposures as long as the results are interpreted with the caveat that the $^{98}\text{Mo}^{16}\text{O}_2^+$ interference

on ^{114}Cd is not eliminated. The anticipated bias on ^{114}Cd due to the $^{98}\text{Mo}^{16}\text{O}_2^+$ interference is described approximately by equation 1 [9].

$$\text{approximate } \mu\text{g/L difference (bias)} = 0.00175[\text{Mo}] - 0.0136 \quad [\text{equation 1}]$$

iii. Reduction of interference ($^{39}\text{K}^{16}\text{O}$) on manganese, ^{55}Mn , using DRC:

The dynamic reaction cell is used in this method to reduce the potassium oxide ($^{39}\text{K}^{16}\text{O}$) interference on manganese at m/z 55. See Section 1.b for an explanation of this process.

iv. Matrix enhancement of arsenic signal:

Matrix induced signal enhancement in ICP-MS analysis from carbon on arsenic has been previously reported in the literature [10, 11]. When arsenic is being determined by this method for normal biomonitoring purposes, ethanol (1.5% v/v) is added in the diluent and rinse solutions to “normalize” the arsenic signal enhancement in all blanks, calibrators, and samples. If arsenic is combined with multi-element analysis for emergency response situations where ethanol is not part of the diluent and rinse solutions, the positive bias potentially resulting is not anticipated to significantly hinder the identification of acute arsenic exposures.

b. Limitations of method (interferences remaining in method)

i. Calcium chloride ($^{40}\text{Ca}^{35}\text{Cl}$) interference on arsenic (^{75}As):

It has been determined that a small interference remains at m/z 75 when the urine matrix contains **both** high chloride **and** high calcium levels [8]. Even at extreme calcium and chloride levels, this interference is has not been found to be significant (approximately 0.4 $\mu\text{g/L}$).

ii. Gallium oxide ($^{71}\text{Ga}^{17}\text{O}$) interference on strontium (^{88}Sr):

Gallium at 10 $\mu\text{g/L}$ in the method diluent or 100 $\mu\text{g/L}$ in a urine specimen (extremely unlikely) will produce approximately a +2 $\mu\text{g/L}$ bias on observed Sr concentration due to the formation of $^{71}\text{Ga}^{17}\text{O}^+$, which occurs at the same m/z as $^{88}\text{Sr}^+$. While this interference is not large in comparison to typical Sr concentrations in urine, ~150 $\mu\text{g/L}$ (see Table 11 in Appendix B), best urine Sr accuracy for biomonitoring purposes will be achieved when Ga is not present in the diluent as an internal standard. If As is combined with the multi-element method for emergency response purposes, the resulting bias on Sr from the $^{71}\text{Ga}^{17}\text{O}^+$ interference is not likely to hinder the ability to hinder the identification of acute Sr exposures.

3) **Procedures for collecting, storing, and handling specimens; criteria for specimen rejection; specimen accountability and tracking**

a. Procedures for collecting, storing, and handling specimens: Guidelines for receiving and shipping packages and an example case study specimen collection protocol are presented in the laboratory Policies and Procedures Manual [12]. Special specimen handling conditions, requirements, and procedures for this method include:

i. No fasting or special diets are required before collection of urine.

- ii. Use sterile, lot-screened collectors for specimen acquisition.
 - iii. Transport urine specimens frozen (packed in dry ice during shipment is preferred when possible).
 - iv. Once received, store long term at ≤ -20 °C until time for analysis. Short-term storage at 2-8 °C is acceptable. Refreeze at ≤ -20 °C portions of the sample that remain after analytical aliquots are withdrawn. Thawing and refreezing samples has not been found to compromise sample results.
 - v. Acceptable containers for analytical aliquots include lot screened polypropylene (PP) cryovials or tubes (i.e. 2 to 5 mL cryogenic vial or 15mL centrifuge tube). Avoid colored plastics and containers containing o-rings when possible due to increased risk of trace element contamination from coloring agents or o-ring materials. Externally threaded containers are preferred because they are less prone to contamination of the specimen and to leaks (internally threaded containers can develop leaks when biological material dries within the threads, compromising resealing).
- b. Criteria for specimen rejection: Specimen characteristics that compromise test results are indicated above. Reasons for rejection of a sample for analysis include the following (in all cases, request a second urine specimen):
- i. Low volume:
 - 1. Method 3018 (15 element): Optimal amount of urine is 1.8+ mL. The volume of urine used for one analysis is 0.5 mL. Less volume is consumed when only subsets of the 15 elements are analyzed.
 - 2. Method 3018A (As): Optimal amount of urine is 1+ mL. The volume of urine used for one analysis is 0.25 mL.
 - ii. Contamination: Improper collection procedures or collection devices can contaminate the urine by contact with dust, dirt, etc.

- c. Transfer or referral of specimens; procedures for specimen accountability and tracking: Location, status, and final disposition of the specimens will be tracked and records are maintained according to the Division's Policies and Procedures Manual [12]. Use only numerical identifiers for samples within the laboratory (e.g., case ID numbers) in order to safeguard confidentiality. Only the medical supervisor (MS) or project coordinator (PC) i.e. non-CDC personnel will have access to the personal identifiers.

4) Safety precautions

a. General safety

- i. Observe all safety regulations as detailed in the Laboratory Safety Manual and the Chemical Hygiene Plan.
- ii. Wear gloves, lab coat, and safety glasses while handling reagents, prepared solutions, or urine specimens.

Stock calibration standard B contains 1% hydrofluoric acid (HF). Latex is not an appropriate protective barrier against HF. When handling these solutions nitrile gloves (or other glove material known to be an adequate protective barrier to HF) are required.

- iii. Observe universal precautions when working with urine.
- iv. Exercise special care when handling and dispensing concentrated nitric or hydrochloric acid. Use additional personal protective equipment, which protects face, neck, and front of body. Add acid to water. Nitric and hydrochloric acids are caustic chemicals that are capable of causing severe eye and skin damage. ***If concentrated acids come in contact with any part of the body, quickly wash the affected area with copious quantities of water for at least 15 minutes.***
- v. Use secondary containment for containers of biological or corrosive liquids.
- vi. The use of the foot pedal on the benchtop automatic pipette is recommended because it reduces analyst contact with work surfaces that have been in contact with urine and also keeps the analyst's hands free to hold the specimen cups and autosampler tubes and to wipe off the tip of benchtop automatic pipette.
- vii. There are many potential hazards on an operating ICP-MS instrument including ultraviolet radiation, high voltages, radio-frequency radiation, and high temperatures. This information is detailed in the ICP-MS System Safety Manual.
- viii. Transport and store compressed gas cylinders with proper securing harnesses. For compressed oxygen gas, use regulators, which are oil-free and are equipped with a flash arrestor. Use of flash arrestors is recommended when working with compressed oxygen and compressed hydrogen at greater than a 7% composition.
- ix. Wipe down all work surfaces at the end of the day with freshly prepared 10% (v/v) sodium-hypochlorite solution or comparable disinfectant.

- b. Radiation safety: Calibration standards used in this method contain $\mu\text{g/L}$ natural uranium. Staff performing this method in the CDC laboratory must maintain the status of Radiation Worker through the Radiation Safety Office and practice appropriate radiation safety when handling these solutions in accordance to the CDC's license with the Nuclear Regulatory Commission (NRC).
- c. Waste disposal:
 - i. Autoclaving: All diluted biological specimens, original biological specimens being disposed, or consumables, which come into contact with biological specimens (even diluted or aerosolized). Use sharps containers or special autoclave pans for broken glass / quartz or items, which puncture autoclave bags (e.g. pipette tips).
 - ii. Other liquid waste
 1. Waste discarded down sink: Only non-corrosive liquid waste (EPA defines as $\text{pH} > 2$ and $\text{pH} < 12.5$, 40CFR §261.22) from the ICP-DRC-MS instrument can be discarded at the sink. Flush the sink with copious amounts of water.
 2. Waste to be picked up by the CDC radiation safety office: Solutions used in the CDC laboratory having uranium concentrations equal to that of the single element standard, intermediate stock standard, or intermediate working standards.
 3. Waste to be picked up by CDC hazardous waste program: Submit request for hazardous waste removal of all other liquid waste generated in the CDC laboratory for this method.

4.

5) Instrument & material sources

- a. Sources for ICP-MS instrumentation
 - i. ICP-MS: Inductively coupled plasma mass spectrometer with dynamic reaction cell technology (ELAN[®] DRC II or NexION) (PerkinElmer Norwalk, CT, www.perkinelmer.com).
 - ii. Recirculating chiller / heat exchanger for ICP-MS: Refrigerated chiller (PolyScience 6105PE) or heat exchanger (PolyScience 3370) (PerkinElmer Norwalk, CT, www.perkinelmer.com).
 - iii. Autosampler: ESI SC4-DX autosampler (Elemental Scientific Inc., Omaha, NE) or equivalent.
 - iv. Computer: Computer controller provided or recommended by ICP-MS manufacturer is recommended to ensure proper communication between computer and ICP-MS. Recommend 1-2 Gb RAM and secondary internal hard disk for nightly backups (if network backups are not possible).

v. FAST sample introduction system (optional): Standard peristaltic pump on ICP-MS replaced by DXi-FAST micro-peristaltic pump / FAST actuator and valve combination unit. For NEXION, like part # DXI-54-P4-F6. If DXi-FAST upgrade on ICP-MS is not used, a separate FAST actuator (built-in option on ESI SC4-DX autosampler or stand-alone FAST actuator) will be necessary to complete the FAST sample introduction system.

b. Sources for ICP-MS parts & consumables

NOTE: The minimum number of spares recommended before reordering (if owning one instrument) are listed as “# Spares =” in the descriptions below.

i. Adapter, PEEK: Securely connects 1.6mm O.D. PFA tubing to 0.03” I.D. peristaltic tubing. Composed of three PEEK parts.

1. Female nut for 1.6mm O.D. (1/16”) tubing. Like part P-420 (Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).
2. PEEK ferrule. Like part P-260x (10pk SuperFlangeless ferrule, Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).
3. Conical Adapter Body. Like part P-692 (Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).

ii. Bottles (for rinse solution): Four liter screw-cap polypropylene container with built-in luer connections (2) designed for use with FAST sample introduction system (like catalog# SC-0305-1, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).

iii. Carboy and cap assembly for waste collection: 10-15 L, polypropylene wide-mouth carboy (100 mm neck size) with handles and no spigot (Like part #7BE-25126, Lab Safety Supply, Janesville, WI, www.lss.com) with cap assembly like part # N0690271 (PerkinElmer, Norwalk, CT, www.perkinelmer.com) with tubing connections built into the cap for addition of liquid waste.

iv. Coolant, for Polyscience chiller or heat exchanger: Only PerkinElmer part # WE01-6558 (PerkinElmer Norwalk, CT, www.perkinelmer.com) is approved for use by PerkinElmer. # Spares = 6.

v. Cone, hyperskimmer (NexION): PerkinElmer part # W1033995 (PerkinElmer Norwalk, CT, www.perkinelmer.com).

1. Screws (for hyper skimmer cone, NexION): PerkinElmer part # 09919737 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 4 screws per instrument.

vi. Cone, sampler (nickel/platinum):

1. ELAN ICP-MS: PerkinElmer part # WE021140 / WE027802 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Part # SC2011-Ni / SC2013-Pt (Testing has also found Spectron, Ventura, CA, www.spectronus.com cones to be comparable). # Spares = 4.
2. NexION ICP-MS: PerkinElmer part # W1033612 / W1033614 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Part # SC4011-Ni / SC4013-Pt (Testing has

also found Spectron, Ventura, CA, www.spectronus.com cones to be comparable).
Spares = 4.

vii. Cone, skimmer (nickel / platinum):

1. ELAN ICP-MS: PerkinElmer part # WE021137 / WE027803 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Part # SC2012-Ni / SC2014-Pt (Testing has also found Spectron, Ventura, CA, www.spectronus.com cones to be comparable) # Spares = 4.
2. NexION ICP-MS: PerkinElmer part # W1026356 / W1026907 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Part # SC4012-Ni / SC4014-Pt (Testing has also found Spectron, Ventura, CA, www.spectronus.com cones to be comparable) # Spares = 4.

viii. Connector (for tubing): Use to connect 1/8" I.D. PVC tubing to 0.125" I.D peristaltic pump tubing. Use part # 3140715 (PerkinElmer Norwalk, CT, www.perkinelmer.com) or equivalent. # Spares = 4.

ix. Detector, electron multiplier: Like part # N8125001 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Available direct from manufacturer (part # 14210, SGE Incorporated, Austin, Texas, <http://www.etpsci.com>) or various distributors. # Spares = 1.

x. FAST / ESI SC4-DX autosampler accessories:

1. Valve: CTFE High-flow valve head for SC-FAST (uses 1/4-28 fittings). Like part # SC-0599-1010 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
2. Stator: CTFE Stator for 6 port SC-FAST high flow valve (1/4-28 fittings). Like part # SC-0599-1010-01 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
3. Rotor: Composite rotor for 6 port SC-FAST high flow valve (1/4-28 fittings). Like part # SC-0599-1010-05 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
4. Sample loop:
 - a. DLS 3018: 3 mL Teflon loop with white connector-nuts for high flow valve head, like part # SC-0315-30 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). Subsets of elements can be analyzed using different loop sizes to minimize sample consumption (e.g. 0.5 mL loop for a single element, 1.0 mL loop for 3 element subset, 2.0mL loop for 8 element subset, etc...).
 - b. DLS 3018A: 0.5mL Teflon sample loop with white nut connectors for high flow valve head of FAST sample introduction system, like part # SC-0315-20 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
5. Probe, autosampler: Teflon, carbon fiber support, 0.8mm i.d., blue marker, 1/4-28 fittings. Like part number SC-5037-3751 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.

6. Probe, carrier solution: Teflon, carbon fiber support, 0.5mm i.d., orange marker, 1/4-28 fittings. Like part number SC-5037-3501 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.
 7. Tubing, carrier solution: 0.5mm i.d. Teflon tubing (orange marker) with red ¼-28 male nut. Connects to high flow FAST valve head, port #2. Like part # SC-0316-0500 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
 8. Tubing, nebulizer: See “Nebulizer, PolyPro-ST micro flow”
 9. Tubing, rinse station: Teflon tubing and adapters (to attach to back of SC autosampler for filling rinse stations and to attach to rinse containers). Like part # SC-0302-0500, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
 10. Tubing, vacuum: Vacuum line for SC-FAST high flow valve, connects to port #6, black nut for connection to valve head, natural brown color nut on other end for connection to SC autosampler vacuum port. Like part # SC-0321 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
 - xi. Hose, for connection to recirculator / chiller: Push on hose. I.D. = ½”, O.D. = ¾”. Use part # PB-8 (per inch, Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. No spares necessary.
 - xii. Hose, for exhaust of ICP-MS: Available as part of ICP-MS installation kit from PerkinElmer (PerkinElmer Norwalk, CT, www.perkinelmer.com), or equivalent. Available direct from manufacturer as part # S-LP-10 air connector (Thermafex, Abbeville, SC, www.thermafex.net). # Spares = 10 feet of 4” diameter (ELAN and NexION) and 10 feet of 6” diameter hose (ELAN).
 - xiii. Injector, quartz with ball joint: I.D. = 2.0 mm. PerkinElmer part # WE023948 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Available direct from manufacturer as part # 400-30 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com) or from various distributors. # Spares = 2.
 - xiv. Ion lens (ELAN): PerkinElmer part # WE018034 (PerkinElmer Norwalk, CT, www.perkinelmer.com). # Spares = 1.
 - xv. Nebulizer: PolyPro-ST micro flow polypropylene nebulizer with external 1/4-28 threaded connector for liquid delivery, low pressure version or equivalent. Like part # ES-4040-7010 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 1. Comparable nebulizers are acceptable for substitution, however, the nebulizer gas flow rate, sample flush time, read delay time, loop fill time, loop size, urine sample dilution preparation volume, and sample-to-sample carry-over must be evaluated and optimized.
1. Gas connection:
 - a. Teflon tubing: 4mm o.d., 2.4mm i.d. Teflon tubing (like part # ES-2502, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 1.
 - b. Adapter kit: Plastic adapters to connect Teflon tubing (2.4mm i.d.) to ¼” male Swagelok (compression) port on ICP-DRC-MS. Parts can be obtained as components in a “gas fittings kit for microflow nebulizer”, kit part # ES-2501-

- 1000, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
Spares = 1.
2. Liquid connection: Connects nebulizer to port #3 of high flow FAST valve head with green, 1/4- 28 fitting. Like part # SC-0317-0250 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.
- xvi. Nut: (for flanged connections of 1.59mm (1/16") o.d. PFA tubing) Flanged, for 1/16" o.d. tubing, 1/4-28 threads. Use part # P-406x (pkg. of 10, Upchurch Scientific, Oak Harbor, WA, www.upchurch.com) or equivalent. Use a Teflon-coated Viton o-ring with this nut instead of the stainless steel washer that comes with part # P-406x). # Spares = 10.
- xvii. Nut and ferrule set, 1/8" Swagelok: Such as part # SS-200-NFSET (stainless steel) or part # B-200-NFSET (brass) (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. For part numbers listed here a quantity of 1 means 1 nut, 1 front ferrule, and 1 back ferrule. Spares = 20.
- xviii. Nut and ferrule set, 1/4" Swagelok: Such as part # SS-400-NFSET (stainless steel) or part # B-400-NFSET (brass) (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. For part numbers listed here a quantity of 1 means 1 nut, 1 front ferrule, and 1 back ferrule. Spares = 20.
- xix. Oil for roughing pumps:
1. Welch Directorr Gold: For roughing pumps. Available direct from manufacturer as part # 8995G-15 (1 gallon, Welch Rietschle Thomas, Skokie, IL, www.welchvacuum.com), or equivalent. # Spares = 4.
 2. Fomblin Y14/5 fluid: PerkinElmer part # N8122265 (1 kg bottle, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 1 per instrument.
- xx. O-ring (for hyper skimmer cone, NexION): PerkinElmer part # 09902123 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 o-rings.
- xxi. O-ring / gasket (for sampler cone):
1. ELAN (o-ring): PerkinElmer part # N8120511 (pkg. of 5, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 o-rings.
 2. NexION (aluminum gasket): PerkinElmer part # WE012989 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 gaskets.
- xxii. O-ring (for skimmer cone):
1. ELAN: PerkinElmer part # N8120512 (pkg. of 5, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 o-rings.
- xxiii. O-ring: (for flanged connections of 1.59mm (1/16") o.d. PFA tubing) Teflon-coated Viton o-ring, i.d. = 1/16", thickness = 1/16", o.d. = 3/16". Such as part # V75-003 (O-rings West, Seattle, WA, www.oringswest.com) or equivalent. # Spares = 20.
- xxiv. O-ring: (for injector support).
1. Internal o-rings: ID = 1/4", OD = 3/8", thickness = 1/16". Need 2 o-rings per injector support setup. PerkinElmer part # N8122008 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent (such as part # V75-010, O-rings West, Seattle, WA, www.oringswest.com). # Spares = 20.

2. External o-rings: ID = 3/8", OD = 1/2", thickness = 1/16". Need 2 o-rings for each injector support setup. PerkinElmer part # N8122009 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent (such as part # V75-012, O-rings West, Seattle, WA, www.oringswest.com). # Spares = 20.
- xxv. O-ring (for inside nebulizer port on standard PerkinElmer cyclonic quartz spray chamber for the ELAN): Such as part # 120-56 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com). Additional o-rings can sometimes be obtained free of charge or at reduced price when acquired while purchasing spray chambers. # Spares = 20.
- xxvi. O-ring (for inside of ELAN bayonet torch mount): Part # WE017284 (PerkinElmer, Shelton, CT, www.perkinelmer.com). Do not substitute. The PerkinElmer o-ring is metal impregnated to minimize RF leakage through the torch mount. # Spares = 2.
- xxvii. Photon Stop (ELAN): PerkinElmer part # WE018278 (PerkinElmer, Shelton, CT, www.perkinelmer.com). # Spares = 1.
- xxviii. Plugs, quick change for roughing pump oil: These plugs will only work on the roughing pumps which come standard on ELAN DRC II and NexION ICP-MS instruments. These plugs will not fit the Leybold pumps which come standard on ELAN DRC Plus instruments. Part # W1011013 (PerkinElmer, Shelton, CT, www.perkinelmer.com). No spares typically needed.
- xxix. RF coil. PerkinElmer part # WE02-1816 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 2.
- xxx. Spray chamber, quartz concentric:
 1. ELAN: PerkinElmer part # WE025221 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. Available direct from manufacturer as part # 400-20 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com) or from various distributors. # Spares = 2.
 2. NexION: PerkinElmer part # N8145013 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 2.
- xxxi. Torch, quartz: PerkinElmer part # N812-2006 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. Available direct from manufacturer as part # 400-10 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com) or various distributors. # New Spares = 2.
- xxxii. Tubing, main argon delivery to instrument: I.D. = 1/8", O.D. = 1/4". Such as part # C-06500-02 (pkg. of 100ft, polypropylene, Fisher Scientific International, Hampton, NH, www.fishersci.com) or equivalent. # Spares = 50ft.
- xxxiii. Tubing, peristaltic, 0.03" i.d. (carrier solution for ESI autosampler): *use either*
 1. Standard PVC, 2-stop (black / black) peristaltic pump tubing, i.d. = 0.03". PerkinElmer part # 09908587 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 6 packs of 12 tubes.
 2. Standard PVC, 3-stop (black/ black/black) peristaltic pump tubing, i.d. 0.76 mm. Spectron part # SC0056 (Spectron, Ventura, CA, www.spectronus.com) or equivalent. #Spares = 6 packs of 12 tubes. *Use this type of tubing with ESI DXi micro-peristaltic pump.*

xxxiv. Tubing, peristaltic, 0.125" i.d. (spray chamber drain): use either

1. Standard PVC, 2-stop (black / white) peristaltic pump tubing, i.d. = 0.125" or equivalent. PerkinElmer part # N812-2012 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 6 packs of 12 tubes.
2. Standard Santoprene, 3-stop (grey/ grey/ grey) peristaltic pump tubing, i.d. 1.30 mm. Spectron part # SC0311 (Spectron, Ventura, CA, www.spectronus.com) or equivalent. #Spares = 6 packs of 12 tubes. Use this type of tubing with ESI DXi micro-peristaltic pump.

xxxv. Tubing, PVC, i.d. = 1/8", o.d. = 3/16". Used to transfer liquid

1. between spray chamber waste port and peristaltic pump
2. between peristaltic pump and liquid waste jug

Like part # 14-169-7A (pkg. of 50ft, Fisher Scientific International, Hampton, NH, www.fishersci.com) or equivalent. # Spares = 20ft.

xxxvi. Tubing, stainless steel, o.d. = 1/8", wall thickness = 0.028": Used to connect DRC gas cylinders to ICP-MS gas ports. Like part # SS-T2-S-028-20 (20ft, Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 20ft.

xxxvii. Tubing, Teflon, corrugated, 1/4" o.d.: Connects to the auxiliary and plasma gas side-arms of the torch. Part # WE015903 (PerkinElmer, Shelton, CT, www.perkinelmer.com). # Spares = 2.

xxxviii. Union elbow, PTFE 1/4" Swagelok (ELAN bayonet mount): Connects argon tubing to torch auxiliary gas sidearm on bayonet mount NEXION ICP-MS instruments. Like part # T-400-9 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 2.

xxxix. Union tee, PTFE, 1/4" Swagelok (ELAN bayonet mount): Connects argon tubing to torch plasma gas sidearm and holds igniter inside torch sidearm on bayonet mount NEXION ICP-MS instruments. Like part # T-400-3 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 2.

c. Sources for ICP-MS maintenance equipment & supplies

- i. Anemometer: Like digital wind-vane anemometer (Model 840032, SPER Scientific LTD., Scottsdale, AZ, www.sperscientific.com) or equivalent. Use to verify adequate exhaust ventilation for ICP-MS (check with hoses fully disconnected).
- ii. Pan, for changing roughing pump oil: Like part # 53216 (United States Plastics Corporation, Lima, OH, www.usplastic.com) or equivalent. # Spares = 1.
- iii. Container, to hold acid baths for glassware: Polypropylene or polyethylene containers with lids (must be large enough for torch, injector, or spray chamber submersion). Available from laboratory or home kitchen supply companies. # Spares = 4.
- iv. Cotton swabs: Any vendor. For cleaning of cones and glassware.
- v. Cutter (for 1/8" o.d. metal tubing): Terry tool with 3 replacement wheels. Like part # TT-1008 (Chrom Tech, Inc., Saint Paul, MN, www.chromtech.com) or equivalent.

- vi. Getter regeneration kit: Part # WE023257 (PerkinElmer, Shelton, CT, www.perkinelmer.com). Use this as needed (at least annually) to clean the getter in the pathway of channel A DRC gas.
 - vii. Magnifying glass: Any 10x+ pocket loupe for inspection of cones and other ICP-MS parts. Plastic body is preferred for non-corrosion characteristics. Like part # 5BC-42813 (Lab Safety Supply, Janesville, WI, www.labsafety.com).
 - viii. Screw driver, for ion lens removal: Screw driver with long, flexible shaft, and 2mm ball-Allen end for removal of ion lens screws part # W1010620. Extra 2mm bits, part # W1010598 (PerkinElmer, Shelton, CT, www.perkinelmer.com). This is not necessary if the lens is mounted in a quick-release mount.
 - ix. Ultrasonic bath: Like ULTRASONIK™ Benchtop Cleaners (NEYTECH, Bloomfield, CT, www.neytech.com) or equivalent.
- d. Sources for general laboratory consumable supplies
- i. Bar code scanner: Like Code Reader 2.0 (Code Corporation, Draper, UT, www.codecorp.com) or equivalent. For scanning sample IDs during analysis setup. Any bar code scanner capable of reading Code 128 encoding at a 3 mil label density can be substituted.
 - ii. Carboy (for preparation of urine quality control pool and waste jug for ICP-MS sample introduction system): Polypropylene 10-L carboy (like catalog # 02-960-20C, Fisher Scientific, Pittsburgh, PA, www.fishersci.com) or equivalent. Carboys with spouts are not advised due to potential for leaking.
 - iii. Containers for diluent and rinse solution: Two liter Teflon™ containers (like catalog# 02-923-30E, Fisher Scientific, Pittsburgh, PA., www.fishersci.com) and 4L polypropylene jugs (like catalog# 02-960-10A, Fisher Scientific, Pittsburgh, PA, www.fishersci.com) have both been used. Acid rinse before use.
 - iv. Cups for urine collection: Like polypropylene 4.5 oz cup, catalog # 354013 (Becton Dickinson Labware, Franklin Lakes, NJ, www.bd.com) or equivalent. Each lot of cups used must be lot screened (tested to be free of trace metal contamination). Colorless plastics tend to have lowest trace metal contamination.
 - v. Gloves: Powder-free, low particulate nitrile (like Best Clean-DEX™ 100% nitrile gloves, any vendor), or equivalent. Use only nitrile for handling calibration stock standard solution B which contains 1% HF.
 - vi. Paper towels: For general lab use, any low-lint paper wipes such as KIMWIPES®EX-L Delicate Task Wipers or KAYDRY®EX-L Delicate Task Wipers (Kimberly-Clark Professional, Atlanta, GA, www.kcprofessional.com). For sensitive applications in cleanrooms, use a wipe designed for cleanrooms such as the Econowipe or Wetwipe (Liberty, East Berlin, CT, www.liberty-ind.com).
 - vii. Pipette, benchtop automatic (for preparation of urine dilutions to be analyzed): Like the Microlab 625 advanced dual syringe diluter (Hamilton, Reno, NV, <http://www.hamilton.com/>) equipped with a 10.0 mL left syringe, a 1.0 mL right syringe, a 12 gauge Concorde CT probe dispense tip, the Microlab cable management system and a foot pedal. Alternatives are acceptable, including the

Micromedic Digiflex™ (Titertek, Huntsville, AL, <http://www.titertek.com/>) equipped with 10.0-mL dispensing syringe, 2 mL sampling syringe, 0.75-mm tip, and foot pedal.

- viii. Pipettes (for preparation of intermediate stock working standards & other reagents): Like Brinkmann Research Pro Electronic pipettes (Brinkmann Instruments, Inc., Westbury, NY, <http://www.brinkmann.com/home/>). 5-100 µL (catalog #4860 000.070), 20-300 µL (catalog #4860 000.089), 50-1000 µL (catalog #4860 000.097), 100-5000 µL (catalog #4860 000.100). Note: pipette catalog numbers are without individual chargers. Can purchase individual chargers (pipette catalog numbers will differ) or a charging stand that will hold four pipettes (catalog #4860 000.860). When purchasing pipette tips (epTips), purchase one or more boxes, then “reloads” for those boxes after that: 5-100 µL (box catalog # 22 49 133-4, reload catalog # 22 49 153-9), 20-300 µL (box catalog # 22 49 134-2, reload catalog # 22 49 154-7), 50-1000 µL (box catalog # 22 49 135-1, reload catalog # 22 49 155-5), 100-5000 µL (box catalog # 22 49 138-5, reload catalog # 22 49 198-9, bulk bag catalog # 22 49 208-0). Equivalent pipettes and tips can be substituted.
- ix. Tubes for sample analysis (for autosampler): Like polypropylene 15-mL conical tubes, BD Falcon model #352097 (Becton Dickinson Labware, Franklin Lakes, NJ, www.bd.com), or equivalent. Clear plastics tend to have lowest trace metal contamination. Blue colored caps have also been used successfully for this method.
- x. Tubes for storage of intermediate working stock standards: Like polypropylene 50-mL conical tubes, BD Falcon model #352098 (Becton Dickinson Labware, Franklin Lakes, NJ, www.bd.com), or equivalent. For use in storage of intermediate working stock standards. Clear plastics tend to have lowest trace metal contamination. Blue colored caps have also been used successfully for this method.
- xi. Vortexer: Like MV-1 Mini Vortexer (VWR, West Chester, PA, www.vwr.com). Used for vortexing urine specimens before removing an aliquot for analysis. Equivalent item can be substituted.
- xii. Water purification system: Like NANOpure Diamond Ultrapure Water System (Barnstead International, Dubuque, Iowa, www.barnstead.com), or equivalent. For ultra-pure water used in reagent and dilution preparations.
- e. Sources of chemicals, gases, and regulators
- i. Acid, hydrochloric acid: Environmental Grade, 30-38% (GFS Chemicals Inc. Columbus, OH, www.gfschemicals.com). This is referred to as “concentrated” hydrochloric acid in this method write-up. For use in preparation of intermediate working stock standards. Equivalent acid products which meet or exceed the trace metals purity are acceptable substitutions.
- ii. Acid, nitric acid: Environmental Grade, 70% (GFS Chemicals Inc. Columbus, OH, www.gfschemicals.com). For use in diluent, rinse solution, intermediate working stock standards, and QC pool preparations. This is referred to as “concentrated” nitric acid in this method write-up. Equivalent acid products which meet or exceed the trace metals purity are acceptable substitutions.
- iii. Ethanol (EtOH): USP dehydrated 200 proof (Pharmco Products, Inc.) or equivalent.

- iv. Argon gas (for plasma & nebulizer) and Regulator: High purity argon (>99.999% purity, Specialty Gases Southeast, Atlanta, GA, www.sgsgas.com) for torch and nebulizer. Minimum tank source is a dewar of liquid argon (180-250L). Bulk tank (1500+L is preferred).
 1. Regulator for argon (at dewar): Stainless steel, single stage, specially cleaned regulator with 3000 psig max inlet, 0-200 outlet pressure range, CGA 580 cylinder connector, and needle valve shutoff on delivery side terminating in a ¼" Swagelok connector. Part number "KPRCGRF415A2/AG10-AR1" (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent. # Spares = 1.
 2. Regulator for argon (between bulk tank and PerkinElmer filter regulator): Single Stage 316SS Regulator, with 0-300 psi Inlet Gauge, 0-200 psi Outlet Gauge, Outlet Spring Range, 0-250 psi, ¼" Swagelok Inlet Connection, ¼ turn Shut off Valve on Outlet with ¼" Swagelok Connection and Teflon Seals. Part number KPR1GRF412A20000-AR1 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent. # Spares = 1.
 3. Regulator for argon (filter regulator on back of ICP-MS):
 - a. ELAN: Argon regulator filter kit. Catalog number N812-0508 (PerkinElmer, Shelton, CT, www.perkinelmer.com).
 - b. NexION: Argon regulator filter kit. Catalog number N8145023 (PerkinElmer, Shelton, CT, www.perkinelmer.com).
- v. Argon / hydrogen: Argon (90%) / hydrogen (10%) for DRC channel A. Initial purity of argon = 99.9997+% ("Research grade 5.7"). Initial purity of hydrogen = 99.9999+% ("Research Grade 6.0"). Mixture is typically purchased in cylinder size 35 (6"x24") (Airgas South, Atlanta, GA, www.airgas.com).
 1. Regulator for argon / hydrogen: Stainless steel, two stage, specially cleaned regulator with 3000 psig max inlet, 0-25 outlet pressure range, CGA 350 cylinder connector, and needle valve shutoff on delivery side terminating in a ¼" Swagelok connector. Like part number KCYADPF412A2AD10 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent. # Spares = 1.
 2. Flash arrestor (stainless steel): Like part # 6104A (Matheson Tri Gas, Montgomeryville, PA, www.mathesontrigas.com) or equivalent.
- vi. Disinfectant, for work surfaces: Diluted bleach (1 part bleach + 9 parts water), but must be re-made daily or equivalent disinfectant.
- vii. Oxygen: Oxygen ("Research Grade Research Grade 5.0", 99.9999% purity) for DRC channel B. Like part # OX R33A (Airgas South, Atlanta, GA, www.airgas.com).
 1. Regulator for oxygen: Stainless steel, two stage regulator for use with high purity oxygen (cleaned to be free of all oils). Maximum inlet pressure 3600-5000 psi. Inlet gauge pressure 0-5000 psi (no oil in gauge). Maximum delivery pressure 50-100 psi with a 0-30 psi outlet gauge (no oil in gauge). CGA 540 cylinder connector on inlet side and an angle pattern (90 degree) stainless steel needle valve on the delivery side terminating in a 1/8" stainless steel Swagelok connector. Like part # GEORG/KCYCFR/ORS2/540 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent.

2. Flash arrestor: Like part # 6104A (Matheson Tri Gas, Montgomeryville, PA, www.mathesontrigas.com), or equivalent. # Spares = 1.
- viii. Standard, dual detector: Like item # SM-2107-052 (High Purity Standards, Charleston, SC, <http://www.hps.net/>).
- ix. Standard, gallium: Like 1,000 mg/L, item # PLGA2-2Y. (SPEX Industries, Inc., Edison, NJ, www.spexcsp.com), or equivalent. Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- x. Standard, iridium: Like 1,000 mg/L iridium, item # PLIR3-2Y (SPEX Industries, Inc., Edison, NJ, www.spexcsp.com), or equivalent. Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- xi. Standard, multi-element stock calibration standard: Item numbers “SM-2107-037 solution A” and “SM-2107-037 solution B” (High Purity Standards, Charleston, SC, <http://www.hps.net/>). This is a set of custom mixes (see Table 3 in Appendix B for concentrations). Both are needed to cover all analytes of methods 3018 and 3018A. These solutions are diluted to prepare the intermediate stock working standards, which are in turn diluted to prepare the working calibrators. This solution can be prepared in-house from NIST traceable single element stock solutions if necessary.
- xii. Standard, rhodium: Like 1,000 mg/L, item # PLRH3-2Y. (SPEX Industries, Inc., Edison, NJ, www.spexcsp.com), or equivalent. Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- xiii. Standard, single element stock standards for preparation of urine quality control pools: National Institute of Standards and Technology (NIST) Standard Reference Materials (SRMs) 3103a (As), 3105a (Be), 3113 (Co), 3132 (Mn), 3134 (Mo), 3108 (Cd), 3102a (Sb), 3111a (Cs), 3104a (Ba), 3163 (W), 3128 (Pb), 3140 (Pt), 3161a (Sn), 3156a (Sr), 3158 (Tl), and 3164 (U) (National Institute of Standards and Technology (NIST), Office of Standard Reference Materials, Gaithersburg, MD, www.nist.gov). Other sources of standards can be used if they are NIST traceable.
- xiv. Triton X-100™ surfactant: Like “Baker Analyzed” TritonX-100™ (J.T. Baker Chemical Co., www.jtbaker.com), or equivalent.
- xv. Standard, gold: Like 10,000 µg/mL, cat # 10M21-2. (High Purity Standards, Charleston, SC, www.highpuritystandards.com), or equivalent. Used in diluent and rinse solution. Standard must be traceable to the National Institute for Standards and Technology.

6) Preparation of reagent and materials

a. Intermediate internal standard mixture

i. Purpose: Preparation of single intermediate solution containing internal standards will simplify the addition of the internal standards into the final diluent solution. This solution can be purchased rather than prepared.

ii. Preparation:

1. For DLS 3018:

To prepare 200 mL of 2% (v/v) HNO₃, 40 µg/L Ir and Rh intermediate internal standard solution:

- a. If not previously dedicated to this purpose, acid wash a 200 mL volumetric flask (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO₃ and ≥18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- b. Partially fill the 200 mL volumetric flask with ≥18 Mega-ohm·cm water (approximately 100-150 mL).
- c. Carefully add 4 mL of concentrated nitric acid. Mix into solution.
- d. Add 8,000 µg of rhodium (e.g. 8 mL of 1,000 µg/mL Rh stock standard).
- e. Add 8,000 µg of iridium (e.g. 8 mL of 1,000 µg/mL Ir stock standard).
- f. Fill to mark (200mL) and mix thoroughly.
- g. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

2. For DLS 3018A:

To prepare 200 mL of 2% (v/v) HNO₃, 40 µg/L Ga intermediate internal standard solution:

- a. If not previously dedicated to this purpose, acid wash a 200 mL volumetric flask (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO₃ and ≥18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- b. Partially fill the 200 mL volumetric flask with ≥18 Mega-ohm·cm water (approximately 100-150 mL).
- c. Carefully add 4 mL of concentrated nitric acid. Mix into solution.
- d. Add 8 µg of Ga (e.g. 8 mL of 1,000 µg/mL Ga standard)
- e. Fill to mark (200mL) and mix thoroughly.
- f. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

3. For emergency response combinations of 3018 and 3018A:

To prepare 200 mL of 2% (v/v) HNO₃, 40 µg/mL Ga, Rh and Ir intermediate internal standard solution:

NOTE: Gallium at 10 µg/L in the method diluent or 100 µg/L in a urine specimen (extremely unlikely) will produce approximately a +2 µg/L bias on observed Sr concentration due to the formation of ⁷¹Ga¹⁷O⁺, which occurs at the same m/z as ⁸⁸Sr⁺. While this interference is not large in comparison to typical Sr concentrations in urine, ~150 µg/L (see Table 11 in Appendix B), best urine Sr accuracy for biomonitoring purposes will be achieved when Ga is not present in the diluent as an internal standard. If As is combined with the multi-element

method for emergency response purposes, the resulting bias on Sr from the $^{71}\text{Ga}^{17}\text{O}^+$ interference is not likely to hinder the ability to identify acute Sr exposures.

- a. If not previously dedicated to this purpose, acid wash a 200 mL volumetric flask (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO_3 and ≥ 18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
 - b. Partially fill the 200 mL volumetric flask with ≥ 18 Mega-ohm·cm water (approximately 100-150 mL).
 - c. Carefully add 4 mL of concentrated nitric acid. Mix into solution.
 - d. Add 8 μg of Rh (e.g. 8 mL of 1,000 $\mu\text{g}/\text{mL}$ Rh standard).
 - e. Add 8 μg of Ir (e.g. 8 mL of 1,000 $\mu\text{g}/\text{mL}$ Ir standard).
 - f. Add 8 μg of Ga (e.g. 8 mL of 1,000 $\mu\text{g}/\text{mL}$ Ga standard).
 - g. Fill to mark (200mL) and mix thoroughly.
 - h. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.
- b. Intermediate Triton X-100™ solution (for DLS 3018 and DLS 3018A):
- i. Purpose: To avoid the time-consuming process of dissolving Triton X-100 on a daily basis for use in rinse solution, prepare an intermediate solution for daily use.
 - ii. Preparation: To prepare 2L of 2% Triton X-100™ in 5% (v/v) HNO_3 :
 1. If not previously dedicated to this purpose, acid wash a 2 L bottle (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO_3 and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
 2. Partially fill the bottle with ≥ 18 Mega-ohm·cm water (approximately 1-1.5 L).
 3. Add 40 mL of Triton X-100™ and stir until completely dissolved. Use a Teflon™ stir bar and stir plate if necessary (acid wash stir bar before use).
 4. Carefully add 100 mL of concentrated nitric acid.
 5. Fill to 2 L and stir thoroughly.
 6. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.
- c. Diluent and carrier
- i. Purpose: All samples (blanks, calibrators, QC, or patient samples) are combined with the diluent during the sample preparation step before analysis. This is where the internal standards are added which during the analysis will compensate for instrumental variations on the analyte signal. If using the FAST sample introduction system, the diluent is also used as the carrier solution.
 - ii. Preparation:

1. For DLS 3018

and emergency response combinations of methods 3018 and 3018A:

To prepare 4 L of an aqueous solution of 10 microgram/L internal standards and 500 µg/L gold in 2% (v/v) nitric acid:

- a. If not previously dedicated to this purpose, acid wash a 4 L container (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO₃ and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- b. Partially fill the 4 L container with ≥18 megaohm·cm water (~2/3 full).
- c. Carefully add 80 mL concentrated nitric acid and mix.
- d. Add spike of internal standard solution (to use other concentrations or volumes, adjust the volumes proportionally).
 - i. If for method 3018, add 1 mL of the 40 µg/mL Rh and Ir internal standard solution.
 - ii. If for emergency response combinations of methods 3018 and 3018A, add 1 mL of the 40 µg/mL Rh, Ir, and Ga internal standard solution.
- e. Add 200 µL of the 10,000 µg/mL gold standard.
- f. Make up to volume (4 L) with ≥18 megaohm·cm water.
- g. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

2. For DLS 3018A

To prepare a 10 µg/L Ga in 2% (v/v) nitric acid and 1.5% (v/v) ethanol:

- i. If not previously dedicated to this purpose, acid wash a 4 L container (PP, PMP, or Teflon™). For example, with 2% (v/v) HNO₃ and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- ii. Partially fill the 4 L container with ≥18 megaohm·cm water (~2/3 full).
- iii. Carefully add 80 mL concentrated nitric acid and mix.
- iv. Carefully add 60 mL dehydrated 200 proof ethanol and mix.
- v. Add 1 mL of the 40 µg/mL Ga internal standard solution. To use other concentrations, adjust the volume proportionally.
- vi. Make up to volume (4 L) with ≥18 megaohm·cm water.
- vii. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

d. ICP-MS rinse solution

i. Purpose: Pump this solution into the sample introduction system between samples to prevent carry-over of the analytes of interest from one sample measurement to the next.

ii. Preparation:

1. For DLS 3018
and emergency response combinations of methods 3018 and 3018A

To Prepare 4 L of 0.002% Triton X-100™, 5% (v/v) nitric acid solution and 500 µg/L gold:

- a. If not previously dedicated to this purpose, acid wash a 4 L container (PP, PMP, or Teflon™). For example, with 5% (v/v) HNO₃ and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- b. Partially fill the bottle with ≥18 Mega-ohm·cm water (approximately 2-3 L).
- c. Add 4 mL of the 2% Triton X-100™ / 5% (v/v) nitric-acid intermediate stock solution and mix well.
- d. Carefully add 200 mL of concentrated nitric acid and mix well.
- e. Add 200 µL of the 10,000 µg/mL gold.
- f. Fill to 4 L using ≥18 Megaohm·cm water.
- g. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

2. For DLS 3018A

To prepare 4 L of 0.002% Triton X-100™, 5% (v/v) HNO₃ solution and 1.5% (v/v) ethanol:

- a. If not previously dedicated to this purpose, acid wash a 4 L container (PP, PMP, or Teflon™). For example, with 5% (v/v) HNO₃ and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- b. Partially fill the bottle with ≥18 Mega-ohm·cm water (approximately 2-3 L).
- c. Add 4 mL of the 2% Triton X-100™ / 5% (v/v) nitric-acid intermediate stock solution and mix well.
- d. Carefully add 200 mL of concentrated nitric acid and mix well.
- e. Carefully add 60 mL dehydrated 200 proof ethanol and mix well.
- f. Fill to 4 L using ≥18 Megaohm·cm water.
- g. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

e. Standards, calibrators, and QC

i. Multi-element stock calibration standard

1. Purpose: These two master solutions will be diluted to prepare intermediate working calibrators.
2. Preparation: Prepared by external vendor or in-house from NIST-traceable standards.
 - a. Solution A: 5% HNO₃ solution containing As, Ba, Be, Cd, Cs, Co, Pb, Mn, Sr, Tl, and U. Concentrations are listed in Table 3 of Appendix B.

- b. Solution B: 5% HNO₃, 1% HF, 0.5% HCl solution containing Sb, Mo, Pt, Sn, and W. Concentrations are listed in Table 3 of Appendix B.
3. Storage: Store at room temperature and label appropriately. Expiration is determined by manufacturer or is 1 year after the container is opened (whichever comes first).

ii. Intermediate multi-element stock calibration standard

1. Purpose: The two stock standards are combined into a single intermediate stock calibration standard preparation.
2. Preparation: To prepare standards in 2% (v/v) nitric acid and 1% (v/v) hydrochloric acid which have final concentrations listed in Table 4 of Appendix B:
 - a. If not previously dedicated to this purpose, acid wash a 100 mL PP, PMP, or Teflon™ volumetric flask. For example, with a 2% (v/v) HNO₃ / 1% (v/v) HCl solution and >18 Mega-ohm·cm water (at least 3 times each) followed by verifying cleanliness through analysis of rinsate. Dedicate to purpose.
 - b. Partially fill the 100 mL volumetric flask with the 2% (v/v) nitric acid and 1% (v/v) hydrochloric acid prepared in Section 6.e.iii.b (50-75% full).
 - c. Using the volume listed in Table 4 of Appendix B, pipette the appropriate volume of the multi-element stock calibration standard solutions (both A and B) into the volumetric flask. Dilute to the volumetric mark with the 2% (v/v) nitric acid and 1% (v/v) hydrochloric acid using a pipette for the final drops. Mix each solution thoroughly. Final concentrations are listed in Table 4 of Appendix B.
 - d. Once mixed, transfer to an acid-cleaned, labeled, 50 mL container (PP, PMP, or Teflon™) for storage.
 - e. Label appropriately and store at room temperature. Expiration is 1 year from the date of preparation.

iii. Multi-element intermediate working calibration standards

1. Purpose: Six multi-element standards (S0 plus 5 spiked standards) used each day of analysis to prepare the final working calibrators.
2. Preparation: To prepare multi-element standards S0-S5 in 2% (v/v) HNO₃, 1% (v/v) HCl according to the volumes and concentrations listed in Table 5 of Appendix B:
 - a. Cleaning flasks: If not previously dedicated to this purpose, acid wash PP, PMP, or Teflon™ volumetric flasks. For example, with a 2% (v/v) HNO₃ / 1% (v/v) HCl solution and >18 Mega-ohm·cm water (at least 3 times each) followed by verifying cleanliness through analysis of rinsate. Dedicate to purpose.
 - b. 2% (v/v) HNO₃ & 1% (v/v) HCl diluent (S0) preparation: In a cleaned 2L volumetric flask, add 1-1.5L ≥18 Megaohm·cm water, 40 mL high purity concentrated HNO₃, and 20 mL high purity concentrated HCl. Fill to the mark and mix thoroughly. Use this diluent to fill the remaining volumetric flasks during preparation of the intermediate working calibration standards.
 - c. Dilutions & storage:

- i. Fill two acid-cleaned 50-mL containers (PP, PMP, or Teflon™) with the HNO₃ & HCl diluent. Label appropriately (including "S0") and store at room temperature. Expiration is 1 year from the date of preparation.
- ii. Partially fill the volumetric flasks with the HNO₃ & HCl diluent (50-75% full).
- iii. Using the volumes listed in Table 5 of Appendix B, pipette the appropriate volume of the multi-element stock calibration standard or the multi-element intermediate stock calibrator solutions (both A and B) into each of the volumetric flasks. Dilute each to the volumetric mark with the HNO₃ & HCl diluent using a pipette for the final drops. Mix each solution thoroughly. Final concentrations are listed in Table 5 of Appendix B.
- iv. Once mixed, transfer to acid-cleaned, labeled, 50-mL containers (PP, PMP, or Teflon™) for storage.
- v. Store at room temperature and label appropriately. Expiration is 1 year from date of preparation.

iv. Working multi-element calibrators

1. Purpose: The working multi-element calibrators will be analyzed in each run to provide a signal-to-concentration response curve for each analyte in the method. The concentration of an analyte in a patient urine sample dilution is determined by comparing the observed signal from the dilution of the patient urine sample to the response curve from the working multi-element calibrators.
2. Preparation: Make dilutions of the intermediate working calibration standards (S0-S5) immediately prior to analysis by combining with base urine (Section 6.e.v) and diluent (Section 6.c) using a benchtop automatic pipette. See Tables 8a and 8b of Appendix B and Section 8.b.ii for details of sample preparation. Expiration of capped dilutions is 3 days from preparation (see Appendix A, test for time between preparation and analysis).

v. Base urine

1. Purpose: This urine pool material will be mixed with the intermediate working calibrators just prior to analysis to matrix-match the calibration curve to the urine matrix of the unknown samples.
2. Contents: A mixture of multiple urine sources collected from anonymous donors are used to approximate an average urine matrix.
3. Preparation & storage:
 - a. Collect urine anonymously by placing screened containers and collection cups in the restrooms with a sign stating the reason the specimens are being collected, the name of the investigator to contact for additional information, and requesting that people provide a urine specimen (see supervisor regarding potential Institutional Review Board, IRB, requirements).
 - b. Once collected, analyze to ensure that concentrations of the analytes in this method are relatively low, so as to not interfere with the proper measurement of calibrators (see Table 2 in Appendix B for suggested maximum base urine concentrations).

- c. Once screened, mix the urine collections together in a larger container (polypropylene (PP), polymethylpentene (PMP), or Teflon™) which has been acid washed. For example, with 2% (v/v) HNO₃ and >18 Mega-ohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate. Add large Teflon™ stir bar and stir for 30+ minutes.
- d. Label appropriately and store long-term in smaller volume tubes (e.g. 50-mL acid-washed or lot screened polypropylene tubes) at ≤ -20 °C. Expiration date is 3 years from the date of preparation.

vi. Internal quality control materials (“Bench” QC)

1. Purpose: Internal (or “bench”) quality control (QC) materials are used to evaluate the accuracy and precision of the analysis process, and to determine if the analytical system is “in control” (is producing results that are acceptably accurate and precise). They are included in the beginning and at the end of each analytical run.
2. Content: The internal (or “bench”) quality control (QC) materials used in this method are pooled human urine, acidified to 1% (v/v) HNO₃, and spiked, if necessary, to reach a desired concentration. The analyte concentrations are in the low-normal concentration range (“low QC”) and high-normal concentration range (“high QC”).
3. Preparation & storage: Quality control materials can be either prepared by and purchased from an external laboratory or prepared within the CDC laboratories. Quality control must always be traceable to the National Institute for Standards and Technology (NIST). The CDC laboratory currently prepares its own bench QC materials using the following procedures:
 - a. Collection of urine: Collect urine anonymously by placing screened containers and / or collection cups in the restrooms with a sign stating the reason the specimens are being collected, the name of the investigator to contact for additional information, and requesting that people provide a urine specimen. Volume of urine to collect is dependent on the desired pool size. This write-up will assume a 10-L pool size for both the low and high bench QC.
 - b. Screening urine: Screen collected samples for metal content before mixing together to make separate pools that will be spiked to low, high, and elevated levels. Samples can be screened individually or after combining several together (reduces number of analyses).
 - i. Keep urine refrigerated whenever possible to minimize microbial growth.
 - ii. Because this is only a quick screen of the metal content, the number of replicates in the urine method can be reduced to one in order to reduce analysis time.
 - iii. Spike analyte concentrations for the low bench QC pool in the low-normal population range. Spike analyte concentrations for the high bench QC pool less than some preselected target concentration values in the high normal population range. See the National Report on Human Exposure to Environmental Chemicals for estimations of the normal population ranges for metals (<http://www.cdc.gov/exposurereport/>).

- c. Combining collected urine: Be attentive not to combine only diluted matrix urine samples into the low pool and only concentrated matrix urine samples into the pool for high and elevated QC. The goal is for combining samples is to approach an 'average' matrix for each pool.
 - i. Graduate four acid-washed 10-L carboys (PP or PMP) in 0.5 L increments (two will be used for decanting into).
 - ii. Combine collected urine samples into separate acid-washed 10-L carboys (PP or PMP), according to their concentrations, for the low high, and elevated bench QC pools.
 - iii. Mix each urine pool using large acid washed, Teflon™ coated stir bars and large stir plates. Keep urine refrigerated whenever possible.
 - iv. Acidify each urine pool to 1% (v/v) HNO₃ by adding the appropriate volume of concentrated HNO₃. Stir for 30+ min on large stir plates.
- d. Settling out of solids:
 - i. Refrigerate the urine (no stirring) for 1-3 days to allow for settling out of solids.
 - ii. For each urine pool, decant the urine into another of the acid-washed 10-L carboys to remove the urine from the solids settled out on the bottom of the carboy.
 - iii. Repeat steps (i) and (ii) until minimal solids are left at the bottom of the carboy after sitting overnight.
- e. Spiking of urine
 - i. Analyze a sample of each urine pool. Record these results for future recovery calculations.
 - ii. Use these results to determine target analyte concentrations possible for the pools
 - iii. Calculate the volume of single element standards needed to spike each pool to the desired concentrations.
 - iv. While stirring the pools on large stir plates, spike each pool with calculated volumes of single element standards (all spiking standards used must be traceable to NIST).
 - v. Continue to stir pools for 30+ minutes after spiking, then reanalyze.
 - vi. Repeat steps 4 and 5 until all analytes reach target concentrations keeping track of the total volume of spiking solution added to each urine pool.
- f. Dispensing and storage of urine
 - i. Container types: Dispense urine into smaller, lot screened containers (e.g. 2 mL cryovials). This allows for one vial of QC to be used in only a small number of analyses, reducing chances of contamination due to long-term use of the same container.
 - ii. Labels: Place labels on vials after dispensing and capping if the vials are originally bagged separately from the caps. This minimizes the chance for

contamination during the process. Include at least the name of QC pool (text and bar code), date of preparation, and a vial number on the labels.

- iii. Dispensing: Dispensing can be accomplished most easily using a benchtop automatic pipette in continuous cycling dispense mode. Complete this process in a clean environment (i.e., a class 100 cleanroom area or hood).
 1. Allow urine pool to reach room temperature before dispensing (to prevent temperature gradients possibly causing concentration gradients across the large number of vials being dispensed and to prevent condensation problems during labeling of vials).
 2. Replace the tubing attached to the dispensing syringe (left when looking at front of benchtop automatic pipette) with a length of clean Teflon™ tubing long enough to reach into the bottom of the 10L carboy while it is sitting on the stir plate.
 3. Check cleanliness of benchtop automatic pipette before use by analyzing 1-2% (v/v) HNO₃ which has been flushed through the benchtop automatic pipette with a portion of the same solution, which has not been through the benchtop automatic pipette.
 4. Approximately one hour before dispensing begins,
 - a. With the large stir plate close to the left side of the benchtop automatic pipette, begin stirring the urine pool to be dispensed.
 - b. Also during this time, flush the benchtop automatic pipette with urine from the pool to be dispensed. Place the ends of the tubing attached to both the sample and dispensing syringes into the carboy of urine so that urine won't be used up during this process. Be sure to secure both ends of tubing in the carboy with Parafilm so they will not come out during the flushing process.
 5. After dispensing the urine into the vials, cap the vials and label them. Placing labels on vials after capping minimizes the chance for contamination during the process.
- iv. Homogeneity testing: After dispensing, check homogeneity of analyte concentrations in pool aliquots by analysis of vials selected from across those dispensed. Seek guidance from a statistician regarding the number of vials needed for homogeneity analysis.
- v. Storage: Store long-term according the same storing and handling criteria described in Section 3.

f. Optimization solutions

i. DRC optimization (cell gas flow rate and RPq):

1. Purpose: For periodic testing of the DRC cell parameters. Procedure requires at a minimum a blank (i), an analyte solution (ii), a blank with interference (iii), and an analyte and interference containing solution (iv).
2. Content:

Diluent in this section refers to sample diluent (10 µg/L internal standards and 500 µg/L gold in 2% (v/v) nitric acid) described in Section 6c.

- a. Solutions for testing elimination of $^{98}\text{Mo}^{16}\text{O}$ interference on ^{114}Cd :
 - i. Base urine in diluent (1 + 9)
 - ii. Base urine in diluent (1 + 9) + 0.24 µg/L Cd
 - iii. Base urine in diluent (1 + 9) + 300 µg/L Mo
 - iv. Base urine in diluent (1 + 9) + 0.24 µg/L Cd + 300 µg/L Mo
 - b. Solutions for testing elimination of $^{39}\text{K}^{16}\text{O}$ interference on ^{55}Mn :
 - i. Base urine in diluent (1 + 9)
 - ii. Base urine in diluent (1 + 9) + 0.3 µg/L Mn
 - iii. Base urine in diluent (1 + 9) + 400 µg/L K
 - iv. Base urine in diluent (1 + 9) + 0.3 µg/L Mn + 400 µg/L K
3. **Preparation & storage:** Prepare different volumes by adding proportionally larger or smaller volumes of solution constituents. Interference concentrations can be prepared higher as needed by adjusting the volume of this spike. Keep interference spike volume small (<0.3 mL) using a high concentration stock solution (i.e. 1000 mg/L). Analyte concentrations can be made higher if needed for sensitivity reasons by preparing a higher concentration calibrator.
- a. Solutions for testing elimination of $^{98}\text{Mo}^{16}\text{O}$ interference on ^{114}Cd :
 - i. Base urine in diluent (1 + 9)
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8a (multiply volumes by 5).
 - ii. Base urine in diluent (1 + 9) + 0.24 µg/L Cd
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8a (multiply volumes by 5).
 - iii. Base urine in diluent (1 + 9) + 300 µg/L Mo
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8a (multiply volumes by 5).
 2. Add 0.015 mL of 1000 mg/L Mo.
 - iv. Base urine in diluent (1 + 9) + 0.24 µg/L Cd + 300 µg/L Mo
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8a (multiply volumes by 5).
 2. Add 0.015 mL of 1000 mg/L Mo.
 - b. Solutions for testing elimination of $^{39}\text{K}^{16}\text{O}$ interference on ^{55}Mn :
 - i. Base urine in diluent (1 + 9)

1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8a (multiply volumes by 5).
 - ii. Base urine in diluent (1 + 9) + 0.3 µg/L Mn
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8a (multiply volumes by 5).
 - iii. Base urine in diluent (1 + 9) + 400 µg/L K
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8a (multiply volumes by 5).
 2. Add 0.02 mL of 1000 mg/L K
 - iv. Base urine in diluent (1 + 9) + 0.3 µg/L Mn + 400 µg/L K
 1. In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8a (multiply volumes by 5).
 2. Add 0.02 mL of 1000 mg/L K.
- c. Store at room temperature and prepare as needed.
- d. Label appropriately, i.e. "Store at room temperature", preparation date, expiration date one year from preparation date, and preparer's initials.

ii. DRC optimization (axial field voltage):

1. Purpose: Use as necessary to verify the elimination of crosstalk phenomenon. The crosstalk phenomenon will present itself as increased sensitivity for an ion of interest in the presence of greater population of spectator ions in the DRC cell. Elimination is verified when the ratio of the intensity for iridium in standard 5 versus standard 0 is the same (1.00 ± 0.05).
2. Content: Working calibration standards 0 and 5 (see Section 6.e.iv).
3. Preparation and Storage: See section 6.e.iv.

iii. For dual detector calibration:

1. Purpose: Use as necessary to perform the dual detector calibration if any element exceeds 1,000,000 cps for calibration standard 5 (typically Sr).
2. Content: Dilutions of single element or special mix stock standards in 2% (v/v) HNO₃. Recommended elements include: As, Ba, Cs, Co, Pb, Mo, Sr, Sn. Other elements can be added as required for optimal instrument performance (esp. if measured intensities approach 500,000 cps in highest working calibrator).
3. Preparation & storage: To prepare elements of interest at 200 µg/L in 2% (v/v) HNO₃:

- a. Partially fill a 50 mL lot screened or acid-washed polypropylene tube with 2% (v/v) HNO₃,
- b. Add a 0.1 mL of 100 ug/mL special mix standard.
- c. Add 0.01 mL of any additional 1,000 ug/mL single element stock standard desired to be added.
- d. Dilute to the 50 mL mark with 2% (v/v) HNO₃.
- e. Label appropriately and store at room temperature. Expiration date is one year from preparation date.

7) Analytical instrumentation setup

(see Section 5 for details on hardware used, including sources)

a. Instrumentation & equipment setup:

i. Configuration for liquid handling

1. FAST valve setup: See Appendix B, Figure 1 for diagram and Section 5.b “FAST / ESI SC4-DX autosampler accessories” for source information.
 - a. Port 1: sample loop (white nut).
 - b. Port 2: 0.5 mm ID probe (red nut) for carrier solution.
 - c. Port 3: nebulizer line (green nut) for transfer of liquid to nebulizer.
 - d. Port 4: sample loop (white nut).
 - e. Port 5: 0.8 mm ID probe (blue nut) for diluted samples.
 - f. Port 6: vacuum line (black nut).
2. Carrier solution uptake: Use peristaltic pump to control uptake flow rate of carrier solution to the SC-FAST valve. Use of a ‘peristaltic to Teflon tubing adapter’ for prevents damage to small i.d. tubing when making connections (see consumables descriptions in Section 5.b).
3. Spray chamber waste removal

Use of a ‘peristaltic to Teflon tubing adapter’ for prevents damage to small i.d. tubing when making connections (see consumables descriptions in Section 5.b).

- a. Between spray chamber and peristaltic tubing:
 - i. Spray chambers with threaded connection: Use vendor-supplied threaded connector on base of chamber, connecting tubing directly to peristaltic pump tubing through a PEEK adapter or directly.
 - ii. Spray chambers without threaded connection: Use of specialized push-on connectors available from various vendors (like UFT-075 from Glass Expansion, Pocasset, MA) are preferred for safety reasons to direct connection of PVC tubing (e.g. 1/8” i.d. x 1/4” o.d.).
- b. Between peristaltic pump tubing and waste container: Connect 1/8” i.d. x 1/4” o.d. PVC tubing to the white / black peristaltic pump tubing using a tubing connector (PerkinElmer item # B3140715). Place the free end of the PVC tubing through the lid of the waste jug (be sure it is secure). Place the waste container in a deep secondary containment tray in case of overflow.

4. Rinse solution for autosampler:

- a. Rinse solution jug: Leave one of the caps on the top of the rinse jug loose to allow air venting into the jug as liquid is removed. Otherwise, the jug will collapse on itself as the liquid is removed and a vacuum is created inside. Use secondary containment tray.
- b. Rinse solution uptake to autosampler rinse station: Use tubing of different lengths and inner diameters between the rinse solution container and the autosampler rinse station to control uptake rate of rinse solution. These can be obtained from the autosampler manufacturer, their distributors, or custom built in the lab. Optimize these factors along with fill time in the software so that waste of rinse solution is minimized and rinse station does not go empty.
- c. Autosampler rinse station waste removal: Gravity drain of waste to the waste container is sufficient. Use minimum drain tubing to make this connection. If this tube is too long, the rinse station will not drain properly.

ii. Gas delivery and regulation

1. ICP-MS modifications:

- a. Plastic tubing between mass flow controllers and dynamic reaction cell have been replaced with stainless steel. Stainless steel tubing is preferred between the reaction gas cylinder / regulator and the back of the ICP-MS instrument.
- b. A second mass flow controller will be needed (channel B) that does not send the DRC gas through a 'getter'.

2. Argon gas: Used for various ICP-MS functions including plasma and nebulizer.

- a. Regulator for argon source (if a dewar): Set delivery pressure of this regulator at least 10 psi higher than the delivery pressure of the step-down regulator to allow for pressure drop across tubing that stretches to the instrument.
- b. Step down regulator (if source of argon is a bulk tank): Place this single stage regulator in the lab so that incoming argon pressure can be monitored and adjusted. Set delivery pressure to 10 psig above the delivery pressure of the filter regulator on the ICP-MS.
- c. Filter Regulator at ICP-MS: Single stage "argon regulator filter kit" supplied with the ICP-DRC-MS. Set the delivery pressure depending on the instrument setup:
 - i. ELAN with a 0-60psi gauge on the filter regulator: 52±1 psi when plasma is running (need 0-150 psi regulator if using a PolyPro or PFA nebulizer made by Elemental Scientific Inc).
 - ii. ELAN or NexION with a 0-150psi gauge on the filter regulator: 90-100 psi when plasma is running.

3. Argon (90%) / hydrogen (10%) gas mixture: Used for dynamic reaction cell interference removal from arsenic isotopes.

- a. Connect to DRC channel A.
- b. Set the delivery pressure of regulator to 5-7 psig when gas is flowing.
- c. Use a flash arrestor is on the outlet side of the regulator.

- d. This gas can be replaced by 100% argon. Argon from the ICP-MS main argon supply can be split off for this purpose by placing a tee in the delivery tubing on the low-pressure side of the filter regulator.
4. Oxygen (99.999±%) gas: Used for dynamic reaction cell interference removal from cadmium and manganese isotopes.
 - a. Connect to DRC channel B.
 - b. Set the delivery pressure of regulator to 5-7 psig when gas is flowing.
 - c. Use a flash arrestor is on the outlet side of the regulator.
- iii. Chiller / Heat Exchanger: If using refrigerated chiller, set temperature control to approximately 18 °C.
- b. Parameters for Instrument and Method: See Tables and Figures in Appendix B for a complete listing of the instrument and method parameters and software screen shots.

8) The run: quality, execution, evaluation, and reporting

a. Bench QC, reference materials and calibration verification:

- i. Bench "QC": Analysis of bench QC permits assessment of methodological imprecision, determination of whether the analytical system is 'in control' during the run, and assessment of time-associated trends. Before QC materials can be used in the QC process, they must be characterized by at least twenty (20) analytical runs to determine appropriate QC parameters.

Bench QC pool analyte concentrations in this method span the analyte concentration range of the calibrators including "low-normal" ('Low QC') and "high-normal" ('High QC') concentrations.

In each analytical run, the analyst will test each of the two bench QC samples two times, subjecting them to the complete analytical process. Bench QC pool samples are analyzed first in the run after the calibration standards but before any patient samples are analyzed. This permits making judgments on calibration linearity and blank levels prior to analysis of patient samples. The second analysis of the bench QC pools is done after analysis of all patient samples in the run (typically 20-30 patient samples total when analyzing for all elements in the method) to ensure analytical performance has not degraded across the time of the run. If more patient samples are analyzed on the same calibration curve after the second run of the bench QC, all bench QC must be reanalyzed before and after the additional samples. For example, the schemes shown in Table 6 in Appendix B are both acceptable ways to analyze multiple consecutive "runs".

- ii. Reference materials: Use standard reference material (SRM, e.g. SRM 2668 levels 1 and 2) from the National Institute of Standards and Technology (NIST) to verify method accuracy. Use previously characterized samples from proficiency testing program or commercially-produced reference materials when NIST SRMs are unavailable.
- iii. Calibration verification: The test system is calibrated as part of each analytical run with NIST-traceable calibration standards. These calibrators, along with the QCs and blanks, are used to verify that the test system is performing properly.

b. Perform, evaluate and report a run

i. Starting the equipment for a run

1. Power on the computer, printer, and autosampler, and instrument computer controller.
2. Peristaltic pump: Set proper tension on peristaltic pump tubing.
3. Software: Start software for the ICP-MS and autosampler control.
4. Daily pre-ignition maintenance checks: Perform and document daily maintenance checks (e.g., Ar supply pressure, interface components cleanliness and positioning, interface pump oil condition, vacuum pressure, etc.).
5. Place probe in adequate volume of carrier or rinse solution: If using an ESI FAST, manually place carrier probe into carrier solution. If not, send the autosampler probe to a rinse solution (e.g. autosampler rinse station).
6. Start the plasma
7. Start the peristaltic pump: Start the pump running slowly, making sure that the rotational direction is correct for the way the tubing is set up.
8. Warm-up time: Allow warm-up time suggested by the manufacturer for the ICP-MS (e.g. RF generator) after igniting the plasma. There will be another warm-up time (or “stability time”) for the DRC later in this procedure.
9. Daily performance check: Perform and document a daily performance check and any optimizations necessary.

Save new parameters to the “default.tun” and “default.dac” files.

10. DRC Stability time: Best analyte-to-internal standard ratio stability is typically observed after 1-1.5 hours of analysis of urine samples using the DRC mode method (~12 measurements of the 15-element panel, or 50 measurements of the total arsenic method can be made in 1 hour). Prepare 50mL⁺ of a calibration standard (e.g. standard 2) to be analyzed repeatedly before the beginning of the run to achieve a stable analyte-to-internal standard ratio. Time to reach stability is instrument-specific and learned from performance of runs. See Table 7 in Appendix B for example of setup in the Samples / Batch window and Tables 8a and 8b in Appendix B for details of making a working standard.
11. Readying the instrument for quick-start analysis: Leave the plasma running to eliminate the need for an initial instrument warm-up period and / or a DRC stabilization period as long as appropriate planning is made for sufficient solution supply and waste collection. Analysis of conditioning samples (diluted urine matrix) can also be scheduled to occur at roughly a predetermined time. Accomplish this by setting up multiple sample analyses with extended rinse times (e.g. one 15-element analysis with a 1400s rinse time will take approximately 30 minutes to complete). Initial samples would be non-matrix, while final samples would be diluted matrix for conditioning. If running a DRC-only method during these scheduled analyses, the ICP-MS will remain in DRC-mode for approximately 45 minutes without depressurizing the cell.

12. Software setup for Analysis:

- a. Workspace (files & folders): Verify & set up the correct files and data directories for your analysis (See Table 1 in Appendix B for defaults).
- b. Samples / Batch Window: Update the software to reflect the current sample set. Use a bar code scanner to input data whenever possible. See Table 1 in Appendix B for times and speeds.

1. Urine vs. Aqueous Method Files:

- a. The difference: There are two method files for this one method (see Table 1 in Appendix B). It is necessary to use both to accomplish each run because the current PerkinElmer software will not allow for more than one blank per method file. The ONLY DIFFERENCE between these two files is on the Sampling tab where one lists the autosampler positions of the urine blank and urine calibrators (the “urblk” method file) and the other lists the autosampler position of the aqueous blank (the “aqblk” method file).
 - b. Use: The ONLY TIME when it matters which of these files is used is when the measurement action *includes* “Run blank” or “Run standards”. When the measurement action is only ‘run sample’, it does not matter whether the “urblk” or “aqblk” method file is used. Analysts typically follow the pattern below, however, for the sake of consistency and as a reminder of which blank must be used for which type of sample. See Table 7 in Appendix B.
 - i. The “urblk” method file: Use to analyze the initial urine blank (blank for the calibration curve), the urine calibrators, and the urine blank checks at the very beginning of the run. The urine blank method defines the autosampler location of the urine blank and the urine calibration standards.
 - ii. The “aqblk” method file must be used to analyze all QC materials and patient samples. The aqueous blank method defines the aqueous blank in autosampler location.
- ii. Preparation of samples for analysis (See Tables 8a and 8b in Appendix B)
1. Thaw urine samples; allow them to reach ambient temperature.
 2. Prepare the following solutions into pre-labeled containers benchtop automatic pipette or other volumetric sample transfer device. See Table 8a or Table 8b in Appendix B for a summary.

Prepare samples in the cleanest environment available to prevent trace element contamination and an environment, which provides personnel protection (e.g. Class II, Type A/B3 biological safety cabinet).

- a. Aqueous Blank: Prepare at least two aqueous blanks. One will be the actual reagent blank for patient and QC samples and the other will be a backup (“Aqueous Blank Check”) in case the original aqueous blank is unusable.

- b. *Calibrators*: Prepare the working calibration standards (S0-S5). Prepare S0 in triplicate. One of these S0 preparations will be the zero standard (urine blank) for the calibration standards; the other two will be analyzed twice after the last calibrator to collect run blank data that can be used in periodically evaluating the method LOD.
- c. *Patient & QC Samples*: Before taking an aliquot for analysis, homogenize the sample.

After preparation, mixed and cover diluted samples. Place prepared dilutions on the autosampler of the ICP-MS in the order corresponding to the sequence setup in the ICP-MS software.

Room temperature storage of original samples for the workday is acceptable.

- iii. Start the analysis using the ICP-MS software.
- iv. Monitor the analysis in real-time as much as possible. If necessary, leave the run to complete itself unattended as long as appropriate planning is made for either overnight operation or Auto Stop (see below).

Monitor the analysis for the following:

1. *Verify proper operation of the instrument* (proper loop filling, sample reaching nebulizer in correct timing, autosampler arm moving properly, etc . . .)
2. *Verify that background signal from instrument and reagents are low*. Helpful checks when diagnosing high background problems include:
 - a. Water to be used in Aq Blank Checks and dilutions.
 - b. Diluent before and after being flushed through the benchtop automatic pipette.
If contamination is observed from the pipette, flush the pipette with ≥ 500 mL of nitric acid solution ($\leq 5\%$ (v/v) HNO_3) and retest.
 - c. Comparison with other instruments.

3. *Verify analyte / internal standard ratio stability*

The net intensity (analyte / internal standard ratio) of the measurements made while stabilizing the UCT can be evaluated to determine the readiness of the system to begin analysis. Continual trending in this ratio indicates that unwanted instrument drift will occur within the run.

4. *Evaluate the Axial Field Voltage (AFV) optimized value*

Monitor the change between S0 and S5 for measured intensities (cps) of the internal standard iridium in DRC mode. If the percent difference between the iridium intensities is greater than 5% (especially if greater in S5), then run the axial field voltage optimization. See Section 6.f.ii for preparation of optimization solutions.

5. *Verify calibration curves* meet R^2 requirements (minimum of 0.98, typically 0.99 to 1.000).
6. *Verify bench QC results are within the acceptable limits*.

If an analyte result for the beginning QC material(s) falls outside of the $\pm 3\text{SD}$ limits, then the following steps are recommended:

- a. Evaluate the blank results.
- b. Evaluate the reproducibility of the 3 replicates within the measurements.
- c. Evaluate the consistency of the internal standard across the measurements (esp. the calibrators).
- d. Evaluate calibration curves. If a particular calibration standard is obviously in error, it can be re-analyzed as a sample (old or new dilution) and incorporated into the curve through data reprocessing as a calibrator. As a last resort, a single calibration point per analyte between and including S2 and S4 can be removed from the curve. Follow up on repeated problems with calibration standards with appropriate corrective actions (e.g. re-preparation of intermediate working standards or troubleshooting instrument parameters).
- e. Prepare a fresh dilution of the failing QC material (same vial) and reanalyze it to see if the QC dilution was not properly made.
- f. Prepare a fresh dilution of the failing QC material (unused vial) and analyze it to see if the QC vial had become compromised.
- g. Prepare and analyze new working calibrators.
- h. Test a different preparation of intermediate working calibration standards.

If these steps do not result in correction of the out-of-control values for QC materials, consult the supervisor for other appropriate corrective actions.

7. *Verify good precision among replicates* of each measurement.
8. *Verify consistent measured intensities of the internal standards.*

Some sample-to-sample variations are to be expected, however, intensities drifting continuously in one direction resulting in failing results for ending QC indicate the instrument needs additional pre-conditioning before the run or environmental conditions are changing too much around the instrument.

9. *Verify elevated patient results.*

Refer to Figure 4 in Appendix B for flowchart.

- a. Confirming an elevated concentration: Repeat for confirmation any sample having a concentration greater than the 1UB threshold. See Table 8 in Appendix B.
- b. Dilution of a sample to within the calibration range: Repeat in duplicate with extra dilution any sample having a concentration greater than the highest calibration standard to bring the observed result within the concentration range of the calibrators. See Table 7 in Appendix B for validated extra dilutions.
- c. Confirming proper washout after an elevated sample: When monitoring the analysis in real-time, if an observed sample concentration is greater than the highest concentration tested for washout (see Table 9 in Appendix B), do the following to verify that the run is still in control for low concentration samples before proceeding with analysis.
 - i. Stop run following elevated sample

Verify that the run is still in control for lower concentration samples before proceeding with analysis. Analyze 2 urine blank checks followed by a low bench QC washout check. If the low bench QC wash check is not in control (within $\pm 3SD$ limits), repeat these 3 check samples until washout is verified before proceeding with analysis.

Example:

3018 UrBlkChk Wash1

3018 UrBlkChk Wash2

LUXXXXX Wash

- ii. If the run is not verified in-control for low concentration samples before the next samples are analyzed, see Section 8.b.vii.2. for directions.
- v. Overnight operation or using auto stop: Ensure sufficient solution supply and waste collection during unattended operation. Turn on the AutoStop feature of the ICP-MS software. Delay the shutdown at least 10 minutes (use peristaltic pump speed approximately that of the method wash) to rinse the sample introduction system of urine matrix before turning off the plasma. It will be necessary to replace the sample peristaltic pump tubing the next day since it will have been clamped shut overnight. Enable "Auto Start/Stop" is on the "AutoStop" tab of the Instrument window.
 1. ELAN specifics: Enable "Auto Start / Stop" is on the "AutoStop" tab of the Instrument window.
 2. NexION specifics: Enable AutoStop in on the Run List window. Select "Batch Completed" for Stop Criteria.
- vi. Records of results: Run results will be documented after each run in both electronic and paper form.
 1. Electronic records: Transfer data electronically to the laboratory information system. When keyboard entry must be used, proofread transcribed data after entry.
 - a. Export data from the ICP-MS software using "original conditions" or files and folders used during the analysis. Use descriptive report filenames (e.g. 2005-0714a_group55.txt). In the ELAN or NexION software under "Report Format" (METHOD window, REPORT tab) choose the "Use Separator" option, and under the "File Write" Section choose "Append."
 - b. Move the generated .TXT data file to the appropriate subdirectory on the network drive where exported data are stored prior to import to the laboratory information management system.
 - c. Import the instrument file into the laboratory information system with appropriate documentation (e.g. instrument ID, analyst, calibration standards lot number, and run or sample specific comments).
 2. Paper records: Run sheets must be documented with
 - i. Analyst initials
 - ii. Instrument ID
 - iii. Date of analysis and run # for the day
- vii. Analyst evaluation of run results:

1. Bench quality control: After completing a run, and importing the results into the laboratory information system, evaluate the run bench QC according to laboratory QC rules [12]. The QC limits are based on the average and standard deviation of the beginning and ending analyses of each of the bench QC pools, so it will not be possible to know if the run is in control until statistically reviewed.
 - a. Rules for bench quality control evaluation: The following are the CDC DLS QC rules for two QC pools per run with two or more QC results per pool.
 - i. If both QC run means are within $2S_m$ limits and individual results are within $2S_i$ limits, then accept the run.
 - ii. If one of the two QC run means is outside a $2S_m$ limit - reject run if:
 1. Extreme Outlier – Run mean is beyond the characterization mean $\pm 4S_m$
 2. 3S Rule - Run mean is outside a $3S_m$ limit
 3. 2S Rule – Two or more of the run means are outside the same $2S_m$ limit
 4. 10 X-bar Rule – Current and previous 9 run means are on same side of the characterization mean
 - iii. If one of the 4 QC individual results is outside a $2S_i$ limit - reject run if:
 1. Extreme Outlier – One individual result is beyond the characterization mean $\pm 4S_i$
 2. R 4S Rule – Within-run ranges for all pools in the same run exceed $4S_w$ (i.e., 95% range limit)
Note: Since runs have multiple results per pool for 2 pools, the R 4S rule is applied within runs only.

Abbreviations:

S_i = Standard deviation of individual results.

S_m = Standard deviation of the run means.

S_w = Within-run standard deviation.

- b. Implications of QC failures: If the DLS SAS program declares the run “out of control” for any analyte, use the following to determine the implications on usability of the data from the run.
 - i. 13 – 15 elements in the run
 1. 1, 2 or 3 analytes “out of control”: ONLY the analytes which were “out of control” are invalid for reporting from the run.
 2. 4 or more analytes “out of control”: All results, regardless of analyte, are invalid for reporting from the run.
 - ii. 4 – 12 elements in the run
 1. 1 or 2 analytes “out of control”: ONLY the analytes which were “out of control” are invalid for reporting from the run.

2. 3 or more analytes “out of control”: All results, regardless of analyte, are invalid for reporting from the run.

iii. 3 elements in the run

1. 1 analyte “out of control”: ONLY the analyte which was “out of control” is invalid for reporting from the run.

2. 2 or more analytes “out of control”: All results, regardless of analyte, are invalid for reporting from the run.

iv. 1 – 2 elements in the run

ONLY the analyte which was “out of control” is invalid for reporting from the run.

2. Patient results:

a. Elevated concentrations: Refer to Figure 5 in Appendix B for flowchart.

i. Boundaries requiring confirmatory measurement:

1. Results greater than the first (1UB) or second (2UB) upper boundaries.

The concentrations assigned to 1UB and 2UB for an element is determined by study protocol but default concentrations are in Table 9 in Appendix B.

a. Results greater than the first upper boundary (1UB): Confirm by repeat analysis of a new sample preparation concentrations observed greater than the “first upper boundary” (defined in the laboratory database as the “1UB”). Report the first analytically valid result, as long as the confirmation is within 10%. Continue repeat analysis until a concentration can be confirmed.

b. Analyst reporting of elevated results: Report any patient results confirmed to be greater than the second upper boundary (2UB) as an “elevated result”.

2. Results greater than highest calibrator: Samples that exceed the high calibrator must be prepared with minimum extra dilution in duplicate to bring the observed result within the calibration range ($\leq S5$). Report the first analytically valid result (i.e. the first one within the calibration range), as long as the confirmation is within 10%. Continue repeat analysis until a concentration can be confirmed.

ii. Concentrations requiring verification of washout: Following a result greater than the highest concentrations tested for washout (see Table 9 of Appendix B) do the following:

1. If the run was determined to be in-control for low concentration samples before the next samples were analyzed, no further action is required.

2. If the run was not determined to be in-control for low concentration samples before the next samples were analyzed confirm by re-analysis the results for the 2 samples immediately following the elevated

sample. Report the results if they confirm the initial results within $\pm 10\%$ or $\pm 3SD$ of the low bench QC, whichever is greater.

- b. Unacceptable reproducibility: If the range of the three replicate readings (maximum replicate concentration value - minimum replicate concentration value) for a single sample analysis is greater than the range maximum criteria listed in Table 9 in Appendix B **and** the range of the three replicate readings is greater than 10% of the observed concentration, do not use the measurement for reporting. Repeat the analysis of the sample.
- viii. Submitting final work for review: All analyses must undergo quality control and quality assurance review. After appropriately documenting the run in the laboratory information system (e.g. sample and run QC, and run and sample comments), inform the first level reviewer of the completed work and submit any printed documentation.

9) Routine equipment maintenance and data backups

Maintenance activities will be documented in the instrument logbook.

- a. Equipment maintenance: Analysts are expected to regularly evaluate the need for, and when necessary perform, cleaning, replacement, or re-positioning of components in ICP-MS the sample introduction system, interface, ion optics region, and equipment required resources (e.g. autosampler, exhaust, compressed gases, and coolant). Frequency of equipment maintenance will be dependent on instrument throughput.
 - i. Parameter optimizations: Analysts are expected to optimize instrument parameters.
 - ii. Dual detector calibration: Perform dual detector calibration regularly (weekly or monthly) for any element exceeding 1,000,000 cps for calibration standard 5. This is typically only Sr. The dual detector calibration solution is described in Section 6.f.iii.
- b. DRC optimizations: DRC conditions (cell gas flow rate and RPq value) can be verified by analyzing the DRC optimization solutions (see Section 6.f.i) as needed to ensure proper reduction of potential ICP-MS interferences.
- c. Data backup: Data on the instrument computer will be backed up via two backup routines. Files used and produced by the ICP-MS in analyzing samples will be backed up and kept a minimum of three years after analysis.
 - i. Daily backups to secondary hard drive: Program automatic backups of the relevant computer files to occur each night onto a secondary hard drive to prevent loss of data from failure of primary hard drive.
 - ii. Weekly backup: Backup relevant computer files weekly either to secondary hard drive which is remote to the laboratory or to removable media which will be placed remote to the laboratory for retrieval in the case of catastrophic data loss elsewhere.

10) Reporting thresholds

- a. Reportable range: Urine multi-element values are reportable in the range between the method LOD and the highest calibrator. Above the highest calibrator, extra dilutions are made of the urine sample to bring the concentration within the reportable range. If extra dilution has been necessary, the reported value will exceed the upper end of the reportable range.

- b. Reference ranges (normal values): In this method the 95% reference ranges (see Appendix, Tables 10 and 11 in Appendix B) for these elements in urine fall within the range of the calibrators.
- c. Action levels: Due to the uncertainty of the health implications of elevated concentrations of many of the elements determined with this method, there is no routine notification for elevated levels of every analyte determined with this method. The present NRC standard for workplace removal is 15 $\mu\text{g/L}$ of U in urine [14]. Other action levels for reporting to supervising physicians are determined on a study-by-study basis.

11) Method calculations

- a. Method limit of detection (LOD): The method detection limits for elements in urine specimens are defined as 3 times s_0 , where s_0 is the estimate of the standard deviation at zero analyte concentration. S_0 is taken as the y-intercept of a linear or 2nd order polynomial regression of standard deviation versus concentration (4 concentration levels of the analytes in urine each measured 60 times across at least a 2-month timeframe). Method LODs are re-evaluated periodically.
- b. Method limit of quantitation (LOQ): The Division of Laboratory Sciences does not currently utilize limits of quantitation in regards to reporting limits [12].
- c. QC limits: Quality control limits are calculated based on concentration results obtained in at least 20 separate runs. It is preferable to perform separate analyses on separate days and using multiple calibrator lot numbers, instruments, and analysts to best mimic real-life variability. The statistical calculations are performed using the SAS program developed for the Division of Laboratory Sciences (DLS_QC_compute_char_stats.sas).

12) Alternate methods for performing test and storing specimens if test system fails:


If the analytical system fails, setup analysis on other ICP-MS instrument, if available. If no other instrument is available, store the specimens at ≤ -20 °C until the analytical system can be restored to functionality.

13) Summary Statistics and QC Graphs

See following pages


Summary Statistics and QC Chart for Urinary arsenic, total (µg/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 110 | 16JAN15 | 08FEB17 | 10.667 | 0.165 | 1.5 |
| 655 | 110 | 16JAN15 | 08FEB17 | 47.192 | 0.689 | 1.5 |




Summary Statistics and QC Chart for Antimony, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 186 | 16JAN15 | 14MAR17 | 0.1600 | 0.0092 | 5.7 |
| 655 | 186 | 16JAN15 | 14MAR17 | 0.5991 | 0.0174 | 2.9 |




Summary Statistics and QC Chart for Barium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|-------|--------------------|--------------------------|
| 654 | 190 | 16JAN15 | 14MAR17 | 1.234 | 0.044 | 3.6 |
| 655 | 190 | 16JAN15 | 14MAR17 | 6.430 | 0.165 | 2.6 |




Summary Statistics and QC Chart for Cadmium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 182 | 16JAN15 | 14MAR17 | 0.1481 | 0.0101 | 6.8 |
| 655 | 182 | 16JAN15 | 14MAR17 | 1.5199 | 0.0535 | 3.5 |




Summary Statistics and QC Chart for Cesium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|---------|--------------------|--------------------------|
| 654 | 195 | 16JAN15 | 14MAR17 | 4.4282 | 0.1020 | 2.3 |
| 655 | 195 | 16JAN15 | 14MAR17 | 12.2242 | 0.2947 | 2.4 |




Summary Statistics and QC Chart for Cobalt, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 187 | 16JAN15 | 14MAR17 | 0.4816 | 0.0329 | 6.8 |
| 655 | 187 | 16JAN15 | 14MAR17 | 2.4142 | 0.0580 | 2.4 |




Summary Statistics and QC Chart for Lead, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|-------|--------------------|--------------------------|
| 654 | 182 | 16JAN15 | 20MAR17 | 0.389 | 0.016 | 4.2 |
| 655 | 182 | 16JAN15 | 20MAR17 | 2.790 | 0.051 | 1.8 |




Summary Statistics and QC Chart for Manganese, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|-------|--------------------|--------------------------|
| 654 | 192 | 16JAN15 | 14MAR17 | 0.806 | 0.101 | 12.5 |
| 655 | 192 | 16JAN15 | 14MAR17 | 2.952 | 0.172 | 5.8 |




Summary Statistics and QC Chart for Molybdenum, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|---------|--------------------|--------------------------|
| 654 | 196 | 16JAN15 | 20MAR17 | 39.142 | 0.697 | 1.8 |
| 655 | 196 | 16JAN15 | 20MAR17 | 131.394 | 2.224 | 1.7 |




Summary Statistics and QC Chart for Strontium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|---------|--------------------|--------------------------|
| 654 | 194 | 16JAN15 | 14MAR17 | 81.384 | 1.980 | 2.4 |
| 655 | 194 | 16JAN15 | 14MAR17 | 167.399 | 4.185 | 2.5 |




Summary Statistics and QC Chart for Thallium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 183 | 16JAN15 | 14MAR17 | 0.1708 | 0.0081 | 4.8 |
| 655 | 183 | 16JAN15 | 14MAR17 | 0.4603 | 0.0128 | 2.8 |




Summary Statistics and QC Chart for Tin, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 183 | 16JAN15 | 14MAR17 | 2.587 | 0.093 | 3.6 |
| 655 | 183 | 16JAN15 | 14MAR17 | 15.255 | 0.513 | 3.4 |




Summary Statistics and QC Chart for Tungsten, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|--------|--------------------|--------------------------|
| 654 | 185 | 16JAN15 | 14MAR17 | 0.1545 | 0.0109 | 7.1 |
| 655 | 185 | 16JAN15 | 14MAR17 | 0.9245 | 0.0221 | 2.4 |



Summary Statistics and QC Chart for Uranium, urine (ug/L)

| Lot | N | Start Date | End Date | Mean | Standard Deviation | Coefficient of Variation |
|-----|-----|------------|----------|---------|--------------------|--------------------------|
| 654 | 185 | 16JAN15 | 14MAR17 | 0.01772 | 0.00094 | 5.3 |
| 655 | 185 | 16JAN15 | 14MAR17 | 0.11695 | 0.00475 | 4.1 |



References

1. Thomas, R., *Practical Guide to ICP-MS (Practical Spectroscopy)*. 2003, New York, New York: Marcel Dekker. 336.
2. Tanner, S.D., Baranov, Vladimir I, *Theory, Design, and Operation of a Dynamic Reaction Cell for ICP-MS*. Atomic Spectroscopy, 1999. **20**(2): p. 45-52.
3. Tanner, S.D., V.I. Baranov, and D.R. Bandura, *Reaction cells and collision cells for ICP-MS: a tutorial review*. Spectrochimica Acta Part B-Atomic Spectroscopy, 2002. **57**(9): p. 1361-1452.
4. PerkinElmer SCIEX Instruments, *ELAN DRC II Hardware Guide*. 2001, Canada.
5. PerkinElmer Sciex, *Service Manual for NexION 300 ICP-MS System*. E9932001A ed. 2010, USA.
6. Bandura, D.R., V.I. Baranov, and S.D. Tanner, *Inductively coupled plasma mass spectrometer with axial field in a quadrupole reaction cell*. Journal of the American Society for Mass Spectrometry, 2002. **13**(10): p. 1176-1185.
7. Mulligan, K.J., T.M. Davidson, and J.A. Caruso, *Feasibility Of The Direct Analysis Of Urine By Inductively Coupled Argon Plasma Mass-Spectrometry For Biological Monitoring Of Exposure To Metals*. Journal Of Analytical Atomic Spectrometry, 1990. **5**(4): p. 301-306.
8. Jarrett, J.M., et al., *Total Urine Arsenic Biomonitoring Using Inductively Coupled Plasma Mass Spectrometry with a Dynamic Reaction Cell*. Atomic Spectroscopy, 2007. **28**(4): p. 113-122.
9. Jarrett, J.M., et al., *Eliminating molybdenum oxide interference in urine cadmium biomonitoring using ICP-DRC-MS*. 2008. **23**: p. 962-967.
10. Larsen, E.H. and S. Sturup, *Carbon-enhanced Inductively Coupled Plasma Mass Spectrometric Detection of Arsenic and Selenium and Its Application to Arsenic Speciation*. Journal Of Analytical Atomic Spectrometry, 1994. **9**: p. 1101-1105.
11. Amarasiriwardena, C.J., et al., *Determination of the total arsenic concentration in human urine by inductively coupled plasma mass spectrometry: a comparison of the accuracy of three analytical methods*. Analyst, 1998. **123**(3): p. 441-445.
12. Division of Laboratory Sciences, *Policies and Procedures Manual*. 2012, Division of Laboratory Sciences (DLS), National Center for Environmental Health, Centers for Disease Control and Prevention, Public Health Service, Department of Health and Human ServicesCenters for Disease Control and Prevention, .
13. Heitland, P. and H.D. Koster, *Biomonitoring of 37 trace elements in blood samples from inhabitants of northern Germany by ICP-MS*. Journal of Trace Elements in Medicine and Biology, 2006. **20**(4): p. 253-262.
14. U.S. Nuclear Regulatory Commission, *Regulatory guide 8.22 (revision 1). Bioassay at uranium mills*. 1988: Atlanta, GA.
15. Centers for Disease Control and Prevention, *Third National Report on Human Exposure to Environmental Chemicals*, <http://www.cdc.gov/exposurereport>. 2005.
16. *Report on Human Biomonitoring of Environmental Chemicals in Canada*. 2010, Health Canada: Ottawa.
17. Heitland, P. and H.D. Koster, *Biomonitoring of 30 trace elements in urine of children and adults by ICP-MS*. Clinica Chimica Acta, 2006. **365**(1-2): p. 310-318.
18. Paschal, D.C., et al., *Trace metals in urine of United States residents: Reference range concentrations*. Environmental Research, 1998. **76**(1): p. 53-59.

19. Agency for Toxic Substances and Disease Registry (ATSDR). 2000. Toxicological profile for Manganese. Atlanta, G.U.S.D.o.H.a.H.S., Public Health Service. , *Toxicological Profile for Manganese*, ATSDR, Editor. 2000. p. 15.
20. Moreno, M.E., et al., *Biomonitoring of metal in children living in a mine tailings zone in Southern Mexico: A pilot study*. International Journal of Hygiene and Environmental Health, 2010. **213**(4): p. 252-258.
21. Gil, F., et al., *Biomonitorization of cadmium, chromium, manganese, nickel and lead in whole blood, urine, axillary hair and saliva in an occupationally exposed population*. Science of the Total Environment, 2011. **409**(6): p. 1172-1180.
22. Wang, H.M., et al., *Urinary heavy metal levels and relevant factors among people exposed to e-waste dismantling*. Environment International, 2011. **37**(1): p. 80-85.
23. Juliao, L., et al., *Exposure of workers in a mineral processing industry in Brazil*. Radiation Protection Dosimetry, 2007. **125**(1-4): p. 513-515.

Appendix A. ruggedness testing results

Parameter Test #1 (15 element panel): Evaluate the impact on analysis results if the RF Power is varied.

Test details:

1. Three different RF power settings were tested in separately prepared, consecutive runs on the instrument without turning off the plasma. At least 15 minutes stabilization time was allowed between each run after the RF power was changed. "Junk urine" samples (20) were analyzed between the beginning and ending QC of each run. All other method parameters were kept per method.
2. Run order was as follows: run #1 (1450W), run #2 (1150W), run #3 (1600W).

| Parameter Test 1 Results (Table 1 of 3 for 15 element panel). All concentrations in µg/L. Test performed 3/24/2010 by Denise Tevis using ELAN DRC-2N. | | | | | |
|--|--|-----------------------|-----------------------|-----------------------|-----------------------|
| ID | RF power tested | Ba | Be | Cd | Co |
| LU-04310_UMP3_e | <i>characterized mean (±2SD range)</i> | 0.76 (0.65 - 0.86) | 0.69 (0.59 - 0.79) | 0.32 (0.28 - 0.36) | 0.42 (0.38 - 0.47) |
| | 1150W | 0.67 | 0.62 | 0.28 | 0.39 |
| | 1450W | 0.70 | 0.66 | 0.31 | 0.32 |
| | 1600W | 0.75 | 0.75 | 0.33 | 0.42 |
| HU-04311_UMP3_e | <i>characterized mean (±2SD range)</i> | 5.01 (4.54 - 5.24) | 5.28 (4.48 - 6.07) | 1.62 (1.47 - 1.78) | 1.88 (1.66 - 2.09) |
| | 1150W | 4.93 | 5.75 | 1.6 | 1.96 |
| | 1450W | 5.55 | 5.28 | 1.75 | 1.67 |
| | 1600W | 4.85 | 5.82 | 1.58 | 1.90 |

Appendix A. Ruggedness testing results. (continued)

| Parameter Test 1 Results (Table 2 of 3 for 15 element panel). All concentrations in µg/L. Test performed 3/23/2011 by Denise Tevis using ELAN DRC-2N. | | | | | |
|---|--|-----------------------|-----------------------|-----------------------|--------------------------|
| ID | RF power tested | Cs | Mo | Pb | Pt |
| LU-04310_UMP_e | <i>characterized mean (±2SD Range)</i> | 2.38 (2.25 - 2.51) | 19.3 (18.6 - 20.0) | 0.42 (0.37 - 0.48) | 0.10 (0.07 - 0.13) |
| | 1150W | 2.13 ‡ | 17.2 * | 0.42 | 0.09 |
| | 1450W | 2.32 | 18.9 | 0.41 | 0.09 |
| | 1600W | 2.42 | 18.8 | 0.43 | 0.10 |
| HU-04311_UMP_e | <i>characterized mean (±2SD Range)</i> | 9.82 (9.03 - 10.6) | 136 (131 - 142) | 2.95 (2.82 - 3.08) | 0.85 (0.71 - 1.00) |
| | 1150W | 9.62 | 136 | 3.03 | 0.93 |
| | 1450W | 10.56 | 133 | 2.89 | 1.02 |
| | 1600W | 9.55 | 135 | 3.05 | 1.08 # |
| ID | RF power tested | Sb | Tl | W | U |
| LU-04310_UMP_e | <i>characterized mean (±2SD Range)</i> | 0.19 (0.17 - 0.21) | 0.18 (0.17 - 0.19) | 0.22 (0.19 - 0.24) | 0.014 (0.011 - 0.016) |
| | 1150W | 0.21 | 0.16 | 0.19 | 0.017 ^α |
| | 1450W | 0.16 | 0.19 | 0.22 | 0.014 |
| | 1600W | 0.20 | 0.18 | 0.22 | 0.017 ^α |
| HU-04311_UMP_e | <i>characterized mean (±2SD Range)</i> | 0.66 (0.60 - 0.71) | 0.58 (0.55 - 0.61) | 0.94 (0.90 - 0.99) | 0.128 (0.115 - 0.141) |
| | 1150W | 0.69 | 0.59 | 0.90 | 0.153 ^α |
| | 1450W | 0.61 | 0.57 | 0.93 | 0.126 |
| | 1600W | 0.66 | 0.60 | 0.91 | 0.150 ^α |

‡Low in one pool only and was within expected precision of method from default setting result.

Appendix A. Ruggedness testing results. (continued)

| Parameter test 1 results (Table 3 of 3 for 15 element panel). All concentrations in µg/L. Test performed 3/23/2011 by Denise Tevis using ELAN DRC-2N. | | | | |
|--|--|-----------------------|-----------------------|-------------------|
| ID | RF power tested | Mn | Sn | Sr |
| NYDOH UE09-05‡ | <i>characterized mean (±2SD Range)</i> | 1.37 (1.55 -1.19) | 2.2 (2.0-2.8) | |
| | 1150W | 1.22 | 2.9 | |
| | 1450W | 0.98 | 2.8 | |
| | 1600W | 1.30 | 3.0 | |
| NYDOH UE09-06‡ | <i>characterized mean (±2SD Range)</i> | 31.1 (26.3 -35.9) | 61 (55.0 - 67.0) | |
| | 1150W | 29.4 | 67.4 | |
| | 1450W | 23.8 | 68.5 | |
| | 1600W | 30.0 | 66.7 | |
| Seronom Trace Elements Urine§ | <i>characterized mean (±2SD Range)</i> | 12.3 (10.9 - 13.7) | 54.6 (51.9 - 57.3) | 110 (104 -116) |
| | 1150W | 10.4 | 62.3 | 113 |
| | 1450W | 8.46 | 62.3 | 111 |
| | 1600W | 10.5 | 61.4 | 111 |

§Purchased from Sero AS, Billingstad, Norway.

‡ Purchased from Wadsworth Center, New York State Department of Health

Conclusion: Results are not compromised by changes in RF power within the range of 1150W to 1600W.

Appendix A. Ruggedness testing results. (continued)

Parameter test #1 (Arsenic): Evaluate the impact on analysis results if the set RF Power is varied.

Test details:

1. Three different RF power settings were tested in separately prepared, consecutive runs on the instrument without turning off the plasma. At least 15 minutes stabilization time was allowed between each run after the RF power was changed. “Junk urine” samples (40) were analyzed between the beginning and ending QC of each run. All other method parameters were kept per method.
2. Run order was as follows: run #1 (1450W), run #2 (1150W), run #3 (1600W).

| Parameter test 1 results (arsenic). All concentrations in µg/L. Test performed 3/26/10 by Graylin Mitchell using ELAN DRC-2G. | | |
|--|--|------------------------------------|
| ID | RF power tested | As |
| LU-04310_UMP_e | characterized mean characterized 2SD range characterized 3SD range | 3.74 3.21 – 4.27 2.95 – 4.53 |
| | 1150W | 3.72 |
| | 1450W | 4.10 |
| | 1600W | 3.66 |
| HU-04311_UMP_e | characterized mean characterized 2SD range characterized 2SD range | 55.8 53.3 – 58.3 52.1 – 59.6 |
| | 1150W | 55.6 |
| | 1450W | 58.8 |
| | 1600W | 52.2 |

Conclusion: Results are not compromised by changes in RF power within the range of 1150W to 1600W.

Appendix A. Ruggedness testing results. (continued)

Parameter test #2 (cadmium and manganese): Evaluate the impact on analysis results if the cell gas flow rate is increased or decreased by 20% for the analytical run.

Test details:

1. Three different cell gas flow rates were tested in separately prepared, consecutive runs on the instrument without turning off the plasma. Samples were prepared with diluent containing the internal standards. At least 15 minutes stabilization time was allowed between each run after the cell gas flow rate was changed. "Junk urine" samples (20) were analyzed between the beginning and ending QC of each run
2. Run #1 (method default = 2.3 mL/min O₂)
3. Run #2 (decreased cell gas flow rate by 20% to 1.84 mL/min O₂).
4. Run #3 (increased cell gas flow rate by 20% to 2.76 mL/min O₂).

| Parameter test 2 results (Table 1 of 2), 15 element, DRC mode (Cd, Mn, and Ir). All concentrations in µg/L. Results collected 12/27/2012 on NexION C. | | | | | | |
|--|-------------------------------------|----------------------------------|-------------|-------------------------------|-------------|-------------|
| ID | cell gas flow rate | Cd | | Mn | | Ir |
| LU10709 | NexION mean (+/- 2SD) | 0.150 (0.114 – 0.186) N=80 | %RSD | 0.79 (0.56 – 1.02) N=80 | %RSD | %RSD |
| | 1.84 mL/min O ₂ | 0.151 | 11.7 | 0.781 | 13.3 | 0.6 |
| | 2.3 mL/min O ₂ (default) | 0.146 | 4.3 | 0.731 | 0.3 | 0.9 |
| | 2.76 mL/min O ₂ | 0.129 | 14.6 | 0.673 | 4.1 | 0.6 |
| HU10710 | NexION mean (+/- 2SD) | 1.52 (1.40 – 1.64) N=80 | %RSD | 2.95 (2.43 – 3.47) N=80 | %RSD | %RSD |
| | 1.84 mL/min O ₂ | 1.47 | 0.6 | 2.89 | 3.5 | 0.9 |
| | 2.3 mL/min O ₂ (default) | 1.51 | 2.5 | 2.76 | 4.0 | 0.7 |
| | 2.76 mL/min O ₂ | 1.39 | 5.7 | 2.48 | 5.5 | 1.0 |
| SRM 2668 Level 2 | NexION mean (+/- 2SD) | 16.0 (14.6 – 17.5) N=48 | %RSD | 51.3 (44.0 – 58.7) N=48 | %RSD | %RSD |
| | target (95% C.I.) | 16.4 (15.3 – 17.5) | - | 47.6 (44.2 – 51.0) | - | - |
| | 1.84 mL/min O ₂ | 16.1 | 1.1 | 49.3 | 1.8 | 0.7 |
| | 2.3 mL/min O ₂ (default) | 16.2 | 0.7 | 48.5 | 1.3 | 0.4 |
| | 2.76 mL/min O ₂ | 14.7 | 1.4 | 43.3 | 1.4 | 1.1 |

Conclusion: The accuracy and intra-measurement precision of Cd and Mn results are not affected by changes in the cell gas flow rate, within the tested range (1.84-2.76 mL/min).

Appendix A. Ruggedness testing results. (continued)

Parameter test #2 (arsenic): Evaluate the impact on analysis results if the cell gas flow rate is increased or decreased by 20% for the analytical run.

Test details:

1. Three different cell gas flow rates were tested in separately prepared, consecutive runs on the instrument without turning off the plasma. At least 15 minutes stabilization time was allowed between each run after the cell gas flow rate was changed. “Junk urine” samples (40) were analyzed between the beginning and ending QC of each run.
2. Run #1 (method default = 0.70 mL/min).
3. Run #2 (decreased cell gas flow rate by 20% to 0.56 mL/min).
4. Run #3 (increased cell gas flow rate by 20% to 0.84 mL/min).

| Parameter Test 2 results, arsenic, DRC mode. All concentrations in µg/L. Results collected on 1/3/2013 NexION C (Ar gas). | | | |
|--|---------------------------|----------------------------|-------------|
| ID | cell gas flow rate | As | |
| LU10709 | NexION mean (+/- 2SD) | 11.0 (9.92 – 12.1) N=74 | %RSD |
| | 0.56 mL/min Ar | 11.0 | 0.4 |
| | 0.70 mL/min Ar (default) | 10.8 | 1.1 |
| | 0.84 mL/min Ar | 10.9 | 0.5 |
| HU10710 | NexION mean (+/- 2SD) | 48.6 (43.5 – 53.6) N=74 | %RSD |
| | 0.56 mL/min Ar | 48.2 | 0.4 |
| | 0.70 mL/min Ar (default) | 47.2 | 0.1 |
| | 0.84 mL/min Ar | 47.9 | 0.2 |
| SRM 2668 Level 2 | NexION mean (+/- 2SD) | 217 (205 – 229) N=34 | %RSD |
| | Target (95% C.I.) | 213 (209 – 218) | - |
| | 0.56 mL/min Ar | 223 | 1.2 |
| | 0.70 mL/min Ar (default) | 219 | 1.7 |
| | 0.84 mL/min Ar | 221 | 0.9 |

Conclusion: Neither accuracy nor within-measurement precision of As results are compromised by changes in cell gas flow rate within the range tested (0.56 – 0.84 mL/min).

Appendix A. Ruggedness testing results. (continued)

Parameter test #3 (DRC elements: cadmium and manganese): Evaluate the impact on analysis results if the RPq is increased or decreased by 20% for the analytical run.

Test details:

1. Three RPq settings were tested for cadmium and manganese in separately prepared, consecutive runs on the instrument without turning off the plasma. Samples were prepared with diluent containing the internal standards. At least 15 minutes stabilization time was allowed between each run after DRC RPq was changed. “Junk urine” samples (20) were analyzed between the beginning and ending QC of each run.
2. Run #1 (instrument default DRC RPq: 0.75).
3. Run #2 (decreased; DRC RPq: 0.65).
4. Run #3 (increased; DRC RPq: 0.85).
5. Run #4 (increased; DRC RPq: 0.80).

| Parameter test 3 results (Table 1 of 2), 15 element, DRC Mode (Cd, Mn, and Ir). All concentrations in µg/L. Results collected 1/3/2013 on NexION D. | | | | | | |
|---|-----------------------|----------------------------------|------|-------------------------------|--------|------|
| ID | DRC RPq | Cd | | Mn | | Ir |
| LU10709 | NexION mean (+/- 2SD) | 0.150 (0.114 – 0.186) N=80 | %RSD | 0.79 (0.56 – 1.02) N=80 | %RSD | %RSD |
| | 0.65 | 0.148 | 8.8 | 0.744 | 0.8 | 0.3 |
| | 0.75 (default) | 0.155 | 10.8 | 0.745 | 5.5 | 0.4 |
| | 0.80 | 0.151 | 5.8 | 1.02 | 6.9 | 0.6 |
| | 0.85 | 0.137 | 3.8 | 0.813 | 27.4 * | 0.7 |
| HU10710 | NexION mean (+/- 2SD) | 1.52 (1.40 – 1.64) N=80 | %RSD | 2.95 (2.43 – 3.47) N=80 | %RSD | %RSD |
| | 0.65 | 1.51 | 3.8 | 2.83 | 1.6 | 0.3 |
| | 0.75 (default) | 1.63 | 1.6 | 3.09 | 3.4 | 0.4 |
| | 0.80 | 1.50 | 3.6 | 3.03 | 2.1 | 0.6 |
| | 0.85 | 1.43 | 1.9 | 2.31 | 19.4 * | 1.1 |
| SRM 2668 Level 2 | NexION mean (+/- 2SD) | 16.0 (14.6 – 17.5) N=48 | %RSD | 51.3 (44.0 – 58.7) N=48 | %RSD | %RSD |
| | target (95% C.I.) | 16.4 (15.3 – 17.5) | - | 47.6 (44.2 – 51.0) | - | - |
| | 0.65 | 16.1 | 0.6 | 50.4 | 1.0 | 0.9 |
| | 0.75 (default) | 15.9 | 0.7 | 50.0 | 1.6 | 0.7 |
| | 0.80 | 15.8 | 0.9 | 50.0 | 0.5 | 0.7 |
| | 0.85 | 15.4 | 2.0 | 46.1 | 0.6 | 0.8 |

Conclusion: Neither accuracy nor within-measurement precision of Cd and Mn results are compromised by changes in RPq within the range of 0.65 – 0.8. However, setting the RPq >0.8 causes problems in with-in measurement precision.

Appendix A. Ruggedness testing results. (continued)

Parameter test #3 (arsenic): Evaluate the impact on analysis results if the RPq is increased or decreased by 20% for the analytical run.

Test details:

1. Three different RPQ settings were tested for Cadmium in separately prepared, consecutive runs on the instrument without turning off the plasma. At least 15 minutes stabilization time was allowed between each run after DRC RPQ was changed. “Junk urine” samples (40) were analyzed between the beginning and ending QC of each run.
2. Run #1 (instrument default DRC RPq: 0.60).
3. Run #2 (decreased; DRC RPq: 0.48).
4. Run #3 (increased; DRC RPq: 0.72).

| Parameter Test 3 Results, As, DRC Mode. All concentrations in µg/L. Results collected 1/4/2013 on NexION C (Ar gas). | | | |
|---|-----------------------|----------------------------|-------------|
| ID | DRC RPq | As | |
| LU10709 | NexION Mean (+/- 2SD) | 11.0 (9.92 – 12.1) N=74 | %RSD |
| | 0.48 | 10.8 | 0.3 |
| | 0.60 (default) | 11.3 | 0.4 |
| | 0.72 | 10.5 | 0.6 |
| HU10710 | NexION Mean (+/- 2SD) | 48.6 (43.5 – 53.6) N=74 | %RSD |
| | 0.48 | 47.1 | 0.3 |
| | 0.60 (default) | 49.4 | 0.8 |
| | 0.72 | 46.8 | 0.5 |
| SRM 2668 Level 2 | NexION Mean (+/- 2SD) | 217 (205 – 229) N=34 | %RSD |
| | Target (95% C.I.) | 213 (209 – 218) | - |
| | 0.48 | 213 | 0.5 |
| | 0.60 (default) | 224 | 0.4 |
| | 0.72 | 210 | 0.3 |

Conclusion: Neither accuracy nor within-measurement precision of Cd and Mn results are compromised by changes in RPq within the range of 0.48 – 0.72.

Appendix A. ruggedness testing results. (continued)

Parameter test #4 (15 element panel): Method descriptions and SOP assume preparation and analysis on same day. Evaluate the impact on analysis results if the analytical run is prepared to analyze but circumstances do not allow for analysis to occur until 24 or 48 hours later.

Test details:

1. Three separate run sets (A, B, and C) were prepared at one sitting from the same starting materials. Set 'A' was analyzed immediately per the assumption of the method. Set's 'B' and 'C' were stored at room temperature for 24 and 48 hours, respectively before analysis. "Junk urine samples (20) were analyzed between the beginning and ending QC of each run, making each a normal length run. All other method parameters were kept per method.
2. On day two, a fresh run set ("D") was prepared and analyzed immediately for comparison to results from set "B" (Run 2 of the day. Results not shown).
3. On day three, another fresh run set ("E") was prepared and analyzed immediately for comparison to results from set "C" (Run 2 of the day. Results not shown).

| Parameter test 4 results (Table 1 of 2 for 15 element panel). All concentrations in µg/L. Test begun 3/23/11 by Denise Tevis using ELAN DRC-2N. | | | | | |
|--|---|-----------------------|-----------------------|-----------------------|-----------------------|
| ID | time from preparation to analysis | Ba | Be | Cd | Co |
| LU-04310_UMP3_e | <i>characterized mean</i> (±2SD range) | 0.76 (0.65 - 0.86) | 0.69 (0.59 - 0.79) | 0.32 (0.28 - 0.36) | 0.42 (0.38 - 0.47) |
| | freshly prepared | 0.70 | 0.66 | 0.32 | 0.32 |
| | 24 hours | 0.70 | 0.69 | 0.31 | 0.42 |
| | 48 hours | 0.86 | 0.67 | 0.31 | 0.43 |
| HU-04311_UMP3_e | <i>characterized mean</i> (±2SD range) | 5.01 (4.54 - 5.24) | 5.28 (4.48 - 6.07) | 1.62 (1.47 - 1.78) | 1.88 (1.66 - 2.09) |
| | freshly prepared | 5.55 | 3.48 | 1.75 | 1.67 |
| | 24 hours | 4.66 | 5.76 | 1.57 | 1.94 |
| | 48 hours | 5.12 | 5.49 | 1.67 | 1.84 |

Appendix A. Ruggedness test #4 results (continued)

| Parameter test 4 results (Table 2 of 3 for 15 element panel). All concentrations in µg/L. Test begun 3/23/11 by Denise Tevis using ELAN DRC-2N. | | | | | |
|--|--|---------------|---------------|---------------|-----------------|
| ID | time from preparation to analysis | Cs | Mo | Pb | Pt |
| LU-04310_UMP_e | <i>characterized mean</i> | 2.38 | 19.3 | 0.42 | 0.10 |
| | <i>(±2SD range)</i> | (2.25 - 2.51) | (18.6 - 20.0) | (0.37 - 0.48) | (0.07 - 0.13) |
| | freshly prepared | 2.32 | 18.9 | 0.41 | 0.09 |
| | 24 hours | 2.35 | 19.2 | 0.44 | 0.11 |
| | 48 hours | 2.30 | 19.0 | 0.46 | 0.10 |
| HU-04311_UMP_e | <i>characterized mean</i> | 9.82 | 136 | 2.95 | 0.85 |
| | <i>(±2SD range)</i> | (9.03 - 10.6) | (131 - 142) | (2.82 - 3.08) | (0.71 - 1.00) |
| | freshly prepared | 10.6 | 133 | 2.89 | 1.02 |
| | 24 hours | 9.36 | 134 | 3.08 | 1.03 |
| | 48 hours | 10.0 | 132 | 3.04 | 1.12 * |
| ID | time from preparation to analysis | Sb | TI | W | U |
| LU-04310_UMP_e | <i>characterized mean</i> | 0.19 | 0.18 | 0.22 | 0.014 |
| | <i>(±2SD range)</i> | (0.17 - 0.21) | (0.17 - 0.19) | (0.19 - 0.24) | (0.011 - 0.016) |
| | freshly prepared | 0.16 | 0.19 | 0.22 | 0.014 |
| | 24 hours | 0.19 | 0.18 | 0.21 | 0.013 |
| | 48 hours | 0.19 | 0.19 | 0.22 | 0.014 |
| HU-04311_UMP_e | <i>characterized mean</i> | 0.61 | 0.58 | 0.94 | 0.128 |
| | <i>(±2SD range)</i> | (0.60 - 0.71) | (0.55 - 0.61) | (0.90 - 0.99) | (0.115 - 0.141) |
| | freshly prepared | 0.61 | 0.57 | 0.93 | 0.126 |
| | 24 hours | 0.65 | 0.60 | 0.90 | 0.128 |
| | 48 hours | 0.69 | 0.59 | 0.90 | 0.126 |

* Results within expected precision of the method of result at fresh preparation.

Appendix A. Ruggedness test #4 results (continued)

| Parameter test 4 results (Table 3 of 3 for 15 element panel). All concentrations in µg/L. Test begun 3/23/11 by Denise Tevis using ELAN DRC-2N. | | | | |
|--|--|----------------------|---------------------|------------|
| ID | time from preparation to analysis | Mn | Sn | Sr |
| NYDOH UE09-05‡ | <i>characterized mean</i> <i>(±2SD range)</i> | 1.37 (1.55 -1.19) | 2.2 (2.0-2.8) | |
| | freshly prepared | 0.98 | 2.8 | |
| | 24 hours | 1.26 | 2.6 | |
| | 48 hours | 1.47 | 2.6 | |
| NYDOH UE09-06‡ | <i>characterized mean</i> <i>(±2SD range)</i> | 31.1 (26.3 -35.9) | 61 (55.0 - 67.0) | |
| | freshly prepared | 23.8 | 68.5 | |
| | 24 hours | 30.6 | 62.8 | |
| | 48 hours | 31.9 | 61.4 | |
| Seronom Trace Elements Urine§ | <i>characterized mean</i> | 12.3 | 54.6 | 110 |
| | <i>(±2SD range)</i> | (10.9 - 13.7) | (51.9 - 57.3) | (104 -116) |
| | freshly prepared | 8.47 | 62.3 | 111 |
| | 24 hours | 10.9 | 57.5 * | 112 |
| | 48 hours | 10.4 | 58.3 * | 114 |

§Purchased from Sero AS, Billingstad, Norway.

‡ Purchased from Wadsworth Center, New York State Department of Health

* Results within expected precision of the method of result at default setting

Conclusions: Results from all times tested were within acceptable boundaries. Therefore, after preparation, capped working samples can be used up to 48 hours after preparation without compromising observed results.

Appendix A. ruggedness test #4 results (continued)

Parameter test #4 (arsenic): Method descriptions and SOP assume preparation and analysis on same day. Evaluate the impact on analysis results if the analytical run is prepared to analyze but circumstances do not allow for analysis to occur until 24 or 48 hours later.

Test details:

1. Three separate run sets (A, B, and C) were prepared at one sitting from the same starting materials. Set 'A' was analyzed immediately per the assumption of the method. Set's 'B' and 'C' were stored at room temperature for 24 and 48 hours, respectively before analysis. "Junk urine samples (20) were analyzed between the beginning and ending QC of each run, making each a normal length run. All other method parameters were kept per method.
2. On day two, a fresh run set ("D") was prepared and analyzed immediately for comparison to results from set "B" (Run 2 of the day. Results not shown).
3. On day three, another fresh run set ("E") was prepared and analyzed immediately for comparison to results from set "C" (Run 2 of the day. Results not shown).

| Parameter test 4 results (arsenic). | | |
|--|--|------------------|
| Test performed 5/26-28/10 by Graylin Mitchell using ELAN DRC2-G. | | |
| ID | time from preparation to analysis | As (µg/L) |
| LU-04310_UMP_e | <i>characterized mean</i> | 3.74 |
| | <i>characterized 2SD range</i> | 3.21 – 4.27 |
| | <i>characterized 3SD range</i> | 2.95 – 4.53 |
| | fresh preparation | 3.40 |
| | after 24 hours | 3.37 |
| | after 48 hours | 3.41 |
| HU-04311_UMP_e | <i>characterized mean</i> | 55.8 |
| | <i>characterized 2SD range</i> | 53.3 – 58.3 |
| | fresh preparation | 54.7 |
| | after 24 hours | 53.8 |
| | after 48 hours | 53.5 |

Conclusions: Results from all times tested were within acceptable boundaries. Therefore, after preparation, capped working samples can be used up to 48 hours after preparation without compromising observed results.

Appendix A. ruggedness test #5 results

Parameter test #5: Evaluate the impact on observed concentration if an extra dilution is performed on the sample relative to the calibration standards.

Test details: A large urine sample spiked to elevated concentrations was prepared for analysis 5 times at various extra dilution levels. *Thallium was re-tested Nov-Dec 2015 using 5 different urine samples (QC and reference materials) to evaluate different urine matrices.

| | Ba | Be | Cd | Co | Cs | Mn |
|---------------------|--|-------------|-------------|-------------|-------------|-------------|
| initial conc (µg/L) | 297 | 98 | 94 | 149 | 252 | 96 |
| | observed concentrations normalized to 'no extra dilution' result | | | | | |
| no extra dilution | 1.00 ± 0.02 | 1.00 ± 0.04 | 1.00 ± 0.01 | 1.00 ± 0.02 | 1.00 ± 0.02 | 1.00 ± 0.03 |
| 2x dilution | 1.00 ± 0.00 | 1.00 ± 0.01 | 1.04 ± 0.00 | 0.99 ± 0.01 | 0.99 ± 0.01 | 1.02 ± 0.01 |
| 5x dilution | 0.99 ± 0.00 | 1.01 ± 0.01 | 1.09 ± 0.00 | 0.99 ± 0.00 | 0.98 ± 0.00 | 1.03 ± 0.01 |
| 10x dilution | 0.98 ± 0.00 | 1.00 ± 0.01 | 1.09 ± 0.00 | 0.98 ± 0.00 | 0.97 ± 0.00 | 1.03 ± 0.00 |
| 20x dilution | 0.97 ± 0.00 | 1.01 ± 0.00 | 1.10 ± 0.00 | 0.98 ± 0.00 | 0.97 ± 0.00 | 1.05 ± 0.00 |

| | Mo | Pb | Pt | Sb | Sn | Sr |
|---------------------|--|-------------|-------------|-------------|-------------|-------------|
| initial conc (µg/L) | 1885 | 243 | 137 | 98 | 325 | 1816 |
| | observed concentrations normalized to 'no extra dilution' result | | | | | |
| no extra dilution | 1.00 ± 0.02 | 1.00 ± 0.01 | 1.00 ± 0.01 | 1.00 ± 0.02 | 1.00 ± 0.05 | 1.00 ± 0.03 |
| 2x dilution | 0.97 ± 0.00 | 1.01 ± 0.00 | 1.00 ± 0.01 | 1.00 ± 0.00 | 1.00 ± 0.02 | 0.92 ± 0.01 |
| 5x dilution | 0.94 ± 0.00 | 1.05 ± 0.00 | 0.99 ± 0.00 | 1.01 ± 0.00 | 1.00 ± 0.01 | 0.90 ± 0.00 |
| 10x dilution | 0.92 ± 0.00 | 1.06 ± 0.00 | 0.94 ± 0.00 | 1.00 ± 0.00 | 0.97 ± 0.01 | 0.88 ± 0.00 |
| 20x dilution | 0.90 ± 0.00 | 1.08 ± 0.00 | 0.89 ± 0.00 | 0.97 ± 0.00 | 0.95 ± 0.00 | 0.87 ± 0.00 |

| | TI* | U | W | As |
|---------------------|--|-------------|-------------|-------------|
| initial conc (µg/L) | 25-115 | 43 | 99 | 2963 |
| | observed concentrations normalized to 'no extra dilution' result | | | |
| no extra dilution | 1.00 | 1.00 ± 0.06 | 1.00 ± 0.01 | 1.00 ± 0.02 |
| 2x dilution | 1.01 ± 0.01 | 1.00 ± 0.03 | 0.98 ± 0.01 | 1.00 ± 0.02 |
| 5x dilution | 1.04 ± 0.01 | 1.01 ± 0.01 | 0.96 ± 0.00 | 1.00 ± 0.02 |
| 10x dilution | 1.07 ± 0.01 | 1.01 ± 0.01 | 0.96 ± 0.00 | 0.99 ± 0.01 |
| 20x dilution | 1.09 ± 0.05 | 1.02 ± 0.00 | 0.95 ± 0.00 | 0.99 ± 0.02 |

Conclusion: Do not use greater than a 2x extra dilution for Sr or a 10x extra dilution for Cd, Mo, Pt, and TI. All others (As with ethanol in reagents, Ba, Be, Co, Cs, Mn, Pb, Sb, Sn, U, and W) can be analyzed at up to a 20x extra dilution without significant effect (> ± 10% error) to the observed concentration.

Appendix B

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|---|--|
| Instrument: PerkinElmer ELAN DRC II ICP-MS or NexION 300D | |
| Autosampler: ESI SC4 autosampler with FAST sample introduction system | |
| Optimization (conditions) window parameters | |
| RF power | 1450W (ELAN), 1600W (NexION) |
| plasma gas flow (Ar) | 15 L/min |
| auxiliary gas flow (Ar) | 1.2 L/min |
| nebulizer gas flow (Ar) | ~0.90 – 1.0 L/min (optimized as needed for sensitivity) |
| ion lens voltage(s) | AutoLens (optimized as needed for sensitivity) |
| AFV, QRO, CRO, CPV, discriminator threshold | Optimized per instrument by service engineer, or advanced user. |
| Parameters of x-y alignment, nebulizer gas flow, AutoLens voltages, mass calibration, dual detector calibration, and detector voltages are optimized regularly. Optimization file name = default.dac. | |
| Configurations window parameters | |
| cell gas changes pause times | Pressurize delay (From Standard to DRC) = 30 Exhaust delay (From DRC to Standard mode) = 30 Flow delay (gas changes while in DRC mode) = 30 Channel delay (channel change in DRC mode) = 30 |
| File names & directories | |
| Method file names | For DLS 3018 CDC_DLS3018_15 element_urblk.mth CDC_DLS3018_15 element_aqblk.mth For DLS 3018A CDC_DLS3018A_UTAS_urblk.mth CDC_DLS3018A_UTAS_aqblk.mth |
| Dataset | Create a new dataset subfolder each day. Name as “2011-0718” for all work done on July 18, 2011 |
| Sample file | Create for each day’s work |
| Report file name | <i>For sample results printouts</i> cdc_quant comprehensive.rop <i>For calibration curve information</i> CDC_Quant Comprehensive (calib curve info).rop |
| Tuning | Default.tun |
| Optimization | Default.dac |
| Calibration | N/A |
| Polyatomic | elan.ply (ELAN), polyatomic.ply (NexION) |
| Report options template (transferring results to the database) | CDC_Database Output.rop <i>Report Format Options: select only “Use Separator”</i> <i>File Write Option: Append</i> <i>Report File name: include date, instrument, and group being analyzed in file name (i.e. 2012-0311b_NexIONA_HM-0364.txt)</i> |
| Method parameters | |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|--|---|
| Method parameters: timing page (see Figures 2a and 3a in the Appendix) | |
| sweeps/reading | 40 |
| readings/replicate | 1 |
| replicates | 3 |
| enable QC checking | On (Figures 2a and 2g for 3018, 3a and 3g for 3018A) |
| isotopes monitored and internal standard associations (exact mass) | <p>For DLS 3018A (use ⁷¹Ga as an internal standard) ⁷¹Ga (70.9249), ⁷⁵As (74.9216)</p> <p>For DLS 3018 <u>Group 1 (use ¹⁰³Rh as an internal standard)</u> ⁹Be (9.0122), ⁵⁹Co (58.9332), ⁸⁸Sr (87.9056), ⁹⁸Mo (97.9055), ¹⁰³Rh (102.905), ¹¹⁸Sn (117.902), ¹²¹Sb (120.904), ¹³³Cs (132.905), ¹³⁸Ba (137.905)</p> <p><u>Group 2 (use ¹⁹³Ir as an internal standard)</u> ¹⁸⁴W (183.951), ¹⁹³Ir (192.963), ¹⁹⁵Pt (194.965), ²⁰⁵Tl (204.975), ²⁰⁸Pb (207.977), ²³⁸U (238.05)</p> <p><u>Group 3 NexION only: (use ¹⁹³Ir as an internal standard)</u> ¹⁹³Ir (192.963), ¹¹⁴Cd (113.904), ⁵⁵Mn (54.9381),</p> <p>3018 can be performed analyzing only a subset of analytes and the appropriate internal standards.</p> |
| dwelt times | <p>30 ms for ⁵⁹Co, ⁸⁸Sr, ⁹⁸Mo, ¹¹⁸Sn, ¹⁰³Rh (vented mode), ¹²¹Sb, ¹³³Cs, ¹³⁸Ba, ¹⁸⁴W, ¹⁹³Ir (vented mode), ²⁰⁵Tl, and ²⁰⁸Pb</p> <p>50 ms for ⁷¹Ga, and ⁷⁵As</p> <p>100 ms for ⁹Be, ⁵⁵Mn, ¹⁰³Rh (DRC mode), ¹¹⁴Cd, ¹⁹³Ir (DRC mode), ¹⁹⁵Pt, and ²³⁸U</p> |
| scan mode | Peak Hopping for all isotopes (1 MCA channel) |
| DRC channel a gas flow rate | <p>For DLS 3018A 10% hydrogen / 90% argon (5-7 psig delivery pressure)</p> <p>Typically 0.7 mL/min (0.56 – 0.84) * *(optimized per instrument, and periodically verified)</p> <p>For 15 DLS 3018 Not used.</p> |
| DRC channel B gas flow rate | <p>For DLS 3018A Not used.</p> <p>For 15 DLS 3018 Oxygen (5-7 psig delivery pressure) Typically 2.3 (1.8 – 2.8) mL/min * * (optimized instrument, and periodically verified)</p> |
| RPa | 0 for all isotopes |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|--|---|
| RPq | <p>For DLS 3018A DRC Mode (As group): Typically* 0.65 - 0.75 for ⁷¹Ga (70.9249), ⁷⁵As (74.9216). Use the same RPQ for each.</p> <p>For DLS 3018 Standard Mode: 0.25 for all standard mode isotopes DRC Mode (Cd and Mn group): default 0.75 (0.65 – 0.80) * for ¹¹⁴Cd (113.904) and ⁵⁵Mn (54.9381), and ¹⁰³Rh (102.905) or ¹⁹³Ir (192.963) in DRC mode. Use the same RPQ for each. (* Optimize per instrument, and periodically verified)</p> |
| Method parameters: processing page (see Figures 2b and 3b in the Appendix) | |
| detector mode | Dual (Pulse for 3018A) |
| process spectral peak | N/A |
| Autolens | On |
| isotope ratio mode | Off |
| enable short settling time | Off |
| blank subtraction | After internal standard |
| measurement units | cps |
| process signal profile | N/A |
| Method parameters: equations page (see Figures 2c and 3c in the Appendix) | |
| equations | <p>On ²⁰⁸Pb, use "+ Pb 206 + Pb 207"</p> <p>On ²³⁸U, use "+ U 235"</p> <p>On 114Cd, use "- 0.027250 * Sn 118"</p> |
| Method parameters: calibration page (see Figures 2d and 3d in the Appendix) | |
| calibration type | external standard |
| curve type | weighted linear |
| sample units | "µg/L" or "ppb" |
| calibration standard concentrations (µg/L) | <p>Be: 0.1, 0.3, 1, 3, 10</p> <p>Co: 0.075, 0.225, 0.75, 2.25, 7.5</p> <p>Sr: 6, 18, 60, 180, 600</p> <p>Mo: 3, 9, 30, 90, 300</p> <p>Sn: 0.3, 0.9, 3, 9, 30</p> <p>Sb: 0.08, 0.24, 0.8, 2.4, 8</p> <p>Cs: 0.2, 0.6, 2, 6, 20</p> <p>Ba: 0.2, 0.6, 2, 6, 20</p> <p>W: 0.06, 0.18, 0.6, 1.8, 6</p> <p>Pt: 0.025, 0.075, 0.25, 0.75, 2.5</p> <p>Tl: 0.04, 0.12, 0.4, 1.2, 4</p> <p>Pb: 0.1, 0.3, 1, 3, 10</p> <p>U: 0.005, 0.015, 0.05, 0.15, 0.5</p> <p>Cd: 0.08, 0.24, 0.8, 2.4, 8</p> <p>Mn: 0.1, 0.3, 1, 3, 10</p> <p>As: 2, 6, 20, 60, 200</p> |
| Method parameters: sampling page (see Figures 2e and 3e in the Appendix) | |
| "peristaltic pump under computer control" | On |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|--|--|
| <p>autosampler tray port sampling device</p> | <p><i>If using ESI autosampler</i> Autosampler Type: AS-93plus Tray Name: esi.try Sampling Device: None</p> <p><i>If using other autosampler</i> Refer to autosampler user guide.</p> |
| <p>sample flush</p> | <p><i>FAST Defaults For DLS 3018A or 3018 single element</i> 2s at 6 rpm (standard ICP-MS peristaltic pump) 2s at 3 rpm (ESI DXi peristaltic pump)</p> <p><i>FAST Defaults For DLS 3018</i> 10s at 6 rpm (standard ICP-MS peristaltic pump) 10s at 3 rpm (ESI DXi peristaltic pump)</p> <p><i>FAST defaults For DLS 3018 subset (~3 elements)</i> 3.5s at 6 rpm (standard ICP-MS peristaltic pump) 3.5s at 3 rpm (ESI DXi peristaltic pump)</p> <p>Can be optimized as needed to adequately fill the FAST loop. As a matter of lab practice, set this time to equal the loop fill time in the ESI FAST program. As long as the combined time of sample flush + read delay is equal to the time required for signal to reach stability, analytical measurement will be good.</p> |
| <p>read delay</p> | <p><i>Default For DLS 3018A or DLS 3018 single element</i> 20s at 6 rpm (standard ICP-MS peristaltic pump) 20s at 3rpm (ESI DXi peristaltic pump)</p> <p><i>Default for DLS 3018 or DLS 3018 subset</i> 30s at 6 rpm (standard ICP-MS peristaltic pump) 30s at 3 rpm (ESI DXi peristaltic pump)</p> <p>Can be optimized as needed to reach signal stability before beginning analysis. As a matter of lab practice, set this time equal to the total time required for the signal to reach stability minus the loop fill time. As long as the combined time of sample flush + read delay is equal to the time required for signal to reach stability, analytical measurement will be good.</p> |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---|----------------------------|----------------|--------------|----------------------------|----|----------|------|----|----------|------|----|------------|------|----|----------|------|----|---------|------|----|---------|------|----|---------|------|----|----------|------|---|---------|------|----|---------|------|----|---------|------|----|----------|------|---|----------|------|----|---------|------|----|---------|------|----|-----------|------|
| wash | <p>Default For DLS 3018A or DLS 3018 single element 20s at 6 rpm (standard ICP-MS peristaltic pump) 20s at 3rpm (ESI DXi peristaltic pump)</p> <p>Can be optimized to allow for changes in FAST loop rinsing (must be greater than total time of steps in FAST program after the initial “on rinse” command).</p> <p>Default For DLS 3018 or DLS 3018 subset 50s at 6 rpm (standard ICP-MS peristaltic pump) 50s at 3 rpm (default for entire panel)</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| extended wash (via ICP-MS software QC checking) | <p>For sample concentrations greater than these, setup the ICP-MS software’s ‘QC checking’ feature to “Wash for X and continue”</p> <table border="1"> <thead> <tr> <th>Analyte</th> <th>Conc.</th> <th>Extended Rinse Time</th> </tr> </thead> <tbody> <tr><td>Be</td><td>300 µg/L</td><td>200s</td></tr> <tr><td>Co</td><td>225 µg/L</td><td>200s</td></tr> <tr><td>Sr</td><td>18000 µg/L</td><td>200s</td></tr> <tr><td>Mo</td><td>900 µg/L</td><td>200s</td></tr> <tr><td>Sn</td><td>90 µg/L</td><td>200s</td></tr> <tr><td>Sb</td><td>24 µg/L</td><td>200s</td></tr> <tr><td>Cs</td><td>60 µg/L</td><td>200s</td></tr> <tr><td>Ba</td><td>600 µg/L</td><td>200s</td></tr> <tr><td>W</td><td>18 µg/L</td><td>200s</td></tr> <tr><td>Pt</td><td>75 µg/L</td><td>200s</td></tr> <tr><td>Tl</td><td>12 µg/L</td><td>200s</td></tr> <tr><td>Pb</td><td>300 µg/L</td><td>200s</td></tr> <tr><td>U</td><td>1.5 µg/L</td><td>200s</td></tr> <tr><td>Cd</td><td>24 µg/L</td><td>200s</td></tr> <tr><td>Mn</td><td>30 µg/L</td><td>200s</td></tr> <tr><td>As</td><td>6000 µg/L</td><td>200s</td></tr> </tbody> </table> | | Analyte | Conc. | Extended Rinse Time | Be | 300 µg/L | 200s | Co | 225 µg/L | 200s | Sr | 18000 µg/L | 200s | Mo | 900 µg/L | 200s | Sn | 90 µg/L | 200s | Sb | 24 µg/L | 200s | Cs | 60 µg/L | 200s | Ba | 600 µg/L | 200s | W | 18 µg/L | 200s | Pt | 75 µg/L | 200s | Tl | 12 µg/L | 200s | Pb | 300 µg/L | 200s | U | 1.5 µg/L | 200s | Cd | 24 µg/L | 200s | Mn | 30 µg/L | 200s | As | 6000 µg/L | 200s |
| Analyte | Conc. | Extended Rinse Time | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Be | 300 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Co | 225 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Sr | 18000 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Mo | 900 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Sn | 90 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Sb | 24 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Cs | 60 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Ba | 600 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| W | 18 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Pt | 75 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tl | 12 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Pb | 300 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| U | 1.5 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Cd | 24 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Mn | 30 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| As | 6000 µg/L | 200s | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|--|--|
| autosampler locations of blanks and standards | <p>For DLS 3018A <i>For calibration curve (points to urine blank)</i> CDC_UMP4_DLS3018A_Urine Arsenic_urblk.mth By default (but can be customized), urine blank, 102; calibration stds, 103-107.</p> <p><i>For QC & patient samples (points to aqueous blank)</i> CDC_UMP4_DLS3018A_Urine Arsenic_aqblk.mth By default (but can be customized), aqueous blank, 149.</p> <p>For 15 DLS 3018 <i>For calibration curve (points to urine blank)</i> CDC_UMP4_DLS3018_15 elem_urblk.mth By default (but can be customized), urine blank, 102; calibration stds, 103-107.</p> <p><i>For QC & patient sample (points to aqueous blank)</i> CDC_UMP4_DLS3018_15 elem_aqblk.mth By default (but can be customized), aqueous blank, 149.</p> |
| FAST parameters: See Figures 4a through 4k in Appendix B for details | |
| configuration file | default.sc (saved at C:\Program Files\ESI\ESI-SC\ OR at C:\Users\Public\ESI\ESI SC) |
| FAST programs | <p>For DLS 3018 e.g. cdc_dls3018_15element_loop3.0ml_scfast.txt others will be needed for different loop sizes and subsets of elements</p> <p>For DLS 3018A cdc_dls3018A_arsenic_loop0.5ml_scfast.txt</p> |
| Potential emergency response modifications: | |
| <u>cadmium:</u> | Analyze cadmium in standard mode with rhodium as the internal standard. Set dwell time to 50ms, DRC gas flow to 0, and RPq to 0.25. |

| Table 1. Instrument and method parameters. Parameters are the same for DLS 3018 and DLS 3018A unless otherwise noted. | |
|--|--|
| <p><u>arsenic:</u></p> | <ul style="list-style-type: none"> • Pure argon can be used in place of 10% hydrogen 90% argon for the DRC gas to remove the $^{40}\text{Ar}^{35}\text{Cl}^+$ interference. Small interferences ($^{40}\text{Ca}^{35}\text{Cl}^+$, $^{59}\text{Co}^{16}\text{O}^+$) are expected, but are not anticipated to be significant for emergency response applications. A tee can be setup on the main argon delivery line for the ICP-MS to provide this argon for the DRC. No modifications of the DRC gas flow rate necessary. • Arsenic can be analyzed along with the 15-element method to create a 16-element panel with no ethanol in the diluent and rinse solutions. This could bias As results 1-5% high and will likely require the use of revised QC limits for arsenic. |
| <p><u>Non-FAST sample introduction system:</u></p> | <p>If the FAST sample introduction system is not available on any instruments, the method can still be implemented, but these changes will need to be made in the ICP-MS software (and ESI software if present).</p> <ul style="list-style-type: none"> • <u>Sample flush:</u> Default is ~90s at 10 rpm. Set so that solution reaches nebulizer. • <u>Read delay:</u> Default is 20s at 10rpm. Set for best reproducibility of replicate measured intensities. • <u>Wash:</u> Default is 120s at 24rpm. Set to prevent significant carry-over from one sample to the next. • If using ESI autosampler without FAST, disable FAST in the ESI software before running analysis. |

Appendix B (continued)

| Table 2. Suggested maximum analyte concentrations for base urine. | |
|--|-----------------------------|
| Analyte | Concentration (µg/L) |
| Be | 0.5 |
| Co | 0.25 |
| Mo | 30 |
| Sb | 0.2 |
| Cs | 3 |
| Ba | 2 |
| W | 0.2 |
| Pt | 0.25 |
| Tl | 0.2 |
| Pb | 0.75 |
| U | 0.03 |
| Cd | 0.25 |
| Mn | 0.1 |
| Sr | 80 |
| Sn | 3 |
| As | 5 |

Appendix B (continued)

| Table 3. Multi-element stock standard concentrations | | |
|---|---|--|
| analyte | stock calibration standard conc. (mg/L) | stock calibration standard conc. (mg/L) |
| | High Purity Standards Item # SM-2107-037 Solution A (5% HNO ₃) | High Purity Standards Item # SM-2107-037 Solution B (5% HNO ₃ , 1% HF, 0.5% HCl) |
| Be | 200 | |
| Co | 150 | |
| Mo | | 6000 |
| Sb | | 160 |
| Cs | 400 | |
| Ba | 400 | |
| W | | 120 |
| Pt | | 50 |
| Tl | 80 | |
| Pb | 200 | |
| U | 10 | |
| As | 4000 | |
| Cd | 160 | |
| Sr | 12,000 | |
| Sn | | 600 |
| Mn | 200 | |

Appendix B (continued)

| Table 4. Preparation of multi-element intermediate stock calibration standard | |
|---|--------------------------------|
| volume of flask (mL) | 100 |
| volume of spike of stock standard A (mL)* | 5 |
| volume of spike of stock standard B (mL)* | 5 |
| | concentrations (mg / L) |
| Be | 10 |
| Co | 7.5 |
| Mo | 300 |
| Sb | 8 |
| Cs | 20 |
| Ba | 20 |
| W | 6 |
| Pt | 2.5 |
| Tl | 4 |
| Pb | 10 |
| U | 0.5 |
| Cd | 8 |
| Sr | 600 |
| Sn | 30 |
| Mn | 10 |
| As | 200 |
| *If preparing from HPS # SM-2107-037, both stock solutions A and B need to be spiked into the flask | |

Appendix B (continued)

| Table 5. Preparation of multi-element intermediate working standards | | | | | |
|---|-------------------|-------------------|-----------------|------------------|-----------------|
| standard # | 1 | 2 | 3 | 4 | 5 |
| flask vol. (mL) | 500 | 200 | 100 | 100 | 100 |
| vol. spike of int. stock std. (mL) | 0.050 | 0.060 | 0.100 | 0.300 | 1.00 |
| concentrations (µg/L) ‡ | | | | | |
| Be | 1 (0.1) ‡ | 3 (0.3) ‡ | 10 (1.0) ‡ | 30 (3.0) ‡ | 100 (10.0) ‡ |
| Co | 0.75 (0.075) ‡ | 2.25 (0.225) ‡ | 7.5 (0.75) ‡ | 22.5 (2.25) ‡ | 75 (7.5) ‡ |
| Mo | 30 (3.0) ‡ | 90 (9.0) ‡ | 300 (30) ‡ | 900 (90) ‡ | 3000 (300) ‡ |
| Sb | 0.8 (0.08) ‡ | 2.4 (0.24) ‡ | 8 (0.8) ‡ | 24 (2.4) ‡ | 80 (8.0) ‡ |
| Cs | 2 (0.2) ‡ | 6 (0.6) ‡ | 20 (2.0) ‡ | 60 (6.0) ‡ | 200 (20) ‡ |
| Ba | 2 (0.2) ‡ | 6 (0.6) ‡ | 20 (2.0) ‡ | 60 (6.0) ‡ | 200 (20) ‡ |
| W | 0.6 (0.06) ‡ | 1.8 (0.18) ‡ | 6 (0.6) ‡ | 18 (1.8) ‡ | 60 (6.0) ‡ |
| Pt | 0.25 (0.025) ‡ | 0.75 (0.075) ‡ | 2.5 (0.25) ‡ | 7.5 (0.75) ‡ | 25 (2.5) ‡ |
| Tl | 0.4 (0.04) ‡ | 1.2 (0.12) ‡ | 4 (0.4) ‡ | 12 (1.2) ‡ | 40 (4.0) ‡ |
| Pb | 1 (0.1) ‡ | 3 (0.3) ‡ | 10 (1.0) ‡ | 30 (3.0) ‡ | 100 (10) ‡ |
| U | 0.05 (0.005) ‡ | 0.15 (0.015) ‡ | 0.5 (0.05) ‡ | 1.5 (0.15) ‡ | 5 (0.5) ‡ |
| Cd | 0.8 (0.08) ‡ | 2.4 (0.24) ‡ | 8 (0.8) ‡ | 24 (2.4) ‡ | 80 (8.0) ‡ |
| Sr | 60 (6.0) ‡ | 180 (18.0) ‡ | 600 (60) ‡ | 1800 (180) ‡ | 6000 (600) ‡ |
| Sn | 3 (0.3) ‡ | 9 (0.9) ‡ | 30 (3.0) ‡ | 90 (9.0) ‡ | 300 (30) ‡ |
| Mn | 1 (0.1) ‡ | 3 (0.3) ‡ | 10 (1.0) ‡ | 30 (3.0) ‡ | 100 (10) ‡ |
| As | 20 (2.0) ‡ | 60 (6.0) ‡ | 200 (20) ‡ | 600 (60) ‡ | 2000 (200) ‡ |
| ‡ A further 1:10 dilution occurs when added to base urine. Enter concentrations in parentheses into the ICP-MS software (method window, calibration page). | | | | | |

Appendix B (continued)

| Table 6. Acceptable ways to perform two consecutive analytical runs, bracketing with bench quality control samples. | |
|--|---|
| setup 1 | setup 2 (typical) |
| <p><i>Run #1</i> calibration standards low bench QC high bench QC patient samples low bench QC high bench QC</p> <p><i>Run #2</i> low bench QC high bench QC patient samples low bench QC high bench QC</p> | <p><i>Run #1</i> calibration standards low bench QC high bench QC patient samples low bench QC high bench QC</p> <p><i>Run #2</i> calibration standards low bench QC high bench QC patient samples low bench QC high bench QC</p> |

Appendix B (continued)

| Table 7. A typical SAMPLE/BATCH window. | | | |
|---|-----------------------------|-------------------------------------|---------------------|
| <u>AS location*</u> | <u>sample id</u> | <u>measurements action</u> | <u>method</u> |
| 236 | DRCstability1 | Run sample | ...15elem_urblk.mth |
| 236 | DRCstability2 | Run sample | ...15elem_urblk.mth |
| 236 | DRCstability3 | Run sample | ...15elem_urblk.mth |
| 236 | DRCstability4 | Run sample | ...15elem_urblk.mth |
| Continue DRC stability samples . . . | | | |
| 236 | DRCstability11 | Run sample | ...15elem_urblk.mth |
| 236 | DRCstability12 [£] | Run sample | ...15elem_urblk.mth |
| 101 | 3018 UrBlkChk Wash1 | Run blank, standards, and sample ** | ...15elem_urblk.mth |
| 113 | 3018 UrBlkChk Wash2 | Run sample | ...15elem_urblk.mth |
| 114 | 3018 UrBlkChk1 | Run sample | ...15elem_urblk.mth |
| 115 | 3018 UrBlkChk2 | Run sample | ...15elem_urblk.mth |
| 150 | 3018 AQBLK | Run blank and sample ‡ | ...15elem_aqblk.mth |
| 136 | L Bench QC | Run sample | ...15elem_aqblk.mth |
| 160 | H Bench QC | Run sample | ...15elem_aqblk.mth |
| 301 | Sample 1 | Run sample | ...15elem_aqblk.mth |
| 302 | Sample 2 | Run sample | ...15elem_aqblk.mth |
| 303 | Sample 3 | Run sample | ...15elem_aqblk.mth |
| 124 | L Bench QC | Run sample | ...15elem_aqblk.mth |
| 148 | H Bench QC | Run sample | ...15elem_aqblk.mth |
| <p>* The exact autosampler positions of QCs and patient samples do not have to be those shown above. QC samples do not have to be run in the order of low, then high, then elevated.</p> <p>** When executing this row, the ICP-MS will first analyze the urine blank (standard 0) at AS position 102, then standards 1-5 at autosampler positions 102-107, <u>then</u> the “UrBlkChk wash1” sample at A/S position 102. The sampling information about AS positions 103-107 are stored in the “urblk” method file and can be customized.</p> <p>‡ When executing this row, the ICP-MS will first analyze the aqueous blank at AS position 149, then the “AQBlkChk” at AS position 150. The sampling information about AS positions 149 is stored in the “aqblk” method file and can be customized.</p> <p>£ A larger number of DRC stability samples will need to be analyzed to make this stability period 1-1.5 hrs when measuring only arsenic (50 measurements ~ 1hour).</p> | | | |

Appendix B (continued)

Table 8a. Preparation of samples, working standards, and QC materials for analysis for DLS 3018A or DLS 3018 single element *

If a different total volume is prepared, adjust the volumes for each component proportionally.

* These directions are written with the expectation of a 10,000 μL syringe on the left side and a 1,000 μL syringe on the right side of the benchtop automatic pipette.

| Description | water (μL) | base urine (μL) | AQ intermediate working standard (μL) | patient or QC urine sample (μL) | Diluent ** (μL) |
|---|-------------------------|------------------------------|--|--|------------------------------|
| working calibration standards (S0-S5) and UrBlkChk (S0) | - | 900 x 1 | 100 x 1 | - | 9,000 (4,500 x 2) |
| AQ Blank | 1,000 x 1 | - | - | - | 9,000 (4,500 x 2) |
| patient urine or urine-based QC | - | - | - | 250 x 1 | 2,250 x 1 |
| patient urine <i>2x dilution</i> ^H | 250 x 1 | - | - | 250 x 1 | 4,500 (2,250 x 2) |
| patient urine <i>5x dilution</i> ^H | 400 x 1 | - | - | 100 x 1 | 4,500 (2,250 x 2) |
| patient urine <i>10x dilution</i> ^H | 900 x 1 | - | - | 100 x 1 | 9,000 (4,500 x 2) |
| patient urine <i>20x dilution</i> ^H | 4,750 (950 x 5) | - | - | 250 x 1 | 45,000 (7,500 x 6) |

** By splitting the dispense step of diluent into two or more portions, liquids pulled up into the right pipette tip are flushed out more completely. For example, when preparing a working calibration standard dilution, do the preparation in two steps: in step 1, dispense 4500 μL diluent + 100 μL ; in step 2, dispense 4500 μL diluent + 900 μL base urine to prepare a 10 mL total volume dilution.

^H Extra dilution is performed on urine samples whose concentration is greater than the highest calibrator listed in Table 5 in the Appendix B.

Maximum extra dilution (see Appendix A, ruggedness test #6 for details)

2x Sr

10x Cd, Mo, Pt, Tl

20x (As with ethanol in reagents, Ba, Be, Co, Cs, Mn, Pb, Sb, Sn, U, and W)

Any extra dilutions within these limits can be prepared as long as the 9:10 ratio of diluent to total dilution volume is maintained. Use of the lowest possible dilution level is preferred to minimize differences between the calibrators and the samples (i.e. 2x dilution is preferred over 10x if 2x is sufficient to dilute analyte into the documented linearity range).

Appendix B (continued)

| Table 8b. Preparation of samples, working standards, and QC materials for analysis for DLS 3018 (or DLS 3018 subset) * | | | | | |
|---|-------------------|------------------------|--|--|------------------------|
| If a different total volume is prepared, adjust the volumes for each component proportionally. | | | | | |
| * These directions are written with the expectation of a 10,000 µL syringe on the left side and a 1,000 µL syringe on the right side of the benchtop automatic pipette. | | | | | |
| dilution ID | water (µL) | base urine (µL) | AQ intermediate working standard (µL) | patient or QC urine sample (µL) | Diluent ** (µL) |
| working calibration standards (S0-S5) and UrBlkChk (S0) | - | 900 x 1 | 100 x 1 | - | 9,000 (4,500 x 2)* |
| AQ blank | 1000 x 1 | - | - | - | 9,000 (4,500 x 2)* |
| patient urine or urine-based QC | - | - | - | 500 x 1 | 4,500 x 1 |
| patient urine 2x dilution ^H | 500 x 1 | - | - | 500 x 1 | 9,000 (4,500 x 2)* |
| patient urine 5x dilution ^H | 800 x 1 | - | - | 200 x 1 | 9,000 (4,500 x 2)* |
| patient urine 10x dilution ^H | 900 x 1 | - | - | 100 x 1 | 9,000 (4,500 x 2)* |
| patient urine 20x dilution ^H | 4,750 (950 x 5) | - | - | 250 x 1 | 45,000 (7,500 x 6)* |
| ** By splitting the dispense step of diluent into two or more portions, liquids pulled up into the right pipette tip are flushed out more completely. For example, when preparing a working calibration standard dilution, do the preparation in two steps: in step 1, dispense 4500 µL diluent + 100 µL; in step 2, dispense 4500 µL diluent + 900 µL base urine to prepare a 10 mL total volume dilution. | | | | | |
| ^H Extra dilution is performed on urine samples whose concentration is greater than the highest calibrator listed in Table 4b in the Appendix B. | | | | | |
| Maximum extra dilution (see Appendix A, ruggedness test #6 for details) | | | | | |
| 2x Sr | | | | | |
| 10x Cd, Mo, Pt, Tl | | | | | |
| 20x (As with ethanol in reagents, Ba, Be, Co, Cs, Mn, Pb, Sb, Sn, U, and W) | | | | | |
| Any extra dilutions within these limits can be prepared as long as the 9:10 ratio of diluent to total dilution volume is maintained. Use of the lowest possible dilution level is preferred to minimize differences between the calibrators and the samples (i.e. 2x dilution is preferred over 10x if 2x is sufficient to dilute analyte into the documented linearity range). | | | | | |

Appendix B (continued)

| Table 9. Boundary concentrations for urine concentrations (µg/L). | | | | |
|--|--|---|--|--|
| analyte | 1st upper boundary ("1UB") * | 2nd upper boundary ("2UB") ** | Range Maximum ("Lim Rep Delta") † | Highest Concentration Validated for Washout |
| Be | 0.2 | 0.4 | 0.3 | 300 |
| Co | 2.83 | 5.66 | 0.3 | 225 |
| Mo | 293.5 | 587 | 4.0 | 9,000 |
| Sb | 0.8 | 1.6 | 0.2 | 240 |
| Cs | 16.5 | 33 | 0.5 | 600 |
| Ba | 17.1 | 34.2 | 0.4 | 600 |
| W | 1.38 | 2.76 | 0.2 | 180 |
| Pt | 0.2 | 0.4 | 0.2 | 75 |
| Tl | 0.62 | 1.24 | 0.2 | 120 |
| Pb | 7.8 | 15.6 | 0.3 | 300 |
| U | 0.277 | 0.554 | 0.03 | 15 |
| Cd | 2.54 | 5.08 | 0.3 | 240 |
| Mn | 4 | 8 | 0.4 | 300 |
| Sr | 400 | 800 | 3 | 18,000 |
| Sn | 25 | 50 | 0.5 | 900 |
| As | 100 | 200 | 10 | 6,000 |

* Typically, the 1UB threshold is based on percentiles of non-weighted, non-creatinine corrected concentration results from NHANES. In the absence of that data, these boundaries can be based on normal ranges reported in the literature. The concentrations assigned to these boundaries is determined by study protocol but default concentrations are listed in this table.

**Typically the 2nd upper boundary (2UB) is set to 2x the 1UB. The concentrations is determined by study protocol but default concentrations are listed in this table.

† Range maximum is the range of the three replicate readings for a single sample analysis. This value is also called the "Lim RepDelta" in the database which handles data for the Inorganic and Radiation Analytical Toxicology Branch.

Appendix B (continued)

| Table 10. Reference ranges for urine concentrations (from the Fourth National Report on Exposure to Environmental Chemicals [15]). All results in µg/L. | | | | | | | |
|---|--|-----------------------|------------------------|------------------------|------------------------|------------------------|----------|
| analyte | survey years | geometric mean | 50th | 75th | 90th | 95th | N |
| Be | 07-08 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | 2627 |
| | 09-10 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | ≤ 0.072 | 2848 |
| Co | 07-08 | 0.369 | 0.380 | 0.580 | 0.920 | 1.32 | 2627 |
| | 09-10 | 0.369 | 0.380 | 0.600 | 0.960 | 1.40 | 2848 |
| Mo | 07-08 | 45.2 | 49.7 | 86.2 | 136 | 163 | 2627 |
| | 09-10 | 42.7 | 45.2 | 79.9 | 121 | 160 | 2848 |
| Sb | 07-08 | 0.610 | 0.600 | 0.100 | 0.170 | 0.240 | 2627 |
| | 09-10 | 0.560 | 0.500 | 0.090 | 0.170 | 0.230 | 2847 |
| Cs | 07-08 | 4.42 | 4.90 | 7.02 | 9.46 | 11.3 | 2627 |
| | 09-10 | 4.06 | 4.47 | 6.54 | 8.98 | 11.0 | 2848 |
| Ba | 07-08 | 1.56 | 1.64 | 3.01 | 4.93 | 7.04 | 2627 |
| | 09-10 | 1.47 | 1.48 | 2.81 | 4.66 | 6.78 | 2848 |
| W | 07-08 | 0.099 | 0.100 | 0.190 | 0.360 | 0.500 | 2627 |
| | 09-10 | 0.081 | 0.080 | 0.160 | 0.310 | 0.460 | 2847 |
| Pt | 07-08 | ≤ 0.009 | ≤ 0.009 | ≤ 0.009 | ≤ 0.009 | 0.014 | 2627 |
| | 09-10 | ≤ 0.009 | ≤ 0.009 | ≤ 0.009 | 0.009 | 0.016 | 2847 |
| Tl | 07-08 | 0.146 | 0.160 | 0.250 | 0.330 | 0.400 | 2627 |
| | 09-10 | 0.144 | 0.160 | 0.240 | 0.340 | 0.410 | 2848 |
| Pb | 07-08 | 0.493 | 0.500 | 0.850 | 1.38 | 1.97 | 2627 |
| | 09-10 | 0.458 | 0.470 | 0.790 | 1.24 | 1.65 | 2848 |
| U | 07-08 | 0.007 | 0.007 | 0.013 | 0.024 | 0.039 | 2627 |
| | 09-10 | 0.007 | 0.007 | 0.013 | 0.022 | 0.036 | 2848 |
| Cd | 07-08 | 0.185 | 0.180 | 0.380 | 0.700 | 1.00 | 2627 |
| | 09-10 | 0.179 | 0.180 | 0.370 | 0.690 | 1.03 | 2848 |
| As | 07-08 | 8.10 | 7.49 | 14.9 | 33.3 | 50.8 | 2605 |
| | 09-10 | 9.28 | 8.15 | 18.0 | 44.6 | 85.6 | 2860 |
| Mn, Sn, and Sr | See Table 11 in Appendix B for reference values. These elements were not part of the analyses reported in the Fourth National Report on Exposure to Environmental Chemicals. | | | | | | |

Appendix B (continued)

| Table 11. Reference concentrations from published literature for urine Mn, Sn and Sr. | | | |
|--|---------------------------|-----------------------------|---|
| analyte | reference | concentration (µg/L) | group type sampled |
| normal (i.e. non-exposed) | | | |
| Mn | Health Canada, 2010[16] | 0.15 | general population |
| Mn | Heitland et al., 2006[17] | 0.1 | German children |
| Mn | Heitland et al., 2006[17] | 0.087 | German adults |
| Mn | Paschal et al., 1998[18] | 1.19 | NHANES III |
| Mn | ASTDR[19] | 1 to 8 | general population |
| Sr | Heitland et al., 2006[17] | 154 | German children |
| Sr | Heitland et al., 2006[17] | 166 | German adults |
| Sr | Usuda et al., 2006 | 143.9 [‡] | Japanese adults |
| Sn | Heitland et al., 2006[17] | 1.2 | German children |
| Sn | Heitland et al., 2006[17] | 8.6 | German adults |
| Sn | Paschal et al., 1998[18] | 6.29 | NHANES III |
| elevated (exposed) | | | |
| Mn | Moreno et al., 2010[20] | 5.2 | children living in a mine tailings zone in Mexico |
| Mn | Gil et al., 2011[21] | 0.43 +/- 4,00 | Iron and steel industry workers in Spain |
| Mn | Wang et al., 2011[22] | 3.15 +/- 3.45 | e-waste dismantling workers in China |
| Sr | Moreno et al., 2010[20] | 49.2 | children living in a mine tailings zone in Mexico |
| Sn | Juliao et al., 2007[23] | 0.45 +/- 0.93 | workers in a niobium mine in Brazil |
| | | | |
| ‡ concentration is a geometric mean and measurements were made via ICP-AES | | | |

Appendix C: Help Sheets

Reagent Preparation 3018, page 1 of 2

NOTE:

mg/L = ppm

µg/L = ppb

µg/mL = ppm

Rinse Solution – 4 L

(5.0% (v/v) HNO₃, 0.002% Triton X-100, 500 µg/L gold)

1. Partially fill a 4 liter bottle with 18 MΩ·cm water.
2. Add 4 mL of the 2% Triton X-100™ / 5% (v/v) nitric-acid intermediate stock solution.
3. Carefully add 200 mL of concentrated HNO₃.
4. Add 200 µL of the 10,000 µg/mL gold internal standard solution.
5. Add enough ≥18 MΩ·cm water to bring to 4 liter mark.
6. Mix well by gently inverting several times.
7. Label appropriately and store at room temperature.

Sample Diluent/Carrier Solution – 2 L

(2.0% (v/v) HNO₃, 10 µg/L Ir and Rh, 500 µg/L Au)

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 40 mL concentrated HNO₃.
3. Add 500 µL of the 40 µg/mL Rh and Ir internal standard solution.
4. Add 100 µL of the 10,000 µg/mL gold standard.
5. Add enough 18 MΩ·cm water to bring to 2 liter mark.
6. Mix well by gently inverting several times.
7. Label appropriately and store at room temperature.

Intermediate Internal Standard Solution – 200 mL

(2.0% (v/v) HNO₃, 40 µg/mL Ir and Rh)

1. Partially fill a 200 mL volumetric flask with 18 MΩ·cm water.
2. Add 4 mL of concentrated HNO₃.
3. Add 8 mL of 1,000 µg/mL Rh standard.
4. Add 8 mL of 1,000 µg/mL Ir standard.
5. Add enough 18 MΩ·cm water to bring to 200 mL mark.
6. Mix well by gently inverting several times.
7. Label appropriately and store at room temperature.

2% Triton X-100 in 5% (v/v) HNO₃ – 2 L

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 100 mL of concentrated HNO₃.
3. Add 40 mL of Triton X-100.
4. Add enough 18 MΩ·cm water to bring to 2 liter mark.
5. Add a clean Teflon magnetic stirring bar and stir on stirrer until dissolved.
6. Label appropriately and store at room temperature.

Appendix C: Help Sheets (continued)

Reagent Preparation 3018, page 2 of 2

0.5% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 10 mL of concentrated HNO₃.
3. Add enough 18 MΩ·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

2% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 40 mL of concentrated HNO₃.
3. Add enough 18 MΩ·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

5% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 100 mL of concentrated HNO₃.
3. Add enough 18 MΩ·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

DRC Stability Solution – 200 mL

1. Add 180 mL of diluent to 250 mL bottle.
2. Add 18 mL of human urine to bottle.
3. Add 2 mL of Standard 2 to bottle.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

Daily Performance Solution (1 μg/L) in 2% (v/v) HNO₃

1. Partially fill a 1 liter volumetric flask with 18 MΩ·cm water.
2. Add 1 mL of High Purity Standard: SM-2107-018.
3. Add 20 mL of concentrated HNO₃.
4. Add enough 18 MΩ·cm water to bring to 1 liter mark.
5. Mix well by gently inverting several times.
6. Label appropriately and store at room temperature.

Dual Detector Solution

1. Partially fill a 50 mL polypropylene tube with 2% (v/v) HNO₃.
2. Add 100 μL of 100 μg/mL High Purity Standard: SM-2107-053.
3. Add 10 μL of any additional 1,000 μg/mL single element stock standard, if required.
4. Dilute to the 50 mL mark with 2% (v/v) (v/v) HNO₃.
5. Label appropriately and store at room temperature.

Appendix C: Help Sheets (continued)

Reagent Preparation 3018A, page 1 of 2

NOTE:

mg/L = ppm

µg/L = ppb

µg/mL = ppm

Rinse Solution – 4 L

(5.0% (v/v) HNO₃, 0.002% Triton X-100, 1.5% (v/v) ethyl alcohol)

1. Partially fill a 4 liter bottle with 18 MΩ·cm water.
2. Add 4 mL of the 2% Triton X-100™ / 5% (v/v) nitric-acid intermediate stock solution.
3. Add 200 mL of concentrated HNO₃.
4. Add 60 mL of ethyl alcohol (200 proof).
5. Add enough 18 MΩ·cm water to bring to 4 liter mark.
6. Mix well by gently inverting several times.
7. Label appropriately and store at room temperature.

Sample Diluent/Carrier Solution – 2 L

(2.0% (v/v) HNO₃, 1.5% (v/v) Ethanol, 10 µg/L Ga)

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 40 mL concentrated HNO₃.
3. Add 30 mL ethyl alcohol (200 proof).
4. Add 500 µL of the 40 µg/mL Ga internal standard solution.
5. Add enough 18 MΩ·cm water to bring to 2 liter mark.
6. Mix well by gently inverting several times.
7. Label appropriately and store at room temperature.

Intermediate Internal Standard Solution – 200 mL

(2.0% (v/v) HNO₃, 40 µg/mL Ga)

1. Partially fill a 200 mL volumetric flask with 18 MΩ·cm water.
2. Add 4 mL of concentrated HNO₃.
3. Add 8 mL of 1,000 µg/mL Ga standard. If initial Ga concentration is different, adjust volume accordingly.
4. Add enough 18 MΩ·cm water to bring to 200 mL mark.
5. Mix well by gently inverting several times.
6. Label appropriately and store at room temperature.

Appendix C: Help Sheets (continued)

Reagent Preparation 3018A, page 2 of 2

DRC Stability Solution – 200 mL

1. Add 180 mL of diluent to 250 mL bottle.
2. Add 18 mL of human urine to bottle.
3. Add 2 mL of Standard 2 to bottle.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

2% Triton X-100 in 5% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 100 mL of concentrated HNO₃.
3. Add 40 mL of Triton X-100.
4. Add enough 18 MΩ·cm water to bring to 2 liter mark.
5. Add a clean Teflon magnetic stirring bar and stir on stirrer until dissolved.
6. Label appropriately and store at room temperature.

1% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 20 mL of conc. HNO₃.
3. Add enough 18 MΩ·cm water to bring to 2 liter mark.
4. Mix well by gently swirling several times.
5. Label appropriately and store at room temperature.

5% (v/v) HNO₃

1. Partially fill a 2 liter bottle with 18 MΩ·cm water.
2. Add 100 mL of conc. HNO₃.
3. Add enough 18 MΩ·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.
5. Label appropriately and store at room temperature.

Daily solution (1ppb) in 2% (v/v) HNO₃

1. Partially fill a 1 liter volumetric flask with 18 MΩ·cm water.
2. Add 1 mL of High Purity Standard: SM-2107-018 (or current lot #)
3. Add 20 mL of HNO₃
4. Add enough 18 MΩ·cm water to bring to 1 liter mark.
5. Mix well by gently inverting several times.
6. Label appropriately and store at room temperature.