



Laboratory Procedure Manual

Analytes: **Cadmium, Lead, Manganese, Mercury, and Selenium**

Matrix: **Whole Blood**

Method: blood multi-element analysis by ICP-DRC-MS

Method No: DLS 3016.8-06

Revised: 3/23/2018

As performed by: Inorganic and Radiation Analytical Toxicology
Branch
Division of Laboratory Sciences
National Center for Environmental Health

Contact: Jeffery M. Jarrett, MS

Phone: 770-488-7906

Fax: 770-488-4097

Email: JJarrett@cdc.gov

James L. Pirkle, M.D., Ph.D.

Director, Division of Laboratory Sciences

Important Information for Users

The Centers for Disease Control and Prevention (CDC) periodically refines these laboratory methods. It is the responsibility of the user to contact the person listed on the title page of each write-up before using the analytical method to find out whether any changes have been made and what revisions, if any, have been incorporated.

Public Release Data Set Information

This document details the Lab Protocol for testing items in the following table:

Data File Name	Variable Name	SAS Label
PBCD_J & PBY_J_R	LBXBCD	Cadmium (µg/L)
	LBDBCDSI	Cadmium (nmol/L)
	LBXBPB	Lead (µg/dL)
	LBDBPBSI	Lead (µmol/L)
	LBXTHG	Mercury, total (µg/L)
	LBPTHGSI	Mercury, total (nmol/L)
	LBXBMN	Manganese (µg/L)
	LBDBMNSI	Manganese (nmol/L)
	LBXBSE	Selenium (ug/L)
	LBDBSESI	Selenium (µmol/L)

1. Clinical relevance and summary of test principle

A. Clinical relevance

Metals ions affect human health in various ways. Some metals (i.e. lead, cadmium, and mercury) show only deleterious effects on human health. Some (i.e. selenium and manganese) play an essential role in the human biological system if within certain concentration ranges, while negative health implications are observed when concentrations in biological systems are in deficit or excess. Determination of a person's level of environmental exposure to chemicals through direct measurement of the substances or their metabolites in human specimens such as blood is called biomonitoring. Biomonitoring reduces the uncertainty of determining levels of exposure over making these determinations through calculations of estimated dose based on analysis of environmental samples and assumptions about exposure pathways[1]. Biomonitoring measurements are the most health-relevant assessments of exposure because they indicate the amount of the chemical that actually gets into people from all environmental sources (e.g., air, soil, water, dust, or food) combined. The laboratory method described here is a multi-element technique for monitoring the concentrations of cadmium (Cd), lead (Pb), manganese (Mn), mercury (Hg), and selenium (Se) in whole human blood for the purpose of biomonitoring.

There is no known biological role of mercury in the human body. The main sources of mercury intake in humans are fish, dental amalgams, and occupational exposures[2]. The main organs affected by mercury are the brain and the kidneys. Exposure of childbearing-aged women is of particular concern because of the potential adverse neurologic effects of Hg in fetuses. The health effects of mercury are diverse and depend on the form of mercury encountered and the severity and length of exposure. The general population is exposed to three forms of mercury: elemental, inorganic, and organic (predominantly methyl). However, this method tests only for the total amount of mercury in the blood without regard to chemical form. In the general population, total blood mercury is due mostly to the dietary intake of organic forms which are formed through microbial action from inorganic mercury that has deposited in aquatic environments and bioaccumulated through the food chain (especially into large predatory fish)[3]. Exposure to inorganic or elemental mercury (e.g. dental amalgams or occupational exposures) is particularly reflected in urine excretion rather than blood. Psychic and emotional disturbances are the initial signs of chronic intoxication by elemental mercury vapors or salts. Those exposed are at increased risk for parasthesia, neuralgias, renal disease, digestive disturbances, and ocular lesions[4]. Massive exposure over a longer period of time results in violent muscular spasms, hallucinations, delirium, and death[5]. Except for methylmercury exposures, blood is considered useful if samples are taken within a few days of exposure. This is because most forms of mercury in the blood decrease by one-half every three days if exposure has been stopped. Thus, mercury levels in the blood provide more useful information after recent exposures than after long-term exposures. Several months after an exposure, mercury levels in the blood and urine are much lower. Blood mercury reference ranges for the U.S. population are listed in Table 11 of Appendix C.

There is no known biological role of lead in the human body. Lead, a naturally occurring metal, has had many different commercial uses from which a person can be exposed either in the occupational / manufacturing process or by the manufactured products such as paint (paint chips, or dust and soil contaminated from deteriorating paint), solder or pipes (only now in older homes), gasoline (now outlawed for all but specialized applications), glazes on pottery, hobby uses (e.g. stained glass),

commercial products (e.g. batteries, lead-containing jewelry), home remedy medicines containing lead compounds and non-Western cosmetics. Soil contains lead naturally, or from man-made uses of lead such as paint (near older homes), gasoline (near roadways), mining, manufacturing, and disposal. The main target for lead toxicity is the nervous system, both in adults and children. The developing biological systems of children are most sensitive to the effects of Pb, where effects are being recognized even at blood lead levels <5 µg/dL [6-10]. Acute, elevated lead exposure is associated with anorexia, dyspepsia, and constipation followed by diffuse paroxysmal abdominal pain. When lead exposure is high, particularly in children, the person is at increased risk for encephalopathy[11]. The alkyl lead species are highly toxic to the central nervous system[12]. The primary screening method for lead exposure is blood lead, which primarily reflects recent exposures (excretory half-life in blood is approximately 30 days)[13]. Lead in blood is primarily (99%) in the red blood cells. Blood lead reference ranges for the U.S. population are listed in Table 11 of Appendix C. The CDC now uses a reference level of 5 µg/dL to identify children with blood lead levels that are much higher than most children's levels. This new level is based on the U.S. population of children ages 1-5 years who are in the highest 2.5% of children when tested for lead in their blood. This reference value is based on the 97.5th percentile of the National Health and Nutrition Examination Survey (NHANES)'s blood lead distribution in children. CDC will update the reference value every four years using the two most recent NHANES surveys[14].

There is no known biological role of cadmium in the human body. The predominant commercial use of cadmium is in battery manufacturing. Other uses include pigment production, coatings and plating, plastic stabilizers, and nonferrous alloys. Since 2001, U.S. cadmium use has declined in response to environmental concerns. In the United States, for nonsmokers the primary source of cadmium exposure is from the food supply. People who regularly consume shellfish and organ meats will have higher exposures. In general, leafy vegetables such as lettuce and spinach, potatoes and grains, peanuts, soybeans, and sunflower seeds contain high levels of cadmium due to bioaccumulation from the soil. Tobacco leaves accumulate high levels of cadmium from the soil, and smoking is the primary non-occupational source of cadmium exposure for smokers. Generally, the critical organ for Cd is the kidney. Kidney dysfunction is one of the most characteristic signs of exposure to Cd. Workers in an environment with high exposure levels have developed proteinuria, renal glucosuria, aminoaciduria, hypercalciuria, phosphaturia, and polyuria. Chronic obstructive lung disease of varying degrees of severities is frequently seen in Cd workers. Concentration of cadmium in blood of healthy unexposed adults are in the range 0.1 – 4 µg/L[15]. Newborn babies are practically free of Cd[16]. Exposure to high concentration of fumes appearing from heated cadmium metal or compounds has led to acute poisoning and in some cases to the death of workers[17]. Principal symptoms reported were respiratory distress due to chemical pneumonitis and edema. It has been estimated that 8 hrs. exposure to 5 g Cd/m³ will be lethal[17]. Ingestion of high amounts of Cd puts a person at increased risk to a rapid onset with severe nausea, vomiting, and abdominal pain. Cadmium levels in blood, urine, feces, liver, kidney, hair, and other tissues have been used as biological indicators of exposure to cadmium. Blood cadmium levels are principally indicative of recent exposure(s) to cadmium rather than whole-body burdens[18-21]. Urine cadmium levels primarily reflect total body burden of cadmium, although urine levels do respond somewhat to recent exposure[22]. Blood cadmium reference ranges for the U.S. population are listed in Table 11 of Appendix C.

Manganese (Mn) is a trace element essential to humans and is associated with the formation of connective and bony tissue, growth and reproductive functions and with carbohydrate and lipid metabolism[23]. Manganese is also a known neurotoxin but little information exists about levels of manganese that cause toxicity. Symptoms of manganese toxicity are similar to Parkinson's Disease and can also include disorientation, memory impairment, anxiety and compulsive behavior[24]. There is much concern for the levels of manganese in humans whom are occupationally exposed to it [25-31]. Recently, there are growing concerns over exposure due to contamination of drinking water with manganese[32-34] and as a result of methylcyclopentadienyl manganese tricarbonyl (MMT) used as an anti-knocking additive in gasoline[35-41]. Populations suffering from iron deficiencies are at an increased risk to manganese toxicity because iron deficiency can result in an accumulation of manganese in the central nervous system[38]. To fully understand the essentiality and toxicity of manganese, further investigations are needed regarding the levels of manganese in biological matrices. Group average levels in blood appear to be related to manganese body burden, while average urinary excretion levels appear to be most indicative of recent exposures [42]. On an individual basis the correlation between the level of workplace exposure and the levels in blood or urine has always been found to be a reliable predictor of exposure [26, 42-44]. Measurements of manganese levels in blood are useful in detecting groups with above-average current exposure, although levels are sometimes related to exposures that have already ceased. In addition to individual variability, another factor that limits the usefulness of measuring manganese in blood, urine, or feces as a measure of excess manganese exposure is the relatively rapid rate of manganese clearance from the body. Excess manganese in blood is rapidly removed by the liver and excreted into the bile, with very little excretion in urine [45, 46]. Thus, levels of manganese in blood or urine are not expected to be the most sensitive indicators of exposure[47]. Blood manganese reference ranges for the U.S. population are listed in Table 11 of Appendix C.

Selenium is an essential element that is required to maintain good health but both selenium deficiency and excessive levels of selenium are associated with several disorders [48, 49]. Selenium is a naturally occurring mineral element that is distributed widely in nature in most rocks and soils. Most processed selenium is used in the electronics industry, but it is also used: as a nutritional supplement; in the glass industry; as a component of pigments in plastics, paints, enamels, inks, and rubber; in the preparation of pharmaceuticals; as a nutritional feed additive for poultry and livestock; in pesticide formulations; in rubber production; as an ingredient in antidandruff shampoos; and as a constituent of fungicides. Radioactive selenium is used in diagnostic medicine. In the body, selenium is incorporated into proteins to make selenoproteins, which are important antioxidant enzymes. The antioxidant properties of selenoproteins help prevent cellular damage from free radicals. Free radicals are natural by-products of oxygen metabolism that increase risk of chronic diseases such as cancer and heart disease [49, 50]. Other selenoproteins help regulate thyroid function and play a role in the immune system [51-54]. Human selenium deficiency is rare in the U.S. but is seen in other countries where soil concentration of selenium is low[55]. There is evidence that selenium deficiency increases the risk of a form of heart disease, hypothyroidism, and a weakened immune system [56, 57]. There is also evidence that selenium deficiency does not usually cause illness by itself. Rather, it can make the body more susceptible to illnesses caused by other nutritional, biochemical or infectious stresses[58]. Symptoms of very high exposure to selenium, a condition called selenosis, include gastrointestinal upsets, hair loss, white blotchy nails, garlic breath odor, fatigue, irritability, and mild nerve damage[48]. Selenium can be detected in the

blood, feces, urine, hair, and nails of exposed individuals, however, field studies have used primarily blood or urine levels to indicate the degree of selenium exposure[48]. Blood selenium reference ranges for the U.S. population are listed in Table 11 of Appendix C.

The laboratory method presented here can be used to achieve rapid and accurate quantification of five elements of toxicological and nutritional interest including cadmium (Cd), lead (Pb), mercury (Hg), manganese (Mn) and selenium (Se) in whole human blood. Use this method to screen blood when people are suspected to be acutely exposed to these elements or to evaluate chronic environmental or other non-occupational exposure.

B. Test principle

This method directly measures the Cd, Mn, Hg, Pb, and Se content of whole blood specimens using mass spectrometry after a simple dilution sample preparation step.

During the sample dilution step, a small volume of whole blood is extracted from a larger whole blood patient specimen after the entire specimen is mixed (vortexed) to create a uniform distribution of cellular components. This mixing step is important because some metals (e.g. Pb) are known to be associated mostly with the red blood cells in the specimen and a uniform distribution of this cellular material must be produced before a small volume extracted from the larger specimen will accurately reflect the average metal concentration of all fractions of the larger specimen. Coagulation is the process in which blood forms solid clots from its cellular components. If steps are not taken to prevent this process from occurring, i.e. addition of anti-coagulant reagents such as EDTA in the blood collection tube prior to blood collection, blood will immediately begin to form clots once leaving the body and entering the tube. These clots prevent the uniform distribution of cellular material in the blood specimen even after rigorous mixing, making a representative sub-sample of the larger specimen unattainable. It is important that prior to or during sample preparation the analyst identify any sample having clots or micro-clots (small clots). Clotted samples are not analyzed by this method due to the inhomogeneity concerns (i.e. all results for the sample are processed as “not reportable”).

Dilution of the blood in the sample preparation step prior to analysis is a simple dilution of 1 part sample + 1 part water + 48 parts diluent. The effects of the chemicals in the diluent are to release metals bound to red blood cells making them available for ionization, reduce ionization suppression by the biological matrix, prevent clogging of the sample introduction system pathways by undissolved biological solids, and allow introduction of internal standards to be utilized in the analysis step. Tetramethylammonium hydroxide (TMAH, 0.4% v/v) and Triton™ X-100 (0.05%) in the sample diluent solubilizes blood components. Triton™ X-100 also helps prevent biological deposits on internal surfaces of the instrument’s sample introduction system and reduce collection of air bubbles in sample transport tubing. Ammonium pyrrolidine dithiocarbamate (APDC) in the sample diluent (0.01%) aids in solubilizing metals released from the biological matrix. Ethanol in the sample diluent (1%) aids solubility of blood components and aids in aerosol generation by reduction of the surface tension of the solution. The internal standards, rhodium, iridium, and tellurium, are at a constant concentration in all blanks, calibrators, QC, and samples. Monitoring the instrument signal ratio of a metal to its internal standard allows correction for instrument noise and drift, and sample-to-sample matrix differences.

Liquid samples are introduced into the mass spectrometer through the inductively coupled plasma (ICP) ionization source. The liquid diluted blood sample is forced through a nebulizer which converts the bulk liquid into small droplets in an argon aerosol. The smaller droplets from the aerosol are selectively passed through the spray chamber by a flowing argon stream into the ICP. By coupling radio-frequency power into flowing argon, plasma is created in which the predominant species are positive argon ions and electrons and has a temperature of 6000-8000 K. The small aerosol droplets pass through a region of the plasma and the thermal energy vaporizes the liquid droplets, atomizes the molecules of the sample and then ionizes the atoms. The ions, along with the argon, enter the mass spectrometer through an interface that separates the ICP (at atmospheric pressure, ~760 torr) from the mass spectrometer (operating at a pressure of 10⁻⁵ torr). The ions first pass through a focusing region, then the dynamic reaction cell (DRC), the quadrupole mass filter, and finally are selectively counted in rapid sequence at the detector allowing individual isotopes of an element to be determined.

Generally, the DRC operates in one of two modes. In 'vented' (or 'standard') mode the cell is not pressurized and ions pass through the cell to the quadrupole mass filter unaffected. In 'DRC' mode, the cell is pressurized with a gas for the purpose of causing collisions and/or reactions between the fill gas and the incoming ions. In general, collisions or reactions with the incoming ions selectively occur to either eliminate an interfering ion, change the ion of interest to a new mass, which is free from interference, or collisions between ions in the beam and the DRC gas can focus the ion beam to the middle of the cell and increase the ion signal. In this method, the instrument is operated in DRC mode when analyzing for manganese, mercury, and selenium. For selenium, the DRC is pressurized with methane gas (CH₄) which reduces the signal from ⁴⁰Ar₂⁺ while allowing the ⁸⁰Se⁺ ions to pass relatively unaffected through the DRC on toward the analytical quadrupole and detector. Manganese and mercury are both measured when the DRC is pressurized with oxygen gas (O₂). They are analyzed at the same flow rate of oxygen to the DRC cell to avoid lengthening analysis time due to pause delays that would be necessary if different gas flows were used for the two analytes. The oxygen reduces the ion signal from several interfering ions (³⁷Cl¹⁸O⁺, ⁴⁰Ar¹⁵N⁺, ³⁸Ar¹⁶O¹H⁺, ⁵⁴Fe¹H⁺) while allowing the ⁵⁵Mn⁺ ion stream to pass relatively unaffected through the DRC on toward the analytical quadrupole and detector. In the case of mercury, collisional focusing of the mercury ions occurs, increasing the observed mercury signal at the detector by approximately a factor of two (2x). The DRC is used to reduce ion signals from polyatomic ions via ion-molecule reaction chemistry [59, 60].

Once ions pass through the DRC cell and electrically selected for passage through the analytical quadrupole, electrical signals resulting from the ions striking the discrete dynode detector are processed into digital information that is used to indicate the intensity of the ions. The intensity of ions detected while aspirating an unknown sample is correlated to an elemental concentration through comparison of the analyte: internal standard signal ratio with that obtained when aspirating calibrators. This method was originally based on the method by Lutz et al [61]. The DRC portions of the method are based on work published by Tanner et al. [59, 60].

2. Limitations of Method; Interfering Substances and Conditions

A. Interferences addressed by this method

i. Reduction of argon dimer ($^{40}\text{Ar}_2^+$) interference on selenium ($^{80}\text{Se}^+$): The reaction cell of the ICP-DRC-MS is used to reduce the $^{40}\text{Ar}_2^+$ polyatomic ion. $^{40}\text{Ar}_2^+$ is formed in the plasma as a result of a reaction between the plasma gas (Ar) and itself. The reaction cell is filled with methane (CH_4) gas which reacts with $^{40}\text{Ar}_2^+$ ions through a charge transfer reaction. The products of the reaction are $^{40}\text{Ar}^+$ (ion at a different mass) and ^{40}Ar (neutral). The background ion signal at m/z 80 is reduced by six orders of magnitude because of this reaction.

ii. Reduction of argon nitride ($^{40}\text{Ar}^{15}\text{N}^+$), argon hydroxide ($^{38}\text{Ar}^{16}\text{O}^1\text{H}^+$) on manganese ($^{55}\text{Mn}^+$): The reaction cell of the ICP-DRC-MS is used to reduce the $^{40}\text{Ar}^{15}\text{N}^+$ and $^{38}\text{Ar}^{16}\text{O}^1\text{H}^+$ polyatomic ions. These ions are formed in the plasma as a result of reactions between the plasma gas (Ar) and atmospheric gases (N_2 , O_2) or the solvent (H_2O). The reaction cell is filled with oxygen (O_2) gas which reacts with $^{40}\text{Ar}^{15}\text{N}^+$ and $^{38}\text{Ar}^{16}\text{O}^1\text{H}^+$ ions through either charge transfer reactions or oxygen transfer reactions. The products of the reactions are either neutral molecules and are not detected (charge transfer), or a new ion with higher mass (oxygen transfer). In either case, attenuation of the background ion signal at m/z 55 occurs.

iii. Reduction of $^{37}\text{Cl}^{18}\text{O}^+$, $^{39}\text{K}^{16}\text{O}^+$, and $^{54}\text{Fe}^1\text{H}^+$ on manganese ($^{55}\text{Mn}^+$): The reaction cell of the ICP-DRC-MS is used to reduce the $^{37}\text{Cl}^{18}\text{O}^+$, $^{39}\text{K}^{16}\text{O}^+$, $^{54}\text{Fe}^1\text{H}^+$ polyatomic ions. These ions are formed in the plasma as a result of reactions between elements present in the blood matrix (Cl, K, and Fe) and the solvent (H_2O). Due to high concentrations of Cl, K, and Fe in the blood matrix, resulting ion signals of $^{37}\text{Cl}^{18}\text{O}^+$, $^{39}\text{K}^{16}\text{O}^+$, $^{54}\text{Fe}^1\text{H}^+$ interfere with the measurement of $^{55}\text{Mn}^+$. The reaction cell is filled with oxygen (O_2) gas which reacts with $^{37}\text{Cl}^{18}\text{O}^+$, $^{39}\text{K}^{16}\text{O}^+$, $^{54}\text{Fe}^1\text{H}^+$ ions through either charge transfer reactions or oxygen transfer reactions. The products of the reactions are either neutral molecules and are not detected (charge transfer), or a new ions with higher mass (oxygen transfer). In either case, attenuation of the background ion signal at m/z 55 occurs.

B. Limitations of method (interferences remaining in method)

i. MoO_2 interference on tellurium (^{130}Te): Molybdenum will combine with oxygen to form a polyatomic ion, $^{98}\text{Mo}^{16}\text{O}_2^+$. At the DRC conditions with oxygen gas used in this method for Hg analysis, $^{98}\text{Mo}^{16}\text{O}_2^+$ will interfere with the measurement of the internal standard for Hg, $^{130}\text{Te}^+$. Increased signal at m/z 130 (due to measuring both $^{130}\text{Te}^+$ and $^{98}\text{Mo}^{16}\text{O}_2^+$) results in an erroneously low net intensity for Hg (net intensity = measured intensity for analyte isotope / measured intensity for internal standard isotope). If this interference occurs during the measurement of the calibrators (i.e. a multi-element calibration stock standard includes high levels of Mo) it can result in a positive bias for observed mercury concentrations as a consequence of a nonlinear calibration curve having an artificially low slope. If this interference occurs during the measurement of an unknown sample, the reduced net intensity observed can result in reporting an erroneously low Hg result. This interference has been verified to be of concern (>5% effect negative bias) at blood molybdenum concentrations greater than 15 $\mu\text{g}/\text{L}$. However, typical levels of molybdenum in whole blood (0.2 – 4.6 $\mu\text{g}/\text{L}$ [62, 63]) are below this. Also, levels of molybdenum in whole blood after acute exposures have been observed to be ≤ 15 $\mu\text{g}/\text{L}$ [63]. Molybdenum concentrations below 5 $\mu\text{g}/\text{mL}$ in stock calibration standard solutions do not produce an observable interference.

3. Procedures for collecting, storing, and handling specimens; criteria for specimen rejection; specimen accountability and tracking

A. Procedures for collecting, storing, and handling specimens

Specimen handling conditions, special requirements, and procedures for collection and transport are discussed in the Division of Laboratory Science's (DLS) Policies and Procedures Manual[64]. In general:

- i. No fasting or special diets are required before collection of blood
- ii. Specimen type – whole blood
- iii. Optimal amount of specimen is 1+ mL. Request a minimum volume of 0.25 mL. Volume for one analytical measurement is 0.05 mL.
- iv. Verify sample collection devices and containers are free of significant contamination (“pre-screened”) before use.
- v. Draw the blood through a stainless steel needle into a pre-screened vacutainer.
- vi. Do not freeze blood collection tubes due to the risk the tubes cracking. Transfer blood to plastic, pre-screened cryovials for freezer storage
- vii. Once received, store blood collection tubes at refrigerated temperatures (2–8 °C). Transfer to plastic, pre-screened cryovials before freezing. Specimen stability has been demonstrated for over 4 years at ≤ -70 °C.

B. Criteria for specimen rejection

The criteria for an unacceptable specimen include

- i. Contamination: Improper collection procedures, collection devices, or sample handling can contaminate the blood through contact with dust, dirt, etc. Manganese is present in the general environment, found often in combination with iron, and is present in many alloys (especially stainless steel).
- ii. Low Volume: Request a minimum volume of 0.25 mL. Volume for one analytical measurement is 0.05 mL.
- iii. In all cases, request a second blood specimen.

C. Transfer or referral of specimens; procedures for specimen accountability and tracking

Location, status, and final disposition of the specimens will be tracked at least by paper document in the “Study Folder” (created before analysts receive the samples). Apart from this specimen tracking form, this folder will also contain the paper print outs of results from analysis of the specimens. Maintain records for a minimum of 3 years. Use only numerical identifiers for samples within the laboratory (e.g., case ID numbers) in order to safeguard confidentiality. Access to personal identifiers for samples will be limited to the medical supervisor or project coordinator (e.g. non-CDC personnel).

4. Safety precautions

A. General safety

- i. Observe all safety regulations as detailed in the Laboratory Safety Manual and the Chemical Hygiene Plan. Participate in training regarding blood-borne pathogens prior to performing this method.
- ii. Observe Universal Precautions when working with blood.
- iii. Wear appropriate gloves, lab coat, and safety glasses while handling all solutions.
- iv. Take special care when handling and dispensing bases and concentrated acids. Use additional personal protective equipment which protects face, neck, and front of body. If 25% w/w TMAH or concentrated hydrochloric acid comes in contact with any part of the body, quickly wash with copious quantities of water for at least 15 minutes.
- v. Use secondary containment for containers holding biological or corrosive liquids.
- vi. The use of the foot pedal on the benchtop automatic pipette is recommended because it reduces analyst contact with work surfaces that have been in contact with blood and also keeps the analyst's hands free to hold the specimen cups and autosampler tubes and to wipe off the tip of benchtop automatic pipette.
- vii. There are many potential hazards on an operating ICP-MS instrument including ultraviolet radiation, high voltages, radio-frequency radiation, and high temperatures. This information is detailed in the ICP-MS System Safety Manual.
- viii. Transport and store compressed gas cylinders with proper securing harnesses. For compressed oxygen gas, use regulators which are oil-free.
- ix. Wipe down all work surfaces at the end of the day with disinfectant. Disinfectant may be either daily remake of diluted bleach (1 part household bleach containing 5.25% sodium hypochlorite + 9 parts water) or an equivalent disinfectant

B. Waste disposal

- i. Autoclaving: All diluted biological specimens, original biological specimens being disposed, or consumables which come into contact with biological specimens (even diluted or aerosolized). Use sharps containers or special autoclave pans for broken glass / quartz or items which puncture autoclave bags (e.g. pipette tips).
- ii. Other liquid waste
 - 1) Waste discarded down sink: Do not discard solutions at the sink having a pH lower than 5.0 or higher than 11.5 (limits defined by Dekalb County, GA). Inactivate biological compounds and cellular constituents in mixed chemical and biological waste, such as the waste carboy of the ICP-MS, by adding an approved disinfectant (e.g. Lysol I.C. or equivalent) prior to drain disposal. Flush the sink with copious amounts of water.
 - 2) Waste to be picked up by CDC hazardous waste program: Submit request for hazardous waste removal of all other liquid waste generated in the CDC laboratory for this method.

5. Instrument & material sources

A. Sources for ICP-MS instrumentation

- i. ICP-MS: Inductively Coupled Plasma Mass Spectrometer with Dynamic Reaction Cell Technology (ELAN® DRC II) (PerkinElmer Norwalk, CT, www.perkinelmer.com).
- ii. Recirculating chiller / heat exchanger for ICP-MS: Refrigerated chiller (PolyScience 6105PE) or heat exchanger (PolyScience 3370) (PerkinElmer Norwalk, CT, www.perkinelmer.com).
- iii. Autosampler: ESI SC4-DX autosampler (Elemental Scientific Inc., Omaha, NE) or equivalent.
- iv. Computer: Computer controller provided or recommended by ICP-MS manufacturer is recommended to ensure proper communication between computer and ICP-MS. Recommend 1-2 Gb RAM and secondary internal hard disk for nightly backups (if network backups are not possible).
- v. FAST sample introduction system (optional): Standard peristaltic pump on ICP-MS replaced by DXi-FAST micro-peristaltic pump / FAST actuator and valve combination unit. Like part # DXI-54-P4-F6. If DXi-FAST upgrade on ICP-MS is not used, a separate FAST actuator (built-in option on ESI SC4-DX autosampler or stand-alone FAST actuator) will be necessary to complete the FAST sample introduction system.

B. Sources for ICP-MS parts & consumables

NOTE: The minimum number of spares recommended before reordering (if owning one instrument) are listed as “# Spares = X amount” in the descriptions below.

- i. Adapter, PEEK: Securely connects 1.6mm O.D. PFA tubing to 0.03” I.D. peristaltic tubing. Composed of three PEEK parts.
 - 1) Female nut: for 1.6mm O.D. (1/16”) tubing. Like part P-420 (Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).
 - 2) PEEK ferrule: Like part P-260x (10pk SuperFlangeless ferrule, Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).
 - 3) Conical Adapter Body: Like part P-692 (Upchurch Scientific, Oak Harbor, WA, www.upchurch.com).
- ii. Bottles (for rinse solution): Four liter screw-cap polypropylene container with built-in luer connections (2) designed for use with FAST sample introduction system (like catalog# SC-0305-1, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- iii. Carboy and cap assembly for waste collection: 10-15 L, polypropylene wide-mouth carboy (100 mm neck size) with handles and no spigot (Like part #7BE-25126, Lab Safety Supply, Janesville, WI, www.lss.com) with cap assembly like part # N0690271 (PerkinElmer, Norwalk, CT, www.perkinelmer.com) with tubing connections built into the cap for addition of liquid waste.
- iv. Coolant, for polyscience chiller or heat exchanger: Only PerkinElmer part # WE01-6558 (PerkinElmer Norwalk, CT, www.perkinelmer.com) is approved for use by PerkinElmer. # Spares = 6.
- v. Cones: Platinum or Nickel cones have been used. Platinum cones are more expensive, but will last longer, can be refurbished, and will frequently yield higher sensitivity.
 - 1) Sampler (nickel/platinum): PerkinElmer part # WE021140/WE027802 (PerkinElmer Norwalk, CT, www.perkinelmer.com) or cross-referenced part number manufactured by Spectron Inc. (Ventura, CA, www.spectronus.com) or Glass Expansion (Pocasset, MA, www.geicp.com). # Spares = 4.

- 2) Skimmer (nickel/platinum): PerkinElmer part # WE021137/WE027803 (PerkinElmer Norwalk, CT, www.perkinelmer.com) or cross-referenced part number manufactured by Spectron Inc. (Ventura, CA, www.spectronus.com) or Glass Expansion (Pocasset, MA, www.geicp.com). # Spares = 4.
- vi. Connector (for tubing): Use to connect 1/8" I.D. PVC tubing to 0.125" I.D peristaltic pump tubing. Use part # 3140715 (PerkinElmer Norwalk, CT, www.perkinelmer.com) or equivalent. # Spares = 4.
- vii. Detector, electron multiplier: Like part # N8125001 (PerkinElmer Norwalk, CT, www.perkinelmer.com). # Spares = 1.
- viii. FAST accessories, Valve: CTFE High-flow valve head for SC-FAST (uses ¼-28 fittings). Like part # SC-0599-1010 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- ix. FAST accessories, Stator: CTFE Stator for 6 port SC-FAST high flow valve (¼-28 fittings). Like part # SC-0599-1010-01 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- x. FAST accessories, Rotor: Composite rotor for 6 port SC-FAST high flow valve (¼-28 fittings). Like part # SC-0599-1010-05 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xi. FAST accessories, Sample Loop: 1 mL Teflon, white connector-nuts for high flow valve head (¼-28 fittings). Like part # SC-0315-10 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xii. FAST accessories, Probe, Autosampler: Teflon, carbon fiber support, 0.8mm i.d., blue marker, 1/4-28 fittings. Like part number SC-5037-3751 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.
- xiii. FAST accessories, Probe, Carrier Solution: Teflon, carbon fiber support, 0.5mm i.d., orange marker, 1/4-28 fittings. Like part number SC-5037-3501 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.
- xiv. FAST accessories, Tubing, FAST vacuum: Vacuum line for SC-FAST high flow valve, connects to port #6, black nut for connection to valve head, natural brown color nut on other end for connection to SC autosampler vacuum port. Like part # SC-0321 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xv. FAST accessories, Tubing and nut for carrier solution: 0.5 mm i.d. Teflon tubing (orange marker) with red ¼-28 male nut. Connects to high flow FAST valve head, port #2. Like part # SC-0316-0500 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xvi. Hose, connection to chiller: Push on hose. I.D. = ½", O.D. = ¾". Use part # PB-8 (per inch, Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Do not normally need spare hose (unless moving instrument into a new location).
- xvii. Hose, exhaust of ICP-MS: Available as part of ICP-MS installation kit from Perkin Elmer (PerkinElmer Norwalk, CT, www.perkinelmer.com). Available direct from manufacturer as part # S-LP-10 air connector (Thermaflex, Abbeville, SC, www.thermaflex.net), or equivalent. # Spares = 10 feet of 4" diameter and 10 feet of 6" diameter hose.
- xviii. Injector, quartz with ball joint: I.D. = 2.0 mm. PerkinElmer part # WE023948 (PerkinElmer Norwalk, CT, www.perkinelmer.com). Available direct from manufacturer as part # 400-30 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com) or from various distributors. # Spares = 2.
- xix. Ion lens: PerkinElmer part # WE018034 (PerkinElmer Norwalk, CT, www.perkinelmer.com). # Spares = 3.

xx. Nebulizer: PolyPro-ST micro flow polypropylene nebulizer with external 1/4-28 threaded connector for liquid delivery, low pressure version or equivalent. Like part # ES-4040-7010 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 1. Different nebulizers are acceptable, however, the nebulizer gas flow rate, sample flush time, read delay time, loop fill time, loop size, blood sample dilution preparation volume, and sample-to-sample carry-over must be evaluated and optimized.

1) Gas connection:

a) Teflon tubing: 4mm o.d., 2.4mm i.d. Teflon tubing (like part # ES-2502, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 1.

b) Adapter kit: Plastic adapters to connect Teflon tubing (2.4 mm i.d) to ¼" male Swagelok (compression) port on ICP-DRC-MS. Parts can be obtained as components in a "gas fittings kit for microflow nebulizer", kit like part # ES-2501-1000 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 1.

2) Liquid connection: Connects nebulizer to port #3 of high flow FAST valve head with green, 1/4- 28 fitting. Like part # SC-0317-0250 (Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com). # Spares = 2.

xxi. Nut: (for flanged connections of 1.59mm (1/16") o.d. PFA tubing) Flanged, for 1/16" o.d. tubing, 1/4-28 threads. Use part # P-406x (pkg. of 10, Upchurch Scientific, Oak Harbor, WA, www.upchurch.com) or equivalent. Use a Teflon-coated Viton o-ring with this nut instead of the stainless steel washer that comes with part # P-406x). # Spares = 10.

xxii. Nut and ferrule set, 1/8" Swagelok: Such as part # SS-200-NFSET (stainless steel) or part # B-200-NFSET (brass) (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. For part numbers listed here a quantity of 1 means 1 nut, 1 front ferrule, and 1 back ferrule. Spares = 20.

xxiii. Nut and ferrule set, 1/4" Swagelok: Such as part # SS-400-NFSET (stainless steel) or part # B-400-NFSET (brass) (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. For part numbers listed here a quantity of 1 means 1 nut, 1 front ferrule, and 1 back ferrule. Spares = 20.

xxiv. Oil, roughing pumps: Welch Directorr Gold Available direct from manufacturer as part # 8995G-15 (1 gallon, Welch Rietschle Thomas, Skokie, IL, www.welchvacuum.com), or equivalent. # Spares = 4. Or Fomblin Y14/5 fluid: PerkinElmer part # N8122265 (1 kg bottle, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares =1 per instrument.

xxv. O-ring, for sampler cone: PerkinElmer part # N8120511 (pkg. of 5, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 o-rings.

xxvi. O-ring, for skimmer cone: PerkinElmer part # N8120512 (pkg. of 5, PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 20 o-rings.

xxvii. O-ring, for flanged connections of 1.59mm (1/16") o.d. PFA tubing: Teflon-coated Viton o-ring, i.d. = 1/16", thickness = 1/16", o.d. = 3/16". Such as part # V75-003 (O-rings West, Seattle, WA, www.oringswest.com) or equivalent. # Spares = 20.

xxviii. O-ring, for injector support: Internal o-rings: ID = ¼", OD = 3/8", thickness = 1/16". Need 2 o-rings per injector support setup. PerkinElmer part # N8122008 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent (such as part # V75-010, O-rings West, Seattle, WA, www.oringswest.com). # Spares = 20.

- xxix. O-ring, for injector support: External o-rings: ID = 3/8", OD = 1/2", thickness = 1/16". Need 2 o-rings for each injector support setup. PerkinElmer part # N8122009 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent (such as part # V75-012, O-rings West, Seattle, WA, www.oringswest.com). # Spares = 20.
- xxx. O-ring, for inside nebulizer port on standard PerkinElmer cyclonic quartz spray chamber for the ELAN: Such as part # 120-56 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com). Additional o-rings can sometimes be obtained free of charge or at reduced price when acquired while purchasing spray chambers. # Spares = 20.
- xxxi. O-ring, for inside of bayonet torch mount: Part # WE017284 (PerkinElmer, Shelton, CT, www.perkinelmer.com). Do not substitute. The PerkinElmer o-ring is specially metal impregnated to minimize RF leakage through the torch mount. # Spares = 2.
- xxxii. Photon stop: PerkinElmer part # WE018278 (PerkinElmer, Shelton, CT, www.perkinelmer.com). # Spares = 1.
- xxxiii. Plugs, for quick change for roughing pump oil: These plugs will only work on the Varian roughing pumps which come standard on ELAN DRC II ICPMS instruments. These plugs will not fit the Leybold pumps which come standard on the ELAN DRC Plus instruments. Part # W1011013 (PerkinElmer, Shelton, CT, www.perkinelmer.com). No spares typically needed.
- xxxiv. Spray chamber, quartz concentric: PerkinElmer part # WE025221 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. Available direct from manufacturer as part # 400-20 (Precision Glass Blowing, Centennial, CO, www.precisionglassblowing.com) or from various distributors. # Spares = 2.
- xxxv. Torch, quartz: PerkinElmer part # N812-2006 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # New Spares = 2.
- xxxvi. Tubing and adapter, for SC autosampler rinse station drain: Tygon tubing and adapter to attach to back of SC autosampler for draining rinse station waste (like part # SC-0303-002, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xxxvii. Tubing and adapters, for SC autosampler rinse station filling: Teflon tubing and adapters (to attach to back of SC autosampler for filling rinse stations and to attach to rinse containers). Like part # SC-0302-0500, Elemental Scientific Inc., Omaha, NE., www.elementalscientific.com).
- xxxviii. Tubing, for main argon delivery to instrument: I.D. = 1/8", O.D. = 1/4". Like part # C-06500-02 (pkg. of 100ft, polypropylene, Fisher Scientific International, Hampton, NH, www.fishersci.com) or equivalent. # Spares = 50 ft.
- xxxix. Tubing, PFA: I.D. = 0.5 mm, O.D. = 1.59 mm (1/16"). Used to transfer liquid between rinse solution jug and peristaltic pump tubing. The Perfluoroalkoxy (PFA) copolymer is a form of Teflon®. Like part # 1548 (20ft length, Upchurch Scientific, Oak Harbor, WA, www.upchurch.com) or equivalent. # Spares = 20ft.
- xl. Tubing, peristaltic, for carrier solution (ESI FAST system): use either
- 1) Standard PVC, 2-stop (black / black) peristaltic pump tubing, i.d. = 0.76 mm. PerkinElmer part # 09908587 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 6 packs of 12 tubes.
 - 2) Standard PVC, 3-stop (black/ black/black) peristaltic pump tubing, i.d. 0.76 mm. Spectron part # SC0056 (Spectron, Ventura, CA, www.spectronus.com) or equivalent. #Spares = 6 packs of 12 tubes. Use this type of tubing with ESI DXi micro-peristaltic pump.
- xli. Tubing, peristaltic, for spray chamber drain: use either

1) Standard PVC, 2-stop (black / white) peristaltic pump tubing, i.d. = 3.18 mm or equivalent. PerkinElmer part # N812-2012 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 6 packs of 12 tubes.

2) Standard Santoprene, 3-stop (grey/ grey/ grey) peristaltic pump tubing, i.d. 1.30 mm. Spectron part # SC0311 (Spectron, Ventura, CA, www.spectronus.com) or equivalent. #Spares = 6 packs of 12 tubes. Use this type of tubing with ESI DXi micro-peristaltic pump.

xlii. Tubing, PVC, i.d. = 1/8", o.d. = 3/16": Used to transfer liquid between spray chamber waste port and peristaltic pump and between peristaltic pump and liquid waste jug. Like part # 14-169-7A (pkg. of 50 ft, Fisher Scientific International, Hampton, NH, www.fishersci.com) or equivalent. # Spares = 20ft.

xliii. Tubing, Stainless Steel, o.d. = 1/8", wall thickness = 0.028": Used to connect gas cylinders to NexIONUCT gas ports. Like part # SS-T2-S-028-20 (20ft, Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 20 ft.

xliv. Tubing, Teflon, corrugated, 1/4" o.d.: Connects to the auxiliary and plasma gas side-arms of the torch. Part # WE015903 (PerkinElmer, Shelton, CT, www.perkinelmer.com) or equivalent. # Spares = 2.

xlv. Tubing, vinyl, for argon delivery to nebulizer: Vinyl Tubing, 1/8" ID x 1/4" OD. Like part # EW-06405-02 (Cole Parmer, Vernon Hills, Illinois, www.coleparmer.com) or equivalent. # Spares = 10 ft.

xlvi. Union elbow, PTFE 1/4" Swagelok (ELAN bayonet mount): Connects argon tubing to torch auxiliary gas sidearm on bayonet mount NEXION ICP-MS instruments. Like part # T-400-9 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 2.

xlvii. Union tee, PTFE, 1/4" Swagelok (ELAN bayonet mount): Connects argon tubing to torch plasma gas sidearm and holds igniter inside torch sidearm on bayonet mount NEXION ICP-MS instruments. Like part # T-400-3 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. Spares = 2.

C. Sources for ICP-MS maintenance equipment & supplies

i. Anemometer: Like digital wind-vane anemometer (Model 840032, SPER Scientific LTD., Scottsdale, AZ, www.sperscientific.com) or equivalent. Use to verify adequate exhaust ventilation for ICP-MS (check with hoses fully disconnected).\

ii. Pan, for changing roughing pump oil: Like part # 53216 (United States Plastics Corporation, Lima, OH, www.usplastic.com) or equivalent.

iii. Container, to hold acid baths for glassware: Polypropylene or polyethylene containers with lids (must be large enough for torch, injector, or spray chamber submersion). Available from laboratory or home kitchen supply companies.

iv. Cotton swabs: Any vendor. For cleaning of cones and glassware.

v. Cutter (for 1/8" o.d. metal tubing): Terry tool with 3 replacement wheels. Like part # TT-1008 (Chrom Tech, Inc., Saint Paul, MN, www.chromtech.com) or equivalent.

vi. Getter regeneration Kit: Part # WE023257 (PerkinElmer, Shelton, CT, www.perkinelmer.com). Use this as needed (at least annually) to clean the getter in the pathway of channel A DRC gas.

vii. Magnifying glass: Any 10x + pocket loupe for inspection of cones and other ICP-MS parts. Plastic body is preferred for non-corrosion characteristics. Like part # 5BC-42813 (Lab Safety Supply, Janesville, WI, www.labsafety.com).

viii. Ultrasonic bath: Like ULTRASONIK™ Benchtop Cleaners (NEYTECH, Bloomfield, CT, www.neytech.com) or equivalent.

D. Sources for general laboratory equipment and consumables

i. Bar code scanner: Like Xenon 1902 cordless area-imaging scanner (Honeywell International Inc., Morristown, NJ, www.honeywellaidc.com). For scanning sample IDs during analysis setup. Any bar code scanner capable of reading Code 128 encoding at a 3 mil label density and 2D bar codes can be substituted.

ii. Carboy (for preparation of blood quality control pool and waste jug for ICPMS sample introduction system): Polypropylene 10-L carboy (like catalog # 02-960-20C, Fisher Scientific, Pittsburgh, PA, www.fishersci.com) or equivalent. Carboys with spouts are not advised due to potential for leaking.

iii. Containers for diluent and rinse solution: Two liter Teflon™ containers (like catalog# 02-923-30E, Fisher Scientific, Pittsburgh, PA., www.fishersci.com, or equivalent) and 4L polypropylene jugs (like catalog# 02-960-10A, Fisher Scientific, Pittsburgh, PA, www.fishersci.com, or equivalent) have both been used. Acid rinse before use.

iv. Flask, volumetric:

1) 50mL volumetric flasks (like catalog # 40000050, Thermo Scientific, Fisher Scientific, Pittsburgh, PA., www.fishersci.com). Plastic or glass is acceptable.

2) 100mL volumetric flasks (like catalog # 40000100, Thermo Scientific, Fisher Scientific, Pittsburgh, PA., www.fishersci.com). Plastic or glass is acceptable.

3) 1L volumetric flask (like catalog # 40001000, Thermo Scientific, Fisher Scientific, Pittsburgh, PA., www.fishersci.com). Plastic or glass is acceptable.

4) 2L volumetric flask (like glass flask catalog # 92812G2000, DWK Life Sciences (Kimble), Fisher Scientific, Pittsburgh, PA., www.fishersci.com). Plastic or glass is acceptable for the making of the S₀ intermediate working calibration standard.

v. Gloves: Powder-free, low particulate nitrile (like Best CleaN-DEX™ 100% nitrile gloves, any vendor).

vi. Paper towels: For general lab use, any low-lint paper wipes such as KIMWIPES®EX-L Delicate Task Wipers or KAYDRY®EX-L Delicate Task Wipers (Kimberly-Clark Professional, Atlanta, GA, www.kcprofessional.com). For sensitive applications in cleanrooms, use a wipe designed for cleanrooms such as the Econowipe or Wetwipe (Liberty, East Berlin, CT, www.liberty-ind.com).

vii. Pipette, benchtop automatic (for preparation of blood dilutions to be analyzed): Like the Microlab 625 advanced dual syringe diluter (Hamilton, Reno, NV, <http://www.hamilton.com/>) equipped with a 5.0 mL left syringe, a 250 µL right syringe, a 12 gauge Concorde CT probe dispense tip, the Microlab cable management system and a foot pedal. PEEK valves like part # 60676-01 (left) and part # 60675-01 (right) may reduce metal (e.g. manganese) background in prepared samples. Alternatives are acceptable, including the Micromedic Digiflex™ (Titertek, Huntsville, AL, <http://www.titertek.com/>) equipped with 10.0-mL dispensing syringe, 200 µL sampling syringe, 0.75-mm tip, and foot pedal.

viii. Pipettes (for preparation of intermediate working calibration standards and other reagents): Like Picus® NxT electronic, single-channel pipettes (Sartorius AG, Göttingen, Germany, www.sartorius.com). 5-120 µL (catalog # LH-745041), 10-300 µL (catalog #LH-745061), 50-1000 µL (catalog #LH-745081), 100-5000 µL (catalog #LH-745101). Equivalent pipettes and tips can be substituted. Tubes for sample analysis (for autosampler): Like polypropylene 15-mL conical tubes,

BD Falcon model #352097 (Becton Dickinson Labware, FranklinLakes, NJ, www.bd.com) or equivalent. Clear plastics tend to have lowest trace metal contamination. Blue colored caps have also been used successfully for this method.

ix. Tubes for storage of intermediate stock calibration standards: Like polypropylene 50-mL conical tubes, BD Falcon model #352098 (Becton Dickinson Labware, FranklinLakes, NJ, www.bd.com) or equivalent. For use in storage of intermediate stock calibration standards. Clear plastics tend to have lowest trace metal contamination. Blue colored caps have also been used successfully for this method.

x. Vortexer: Like MV-1 Mini Vortexer (VWR, West Chester, PA, www.vwr.com). Used for vortexing blood specimens before removing an aliquot for analysis. Equivalent item can be substituted.

E. Sources of chemicals, gases, and regulators

i. Acid, hydrochloric acid: Veritas™ double-distilled grade, 30–35% (GFS Chemicals Inc. Columbus, OH, www.gfschemicals.com) or equivalent. This is referred to as “concentrated” hydrochloric acid in this method write-up. For use in preparation of intermediate stock and working calibration standards.

ii. Acid, nitric acid: Veritas™ double-distilled grade, 68-70% (GFS Chemicals Inc. Columbus, OH, www.gfschemicals.com). For use in cleaning any bottles, vials, tubes, and flasks. This is referred to as “concentrated” nitric acid in this method write-up.

iii. Blood, whole (human or bovine): Bags of human blood can be purchased from various sources such as American Red Cross (<http://www.redcross.org>) or Tennessee Blood services (Memphis, TN, <http://tennesseebloodservices.com/>). Request that human blood be screened for infectious diseases such as Hepatitis B and HIV. Source for bovine blood includes the Wisconsin State Laboratory of Hygiene (WSLH, Madison, WI, <http://www.slh.wisc.edu>).

iv. Ethanol (EtOH): USP dehydrated 200 proof (Pharmco Products, Inc.) or equivalent.

v. Ammonium pyrrolidine dithiocarbamate (APDC), laboratory grade (Fisher Scientific, Fairlawn, NJ) or equivalent.

vi. Argon gas (for plasma & nebulizer) and regulator: High purity argon (99.999+% purity, Specialty Gases Southeast, Atlanta, GA, www.sgs-gas.com) for torch and nebulizer. Minimum tank source is a dewar of liquid argon (180-250 L). Bulk tank (1500+L is preferred).

vii. Regulator for argon (at dewar): Stainless steel, single stage, specially cleaned regulator with 3000 psig max inlet, 0–200 outlet pressure range, CGA 580 cylinder connector, and needle valve shutoff on delivery side terminating in a ¼” Swagelok connector. Part number “KPRCGRF415A2/AG10-AR1” (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. # Spares = 1.

viii. Regulator for argon (between bulk tank and PerkinElmer filter regulator): Single Stage 316SS Regulator, with 0-300 psi Inlet Gauge, 0-200 psi Outlet Gauge, Outlet Spring Range, 0-250 psi, ¼” Swagelok Inlet Connection, ¼ turn Shut off Valve on Outlet with ¼” Swagelok Connection and Teflon Seals. Part number KPR1GRF412A20000-AR1 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com) or equivalent. # Spares = 1.

ix. Regulator for argon (filter regulator on back of ICP-MS): Argon regulator filter kit. Catalog number N812-0508 (PerkinElmer, Shelton, CT, www.perkinelmer.com).

x. Disinfectant, for work surfaces: Diluted bleach (1 part household bleach containing 5.25% sodium hypochlorite + 9 parts water), remade daily, or equivalent disinfectant.

- xi. Methane: Methane (Research Grade 5.0, 99.99% or higher purity), for DRC channel A. Typically purchased in cylinder size 200 (part # ME R200, Airgas South, Atlanta, GA, www.airgas.com).
- xii. Regulator for methane: Stainless steel, two stage, specially cleaned regulator with 3000 psig max inlet, 0-25 outlet pressure range, CGA 350 cylinder connector, and needle valve shutoff on delivery side terminating in a ¼" Swagelok connector. Like part number KCYADPF412A2AD10 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent. # Spares = 1.
- xiii. Oxygen: Oxygen ("Research Grade 5.0", 99.999% purity, equivalent, or higher purity) for DRC channel B. Like part # OX R33A (Airgas South, Atlanta, GA, www.airgas.com).
- xiv. Regulator for oxygen: Stainless steel, two stage regulator for use with high purity oxygen (cleaned to be free of all oils). Maximum inlet pressure 3600-5000 psi. Inlet gauge pressure 0-5000 psi (no oil in gauge). Maximum delivery pressure 50–100 psi with a 0-30 psi outlet gauge (no oil in gauge). CGA 540 cylinder connector on inlet side and an angle pattern (90 degree) stainless steel needle valve on the delivery side terminating in a 1/8" stainless steel Swagelok connector. Like part # GEORG/KCYCFR/ORS2/540 (Georgia Valve and Fitting, Atlanta, GA, www.swagelok.com), or equivalent.
- xv. Standard, iridium: Like 1,000 µg/mL, item #CGIR1-1 (Inorganic Ventures, Christiansburg, VA <http://www.inorganicventures.com>). Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- xvi. Standard, multi-element stock calibration standard: Item number SM-2107-042 (High Purity Standards, Charleston, SC, <http://www.hps.net/>). Standard must be traceable to the National Institute for Standards and Technology.
- xvii. Standard, rhodium: Like 1,000 mg/L, item # PLRH3-2Y. (SPEX Industries, Inc., Edison, NJ, www.spexcsp.com). Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- xviii. Standard, single element stock standards for preparation of calibrators and blood quality control pools: National Institute of Standards and Technology (NIST) Standard Reference Materials (SRMs): 3108 (Cd), 3132 (Mn), 3128 (Pb), 3133 (Hg), 3149 (Se). (Gaithersburg, MD, www.nist.gov). Standard must be traceable to the National Institute for Standards and Technology.
- xix. Standard, tellurium: Like 1,000 mg/L, item #CGTE1-1 (Inorganic Ventures, Christiansburg, VA <http://www.inorganicventures.com>). Used as an internal standard in diluent. Standard must be traceable to the National Institute for Standards and Technology.
- xx. Tetramethylammonium hydroxide, 25% w/w, or equivalent: (AlfaAesar, 30 Bond St., Ward Hill, MA 01835).
- xxi. Triton™ X-100 surfactant: Like "Baker Analyzed" Triton™ X-100 (J.T. Baker Chemical Co., www.jtbaker.com).

6. Preparation of reagents and materials

A. Internal standard intermediate mixture:

i. Purpose: Preparation of single intermediate solution containing all internal standards simplifies the addition of the internal standard(s) into the final diluent solution. This solution can be purchased rather than prepared.

ii. Preparation and storage: To prepare 50 mL of 20 mg/L Rh, Ir, Te in 1% v/v HNO₃:

- 1) If not previously dedicated to this purpose, acid wash a 50 mL volumetric flask (PP, PMP, or Teflon™). For example, with 1% v/v HNO₃ and ≥18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- 2) Partially fill the 50 mL volumetric flask with 1% v/v HNO₃ (approximately 25-40 mL).
- 3) Add 1 mL of 1,000 µg/mL Rh standard, 1 mL of 1,000 µg/mL Ir standard, and 1 mL of 1,000 µg/mL Te standard. If initial Rh, Ir, or Te standard concentration is different, adjust volume proportionally.
- 4) Fill to mark (50 mL) with 1% v/v HNO₃ and mix thoroughly.
- 5) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

B. Intermediate Triton™ X-100 solution:

i. Purpose: To ease daily preparation of the diluent and rinse solutions by first preparing an intermediate Triton™ X-100 solution.

ii. Preparation and storage: To prepare 1 L of 20% Triton™ x-100

- 1) If not previously dedicated to this purpose, acid wash a 1 L volumetric flask (PP, PMP, or Teflon™). For example, with 1% v/v HNO₃ and ≥18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- 2) Add 200 mL of Triton™ X-100 to the 1L container that is partially filled with ≥18 Mohm·cm water.
- 3) Fill to 1 L with ≥18 Mohm·cm water and mix until the Triton™ X-100 has completely dissolved into solution (overnight). A magnetic stirring plate can be used to assist mixing by adding an acid-washed Teflon™ coated stirring bar to the bottle.
- 4) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

C. Sample diluent and carrier

i. Purpose: This solution will be used in the preparation of all samples and calibrators during the dilution process prior to analysis. Make all samples, standards, blanks, QC, etc. . . in a run from the same diluent solution so that the concentration of the internal standards will be the same among all calibrators and samples in the run. When using a flow-injection component in the sample introduction system (i.e. the Elemental Scientific SC4-FAST autosampler), use the same solution for the the 'carrier' and sample diluent. The diluent is an aqueous solution of 5 µg/L internal standard mixture (Rh, Ir, Te), in 0.4% v/v tetramethyl ammonia hydroxide (TMAH), 1% ethanol, 0.01% APDC, and 0.05% v/v Triton™ X-100. Larger volumes of these solutions can be prepared by adjusting component volumes proportionally.

ii. Preparation and storage: To prepare 2L of 5 µg/L Rh, Ir and Te, 0.01% APDC in 0.4% v/v TMAH, 1% ethanol, and 0.05% v/v Triton™ X-100:

- 1) If not previously dedicated to this purpose, acid wash a 2L container (PP, PMP, or Teflon™). For example, with 1% v/v HNO₃ and ≥18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.

- 2) Partially fill the 2 L container with ≥ 18 Mohm·cm water.
- 3) Add 0.2 g of APDC, 8 mL of 25% v/v TMAH, 20 mL of ethanol, and 5 mL of 20% Triton™ X-100.
- 4) Dilute to volume (2L) with ≥ 18 Mohm·cm water.
- 5) Spike 500 μ L of 20 mg/L Rh, Ir, Te to the final diluent.
- 6) Invert bottle a few times to insure thorough mixing. Allow to sit for several hours or overnight before using.
- 7) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

D. ICP-MS rinse solution

i. Purpose: The rinse solution used in this method is an aqueous solution of 0.01% APDC in 0.4% v/v TMAH, 1% ethanol, and 0.05% v/v Triton™ X-100. This solution will be pumped through the autosampler rinse station, probe, and sample loop between sample analyses to prevent carry-over of analytes from one sample measurement to the next.

ii. Preparation and storage: To Prepare 4 L of 0.01% APDC in 0.4% v/v TMAH, 1% ethanol, and 0.05% v/v Triton™ X-100:

- 1) If not previously dedicated to this purpose, acid wash a 4L container (PP, PMP, or Teflon™). For example, with 1% v/v HNO₃ and ≥ 18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
- 2) Partially fill the 4 L bottle with ≥ 18 Mohm·cm water (approximately 2-3 L). Use of volumetric flask is not required.
- 3) Add 0.4 g of APDC
- 4) Add 16 mL of 25% v/v TMAH
- 5) Add 40 mL of ethanol,
- 6) Add 10mL of 20% Triton™ X-100, (See Section 6.b for details on preparation)
- 7) Fill to 4 L using ≥ 18 Mohm·cm water.
- 8) Invert bottle a few times to ensure thorough mixing. Allow to sit for several hours or overnight before using.
- 9) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

E. Base blood

i. Purpose: This blood pool material will be mixed with the intermediate working calibration standards just prior to analysis to matrix-match the calibration curve to the blood matrix of the unknown samples.

ii. Preparation and storage: To prepare a mixture of multiple blood sources collected from anonymous donors to approximate an average blood matrix:

- 1) Purchase several bags of whole blood.
- 2) Screen each individual bag of blood for concentration of analytes of interest. See Table 2 of Appendix C for minimum acceptable values
- 3) Once screened, mix the acceptable blood together in a larger container (i.e. acid washed polypropylene (PP), polymethylpentene (PMP), or Teflon™) and stir for 30+ minutes on a large stir plate (acid wash large Teflon™ stir bar before use).
- 4) Store long-term as smaller portions for daily use (e.g. 2 mL cryovials) according to the same storing and handling criteria described in Section 3.

F. Multi-element stock calibration standards

- i. Purpose: This multi-element stock standard will be used to prepare the intermediate working calibration standards.
- ii. Purchase and storage: Whether purchased or prepared in-house, the starting materials must be NIST-traceable. Matrix and concentrations of Pb, Cd, Hg, Mn and Se are listed in Table 3 of Appendix C. Store at ambient temperature and label appropriately. Expiration is as defined by the manufacturer or 1 year from date of opening, whichever comes first.

G. Intermediate calibration standard (S0):

- i. Purpose: This diluent is used to dilute stock and intermediate stock calibration standards, not to prepare working calibrators or blood samples for analysis.
- ii. Preparation and storage: To prepare 2L of 3% v/v HCl:
 - 1) If not previously dedicated to this purpose, acid wash a 2L volumetric flask (PP, PMP, or Teflon™). For example, with 3% v/v HCl and ≥ 18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate.
 - 2) In the 2 L flask, add 1-1.5L ≥ 18 Mohm·cm water.
 - 3) Add 60 mL high purity concentrated HCl.
 - 4) Fill to the mark and mix thoroughly.
 - 5) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

H. Multi-element intermediate stock calibration standard

- i. Purpose: This multi-element intermediate stock standard will be used to prepare the intermediate working calibration standards.
- ii. Preparation and storage: To prepare 3% v/v HCl solutions containing Cd, Pb, Hg, Se, and Mn with concentrations listed in Table 4 of Appendix C:
 - 1) Acid-rinse one 100 mL, PP (or PMP) volumetric flask. For example, with 3% v/v HCl and ≥ 18 Mohm·cm water (at least 3 times each) and verify cleanliness through analysis of rinsate. Mark flask according to intended use. Dedicate to purpose.
 - 2) Partially fill (50-75% full) the 100 mL flask with the 3% v/v HCl diluent prepared in Section 6.e.ii.
 - 3) Using the volume listed in Table 4 of Appendix C, pipette the appropriate volume of the multi-element stock calibration standard solution into the volumetric flask. Dilute to the volumetric mark with the 3% v/v HCl diluent using a pipette for the final drops. Mix each solution thoroughly. Final concentrations are listed in Table 4 of Appendix C.
 - 4) Once mixed, transfer to acid-cleaned, labeled, 50 mL containers (PP, PMP, or Teflon™) for storage.
 - 5) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.

I. Intermediate working calibration standards

i. Purpose: Used each day of analysis to prepare the final working calibrators that will be placed on the autosampler.

ii. Preparation and storage: To prepare 3% v/v HCl solutions containing Cd, Pb, Hg, Se, and Mn with concentrations listed in Table 5 of Appendix C:

- 1) Acid-rinse eight 100 mL, PP (or PMP) volumetric flasks. For example, with 3% v/v HCl and ≥ 18 Mohm-cm water (at least 3 times each) and verify cleanliness through analysis of rinsate. Mark each flask according to intended use. Dedicate to purpose.
- 2) Fill each 100 mL flask 50-75% with the 3% v/v HCl diluent prepared in Section 6.e.ii.
- 3) Using the volumes listed in Table 5 of Appendix C, pipette the appropriate volume of the multi-element intermediate stock calibration standard solutions into each of the volumetric flasks. Dilute each to the volumetric mark with the 3% v/v HCl diluent using a pipette for the final drops. Mix each solution thoroughly. Final concentrations are listed in Table 5 of Appendix C.
- 4) Once mixed, transfer to acid-cleaned, labeled, 50 mL containers (PP, PMP, or Teflon™) for storage.
- 5) Store at ambient temperature and label appropriately. Expiration is 1 year from date of preparation.
- 6) Pour aliquots of each standard into clean 15mL polypropylene tubes and label for daily use.

J. Working calibrators

i. Purpose: The working calibrators will be analyzed in each run to provide a signal-to-concentration response curve for each analyte in the method. The concentration of the analyte of interest in a patient blood sample dilution is determined by comparing the observed signal ratio (element/internal standard) from the dilution of the patient blood sample to the signal ratio response curve from the working calibrators.

ii. Preparation and storage: Dilutions (1:50) of the corresponding eight intermediate working calibration standards with base blood and sample diluent. Mix with base blood and diluent (Section 6.c) using a benchtop automatic pipette to make 1:50 dilutions of the corresponding eight intermediate working calibrators immediately prior to analysis (see Table 8 of Appendix C).

K. Internal quality control materials (“bench” QC)

i. Purpose: Internal (or “bench”) quality control (QC) materials are used to evaluate the accuracy and precision of the analysis process, and to determine if the analytical system is “in control” (is producing results that are acceptably accurate and precise). They are included in the beginning and at the end of each analytical run.

ii. Preparation and storage: To prepare pooled animal or human blood at low-normal and high-normal concentrations: Both purchased or in-house prepared quality control materials are suitable for this purpose if volumes, concentrations meet method requirements and any spikes of elemental levels are traceable to the National Institute for Standards and Technology (NIST).

- 1) Screening blood: Screen bags of blood for analytes of interest concentration before mixing together to make separate base blood pools for each QC material that will be created.
- 2) Keep blood refrigerated whenever possible to minimize microbial growth.
- 3) Because this is only a quick screen of the analyte of interest concentration, the number of replicates in the blood method can be reduced to one in order to reduce analysis time.

- 4) Select blood for the low bench QC pool which has analyte concentrations in the low-normal population range. Select blood for the high and elevated bench QC pools which has analyte concentrations less than some pre-selected target concentration values in the high normal population range. See Table 11 of Appendix C for normal population reference ranges.
- 5) Combining collected blood: The goal for combining samples is to approach an 'average' matrix for each pool.
 - a) Graduate four acid-washed 10 L carboys (PP or PMP) in 0.5 L increments (two will be used for decanting into).
 - b) Combine collected blood samples into two separate acid-washed 10 L carboys (PP or PMP), according to their concentrations, for the low bench and high bench QC pools.
 - c) Mix each blood pool using carboy stirrers and large stir plates. Keep blood refrigerated whenever possible.
- 6) Spiking of blood
 - a) Analyze three samples of each blood pool. Record these results for future recovery calculations.
 - b) Use these results to determine target analyte concentrations possible for the pools
 - c) Calculate the volume of single element standards needed to spike each pool to the desired concentrations. See Table 11 of Appendix C for normal population reference ranges.
 - d) While stirring the pools on large stir plates, spike each pool with calculated volumes of single element standards (all spiking standards used must be traceable to NIST).
 - e) Continue to stir pools overnight after spiking, then reanalyze.
 - f) Repeat steps 4 and 5 until all analytes reach target concentrations keeping track of the total volume of spiking solution added to each blood pool.
- 7) Dispensing and storage of blood
 - a) Container types: Dispense blood into lot screened containers (i.e. – 2 mL polypropylene tubes). If possible, prepare tubes of QC which have only enough volume for one typical run + 1 repeat analysis. This allows for one vial of QC to be used per day of analysis, reducing chances of contamination of QC materials due to multi-day use.
 - b) Labels: Place labels on vials after dispensing and capping if the vials are originally bagged separately from the caps. This minimizes the chance for contamination during the process. Include at least the name of QC pool (text and bar code), date of preparation, and a vial number on the labels.
 - c) Dispensing: Dispensing can be accomplished most easily using a benchtop automatic pipette in continuous cycling dispense mode. Dispense the pools in a clean environment (i.e. a class 100 cleanroom area or hood).
 - d) Allow blood to reach ambient temperature before dispensing (to prevent temperature gradients possibly causing concentration gradients across the large number of vials being dispensed and to prevent condensation problems during labeling of vials).
 - e) Replace the tubing attached to the dispensing syringe (left when looking at front of the benchtop automatic pipette) with a length of clean Teflon™ tubing long enough to reach into the bottom of the 10 L carboy while it is sitting on the stir plate.
 - f) Check cleanliness of the benchtop automatic pipette before use by analyzing 1-2% v/v HNO₃ which has been flushed through the benchtop automatic pipette with a portion of the same solution which has not been through the benchtop automatic pipette.

- g) Approximately one hour before dispensing begins,
- h) With the large stir plate close to the left side of the benchtop automatic pipette, begin stirring the blood pool to be dispensed.
- i) Also during this time, flush the benchtop automatic pipette with blood from the pool to be dispensed. Place the ends of the tubing attached to both the sample and dispensing syringes into the carboy of blood so that blood won't be used up during this process. Be sure to secure both ends of tubing in the carboy with Parafilm so they will not come out during the flushing process.
- j) After dispensing the blood into the vials, cap the vials and label them. Placing labels on vials after capping minimizes the chance for contamination during the process.
- k) Homogeneity test: Check homogeneity of analyte concentrations in pool aliquots.
- l) Storage: Store long-term as smaller portions for daily use (e.g. 2 mL cryovials) according to the same storing and handling criteria described in Section 3.

L. DRC optimization solutions:

i. Purpose: For periodic testing of the DRC cell parameters. Procedure requires at a minimum a blank (i), an analyte solution (ii), a blank with interference (iii), and an analyte and interference containing solution (iv). For Se, only the blank (i), an analyte solution (ii) are needed because the interference on Se is plasma based. Interferences are discussed in Sections 1.B. and 2.A. Interference concentrations can be prepared higher as needed by adjusting the volume of the spikes. Keep interference spike volume small (<0.3 mL) using a high concentration stock solution (i.e. 1000 µg/mL). Analyte concentrations can be made higher if needed for sensitivity reasons by preparing a higher concentration calibrator. If elimination of the interference is difficult to verify, replace the use of blood in these preparations with ultrapure water to minimize trace amounts of the analyte in the preparation. Diluent in this section refers to sample diluent (5 µg/L internal standard mixture (Rh, Ir, Te), 0.4% v/v tetramethyl ammonia hydroxide (TMAH), 1% ethanol, 0.01% APDC, and 0.05% v/v Triton™ X-100 as described in Section 6c.

ii. Preparation and storage: (⁵⁴Fe¹H interference on ⁵⁵Mn):

- 1) Base blood in diluent (1 + 49): In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8 of Appendix C (multiply volumes by 5).
- 2) Base blood in diluent (1 + 49) + 4.5 µg/L Mn: In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8 of Appendix C (multiply volumes by 5).
- 3) Base blood in diluent (1 + 49) + 500 µg/L Fe: In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8 of Appendix C (multiply volumes by 5). Add 0.025 mL of 1000 µg/mL Fe.
- 4) Base blood in diluent (1 + 49) + 4.5 µg/L Mn + 500 µg/L Fe: In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8 of Appendix C (multiply volumes by 5). Add 0.025 mL of 1000 µg/mL Fe.
- 5) Store at ambient temperature and prepare as needed.

iii. Label appropriately (see Section 6.f.i.2), "Store at ambient temperature", preparation date, expiration date one year from preparation date, and preparer's initials.

iv. Preparation and storage: (⁴⁰Ar₂ interference on ⁸⁰Se):

- 1) Base blood in diluent (1 + 49): In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 0 as described in Table 8 of Appendix C (multiply volumes by 5).
- 2) Base blood in diluent (1 + 49) + 90 µg/L Se: In a 50 mL lot screened or acid-washed polypropylene tube, prepare a 50 mL portion of working calibrator 2 as described in Table 8 of Appendix C (multiply volumes by 5).
- 3) Store at ambient temperature and prepare as needed.
- 4) Label appropriately (see Section 6.f.i.2), “Store at ambient temperature”, preparation date, expiration date one year from preparation date, and preparer’s initials.

M. Dual detector calibration (aqueous dilutions of stock standard(s) in 2% v/v nitric acid):

i. Purpose: Use as necessary to perform the dual detector calibration procedure. Lead is typically the only element requiring dual detector calibration (exceeding 1,000,000 cps for the highest calibrators), but add others as is recommended by service representatives.

ii. Preparation & storage: Prepare different volumes, if needed, by adding proportionally larger or smaller volumes of solution constituents. To prepare a total of 50 mL of 200 µg/L Pb in 2% v/v nitric acid:

- 1) Partially fill a 50 mL lot screened or pre-washed polypropylene tube with 10-40 mL 2% v/v nitric acid.
- 2) Add 0.1 mL of 100 µg/mL stock standard for each element desired in the final solution.
- 3) Dilute to the 50 mL mark with 2% v/v nitric acid.
- 4) Label appropriately and store at ambient temperature. Expiration date one year from preparation date.

7. Analytical instrumentation setup

See Section 5 for details on hardware used, including sources. See Table 1 and the Figures in Appendix C for a complete listing of the instrument and method parameters, and software screen shots.

A. Configuration for liquid handling

i. FAST valve setup: See Figure 1 of Appendix C for diagram and Section 5.b “FAST / ESI SC4-DX autosampler accessories” for source information.

- 1) Port 1: sample loop (white nut).
- 2) Port 2: 0.5 mm ID probe (red nut) for carrier solution.
- 3) Port 3: nebulizer line (green nut) for transfer of liquid to nebulizer.
- 4) Port 4: sample loop (white nut).
- 5) Port 5: 0.8 mm ID probe (blue nut) for diluted samples.
- 6) Port 6: vacuum line (black nut).

ii. Carrier solution uptake: Use peristaltic pump to control uptake flow rate of carrier solution to the SC-FAST valve. The carrier probe tubing can be connected directly to the peristaltic pump tubing. The other side of the peristaltic pump tubing connects directly to “carrier in” line with the red nut (see consumables descriptions in Section 5.B).

iii. Spray chamber waste removal: Use the peristaltic pump to control the removal of liquid waste from the spray chamber. The spray chamber drain tubing connects directly to the Santoprene™ peristaltic pump tubing. Connect the other end of the peristaltic pump tubing to 0.5 mm i.d. PFA

tubing. Place the free end of the PFA tubing through the lid of the waste jug (be sure it is secure). Between peristaltic pump tubing and waste container: Connect 1/8" i.d. x 1/4" o.d. PVC tubing to the white / black peristaltic pump tubing using a tubing connector (PerkinElmer item # B3140715). Place the free end of the PVC tubing through the lid of the waste jug (be sure it is secure). Place waste container in a deep secondary containment tray in case of overflow.

iv. Between spray chamber and peristaltic tubing:

- 1) Spray chambers with threaded connection: Use vendor-supplied threaded connector on base of chamber, connecting tubing directly to peristaltic pump tubing through a PEEK adapter or directly.
- 2) Spray chambers without threaded connection: Use of specialized push-on connectors available from various vendors (like UFT-075 from Glass Expansion, Pocasset, MA) are preferred for safety reasons to direct connection of PVC tubing (e.g. 1/8" i.d. x 1/4" o.d.).

v. Rinse solution for autosampler:

- 1) Rinse solution jug: Leave one of the caps on the top of the rinse jug loose to allow air venting into the jug as liquid is removed. Otherwise the jug will collapse on itself as the liquid is removed and a vacuum is created inside. Use secondary containment tray.
- 2) Rinse solution uptake to autosampler rinse station: Use tubing of different lengths and inner diameters between the rinse solution container and the autosampler rinse station to control uptake rate of rinse solution. These can be obtained from the autosampler manufacturer, their distributors, or custom built in the lab. Optimize these factors along with fill time in the software so that waste of rinse solution is minimized and rinse station does not go empty.
- 3) Autosampler rinse station waste removal: Gravity drain of waste to the waste container is sufficient. Use minimum drain tubing to make this connection. If this tube is too long, the rinse station will not drain properly.

B. Gas delivery and regulation

i. ICP-MS modifications:

- 1) Plastic tubing between mass flow controllers and dynamic reaction cell have been replaced with stainless steel. Stainless steel tubing is preferred between the reaction gas cylinder / regulator and the back of the ICP-MS instrument.
- 2) A second mass flow controller is used on channel B that does not send the DRC gas through a 'getter'.

ii. Argon gas: Used for various ICP-MS functions including plasma and nebulizer.

- 1) Regulator for argon source (if a dewar): Set delivery pressure of this regulator at least 10 psi higher than the delivery pressure of the step-down regulator to allow for pressure drop across tubing that stretches to the instrument.
- 2) Step down regulator (if source of argon is a bulk tank): Place this single stage regulator in the lab so that incoming argon pressure can be monitored and adjusted. Set delivery pressure to 10 psig above the delivery pressure of the filter regulator on the ICP-MS.
- 3) Filter regulator at ICP-MS: Single stage "argon regulator filter kit" supplied with the ICP-DRC-MS. Set the delivery pressure depending on the instrument setup:
- 4) ELAN with a 0-60psi gauge on the filter regulator: 52±1 psi when plasma is running (need 0-150 psi regulator if using a PolyPro or PFA nebulizer made by Elemental Scientific Inc).
- 5) ELAN with a 0-150psi gauge on the filter regulator: 90-100 psi when plasma is running.

iii. Methane gas: Used for dynamic reaction cell interference removal from selenium isotopes.

- 1) Connect to DRC channel A
- 2) Set the delivery pressure of regulator to 5-7 psig when gas is flowing. See section 5.e for part numbers and details.

iv. Oxygen gas: Used for dynamic reaction cell interference removal from manganese isotopes.

- 1) Connect to DRC channel B.
- 2) Set the delivery pressure of regulator to 5-7 psig when gas is flowing. See Section 5.e for part numbers and details.

v. Chiller / heat exchanger: If using refrigerated chiller, set temperature control to approximately 18 °C.

8. The run: quality, execution, evaluation, and reporting

A. Quality: bench QC, reference materials, and calibration verification:

i. Bench "QC": Analysis of bench QC permits assessment of methodological imprecision, determination of whether the analytical system is 'in control' during the run, and assessment of time-associated trends. Before QC materials can be used in the QC process, they must be characterized by at least twenty (20) analytical runs to determine appropriate QC parameters. Bench QC pool analyte concentrations in this method span the analyte concentration range of the calibrators including "low-normal" ('Low QC'), "high-normal" ('High QC'), and "above-normal" ('Elevated QC') concentrations.

In each analytical run the analyst will test each of the three bench QC samples two times, subjecting them to the complete analytical process. Bench QC pool samples are analyzed first in the run after the calibrators but before any patient samples are analyzed. This permits making judgments on calibration linearity and blank levels prior to analysis of patient samples. The second analysis of the bench QC pools is done after analysis of all patient samples in the run (typically 40-50 patient samples total when analyzing for all elements in the method) to ensure analytical performance has not degraded across the time of the run. If more patient samples are analyzed on the same calibration curve after the second run of the bench QC, all bench QC must be reanalyzed before and after the additional samples. For example, the schemes shown in Table 6 of Appendix C are both acceptable ways to analyze multiple consecutive "runs".

ii. Reference materials: Use standard reference materials (SRM) from the National Institute of Standards and Technology (NIST) (i.e. SRM 955c Levels 1-4) to verify method accuracy. Use previously characterized samples from proficiency testing program or commercially-produced reference materials when NIST SRMs are unavailable.

iii. Calibration verification: The test system is calibrated as part of each analytical run with NIST-traceable calibrators. These calibrators, along with the QCs and blanks, are used to verify that the test system is performing properly.

B. Execution: perform, evaluate and report a run

i. Starting the equipment for a run:

- 1) Power on the computer, printer, and autosampler, and instrument computer controller.
- 2) Peristaltic pump: Set proper tension on peristaltic pump tubing.
- 3) Software: Start software for the ICP-MS and autosampler control.

- 4) Daily pre-ignition maintenance checks: Perform and document daily maintenance checks (e.g., Ar supply pressure, interface components cleanliness and positioning, interface pump oil condition, vacuum pressure, etc.).
- 5) Place probe in adequate volume of carrier or rinse solution: If using an ESI FAST, manually place carrier probe into carrier solution. If not, send the autosampler probe to a rinse solution (e.g. autosampler rinse station).
- 6) Start the plasma
- 7) Start the peristaltic pump: Start the pump running slowly, making sure that the rotational direction is correct for the way the tubing is set up.
- 8) Warm-up time: Allow warm-up time suggested by the manufacturer for the ICP-MS (e.g. RF generator) after igniting the plasma. There will be another warm-up time (or “stability time”) for the DRC later in this procedure.
- 9) Daily performance check: Perform and document a daily performance check and any optimizations necessary. Save new parameters to the “default.tun” and “default.dac” files.
- 10) DRC stability time: Best analyte-to-internal standard ratio stability is typically observed after 1-1.5 hours of analysis of diluted blood samples using the DRC mode method (~15 measurements of the 5 element panel can be made in 1 hour). Prepare 50mL+ of a calibration standard (e.g. standard 2) to be analyzed repeatedly before the beginning of the run to achieve a stable analyte-to-internal standard ratio. Time to reach stability is instrument-specific and learned from performance of runs. See Table 7 of Appendix C for example of setup in the Samples / Batch window and Table 8 of Appendix C for details of making a working calibrator.
- 11) (Optional) Ready the instrument for priority samples: If priority samples are expected for analysis, the plasma can be started well in advance and left running to eliminate the need for an initial instrument warm-up period and / or a DRC stabilization period as long as appropriate planning is made for sufficient solution supply and waste collection. Analysis of conditioning samples (diluted blood matrix) can also be scheduled to occur at roughly a predetermined time. Accomplish this by setting up multiple sample analyses with extended rinse times (e.g. one 5 element analysis with a 1500s rinse time will take approximately 30 minutes to complete). Initial samples would be non-matrix, while final samples would be diluted matrix for conditioning. If running a DRC-only method during these scheduled analyses, the ICP-MS will remain in DRC-mode for approximately 45 minutes without depressurizing the cell.

ii. Software setup for analysis:

- 1) Workspace (files & folders): Verify & set up the correct files and data directories for your analysis (See Table 1 of Appendix C for defaults).
- 2) Samples / batch window: Update the software to reflect the current sample set. Use a bar code scanner to input data whenever possible. See Table 1 of Appendix C for times and speeds.
- 3) Blood vs. aqueous method files:
 - a) The difference: There are two method files for this one method (see Table 1 of Appendix C). It is necessary to use both to accomplish each run because the current PerkinElmer software will not allow for more than one blank per method file. The ONLY DIFFERENCE between these two files is on the Sampling tab where one lists the autosampler positions of the blood blank and blood-based calibrators (the “bldblk” method file) and the other lists the autosampler position of the aqueous blank (the “aqblk” method file).

b) Use: The ONLY TIME when it matters which of these files is used is when the measurement action includes “Run blank” or “Run standards”. When the measurement action is only ‘run sample’, it does not matter whether the “bldblk” or “aqblk” method file is used. Analysts typically follow the pattern below, however, for the sake of consistency and as a reminder of which blank must be used for which type of sample. See Table 7 of Appendix C.

(i) The “bldblk” method file: Use to analyze the initial blood blank (blank for the calibration curve), the blood calibrators, and the blood blank checks at the very beginning of the run. The blood blank method defines the autosampler location of the blood blank and the blood calibration standards.

(ii) The “aqblk” method file must be used to analyze all QC materials and patient samples. The aqueous blank method defines the aqueous blank in autosampler location.

iii. Preparation of samples for analysis (See Table 8 of Appendix C)

- 1) Thaw blood samples; allow them to reach ambient temperature.
- 2) Prepare the following solutions into pre-labeled containers using the benchtop automatic pipette or other volumetric sample transfer device. See Table 8 of Appendix C for a summary.
 - a) Aqueous Blank: Prepare a minimum of two aqueous blanks. One will be the actual aqueous blank and the other will be a backup (“Aqueous Blank Check”) in case the original aqueous blank is unusable.
 - b) Calibrators: Prepare the working calibrators (S0-S8). Prepare S0 in triplicate to use for both the zero calibrator and blood blank checks after the calibrators.
 - c) Patient & QC Samples: Before taking an aliquot for analysis, homogenize the sample thoroughly.
- 3) After preparation, mix and cover. Place prepared dilutions on the autosampler of the ICP-MS in the order corresponding to the sequence setup in the ICP-MS software.
- 4) Ambient temperature storage is acceptable for the original samples during the work day.
- 5) Samples must be analyzed within 24 hours of preparation to obtain valid results for selenium. The method has been validated to produce valid results for other Pb, Cd, Hg, and Mn even 48 hrs after sample preparation. See critical parameter test results in Appendices A and B for details.

iv. Start the analysis using the ICP-MS software

v. Monitor the analysis in real-time as much as possible. If necessary, leave the run to complete itself unattended as long as appropriate planning is made for either overnight operation or Auto Stop (see below).

vi. Monitor the analysis for the following:

- 1) Verify proper operation of the instrument (proper loop filling, sample reaching nebulizer in correct timing, autosampler arm moving properly, etc...).
- 2) Verify that background signal from instrument and reagents are low. Helpful checks when diagnosing high background problems include:
 - a) Water to be used in Aq Blank Checks and dilutions.
 - b) Diluent before and after being flushed through the benchtop automatic pipette.
 - c) If contamination is observed from the pipette, flush the pipette with ≥ 500 mL of nitric acid solution ($\leq 5\%$ v/v HNO_3) and retest.
 - d) Comparison with other instruments.

- 3) Verify analyte / internal standard ratio stability. The net intensity (analyte / internal standard ratio) of the measurements made while stabilizing the DRC can be evaluated to determine the readiness of the system to begin analysis. Continual trending in this ratio indicates that unwanted instrument drift will occur within the run.
- 4) Verify calibration curves meet R2 requirements (minimum of 0.98, typically 0.99 to 1.000).
- 5) Verify bench QC results are within the acceptable limits. If an analyte result for the beginning QC material(s) falls outside of the $\pm 3SD$ limits, then the following steps are recommended:
 - a) Evaluate the blank results.
 - b) Evaluate the reproducibility of the 3 replicates within the measurements.
 - c) Evaluate the consistency of the internal standard across the measurements (esp. the calibrators).
 - d) Evaluate calibration curves. If a particular calibration standard is obviously in error, it can be re-analyzed as a sample (old or new dilution) and incorporated into the curve through data reprocessing as a calibrator. As a last resort, a single calibration point per analyte between or including S2 and S7 can be removed from the curve (Do not drop S1 or S8). Follow up problems with calibrators with appropriate corrective actions (e.g. re-preparation of intermediate working calibration standards or troubleshooting instrument parameters).
 - e) Prepare a fresh dilution of the failing QC material (same vial) and reanalyze it to see if the QC dilution was not properly made.
 - f) Prepare a fresh dilution of the failing QC material (unused vial) and analyze it to see if the QC vial had become compromised.
 - g) Prepare and analyze new working calibrators.
 - h) Test a different preparation of intermediate working calibration standards.
 - i) If these steps do not result in correction of the out-of-control values for QC materials, consult the supervisor for other appropriate corrective actions.
- 6) Verify good precision among replicates of each measurement.
- 7) Verify consistent measured intensities of the internal standards. Some sample-to-sample variations are to be expected, however, intensities drifting continuously in one direction resulting in failing results for ending QC indicate the instrument needs additional pre-conditioning before the run or environmental conditions are changing too much around the instrument.
- 8) Verify elevated patient results. Refer to Figure 17 of Appendix C for flowchart.
 - a) Confirming an elevated concentration: Repeat for confirmation any sample having a concentration greater than the 1UB threshold. See Table 10 of Appendix C.
 - b) Dilution of a sample to within the calibration range: Repeat in duplicate with extra dilution any sample having a concentration greater than the highest calibrator to bring the observed result within the concentration range of the calibrators. See Table 9 of Appendix B for high calibrator concentrations and validated extra dilutions.
 - c) Confirming proper washout after an elevated sample: When monitoring the analysis in real-time, if a sample concentration is greater than the highest concentration validated for washout (see Table 10 of Appendix C), do the following to verify that the run is still in control for low concentration samples before proceeding with analysis.
 - (i) Stop run following elevated sample

(ii) Verify that the run is still in control for lower concentration samples before proceeding with analysis. Analyze 2 blood blank checks followed by a low bench QC washout check. If the low bench QC wash check is not in control (within $\pm 3SD$ limits), repeat these 3 check samples until washout is verified before proceeding with analysis. Example:

3016 BldBlkChk Wash1

3016 BldBlkChk Wash2

LBXXXXX Wash

(iii) If the run is not verified in-control for low concentration samples before the next samples are analyzed, see Section 8.b.vii.2. for directions.

vii. Overnight operation or using auto stop: The run may be left to complete itself unattended as long as appropriate planning is made (e.g. sufficient solution supply and waste collection). Turn on the AutoStop feature of the ICP-MS software. Delay the shutdown at least 10 minutes (use peristaltic pump speed approximately that of the method wash) to rinse the sample introduction system of blood matrix before turning off the plasma. It will be necessary to replace the sample peristaltic pump tubing the next day since it will have been clamped shut overnight. Enable "Auto Start / Stop" is on the "AutoStop" tab of the Instrument window.

viii. Records of results: Run results will be documented after each run in both electronic and paper form.

1) Electronic records: Transfer data electronically to the laboratory information system. When keyboard entry must be used, proofread transcribed data after entry.

a) Export data from the ICP-MS software using "original conditions" or files and folders used during the analysis. Use descriptive report filenames (e.g. 2005-0714a_group55.txt). In the ICP-MS software under "Report Format" (METHOD window, REPORT tab) choose the "Use Separator" option, and under the "File Write" Section choose "Append."

b) Move the generated .TXT data file to the appropriate subdirectory on the network drive where exported data are stored prior to import to the laboratory information management system.

c) Import the instrument file into the laboratory information system with appropriate documentation (e.g. instrument ID, analyst, calibration standards lot number, and run or sample specific comments).

2) Run summary records: Printed run sheets, or PDF equivalent, must be documented with

a) Analyst initials

b) Instrument ID

c) Date of analysis and run # for the day

C. Evaluation: Analyst evaluation of run results:

i. Bench quality control: After completing a run, and importing the results into the laboratory information system, evaluate the run bench QC according to laboratory QC rules. The QC limits are based on the average and standard deviation of the beginning and ending analyses of each of the bench QC pools, so it will not be possible to know if the run is in control until statistically reviewed.

1) Rules for bench quality control evaluation: The following are the CDC DLS QC rules for three QC pools per run with two or more QC results per pool

- a) If all three QC run means are within $2S_m$ limits and individual results are within $2S_i$ limits, then accept the run.
- b) If one of the three QC run means is outside a $2S_m$ limit - reject run if
 - (i) Extreme Outlier – Run mean is beyond the characterization mean $\pm 4S_m$
 - (ii) 3S Rule - Run mean is outside a $3S_m$ limit
 - (iii) 2S Rule – Two or more of the run means are outside the same $2S_m$ limit
 - (iv) 10 X-bar Rule – Current and previous 9 run means are on same side of the characterization mean
- c) If one of the QC individual results is outside a $2S_i$ limit - reject run if:
 - (i) Extreme Outlier – One individual result is beyond the characterization mean $\pm 4S_m$
 - (ii) R 4S Rule – 2 or more of the within-run ranges in the same run exceed $4S_w$ (i.e., 95% range limit). Note: Since runs have multiple results per pool for 3 pools, the R 4S rule is applied within runs only.

2) Abbreviations:

- a) S_i = Standard deviation of individual results.
- b) S_m = Standard deviation of the run means.
- c) S_w = Within-run standard deviation.

3) Implications of QC failures: If the DLS SAS program declares the run “out of control” for any analyte, use the following to determine the implications on usability of the data from the run.

- a) 4–5 elements in the run
 - (i) 1 or 2 analytes “out of control”: ONLY the analytes which were “out of control” are invalid for reporting from the run.
 - (ii) 3 or more analytes “out of control”: All results, regardless of analyte, are invalid for reporting from the run.
- b) 1–3 elements in the run
 - (i) 1 analyte “out of control”: ONLY the analyte which is “out of control” is invalid for reporting from the run.
 - (ii) 2 or more analytes “out of control”: All results, regardless of analyte, are invalid for reporting from the run.

ii. Patient results:

1) Elevated concentrations: Refer to Figure 17 of Appendix C for flowchart.

a) Boundaries requiring confirmatory measurement:

(i) Results greater than the first (1UB) or second (2UB) upper boundaries. The concentrations assigned to 1UB and 2UB for an element is determined by study protocol but default concentrations are in Table 10 of Appendix C.

1. Results greater than the first upper boundary (1UB): Confirm by repeat analysis of a new sample preparation concentrations observed greater than the “first upper boundary” (defined in the laboratory database as the “1UB”). Report the first analytically valid result, as long as the confirmation is within 10%. Continue repeat analysis until a concentration can be confirmed.

2. Analyst reporting of elevated results: Report any patient results confirmed to be greater than the second upper boundary (2UB) as an “elevated result”.

(ii) Results greater than highest calibrator: Samples that exceed the high calibrator must be prepared with minimum extra dilution in duplicate to bring the observed result within the

calibration range ($\leq S8$). Report the first analytically valid result (i.e. the first one within the calibration range), as long as the confirmation is within 10%. Continue repeat analysis until a concentration can be confirmed.

- b) Concentrations requiring verification of washout: Following a result greater than the highest concentrations validated for washout (see Table 10 of Appendix C) do the following:
- (i) If the run was determined to be in-control for low concentration samples before the next samples were analyzed, no further action is required.
 - (ii) If the run was not determined to be in-control for low concentration samples before the next samples were analyzed confirm by re-analysis the results for the 2 samples immediately following the elevated sample. Report the results if they confirm the initial results within $\pm 10\%$ or $\pm 3SD$ of the low bench QC, whichever is greater.
- 2) Unacceptable reproducibility: If the range of the three replicate readings (maximum replicate concentration value - minimum replicate concentration value) for a single sample analysis is greater than the range maximum criteria listed in Table 10 of Appendix C and the range of the three replicate readings is greater than 10% of the observed concentration, do not use the measurement for reporting. Repeat the analysis of the sample.

D. Reporting: Submitting final work for review:

All analyses must undergo quality control and quality assurance review. After appropriately documenting the run in the laboratory information system (e.g. sample and run QC, and run and sample comments), inform the first level reviewer of the completed work and submit any printed documentation.

9. Routine equipment maintenance and data backups

Maintenance activities will be documented in the instrument logbook.

A. Equipment maintenance:

Analysts are expected to regularly evaluate the need for, and when necessary perform, cleaning, replacement, or re-positioning of components in ICP-MS the sample introduction system, interface, ion optics region, and equipment required resources (e.g. autosampler, exhaust, compressed gases, and coolant). Frequency of equipment maintenance will be dependent on instrument throughput.

- i. Parameter optimizations: Analysts are expected to optimize instrument parameters.

B. Dual detector calibration:

Perform dual detector calibration regularly for any element exceeding 1,000,000 cps for calibration standard 8. This is typically only Pb. Dual detector calibration solution is described in Section 6.f.ii. Frequency of dual detector calibration is typically monthly when throughput requires multiple analytical runs per week, or as needed for optimized linearity.

C. DRC optimizations:

DRC conditions (cell gas flow rate and RPq value) can be verified by analyzing the DRC optimization solutions (see Section 6.f.i) as needed to ensure proper reduction of potential ICP-MS interferences.

D. Data backup:

Data on the instrument computer will be backed up via two backup routines. Files used and produced by the ICP-MS in analyzing samples will be backed up and kept a minimum of two years after analysis.

- i. Daily backups to secondary hard drive: Program automatic backups of the relevant computer files to occur each night onto a secondary hard drive to prevent loss of data from failure of primary hard drive.
- ii. Weekly backup: Backup relevant computer files weekly either to secondary hard drive which is remote to the laboratory or to removable media which will be placed remote to the laboratory for retrieval in the case of catastrophic data loss elsewhere.

10. Reporting thresholds

A. Reportable range:

Blood element concentrations are reportable in the range between the method LOD and the high calibrator times the maximum permitted extra dilution (see Table 9 of Appendix C). Above the high calibrator, extra dilutions are made of the blood sample to bring the observed concentration within the calibration range.

B. Reference ranges (normal values):

In this method the 95% reference ranges (see Table 11 of Appendix C) for these elements in blood fall within the range of the calibrators.

C. Action levels:

Report concentrations observed greater than the “second upper boundary” (defined in the laboratory database as the “2UB”) to the QC reviewer as an “elevated result”. The concentration assigned to the 2UB for an element is determined by study protocol but default concentrations are listed in Table 10 of Appendix C. The protocol for supervisors reporting elevated results to medical personnel is defined according to the study protocol. But typically,

- i. Lead: Experts now use a reference level of 5 micrograms per deciliter to identify children with blood lead levels that are much higher than most children’s levels. This reference value is based on the 97.5th percentile of the National Health and Nutrition Examination Survey (NHANES)’s blood lead distribution in children. The current reference value is based on NHANES data from 2007-2008 and 2009-2010. CDC will assess the reference value every 4 years using the two most recent NHANES surveys. Chelation treatment is recommended at blood lead levels ≥ 45 $\mu\text{g}/\text{dL}$ [14]. The Occupational Safety and Health Administration regulations use a blood lead level of 40 $\mu\text{g}/\text{dL}$ as cause for written notification and a medical exam, and a blood lead level of 60 $\mu\text{g}/\text{dL}$ as cause for medical removal from exposure[65].
- ii. Cadmium: Levels of concern for cadmium in blood is >5 $\mu\text{g}/\text{L}$ [66, 67].
- iii. Mercury: The American Conference of Governmental Industrial Hygienists has a biological exposure index (BEI) of 15 $\mu\text{g}/\text{L}$ for inorganic mercury in blood (end of shift at end of work week)[67].
- iv. Manganese: Insufficient data to establish an action level.
- v. Selenium: >500 $\mu\text{g}/\text{L}$ [68, 69]

11. Method Calculations

A. Method limit of detection (LODs):

The method detection limits for elements in blood specimens are defined as 3 times s_0 , where s_0 is the estimate of the standard deviation at zero analyte concentration. S_0 is taken as the y-intercept of a linear or 2nd order polynomial regression of standard deviation versus concentration (4 concentration levels of the analytes in blood each measured 60 times across at least a 2-month timeframe). Method LODs are re-evaluated periodically.

B. Method limit of quantitation (LOQ):

The Division of Laboratory Sciences does not currently utilize limits of quantitation in regards to reporting limits[64].

C. QC Limits:

Quality control limits are calculated based on concentration results obtained in at least 20 separate runs. It is preferable to perform separate analyses on separate days and using multiple calibrator lot numbers, instruments, and analysts to best mimic real-life variability. The statistical calculations are performed using the SAS program developed for the Division of Laboratory Sciences.

12. Alternate methods for performing test and storing specimens if test system fails:

If the analytical system fails, the analysis may be setup on other ICP-MS instruments in the laboratory. If no other instrument is available, store the specimens at $\sim 4^\circ\text{C}$ until the analytical system can be restored to functionality. If interruption longer than 4 weeks is anticipated, then store blood specimens at $\leq -20^\circ\text{C}$.

13. Method performance documentation


Method performance documentation for this method including accuracy, precision, sensitivity, specificity, and stability is provided in Appendix A of this method documentation. **The signatures of the branch chief and director of the Division of Laboratory Sciences on the first page of this procedure denote that the method performance is fit for the intended use of the method.**

14. Summary Statistics and QC Graphs

See following pages


2017-2018 Summary Statistics and QC Chart for Blood cadmium ($\mu\text{g/L}$)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
1094	358	20JAN17	31JAN19	3.0057	0.0837	2.8
1156	358	20JAN17	31JAN19	0.4482	0.0425	9.5
1507	358	20JAN17	31JAN19	44.4582	1.1024	2.5




2017-2018 Summary Statistics and QC Chart for Blood lead ($\mu\text{g/dL}$)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
1094	362	20JAN17	31JAN19	10.0749	0.1570	1.6
1156	362	20JAN17	31JAN19	2.1053	0.0583	2.8
1507	362	20JAN17	31JAN19	88.5323	1.4580	1.6




2017-2018 Summary Statistics and QC Chart for Blood manganese ($\mu\text{g/L}$)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
1094	373	20JAN17	31JAN19	14.346	0.529	3.7
1156	373	20JAN17	31JAN19	7.774	0.385	4.9
1507	373	20JAN17	31JAN19	42.382	1.398	3.3




2017-2018 Summary Statistics and QC Chart for Blood mercury,total (µg/L)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
1094	358	20JAN17	31JAN19	5.921	0.159	2.7
1156	358	20JAN17	31JAN19	0.587	0.067	11.4
1507	358	20JAN17	31JAN19	43.258	4.274	9.9



2017-2018 Summary Statistics and QC Chart for Blood selenium ($\mu\text{g/L}$)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
1094	356	24JAN17	07FEB19	256.688	8.514	3.3
1156	356	24JAN17	07FEB19	190.201	7.421	3.9
1507	356	24JAN17	07FEB19	2697.385	84.010	3.1



15. Appendix A. Method performance documentation

A. Accuracy

i. Cadmium

Accuracy compared to Reference Material
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3*LOD, where it should be within $\pm 20\%$

Method name: Blood multi-element analysis by ICP-DRC-MS
 Method #: 3016
 Matrix: Blood
 Units: $\mu\text{g/L}$
 Reference material: NIST SRM 955c levels 2, 3, and 4
 Analyte: cadmium

Reference material	Replicate	Nominal value	Measured concentration					Mean	SD	CV (%)	Difference from nominal value (%)
			Day 1	Day 2	Day 3	Day 4	Day 5				
Material 1	1	2.14	2.0	2.1	2.1	2.1	2.1	2.11	0.07	3.30	-1.3
	2		2.0	2.1	2.3	2.2	2.1				
Material 2	1	5.20	5.1	5.1	5.0	5.2	4.9	5.05	0.11	2.16	-2.9
	2		5.1	5.1	4.8	5.2	5.0				
Material 3	1	9.85	10	10	9.9	10	10	10.18	0.28	2.76	3.4
	2		10	10	10	11	9.9				

ii. Lead

Accuracy compared to Reference Material
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3*LOD, where it should be within $\pm 20\%$

Method name: Blood multi-element analysis by ICP-DRC-MS
 Method #: 3016
 Matrix: Blood
 Units: $\mu\text{g/dL}$
 Reference material: NIST SRM 955c levels 2, 3, and 4
 Analyte: lead

Reference material	Replicate	Nominal value	Measured concentration					Mean	SD	CV (%)	Difference from nominal value (%)
			Day 1	Day 2	Day 3	Day 4	Day 5				
Material 1	1	13.95	14	14	14	14	14	13.91	0.16	1.15	-0.3
	2		14	14	14	14	14				
Material 2	1	27.8	27	28	28	28	27	27.64	0.24	0.85	-0.4
	2		27	28	28	28	28				
Material 3	1	45.5	46	45	45	46	45	45.57	0.48	1.04	0.1
	2		46	45	45	46	45				

iii. Manganese

Accuracy compared to Reference Material														
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3*LOD, where it should be within $\pm 20\%$														
Method name:		Blood multi-element analysis by ICP-DRC-MS												
Method #:		3016												
Matrix:		Blood												
Units:		$\mu\text{g/L}$												
Reference material:		ClinChek I, Seronorm L1 (lot # 1406263), Seronorm L3 (lot # 1509408)												
Analyte:		manganese												
Reference material	Replicate	Nominal value	Measured concentration								Mean	SD	CV (%)	Difference from nominal value (%)
			Day 1	Day 2	Day 3	Day 4	Day 5							
Level 1	1	8.01	10	7.4	8.3	7.8	7.9	8.31	1.05	12.69	3.7			
	2		10	7.5	8.2	7.7	7.8							
Level 2	1	18.4	21	17	18	18	18	18.73	1.37	7.30	1.8			
	2		21	18	18	18	19							
Level 3	1	33.3	39	35	37	31	39	36.05	2.95	8.19	8.3			
	2		39	35	37	31	38							

iv. Mercury

Accuracy compared to Reference Material														
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3*LOD, where it should be within $\pm 20\%$														
Method name:		Blood multi-element analysis by ICP-DRC-MS												
Method #:		3016												
Matrix:		Blood												
Units:		$\mu\text{g/L}$												
Reference material:		NIST SRM 955c levels 2, 3, and 4												
Analyte:		mercury (total)												
Reference material	Replicate	Nominal value	Measured concentration								Mean	SD	CV (%)	Difference from nominal value (%)
			Day 1	Day 2	Day 3	Day 4	Day 5							
Level 1	1	4.95	5.5	5.7	5.6	5.6	5.6	5.59	0.13	2.33	13.0			
	2		5.6	5.6	5.8	5.6	5.3							
Level 2	1	17.8	19	19	19	19	19	19.30	0.33	1.70	8.4			
	2		19	19	20	20	19							
Level 3	1	33.9	36	35	36	35	34	35.30	0.73	2.06	4.1			
	2		36	35	36	36	34							

v. Selenium

Accuracy compared to Reference Material											
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3*LOD, where it should be within $\pm 20\%$											
Method name:		Blood multi-element analysis by ICP-DRC-MS									
Method #:		3016									
Matrix:		Blood									
Units:		$\mu\text{g/L}$									
Reference material:		Seronorm L1 (lot # 1406263), Seronorm L2 (lot # 1406264), Seronorm L3 (lot # 1509408)									
Analyte:		selenium									
Reference material	Replicate	Nominal value	Measured concentration								Difference from nominal value (%)
			Day 1	Day 2	Day 3	Day 4	Day 5	Mean	SD	CV (%)	
Level 1	1	60	76	77	75	65	52	69.11	9.81	14.19	15.2
	2		76	76	75	65	53				
Level 2	1	161	182	172	171	161	155	169.58	9.03	5.32	5.3
	2		182	172	174	162	163				
Level 3	1	198	236	225	228	207	219	223.43	10.51	4.70	12.8
	2		238	226	231	209	216				

B. Precision

i. Cadmium

Precision						
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)						
Method name:	Blood multi-element analysis by ICP-DRC-MS					
Method #:	3016					
Matrix:	Blood					
Units:	$\mu\text{g/L}$					
Analyte:	cadmium					
Quality material 1						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	2.8	2.8	2.83	2.97025E-05	2.97025E-05	16.06821361
2	3.1	3.1	3.10	0.00140625	0.00140625	19.25349458
3	2.9	3.0	2.98	0.002261003	0.002261002	17.71136645
4	2.9	3.1	2.99	0.005133723	0.005133723	17.82343513
5	3.0	3.1	3.06	0.00470596	0.00470596	18.6843845
6	3.0	3.0	2.98	0.00037636	0.00037636	17.70363008
7	3.0	3.0	2.97	0.00021609	0.00021609	17.68102578
8	3.0	2.9	2.95	0.00210681	0.00210681	17.45341362
9	2.9	2.8	2.87	0.004362602	0.004362603	16.48585621
10	3.3	3.2	3.23	0.001444	0.001444	20.85159042
Grand sum	59.9146	Grand mean	2.99573			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.044085	0.0044085	0.066396536	2.22		
Between Run	0.228445702	0.025382856	0.102406923	3.42		
Total	0.272530702		0.122047851	4.07		
Quality material 2						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	45	43	43.99	0.309859222	0.309859222	3869.703541
2	46	45	45.42	0.34644996	0.34644996	4125.734787
3	43	43	42.99	0.054639063	0.054639063	3696.51235
4	45	41	42.66	4.62078016	4.62078016	3639.171047
5	46	47	46.18	0.140212803	0.140212802	4264.76919
6	45	44	44.75	0.457314063	0.457314062	4004.256897
7	45	43	43.75	0.64657681	0.64657681	3827.337541
8	45	43	44.13	1.008317223	1.008317222	3894.234228
9	43	42	42.54	0.725478063	0.725478063	3620.111505
10	45	46	45.34	0.04264225	0.04264225	4112.265498
Grand sum	883.4752	Grand mean	44.17376			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev		
Within Run	16.70453923	1.670453923	1.292460414	2.93		
Between Run	27.67513226	3.075014695	0.83802171	1.90		
Total	44.37967149		1.540368238	3.49		

ii. Lead

Precision						
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)						
Method name:	Blood multi-element analysis by ICP-DRC-MS					
Method #:	3016					
Matrix:	Blood					
Units:	$\mu\text{g/dL}$					
Analyte:	lead					
Quality material 1						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	10.1	10.1	10.10	0.00076729	0.00076729	204.1371768
2	10.1	10.2	10.14	0.00013456	0.00013456	205.5783645
3	9.8	9.9	9.87	0.00264196	0.00264196	194.8338
4	10.3	10.2	10.24	0.00084681	0.00084681	209.6496691
5	9.7	9.8	9.73	0.00131769	0.00131769	189.4664712
6	10.0	10.0	10.01	0.000504003	0.000504002	200.594441
7	10.1	10.3	10.19	0.00425104	0.00425104	207.7781895
8	10.1	10.1	10.06	4.624E-05	4.624E-05	202.4273205
9	10.0	10.1	10.04	0.000315062	0.000315062	201.6252886
10	10.1	10.0	10.05	8.19025E-05	8.19025E-05	202.111544
Grand sum	200.8881	Grand mean	10.044405			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.021813115	0.002181311	0.046704513	0.46		
Between Run	0.400829295	0.044536588	0.145525387	1.45		
Total	0.42264241		0.15283635	1.52		
Quality material 2						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	88	88	87.81	0.050176	0.050176	15419.50629
2	88	88	88.02	0.056620202	0.056620202	15496.78365
3	89	89	89.02	0.00038025	0.00038025	15848.01697
4	90	91	90.20	0.36687249	0.36687249	16270.8533
5	83	84	83.52	0.08323225	0.08323225	13950.24539
6	89	89	89.16	0.003463323	0.003463322	15899.02903
7	90	91	90.10	0.29964676	0.29964676	16235.69564
8	89	89	88.94	0.00488601	0.00488601	15821.42988
9	89	90	89.36	0.083781303	0.083781302	15969.57923
10	86	86	86.16	0.001235522	0.001235523	14846.90165
Grand sum	1764.5586	Grand mean	88.22793			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev		
Within Run	1.90058822	0.190058822	0.435957363	0.49		
Between Run	74.68839436	8.298710485	2.013535654	2.28		
Total	76.58898258		2.060190441	2.34		

iii. Manganese

Precision						
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)						
Method name:	Blood multi-element analysis by ICP-DRC-MS					
Method #:	3016					
Matrix:	Blood					
Units:	$\mu\text{g/L}$					
Analyte:	manganese					
Quality material 1						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean ²
1	14	14	13.86	0.002505003	0.002505003	383.9691584
2	15	15	14.98	0.021025	0.021025	448.5671424
3	15	15	15.08	0.130718403	0.130718403	454.966629
4	14	15	14.42	0.16842816	0.16842816	415.642112
5	15	15	15.07	0.00101761	0.00101761	454.1977441
6	15	14	14.14	0.310193303	0.310193302	399.7067106
7	15	16	15.64	0.03598609	0.03598609	489.3505848
8	13	14	13.79	0.380133903	0.380133903	380.5461132
9	14	13	13.81	0.126273623	0.126273623	381.3631531
10	15	15	14.62	0.000382203	0.000382203	427.2870678
Grand sum	290.7972	Grand mean	14.53986			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	2.35332659	0.235332659	0.485110976	3.34		
Between Run	7.445839118	0.827315458	0.544050916	3.74		
Total	9.799165708		0.728919789	5.01		
Quality material 2						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean ²
1	45	43	43.85	0.43112356	0.43112356	3845.153896
2	45	43	44.02	1.526089622	1.526089623	3875.177452
3	43	41	42.36	0.984560062	0.984560063	3588.036058
4	43	40	41.53	3.296221803	3.296221803	3449.124651
5	45	44	44.66	0.21566736	0.21566736	3989.495678
6	44	40	42.01	3.944196	3.944196	3529.713808
7	45	47	45.75	0.70408881	0.70408881	4186.1433
8	43	44	43.17	0.167977022	0.167977022	3727.358238
9	41	40	40.57	0.166831402	0.166831403	3291.760547
10	43	44	43.27	0.233627223	0.233627223	3745.061785
Grand sum	862.3688	Grand mean	43.11844			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev		
Within Run	23.34076573	2.334076573	1.527768495	3.54		
Between Run	43.0280519	4.780894655	1.106078226	2.57		
Total	66.36881763		1.886129798	4.37		

iv. Mercury

Precision						
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)						
Method name:	Blood multi-element analysis by ICP-DRC-MS					
Method #:	3016					
Matrix:	Blood					
Units:	$\mu\text{g/L}$					
Analyte:	mercury (total)					
Quality material 1						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean ²
1	6.0	6.1	6.07	0.00201601	0.00201601	73.76994578
2	5.6	5.7	5.61	0.00251001	0.00251001	62.97113088
3	5.7	5.8	5.71	0.003666303	0.003666303	65.25274561
4	5.8	5.9	5.89	0.002943063	0.002943062	69.36888685
5	5.6	6.0	5.77	0.033948063	0.033948062	66.61465313
6	5.7	5.6	5.63	0.00087616	0.00087616	63.30600242
7	5.8	5.9	5.83	0.001447803	0.001447802	68.08161361
8	5.8	6.1	5.97	0.02033476	0.02033476	71.27941202
9	5.7	5.9	5.83	0.010110303	0.010110303	67.95797945
10	6.0	5.9	5.95	0.001922823	0.001922823	70.80381001
Grand sum	116.5332	Grand mean	5.82666			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.15955059	0.015955059	0.126313337	2.17		
Between Run	0.406844618	0.045204958	0.120933657	2.08		
Total	0.566395208		0.174871405	3.00		
Quality material 2						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean ²
1	36	37	36.72	0.41024025	0.41024025	2697.12808
2	45	34	39.53	32.40171006	32.40171006	3125.676645
3	38	37	37.73	0.166586422	0.166586423	2847.234083
4	39	39	38.95	0.16589329	0.16589329	3034.719162
5	47	44	45.42	1.883893502	1.883893503	4125.743871
6	37	44	40.16	12.1874301	12.1874301	3225.916261
7	45	44	44.42	0.015388402	0.015388403	3945.677594
8	43	46	44.51	1.149076802	1.149076803	3963.072518
9	34	35	34.08	0.312425102	0.312425103	2322.599721
10	31	30	30.36	0.045903062	0.045903062	1843.914628
Grand sum	783.7858	Grand mean	39.18929			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev		
Within Run	97.477094	9.7477094	3.122132188	7.97		
Between Run	415.6735495	46.18594994	4.268386143	10.89		
Total	513.1506435		5.288367392	13.49		

v. Selenium

Precision						
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)						
Method name:	Blood multi-element analysis by ICP-DRC-MS					
Method #:	3016					
Matrix:	Blood					
Units:	$\mu\text{g/L}$					
Analyte:	selenium					
Quality material 1						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	247	238	242.59	21.42902972	21.42902972	117696.66255
2	265	258	261.53	13.9760084	13.9760084	136793.63265
3	246	250	248.10	3.148317923	3.148317922	123106.57494
4	251	250	250.51	0.757857303	0.757857302	125508.26562
5	245	243	244.01	0.58660281	0.58660281	119081.85780
6	252	243	247.72	21.05157924	21.05157924	122729.50501
7	244	242	243.04	0.535824	0.535824	118137.66093
8	226	236	231.15	26.60909056	26.60909056	106864.43589
9	234	238	235.82	4.703910322	4.703910323	111223.88987
10	263	253	257.54	25.1597544	25.1597544	132657.77236
Grand sum	4924.0232	Grand mean	246.20116			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	235.9159494	23.59159494	4.857117966	1.97		
Between Run	1500.033931	166.6704368	8.458097949	3.44		
Total	1735.94988		9.753513001	3.96		
Quality material 2						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	2826	2740	2782.65	1856.200897	1856.200897	15486285.94
2	2858	2770	2813.73	1925.098988	1925.098988	15834170.47
3	2489	2478	2483.40	27.83090025	27.83090025	12334595.82
4	2715	2336	2525.41	35763.29181	35763.29181	12755374.67
5	2754	2715	2734.86	378.1702516	378.1702516	14958916.25
6	2771	2722	2746.54	605.9277634	605.9277634	15086995.8
7	2774	2696	2735.22	1536.369532	1536.369532	14962890.27
8	2771	2635	2703.07	4666.898236	4666.898236	14613200.8
9	2408	2524	2465.78	3356.863988	3356.863988	12160189.85
10	2402	2316	2359.07	1862.000171	1862.000171	11130384.78
Grand sum	52699.4875	Grand mean	2634.974375			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	103957.3051	10395.73051	101.9594552	3.87		
Between Run	461205.5211	51245.0579	142.9148827	5.42		
Total	565162.8262		175.5573815	6.66		

C. Stability

i. Cadmium

Stability								
The initial measurement can be from the same day for all stability experiments.								
Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions								
Describe condition: Three times frozen at -20°C and then thawed (3 freeze-thaw cycles)								
Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)								
Describe condition: Original samples in cryovial stored at room temperature for 1 day								
Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler								
Describe condition: Processed samples (ready for instrument analysis) stored at room temperature for 1 day								
Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis								
Describe condition: Samples from 2 characterized pools stored at -70°C for 4 years								
All stability sample results should be within ±15% of nominal concentration								
Method name:	Blood multi-element analysis by ICP-DRC-MS							
Method #:	3016							
Matrix:	Blood							
Units:	µg/L							
Analyte:	cadmium							
Quality material 1								
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	3.0	2.9	3.1	3.1	3.0	3.0	2.9	3.0
Replicate 2	3.1	2.9	3.2	3.2	2.8	2.9	2.8	3.1
Replicate 3	3.2	3.0	3.0	3.1	3.1	3.1	3.2	3.2
Mean	3.101666667	2.952566667	3.108066667	3.1	2.9629	2.9847	2.968033333	3.1
% difference from initial measurement	--	-4.8	--	0.3	--	0.7	--	4.5
Quality material 2								
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	45	44	45	45	44	45	43	45
Replicate 2	43	45	45	44	42	45	42	43
Replicate 3	44	46	45	45	45	44	46	44
Mean	44.31223333	44.9	45.1293	44.8	43.64023333	44.75216667	43.45536667	44.3
% difference from initial measurement	--	1.4	--	-0.8	--	2.5	--	2.0

ii. Lead

Stability									
The initial measurement can be from the same day for all stability experiments.									
Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions									
Describe condition: Three times frozen at -20°C and then thawed (3 freeze-thaw cycles)									
Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)									
Describe condition: Original samples in cryovial stored at room temperature for 1 day									
Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler									
Describe condition: Processed samples (ready for instrument analysis) stored at room temperature for 1 day									
Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis									
Describe condition: Samples from 2 characterized pools stored at -70°C for 4 years									
All stability sample results should be within ±15% of nominal concentration									
Method name:	Blood multi-element analysis by ICP-DRC-MS								
Method #:	3016								
Matrix:	Blood								
Units:	µg/dL								
Analyte:	lead								
Quality material 1									
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability	
Replicate 1	10.0	9.6	10.1	10.2	10.0	9.5	10.2	10.0	
Replicate 2	10.0	10.0	10.1	10.4	9.7	10.0	10.1	10.0	
Replicate 3	10.5	9.8	10.0	10.3	10.7	10.7	10.0	10.5	
Mean	10.1659	9.8336	10.05456667	10.3	10.13676667	10.05676667	10.10496667	10.2	
% difference from initial measurement	--	-3.3	--	2.4	--	-0.8	--	0.6	
Quality material 2									
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability	
Replicate 1	89	86	88	89	89	89	88	89	
Replicate 2	86	89	89	90	85	87	90	86	
Replicate 3	87	89	88	89	90	87	86	87	
Mean	87.4425	88.35676667	88.34143333	89.6	88.04316667	88.0469	87.87046667	87.4	
% difference from initial measurement	--	1.0	--	1.4	--	0.0	--	-0.5	

iii. Manganese

Stability

The initial measurement can be from the same day for all stability experiments.

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: Three times frozen at -20°C and then thawed (3 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: Original samples in cryovial stored at room temperature for 1 day

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: Processed samples (ready for instrument analysis) stored at room temperature for 1 day

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: Samples from 2 characterized pools stored at -70°C for 4 years

All stability sample results should be within ±15% of nominal concentration

Method name: Blood multi-element analysis by ICP-DRC-MS

Method #: 3016

Matrix: Blood

Units: µg/L

Analyte: manganese

Quality material 1								
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	14	14	15	15	15	14	14	14
Replicate 2	14	14	15	15	14	14	13	14
Replicate 3	15	14	15	15	15	15	15	15
Mean	14.46053333	14.06546667	14.78956667	15.1	14.70273333	14.33466667	14.15363333	14.5
% difference from initial measurement	--	-2.7	--	2.2	--	-2.5	--	2.2
Quality material 2								
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	43	42	43	43	43	42	44	43
Replicate 2	42	43	43	44	43	42	40	42
Replicate 3	43	43	43	44	44	42	44	43
Mean	42.4561	42.7184	42.9118	43.6	43.30506667	42.3925	42.4991	42.5
% difference from initial measurement	--	0.6	--	1.5	--	-2.1	--	-0.1

iv. Mercury

Stability									
The initial measurement can be from the same day for all stability experiments.									
Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions									
Describe condition: Three times frozen at -20°C and then thawed (3 freeze-thaw cycles)									
Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)									
Describe condition: Original samples in cryovial stored at room temperature for 1 day									
Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler									
Describe condition: Processed samples (ready for instrument analysis) stored at room temperature for 1 day									
Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis									
Describe condition: Samples from 2 characterized pools stored at -70°C for 4 years									
All stability sample results should be within ±15% of nominal concentration									
Method name:	Blood multi-element analysis by ICP-DRC-MS								
Method #:	3016								
Matrix:	Blood								
Units:	µg/L								
Analyte:	mercury (total)								
Quality material 1									
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability	
Replicate 1	6.0	5.8	6.1	6.1	5.7	5.7	5.8	6.0	
Replicate 2	5.9	6.0	6.0	6.1	5.5	5.9	6.1	5.9	
Replicate 3	6.1	6.0	5.9	6.2	6.2	6.2	5.5	6.1	
Mean	5.999766667	5.929033333	6.022666667	6.2	5.781533333	5.923666667	5.792633333	6.0	
% difference from initial measurement	--	-1.2	--	2.3	--	2.5	--	3.6	
Quality material 2									
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability	
Replicate 1	50	42	48	48	36	36	38	50	
Replicate 2	36	43	38	37	36	38	39	36	
Replicate 3	35	41	44	41	38	37	39	35	
Mean	40.47766667	42.03253333	43.4925	42.0	36.6215	36.88343333	38.71813333	40.5	
% difference from initial measurement	--	3.8	--	-3.4	--	0.7	--	4.5	

v. Selenium

Stability

The initial measurement can be from the same day for all stability experiments.

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: Three times frozen at -20°C and then thawed (3 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: Original samples in cryovial stored at room temperature for 1 day

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: Processed samples (ready for instrument analysis) stored at room temperature for 1 day

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: Samples from 2 characterized pools stored at -70°C for 4 years

All stability sample results should be within ±15% of nominal concentration

Method name: Blood multi-element analysis by ICP-DRC-MS

Method #: 3016

Matrix: Blood

Units: µg/L

Analyte: selenium

Quality material 1		Quality material 1		Quality material 1		Quality material 1		
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	270	264	247	245	275	254	246	270
Replicate 2	262	264	252	244	263	262	250	262
Replicate 3	268	266	255	235	284	268	251	268
Mean	266.4195333	264.7295667	251.0805333	241.4	274.0861667	261.3381333	249.1923333	266.4
% difference from initial measurement	--	-0.6	--	-3.8	--	-4.7	--	6.9
Quality material 2		Quality material 2		Quality material 2		Quality material 2		
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	2765	2750	2688	2485	2872	2755	2489	2765
Replicate 2	2660	2805	2703	2226	2771	2682	2478	2660
Replicate 3	2688	2771	2712	2501	2934	2684	2715	2688
Mean	2704.1685	2775.0806	2700.885533	2403.9	2858.982	2706.9993	2560.443067	2704.2
% difference from initial measurement	--	2.6	--	-11.0	--	-5.3	--	5.6

D. Analytical Sensitivity and Specificity

LOD, specificity and fit for intended use			
Method name:	Blood multi-element analysis by ICP-DRC-MS		
Method #:	3016		
Matrix:	Blood		
Units:	µg/L (lead in µg/dL)		
Analytes	Limit of Detection (LOD)	Interferences successfully checked in at least 50 human samples	Accuracy, precision, LOD, specificity and stability meet performance specifications for intended use
cadmium (BCD)	0.10	Yes	Yes
lead (BPB)	0.07	Yes	Yes
manganese (BMN)	0.99	Yes	Yes
mercury (THG)	0.28	Yes	Yes
selenium (BSE)	24.48	Yes	Yes

16. Appendix B. Ruggedness testing results

A. Ruggedness Parameter #1: Stability of sample preparations

This test evaluated the stability of prepared samples for up to 48 hours after preparation, simulating a situation where immediate sample analysis was not possible.

i. Test Details:

- 1) Day 1: Prepared a set of dilutions (calibrators, blanks, reference material, fake samples) for analysis in triplicate. Analyzed set 1 immediately (normal practice). Cap sets 2 and 3 and leave at ambient temperature for later analysis.
- 2) Day 2: Prepared run set 4 and analyzed it sequentially with run set 2
- 3) Day 3: Prepared run set 5 and analyzed it sequentially with run set 3

ii. Results: See Ruggedness Table 1.

iii. Conclusions: Samples which have been diluted 1+1+48 for analysis up to one (1) day previously can still be analyzed.

Ruggedness Table 1. Measured concentrations of Hg, Pb, Cd, Mn, and Se in bench QC or reference samples analyzed 0, 24, and 48 hours after preparation

ID	Time, prep to analysis	Hg** (µg/L)	Pb** (µg /dL)	Cd** (µg /L)	Mn** (µg /L)	Se** (µg /L)
LB08707 WB2	<i>target mean and 3SD range</i>	0.585 0.318 – 0.852	2.12 1.99 – 2.25	0.488 0.353 – 0.623	7.98 6.38 – 9.59	
	0 hr	0.418	2.03	0.399	6.09	
	24 hr (fresh)	0.504 (0.522)	1.99 (2.18)	0.419 (0.47)	7.06 (7.88)	
	48 hr	0.396 (0.418)	2.04 (2.03)	0.509 (0.40)	7.82 (6.09)	
HB08708 WB2	<i>target mean and 3SD range</i>	6.19 5.74 – 6.63	10.1 9.73 – 10.4	3.14 2.84 – 3.44	14.9 12.8 – 17.1	
	0 hr	5.86	10.0	3.03	12.5	
	24 hr	5.46 (5.7)	9.5 (10.7)	2.85 (3.17)	13.6 (14.7)	
	48 hr	2.64 (5.9)	9.2 (10.0)	2.79 (3.03)	13.5 (12.5)	
QMEQAS 07B-03*	<i>target mean and 3SD range</i>					228 206 – 251
	0 hr					192
	24 hr					202 (217)
	48 hr					56 (192)
QMEQAS 10B-06*	<i>target mean and 2SD range</i>					239 215 – 253
	0 hr					212
	24 hr					221 (238)
	48 hr					62 (212)

*samples purchase from Le centre de toxicology du Quebec (Quebec, Canada)

** Test performed 12/6-8/10 by Deanna Jones. Results are the average of the beginning and ending QC results for each analytical run.

B. Ruggedness Parameter Test #2: RF Power

This test evaluated the significance of the RF Power setting of the ICP when analyzing blood samples for whole blood metals.

i. Test Details:

- 1) Prepare a set of dilutions (calibrators, blanks, reference material, dummy samples) for analysis in triplicate (three separate sets of tubes).
- 2) Analyze them in three separate runs on the same day, same instrument.
- 3) Change the RF Power across the runs.
- 4) Allow 15 minutes equilibration time between runs for RF Power to stabilize.

ii. Results: See Ruggedness Table 2.

iii. Conclusions: Results are not compromised by changes in RF power within the range of 1150W to 1600W.

Ruggedness Table 2. Measured concentrations of Hg, Pb, Cd, Mn, and Se in bench QC or reference samples at RF powers of 1150W, 1450W, and 1600W

ID	RF power (W)	Hg** (µg /L)	Pb** (µg /dL)	Cd** (µg /L)	Mn** (µg /L)	Se** (µg /L)
LB08707_WB2	<i>target mean and 2SD range</i>	0.585 0.407 – 0.763	2.12 2.03 – 2.21	0.488 0.398 – 0.578	7.98 6.91 – 9.05	
	1150 W	0.517	2.09	0.432	7.35	
	1450 W (default)	0.512	2.03	0.369	6.76	
	1600 W	0.529	2.02	0.418	7.17	
HB08708_WB2	<i>target mean and 2SD range</i>	6.19 5.89 – 6.48	10.1 9.84 – 10.3	3.14 2.94 – 3.34	14.9 13.5 – 16.4	
	1150 W	5.90	10.0	2.93	13.7	
	1450 W (default)	6.23	10.2	2.90	12.8	
	1600 W	5.99	10.1	3.07	13.3	
QMEQAS08B-02*	<i>target mean and 2SD range</i>					293 273 - 313
	1150 W					269
	1450 W (default)					288
	1600 W					314
QMEQAS08B-08*	<i>target mean and 2SD range</i>					165 154 - 176
	1150 W					179
	1450 W (default)					147
	1600 W					146

*samples purchase from Le centre de toxicology du Quebec (Quebec, Canada)

** Test performed on December 6 and December 10, 2010 by Deanna Jones. Results below are the average of the beginning and ending QC results for each analytical run.

C. Ruggedness parameter test #3: Dynamic reaction cell (DRC) gas flow rate

This test evaluated the significance of the DRC gas flow rate (e.g. oxygen or methane) while analyzing blood samples for elements analyzed in DRC mode (Hg, Mn, and Se). The cell gas for Mn and Hg is oxygen (O₂) and the typical flow rate is 1.2 mL/min. The cell gas for Se is methane (CH₄) and the typical flow rate is 0.84 mL/min.

i. Test Details:

- 1) Prepare three sets of dilutions in three separate sets of tubes (calibrators, blanks, reference materials, dummy samples) for analysis.
- 2) Analyze each set of dilutions in separate runs on the same day using the same instrument.
- 3) Change the gas flow rate by 20% of the method default for each cell gas.
 - a) Run 1: decreased cell gas flow rates by 20% of method default O₂ = 0.96 mL/min; CH₄ = 0.7 mL/min.
 - b) Run 2: method default O₂ = 1.2 mL/min; CH₄ = 0.84 mL/min
 - c) Run 3: increased cell gas flow rate by 20% of the method default O₂ = 1.44 mL/min; CH₄ = 1.0 mL/min.

ii. Results: See Ruggedness Table 3 and Ruggedness Table 4.

iii. Conclusion: Accuracy of Mn and Hg results are not compromised by changes in cell gas flow rate within the range tested (0.96 – 1.44 mL/min).

Ruggedness Table 3. Measured concentrations of Hg and Mn with oxygen gas at 0.96 mL/min, 1.2 mL/min, and 1.44 mL/min, and Se with methane gas at 0.7 mL/min, 0.84 mL/min, and 1.0 mL/min, in bench QC or reference samples

ID	cell gas flow rate	Hg* (µg /L)	Mn* (µg /L)	ID	cell gas flow rate	Se** (µg /L)
LB08707_WB2	<i>target mean and 2SD range</i>	0.585 0.407 – 0.763	7.98 6.91 – 9.05	QMEQAS07B-09*	<i>target mean and 2SD range</i>	157 146 - 168
	0.96 mL/min O ₂	0.457	8.49		0.7 mL/min CH ₄	187
	1.2 mL/min O ₂	0.479	8.15		0.84 mL/min CH ₄	186
	1.44 mL/min O ₂	0.555	8.12		1.0 mL/min CH ₄	191
HB08708_WB2	<i>Target Mean and 2SD Range</i>	6.19 5.89 – 6.48	14.9 13.5 – 16.4	QMEQAS08B-02*	<i>target mean and 2SD range</i>	293 273 - 313
	0.96 mL/min O ₂	4.71	14.4		0.7 mL/min CH ₄	328
	1.2 mL/min O ₂	5.45	15.2		0.84 mL/min CH ₄	334
	1.44 mL/min O ₂	5.34	14.6		1.0 mL/min CH ₄	339
* Test performed on December 6, 2010 and January 4, 2010 by Deanna Jones. Results for LB08707_WB2 and HB08708_WB2 are the average of the beginning and ending results for each analytical run.						

D. Ruggedness parameter test #4: RPq value

This test evaluated the significance of the RPq value while analyzing blood samples for Se, Mn and Hg.

i. Test Details:

- 1) Prepare three sets of dilutions in three separate sets of tubes (calibrators, blanks, reference materials, dummy samples) for analysis.
- 2) Analyze each set of dilutions in separate runs on the same day using the same instrument.
- 3) Change the RPq value by 20% of the method default for each analyte.
 - a) Run 1: decreased RPq to 0.48 for Mn and Hg; and 0.52 for Se.
 - b) Run 2: method default RPq of 0.6 for Mn and Hg; and 0.65 for Se.
 - c) Run 3: increased RPq to 0.72 for Mn and Hg; 0.78 for Se.

ii. Results: See Ruggedness Table 5 and Ruggedness Table 6.

iii. Conclusion: Accuracy of Mn and Hg results are not compromised by changes in RPq settings within the range tested (0.48 – 0.72). Accuracy of Se results are not compromised by changes in RPq settings within the range tested (0.52 – 0.78 for Se).

Ruggedness Table 4. Measured concentrations of Hg, Mn at RPq values of 0.48, 0.6, and 0.72, and Se at RPq values of 0.52, 0.7, 0.78, in bench QC or reference samples

ID	RPq*	Hg* (µg /L)	Mn* (µg /L)	ID	RPq*	Se** (µg /L)
LB08707_WB2	<i>target mean and 2SD range</i>	0.585 0.407 – 0.763	7.98 6.91 – 9.05	QMEQAS07B-09**	<i>target mean and 2SD range</i>	293 273 – 313
	RPq = 0.48	0.455	7.86		RPq = 0.52	262
	RPq = 0.6	0.418	6.09		RPq = 0.7	250
	RPq = 0.72	0.402	7.99		RPq = 0.78	277
HB08708_WB2	<i>Target Mean and 2SD Range</i>	6.19 5.89 – 6.48	14.9 13.5 – 16.4	QMEQAS08B-02**	<i>target mean and 2SD range</i>	361 337 - 385
	RPq = 0.48	5.54	14.4		RPq = 0.52	347
	RPq = 0.6	5.86	12.5		RPq = 0.7	349
	RPq = 0.72	5.53	14.9		RPq = 0.78	364
* Test performed on December 21, 2010 by Deanna Jones. Results for LB08707_WB2 and HB08708_WB2 are the average of the beginning and ending results for each analytical run. **samples purchase from Le centre de toxicology du Quebec (Quebec, Canada). Test performed on December 21, 2010 by Deanna Jones.						

E. Ruggedness parameter test #5: Axial Field Voltage (AFT)

This test evaluated the significance of the Axial Field Voltage (AFT) while analyzing blood samples for whole blood metals. The Axial Field Voltage may vary on each instrument.

i. Test Details:

- 1) Prepare three sets of dilutions in three separate sets of tubes (calibrators, blanks, reference materials, dummy samples) for analysis.
- 2) Analyze each set of dilutions in separate runs on the same day using the same instrument.
- 3) Change the AFT value by ± 100 V of the optimized AFT value for the instrument.
 - a) Run 1: decrease the AFT from the optimized value by 100 V.
 - b) Run 2: method default AFT.
 - c) Run 3: increase the AFT from the optimized value by 100 V.

ii. Results: See Ruggedness Table 7.

iii. Conclusion: Accuracy of Mn, Hg and Se results are not compromised by changes in AFV settings within the range tested (optimized setting ± 100 V).

Ruggedness Table 5. Evaluating the significance of Axial Field Voltage on sample stability.

ID	axial field voltage	Hg* ($\mu\text{g/L}$)	Pb* ($\mu\text{g/dL}$)	Cd* ($\mu\text{g/L}$)	Mn* ($\mu\text{g/L}$)	Se ($\mu\text{g/L}$)
LB08707_W B2	Target Mean and 2SD Range	0.585 0.407 – 0.763	2.12 2.03 – 2.21	0.488 0.398 – 0.578	7.98 6.91 – 9.05	
	(optimized - 100V)	0.511	2.00	40.415	7.77	
	(optimized)	0.461	2.04	0.394	6.36	
	(optimized + 100V)	0.414	2.01	0.376	6.95	
HB08708_W B2	Target Mean and 2SD Range	6.19 5.89 – 6.48	10.1 9.84 – 10.3	3.14 2.94 – 3.34	14.9 13.5 – 16.4	
	(optimized - 100V)	5.50	9.8	2.91	14.3	
	(optimized)	5.62	9.8	2.84	12.0	
	(optimized + 100V)	5.75	10.1	2.99	12.8	
QMEQAS07B -09**	Target Mean and 2SD Range					157 146 – 168
	(optimized - 100V)					139
	(optimized)					147
	(optimized + 100V)					138
QMEQAS09B -08**	Target Mean and 2SD Range					548 511 - 585
	(optimized - 100V)					501
	(optimized)					556
	(optimized + 100V)					532
**Test performed on December 20, 2010 by Deanna Jones. Results are the average of the beginning and ending QC results for each analytical run.						
**samples purchase from Le centre de toxicology du Quebec (Quebec, Canada)						

F. Ruggedness parameter test #6: Extra dilution of samples

Evaluate the impact on observed concentration if an extra dilution is performed on the sample relative to the calibration standards.

i. Test details: A large blood sample was spiked to elevated concentrations, and mixed well. The spiked sample was then prepared for analysis at various extra dilution levels and the observed results compared to results obtained with no extra dilution performed.

ii. Results: See Ruggedness Table 8.

iii. Conclusions: Results show that all analytes of the method (Pb, Cd, Hg, Mn, and Se) can be analyzed at up to a 20x extra dilution without significant effect ($> \pm 10\%$ error) to the observed concentration.

Ruggedness Table 6. Evaluating the impact on observed concentration of an extra dilution performed on the sample relative to the calibration standards.

Dilution level	Mn	Hg	Se	Cd	Pb
No Extra (N=8)	1.00	1.00	1.00	1.00	1.00
2x dilution (N=8)	1.00 ± 0.01	1.03 ± 0.05	1.02 ± 0.03	1.00 ± 0.01	1.01 ± 0.01
5x dilution (N=6)	1.01 ± 0.01	1.06 ± 0.06	1.01 ± 0.02	1.01 ± 0.01	1.02 ± 0.01
10x dilution (N=8)	1.01 ± 0.03	1.04 ± 0.06	1.04 ± 0.06	1.00 ± 0.02	1.02 ± 0.02
20x dilution (N=8)	1.02 ± 0.04	1.09 ± 0.05	1.06 ± 0.08	1.01 ± 0.03	1.02 ± 0.02
Results are the average normalized concentration of N measurements ± 1RSD.					

17. Appendix C. Tables and Figures

Table 1. Instrument and Method Parameters

Instrument: PerkinElmer ELAN DRC II ICP-MS ESI SC4 autosampler with (optional) PC3 Peltier cooled spray chamber	
Optimization window parameters	
RF power	1450 W
Plasma Gas Flow (Ar)	15 L/min
Auxiliary Gas Flow (Ar)	1.2 L/min
Nebulizer Gas Flow (Ar)	~0.90 – 1.0 L/min (optimized as needed for sensitivity)
Ion Lens Voltage(s)	AutoLens (optimized as needed for sensitivity)
AFV, QRO, CRO, CPV, Discriminator Threshold	Optimized per instrument by service engineer, or advanced user.
Parameters of x-y alignment, nebulizer gas flow, AutoLens voltages, mass calibration, dual detector calibration and detector voltages are optimized regularly. Optimization file name = default.dac.	
Configurations window parameters	
cell gas changes pause times	Pressurize Delay (From Standard to DRC mode) = 60 Exhaust Delay (From DRC to Standard mode) = 30 Flow Delay (Gas changes while in DRC mode) = 30 Channel Delay (Gas channel change in DRC mode) = 30
File names & directories	
method file names	<i>calibration curve (programmed for blood blank)</i> CDC_DLS3016_bldblk.mth <i>For QC & patient sample analysis (programmed for aqueous blank)</i> CDC_DLS3016_aqblk.mth
dataset	Create a new dataset subfolder each day. Name as “2011-0820” for all work done on August 20, 2011
sample file	Create for each day’s work
report file name (See Figure 8 in Appendix C)	<i>For sample results printouts</i> cdc_quant comprehensive.rop <i>For calibration curve information</i> CDC_Quant Comprehensive (calib curve info).rop
tuning	Default.tun
optimization	Default.dac
calibration	N/A
polyatomic	elan.ply

report options template (transferring results to the database)	CDC_Database Output.rop <i>Report Format Options: select only "Use Separator"</i> <i>File Write Option: Append</i> <i>Report File name: make descriptive including date</i> <i>(e.g. 2005-0311b_DRC2A_HM-0364.txt)</i>
Method Parameters	
Method Parameters: Timing Page (see Figure 2 of Appendix C)	
sweeps/reading	30
readings/replicate	1
replicates	3
enable qc checking	On
isotopes monitored and internal standard associations (exact mass)	use ¹⁰³ Rh, ¹³⁰ Te, ¹⁹³ Ir as internal standards ¹⁰³ Rh (102.905): ⁵⁵ Mn (54.93805) ¹³⁰ Te(129.907): ²⁰² Hg (201.971), ⁸⁰ Se(79.9165) ¹⁹³ Ir(192.963): ²⁰⁸ Pb(207.977), ¹¹⁴ Cd(113.904)
dwelt times	100 ms for ⁵⁵ Mn, ²⁰² Hg, ⁸⁰ Se, ²⁰⁸ Pb, and ¹¹⁴ Cd 50 ms for ¹³⁰ Te, ¹⁰³ Rh, and ¹⁹³ Ir
scan mode	Peak Hopping for all isotopes (1 MCA channel)
DRC channel A gas flow rate	methane (5-7 psig delivery pressure) typically 0.84 L/min (0.7 – 1.0) * *optimized per instrument, and periodically verified
DRC channel B gas flow rate	oxygen (5-7 psig delivery pressure) typically 1.2 L/min (0.96 – 1.44) * *optimized per instrument, and periodically verified
RPa	0 for all isotopes
RPq	Typically* 0.6 (0.48 – 0.72) for ¹⁰³ Rh, ⁵⁵ Mn, ¹³⁰ Te, and ²⁰² Hg. 0.65 (0.52 – 0.78) for ¹³⁰ Te and ⁸⁰ Se. 0.25 for ¹⁹³ Ir, ²⁰⁸ Pb, and ¹¹⁴ Cd Use the same RPQ for each analyte and its IS. (* Optimize per instrument, and periodically verified)
Method parameters: processing page (see Figure 3 of Appendix C)	
detector mode	Dual
process spectral peak	N/A
autolens	On
isotope ratio mode	Off
enable short settling time	Off
blank subtraction	After internal standard
measurement units	cps
process signal profile	N/A

Method parameters: equations page (see Figure 4 of Appendix C)	
equations	+Hg 200 -0.027250 * Sn118 +Pb 206 +Pb 207
Method parameters: calibration page (see Figure 5 of Appendix C)	
calibration type	external std.
curve type	weighted linear
sample units	"µg/L" or "ppb"
calibrator concentrations (µg/L)	Mn (µg /L): 1.5, 4.5, 10.5, 15, 30, 75, 225, 600 Cd and Hg (µg /L): 0.5, 1.5, 3.5, 5, 10, 25, 75, 200 Pb (µg /dL): 1, 3, 7, 10, 20, 50, 150, 400 Se (µg /L): 30, 90, 210, 300, 600, 1500, 4500, 12000
Method parameters: sampling page (see Figures 6 and 7 of Appendix C)	
"peristaltic pump under computer control"	On
autosampler tray port sampling device	<i>If using ESI autosampler</i> Autosampler Type: AS-93plus Tray Name: esi.try Sampling Device: None If using other autosampler, refer to user guide.
sample flush	default is 4s at 1.5 rpm (~160 uL/min, ESI DXi peristaltic pump, FAST sample introduction system) Time can be optimized as needed to adequately fill the FAST loop. Time and rpm can be optimized as needed to using a different style peristaltic pump (maintaining approximate liquid flow rate). As a matter of lab practice, set this time to equal the loop fill time in the ESI FAST program. As long as the combined time of sample flush + read delay is equal to the time required for signal to reach stability, analytical measurement will be good.

read delay	<p>60s at 1.5 rpm (~160 uL/min, ESI DXi peristaltic pump, FAST sample introduction system)</p> <p>Time can be optimized as needed to reach signal stability before beginning analysis. Time and rpm can be optimized as needed to using a different style peristaltic pump (maintaining approximate liquid flow rate). As a matter of lab practice, set this time equal to the total time required for the signal to reach stability minus the loop fill time. As long as the combined time of sample flush + read delay is equal to the time required for signal to reach stability, analytical measurement will be good.</p>																		
wash	<p>30s at 1.5 rpm (~160 uL/min, ESI DXi peristaltic pump, FAST sample introduction system)</p> <p>Time can be optimized to allow for changes in FAST loop rinsing (must be greater than total time of steps in FAST program after the initial “on rinse” command). Time and rpm can be optimized as needed to using a different style peristaltic pump (maintaining approximate liquid flow rate).</p>																		
QC / Sample : See Figures 6, 7, and 9 of Appendix C for details																			
<p>extra wash (via ICP-MS software QC checking)</p>	<p>For sample concentrations greater than these, setup the ICP-MS software’s ‘QC checking’ feature to “Wash for X and continue”</p> <table border="1" data-bbox="597 1119 1252 1346"> <thead> <tr> <th>Analyte</th> <th>Concentration</th> <th>Extra Rinse Time</th> </tr> </thead> <tbody> <tr> <td>Cd</td> <td>200 µg/L</td> <td>200s</td> </tr> <tr> <td>Hg</td> <td>200 µg/L</td> <td>200s</td> </tr> <tr> <td>Mn</td> <td>600 µg/L</td> <td>200s</td> </tr> <tr> <td>Pb</td> <td>400 µg/dL</td> <td>200s</td> </tr> <tr> <td>Se</td> <td>1200 µg/L</td> <td>200s</td> </tr> </tbody> </table>	Analyte	Concentration	Extra Rinse Time	Cd	200 µg/L	200s	Hg	200 µg/L	200s	Mn	600 µg/L	200s	Pb	400 µg/dL	200s	Se	1200 µg/L	200s
Analyte	Concentration	Extra Rinse Time																	
Cd	200 µg/L	200s																	
Hg	200 µg/L	200s																	
Mn	600 µg/L	200s																	
Pb	400 µg/dL	200s																	
Se	1200 µg/L	200s																	
<p>autosampler locations of blanks and standards</p>	<p><i>For calibration curve (points to blood blank)</i> CDC_DLS3016_bldblk.mth Calibration Stds 0 – 8 in autosampler positions 105 – 113 by default, but can be customized.</p> <p><i>For QC & patient sample analysis (points to aqueous blank)</i> CDC_DLS3016_aqblk.mth Aqueous Blank in autosampler position 117 by default, but can be customized.</p>																		
FAST parameters: See Figures 10, 11, 12, 13, 14, 15, and 16 of Appendix C for details																			
configuration file	<p>default.sc (saved at C:\Program Files\ESI\ESI-SC\)</p>																		
FAST program	cdc_dls3016_5element_loop1ml_scfast.txt																		

Potential Emergency Response Modifications:	
<u>mercury:</u>	Analyze mercury in standard mode with tellurium as the internal standard. Set dwell time to 100ms, DRC gas flow to 0, and RPq to 0.25.
<u>Non-FAST sample introduction system:</u>	<p>If the FAST sample introduction system is not available on any instruments, the method can still be implemented, but these changes will need to be made in the ELAN (and ESI software if present). Peristaltic pump speeds are for DXi pump; adjust accordingly if another pump is installed.</p> <ul style="list-style-type: none"> • <u>sample flush</u>: Default is ~30s at -16 rpm. Set so that solution reaches nebulizer. • <u>read delay</u>: Default is 45s at -5 rpm. Set for best reproducibility of replicate measured intensities. • <u>wash</u>: Default is 60s at -11rpm. Set to prevent significant carry-over from one sample to the next. • If using ESI autosampler without FAST, disable FAST in the ESI software before running analysis.

Table 2. Suggested concentrations for base blood

analyte (units)	suggested concentration
Cd (µg/L)	<0.5
Hg (µg/L)	<0.5
Mn (µg/L)	< 8
Pb (µg/dL)	<2
Se (µg/L)	<200

Table 3. Stock calibration standard concentrations.

Analyte	Stock calibration concentration (mg/L) High Purity Standards Item # SM-2107-042 10% v/v HCl
Cd	50
Hg	50
Mn	150
Pb	1000
Se	3000

Table 4. Preparation of intermediate stock calibration standard.

volume of flask (mL)	100
volume of spike of stock standard solution (mL)	2
	concentrations (mg /L)
Cd	1
Hg	1
Mn	3
Pb (mg /dL)	20
Se	60

Table 5. Preparation of intermediate working calibration standards

Standard #	1	2	3	4	5	6	7	8
volume of flask (mL)	100	100	100	100	100	100	100	100
volume spike of stock std. (mL)						0.05	0.15	0.4
volume spike of int. stock Std. (mL)	0.05	0.15	0.35	0.50	1.00			
concentrations (µg /L) *								
Cd	0.5	1.5	3.5	5	10	25	75	200
Hg	0.5	1.5	3.5	5	10	25	75	200
Mn	1.5	4.5	10.5	15	30	75	225	600
Pb (µg /dL)	1	3	7	10	20	50	150	400
Se	30	90	210	300	600	1500	4500	12000

* These same concentrations are entered in the ICP-MS software's calibration page to describe the concentrations of the working calibrators (preparations analyzed during a run). This eliminates the need to multiply ICP-MS observed results by a dilution factor except for the case of extra dilutions (see Table 8 of Appendix C).

Table 6. Acceptable ways to perform two consecutive analytical runs, bracketing with bench quality control samples.

setup 1	setup 2
<p><i>Run #1</i> calibration standards low bench QC high bench QC elevated bench QC patient samples low bench QC high bench QC elevated bench QC</p> <p><i>Run #2</i> low bench QC high bench QC elevated bench QC patient samples low bench QC high bench QC elevated bench QC</p>	<p><i>Run #1</i> calibration standards low bench QC high bench QC elevated bench QC patient samples low bench QC high bench QC elevated bench QC</p> <p><i>Run #2</i> <i>calibration standards</i> low bench QC high bench QC elevated bench QC patient samples low bench QC high bench QC elevated bench QC</p>

Table 7. A typical SAMPLE/BATCH window.

<u>AS Location*</u>	<u>Sample ID</u>	<u>Measurements Action</u>	<u>Method</u>
233	DRCstability1	Run sample	...DLS3016_bldblk.mth
233	DRCstability2	Run sample	...DLS3016_bldblk.mth
233	DRCstability3	Run sample	...DLS3016_bldblk.mth
233	DRCstability4	Run sample	...DLS3016_bldblk.mth
Continue DRC stability samples . . .			
233	DRCstability9	Run sample	...DLS3016_bldblk.mth
233	DRCstability10	Run sample	...DLS3016_bldblk.mth
114	3016 BldBlkChk1	Run blank, standards, and sample **	...DLS3016_bldblk.mth
115	3016 BldBlkChk2	Run sample	...DLS3016_bldblk.mth
116	3016 AQBLK	Run blank and sample †	...DLS3016_aqblk.mth
125	L Bench QC	Run sample	...DLS3016_aqblk.mth
126	H Bench QC	Run sample	...DLS3016_aqblk.mth
127	E Bench QC	Run sample	...DLS3016_aqblk.mth
137	Sample 1	Run sample	...DLS3016_aqblk.mth
138	Sample 2	Run sample	...DLS3016_aqblk.mth
125	L Bench QC	Run sample	...DLS3016_aqblk.mth
126	H Bench QC	Run sample	...DLS3016_aqblk.mth
127	E Bench QC	Run sample	...DLS3016_aqblk.mth
<p>* The exact autosampler positions of QCs and patient samples do not have to be those shown above. QC samples do not have to be run in the order of low, then high, then elevated.</p> <p>** When executing this row, the ELAN will first analyze the standard 0 (blood blank) at AS position 105, then standards 1-8 at autosampler positions 106-113, <u>then</u> the “3016 BldBlkChk1” sample at A/S position 114. The sampling information about AS positions 105-113 are stored in the “bldblk” method file.</p> <p>† When executing this row, the ELAN will first analyze the aqueous blank at AS position 117, then the “Aq blank ” at AS position 103. The sampling information about AS positions 117 is stored in the “aqblk” method file.</p>			

Table 8. Preparation of samples, working calibrators, and QC materials for analysis *.

Description	Water (µL)	Base Blood (µL)	AQ Intermediate Working Calibration Standard (µL)	Patient or QC blood sample (µL)	Diluent (µL)**	Total volume (µL)
Working Calibrators (S0-S8) and Bldblkchk (S0)	-	50 x 1	50 x 1	-	2,400	2,500
AQ Blank	100 x 1	-	-	-	2,400	2,500
Patient blood or Blood-Based QC	50 x 1	-	-	50 x 1	2,400	2,500
Patient Blood <i>2x Extra Dilution</i> ^H	150 x 1	-	-	50 x 1	4,800	5,000
Patient Blood <i>5x Extra Dilution</i> ^H	450 (225 x 2)			50 x 1	12,000	12,500
Patient Blood <i>10x Extra Dilution</i> ^H	950 (190 x 5)			50 x 1	24,000	25,000
Patient Blood <i>20x Extra Dilution</i> ^H	1950 (195 x 10)			50 x 1	48,000	50,000
<p>If a different total volume is prepared, adjust the volumes for each component proportionally.</p> <p>* These directions are written with the expectation of a 5,000 µL syringe on the left side and a 250 µL syringe on the right side of the benchtop automatic pipette.</p> <p>** By splitting the dispense step of diluent into two or more portions, liquids pulled up into the right pipette tip are flushed out more completely. For example, when preparing a working calibrator, do the preparation in two steps: in step 1, dispense 2400 µL diluent + 50 µL; in step 2, dispense 2400 µL diluent + 50 µL base blood to prepare a 2.5 mL total volume dilution.</p> <p>^H Extra dilution is performed on blood samples whose concentration is greater than the concentration of the highest calibrator listed in Table 9 of Appendix C. Any extra dilution within these limits can be prepared as long as the 4.8:5 ratio of diluent to total dilution volume is maintained. Use of the lowest possible dilution level is preferred to minimize differences between the calibrators and the samples (i.e. 2x dilution is preferred over 5x if 2x is sufficient to dilute analyte into the documented linearity range).</p>						

Table 9. Limit of Detection, highest calibrator concentration, maximum allowable extra dilutions, and upper reportable range

Analyte	Limit of Detection (LOD)*	High Calibrator	Maximum Extra Dilution**	Reportable Range Upper Boundary
Mn (µg/L)	0.99	600	20	12,000
Pb (µg/dL)	0.07	400	20	8,000
Cd (µg/L)	0.10	200	20	4,000
Hg (µg/L)	0.28	200	20	4,000
Se (µg/L)	24.48	12,000	20	240,000

*Re-evaluated periodically (2+ years) or at significant method changes. LODs shown were calculated 1/24/2013.
**See ruggedness test 6 in Appendix B for supporting validation data.

Table 10. Boundary concentrations for whole blood concentrations.

Analyte (units)	1 st upper boundary ("1UB") *	2 nd upper boundary ("2UB") **	Range maximum ("Lim Rep Delta") †	Highest Concentration Validated for Washout
Mn (µg/L)	20	35	2.0 for values <20 10% of value at ≥20	600
Pb (µg/dL)	5.0	5.0	1.0 for values < 10 10% of value at ≥10	400
Cd (µg/L)	5.0	5.0	1.0 for values <10 10% of value at ≥10	200
Hg (µg/L)	10.0	10.0	1.0 for values < 10 10% of value at ≥10	200
Se (µg/L)	400	400	20 for values <200 10% of value at ≥200	12,000

* Typically, the 1st upper boundary (1UB) is the 99th percentile of non-weighted concentration results from the NHANES subset groups; a concentration significant to public health; or a concentration defined by study protocol. The default 1UB concentrations are listed in this table.
** The 2nd upper boundary (2UB) may be 2x the 1UB; a concentration significant to public health; or defined by study protocol.
† Range maximum is the range of the three replicate readings for a single sample analysis also called the "Lim Rep Delta" or "Rep Delta" in the division LIMS.

Table 11. Reference ranges for blood concentrations [70].

analyte (units)	survey years	geometric mean	50th	75th	90th	95th	N
Cd (µg/L)	07-08	0.315	0.270	0.500	1.00	1.52	8266
	09-10	0.302	0.260	0.480	0.960	1.40	8793
	11-12	0.279	0.250	0.460	0.960	1.50	7920
	13-14	0.235	0.210	0.410	0.840	1.22	5215
Hg (µg/L)	07-08	0.769	0.740	1.48	2.95	4.64	8266
	09-10	0.863	0.790	1.68	3.43	5.13	8793
	11-12	0.703	0.640	1.38	2.87	4.40	7920
	13-14	0.683	0.620	1.29	2.65	4.36	5215
Pb (µg/dL)	07-08	1.27	1.22	1.90	2.80	3.70	8266
	09-10	1.12	1.07	1.70	2.58	3.34	8793
	11-12	0.973	0.930	1.52	2.38	3.16	7920
	13-14	0.858	0.830	1.32	2.10	2.81	5215
Mn (µg/L)	11-12	9.35	9.22	11.5	14.4	16.7	7920
	13-14	9.52	9.41	11.8	14.6	16.7	5215
Se (µg/L)	11-12	190	190	206	223	236	7920
	13-14	193	193	208	223	235	5215

Figure 16. ESI SC4 autosampler screen shots (“FAST” page)


Timer A can be optimized to achieve proper filling of loop with diluted sample digestate. Timers B, C, D, E, and F control rinsing the loop after analysis and can be optimized for eliminating carry-over from one sample to the next while using the minimum amount of rinse solution. Save the file with the name “DLS 3016.8 FAST parameters.txt”. It can be found in the directory C:\Program Files\ESI\ESI-SC\.

Manually clicking the “Load” button prior to starting analysis will ensure the position of the actuator is always the same at the beginning of the analysis.

Manually clicking the “Vacuum On” button prior to starting the analysis will help initial sample uptake to be consistent (the vacuum pump may be slow to start for the first sample if this is not done, possibly resulting in loop filling inconsistencies).

Event	Action	Parameters	Parameter Units	Event Parameter
On Probe Down	Vacuum1 On			
On Probe Down	Load1			
Probe In Sample	Timer A	4	seconds	
Timer A Expires	Inject1			
Timer A Expires	Move Rinse			
Rinse Completed	Probe Up			
On Rinse	Load1			
On Rinse	Probe Down			
On Rinse	A2 On			
On Rinse	Timer B	2	seconds	
Timer B Expires	Probe Up			
Timer B Expires	Timer C	2	seconds	
Timer C Expires	Probe Down			
Timer C Expires	Timer D	2	seconds	
Timer D Expires	Probe Up			
Timer D Expires	Timer E	2	seconds	
Timer E Expires	Probe Down			
Timer E Expires	Timer F	2	seconds	
Timer F Expires	Probe Up			
Timer F Expires	A2 Off			
Timer F Expires	Move Next			

Figure 1. Chart for handling an elevated result.



18. Appendix D. Help Sheets

Reagent Preparation

NOTE:

mg/L = ppm

µg/L = ppb

µg/mL = ppm

Rinse solution

(0.4% (v/v) TMAH, 0.05% Triton™ X-100, 1% ethanol, 0.01% APDC)

1. Partially fill a 4 liter bottle with ≥ 18 Mohm·cm water..
2. Add 0.4 grams of APDC.
3. Add 16 mL of 25% v/v TMAH (Tetramethylammonium hydroxide).
4. Add 40 mL of ethanol (200 proof).
5. Add 200 mL of 1% Triton™ X-100 (OR add 10mL of 20% Triton™ X-100).
6. Add enough ≥ 18 Mohm·cm water to bring to 4 liter mark.
7. Mix well by gently inverting several times.

Sample diluent

(0.4% (v/v) TMAH, 0.01% APDC, 0.05% Triton™ X-100, 1% ethanol, 5 ppb Te, Rh, Ir)

1. Partially fill a 2 liter bottle with ≥ 18 Mohm·cm water.
2. Add 0.2 grams of APDC.
3. Add 8 mL of 25% v/v TMAH.
4. Add 20 mL of ethanol.
5. Add 500 µL of a 20 mg/L stock solution of Te, Rh, and Ir.
6. Add 100 mL of 1% Triton™ X-100 (OR, if using a 20% Triton™ X-100 solution, add 5mL)
7. Add enough ≥ 18 Mohm·cm water to bring to 2 liter mark.
8. Mix well by gently inverting several times.

0.5% v/v HNO₃

(Carrier solution for optimization)

1. Partially fill a 2 liter bottle with ≥ 18 Mohm·cm water.
2. Add 10 mL of concentrated HNO₃.
3. Add enough ≥ 18 Mohm·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.

Reagent Preparation (continued)

1% v/v HNO₃

1. Partially fill a 4 liter bottle with ≥ 18 Mohm·cm water.
2. Add 40 mL of concentrated HNO₃.
3. Add enough ≥ 18 Mohm·cm water to bring to 4 liter mark.
4. Mix well by gently swirling several times.

5% v/v HNO₃

1. Partially fill a 2 liter bottle with ≥ 18 Mohm·cm water.
2. Add 100 mL of concentrated HNO₃.
3. Add enough ≥ 18 Mohm·cm water to bring to 2 liter mark.
4. Mix well by gently inverting several times.

20% Triton™ X-100

1. Partially fill a 1 liter bottle with ≥ 18 Mohm·cm water.
2. Add 200 mL of Triton™ X-100.
3. Add enough ≥ 18 Mohm·cm water to bring to 1 liter mark.
4. Allow to dissolve overnight (or add a Teflon magnetic stirring bar and stir on stirrer until dissolved).
5. Mix well by gently inverting several times.

1% Triton™ X-100

9. Partially fill a 1 liter bottle with ≥ 18 Mohm·cm water.
10. Add 10 mL of Triton™ X-100.
11. Add enough ≥ 18 Mohm·cm water to bring to 1 liter mark.
12. Allow to dissolve overnight (or add a Teflon magnetic stirring bar and stir on stirrer until dissolved).
13. Mix well by gently inverting several times.

Reagent Preparation (continued)

20 ppm Rh, Te and Ir internal standard solution

1. Partially fill an acid rinsed, 50 mL flask with 1% v/v HNO₃.
2. Add 1 mL of Rh from 1000 ppm stock standard.
3. Add 1 mL of Te from 1000 ppm stock standard.
4. Add 1 mL of Ir from 1000 ppm stock standard.
5. Add enough 1% v/v HNO₃ to fill to 50 mL mark.
6. Mix well by gently inverting several times.
7. Pour the standard solution over into an appropriately labeled 50 mL polypropylene tube.

Daily solution (1ppb) in 2% v/v HNO₃

1. Partially fill a 1 liter volumetric flask with ≥ 18 Mohm·cm water.
2. Add 1 mL of High Purity Standard: SM-2107-018 (or current lot #)
3. Add 20 mL of concentrated HNO₃.
4. Add enough ≥ 18 Mohm·cm water to bring to 1 liter mark.
5. Mix well by gently inverting several times.

Stability test solution (1 liter bulk prep)

1. Use a 1 liter bottle dedicated to stability test solution preparation
2. Add 960 mL of Sample Diluent
3. Add 20 mL of “junk” whole blood
4. Add 20 mL of Intermediate Working Calibration Standard (may use S1 or S2)
OR add 1.5mL of Intermediate Stock Calibration Standard.
5. Mix well by gently inverting several times.
6. Store in the refrigerator (when not using).

Standard Preparation (from single element stock standards)

Prepare 3% v/v HCl solution:

1. Store in the refrigerator (when not using).
2. Partially fill a clean 2 liter bottle with ≥ 18 Mohm·cm water.
3. Using a clean 50 mL polypropylene tube to measure, add 60 mL of high purity concentrated HCl.
4. Add enough ≥ 18 Mohm·cm water to bring to 2 liter mark.
5. Gently invert to mix.

Prepare intermediate stock standard (see Table 4 of Appendix C):

1. Partially fill a 100 mL volumetric flask with 3% v/v HCl solution.
2. Label as: "HgPbCdMnSe Intermediate Stock Std"
3. Add 2 mL of HgPbCdMnSe multi-element stock solution.
4. Add enough 3% v/v HCl to bring to 100 mL mark.
5. Mix well by gently inverting several times.

Prepare intermediate working calibration standards (see Table 5 of Appendix C):

1. Partially fill each of eight, 100 mL volumetric flasks with 3% v/v HCl solution.
2. Label as: Intermediate Working Std "S1", "S2", "S3" and "S4", "S5", "S6", "S7" and "S8".
3. For "S1 Intermediate Working Std": add 50 μ L of the Intermediate Stock Std.
4. For "S2 Intermediate Working Std": add 150 μ L of the Intermediate Stock Std.
5. For "S3 Intermediate Working Std": add 350 μ L of the Intermediate Stock Std.
6. For "S4 Intermediate Working Std": add 500 μ L of the Intermediate Stock Std.
7. For "S5 Intermediate Working Std": add 1 mL of the Intermediate Stock Std.
8. For "S6 Intermediate Working Std": add 50 μ L of the Multi-Element Stock Std.
9. For "S7 Intermediate Working Std": add 150 μ L of the Multi-Element Stock Std.
10. For "S8 Intermediate Working Std": add 400 μ L of the Multi-Element Stock Std.
11. Add enough 3% v/v HCl solution to bring to 100 mL mark.
12. Mix well by gently inverting several times.
13. These intermediate working calibration standards may be poured over into clean 15 mL Falcon tubes for daily use (NOTE: "S0 Intermediate Working Std" is 3% v/v HCl only).

19. References

1. Pirkle, J.L., et al., *National exposure measurements for decisions to protect public health from environmental exposures*. International Journal of Hygiene and Environmental Health, 2005. **208**(1-2): p. 1-5.
2. *Agency for Toxic Substances and Disease Registry (ATSDR). Toxicological profile for Mercury*. 1999, U.S. Department of Health and Human Services, Public Health Service: Atlanta, GA.
3. Mahaffey, K.R. *NHANES 1999 - 2002 Update on Mercury*. in *Northeast Regional Mercury Conference*. 2005.
4. Drasch, G.A., *Mercury*, in *Handbook on metals in clinical and analytical chemistry*, H.G. Seiler, H. Sigel, and A. Sigel, Editors. 1994, Marcel Dekker: New York. p. 479-493.
5. World Health Organization, *Environmental Health Criteria 118: Inorganic Mercury*. 1991, Geneva.
6. Centers for Disease Control and Prevention, *Preventing Lead Poisoning in Young Children*. 2005.
7. Needleman, H., et al., *Bone lead levels in adjudicated delinquents. A case control study*. Neurotoxicology and teratology, 2002. **24**(6): p. 711-7.
8. Dietrich, K., et al., *Early exposure to lead and juvenile delinquency*. Neurotoxicology and teratology, 2001. **23**(6): p. 511-518.
9. Bellinger, D.C., *Low-level lead exposure, intelligence and academic achievement: A long-term follow-up study*. Pediatrics, 1992. **90**(6): p. 855-861.
10. Bellinger, D.C., *Intellectual Impairment and Blood Lead Levels*. The New England Journal of Medicine, 2003. **349**(5): p. 500-502.
11. Goyer, R.A., *Lead*, in *Handbook on toxicity of inorganic compounds*, H.G. Seiler, S. Helmut, and S. Astrid, Editors. 1994, Marcel Dekker: New York. p. 359-382.
12. Batley, G.E., *Handbook of Trace Element Speciation: Analytical Methods*. 1991, Boca Raton: CDC Press.
13. *Agency for Toxic Substances and Disease Registry (ATSDR). Toxicological profile for Lead*. 2007, U.S. Department of Health and Human Services, Public Health Service: Atlanta, GA.
14. *Centers for Disease Control and Prevention. CDC Response to Advisory Committee on Childhood Lead Poisoning Prevention. Recommendations in "Low Level Lead Exposure Harms Children: A Renewed Call of Primary Prevention", 2012. Atlanta, GA: U.S. Department of Health and Human Services, Centers for Disease Control and Prevention.*
http://www.cdc.gov/nceh/lead/acclpp/cdc_response_lead_exposure_recs.pdf Accessed 30 June 2016.
15. World Health Organization, *Environmental Health Criteria 134: Cadmium*. 1992.
16. Elinder, C.G., International Journal of Environmental Studies, 1982. **19**(3-4): p. 187-193.
17. Herber, R.F.M., *Cadmium*, in *Handbook on metals in clinical and analytical chemistry*, H.G. Seiler, H. Sigel, and A. Sigel, Editors. 1994, Marcel Dekker: New York. p. 283-297.

18. Ghezzi, I., et al., *Behavior of biological indicators of cadmium in relation to occupational exposure*. International archives of occupational and environmental health, 1985. **55**(2): p. 133-140.
19. Jarup, L., C. Elinder, and G. Spang, *Cumulative blood-cadmium and tubular proteinuria: a dose-response relationship*. International archives of occupational and environmental health, 1988. **60**(3): p. 223-229.
20. Lauwerys, R., et al., *Cadmium - Exposure Markers as Predictors of Nephrotoxic Effects*. Clinical Chemistry, 1994. **40**(7B): p. 1391-1394.
21. Roels, H., et al., *Health significance of cadmium induced renal dysfunction: a five year follow up*. British journal of industrial medicine, 1989. **46**(11): p. 755-764.
22. Bernard, A. and R. Lauwerys, *Cadmium in human population, in Cadmium in the Environment*. 1986, Springer. p. 114-123.
23. Milne, D.B., *Trace Elements*, in *Tietz textbook of clinical chemistry*, C.A. Burtis, Ashwood, Edward R., Editor. 1999, W. B. Saunders Company: Philadelphia. p. 1029-1055.
24. Chiswell, B. and D. Johnson, *Manganese*, in *Handbook on Metals in Clinical and Analytical Chemistry*, A.S. Hans G. Seiler, Helmut Sigel, Editor. 1994, Marcel Dekker: New York. p. 467-478.
25. Smargiassi, A., et al., *Peripheral Markers of Catecholamine Metabolism among Workers Occupationally Exposed to Manganese (Mn)*. Toxicology Letters, 1995. **77**(1-3): p. 329-333.
26. Roels, H.A., et al., *Assessment of the permissible exposure level to manganese in workers exposed to manganese dioxide dust*. Occupational and Environmental Medicine, 1992. **49**(1): p. 25-34.
27. Cowan, D.M., et al., *Manganese exposure among smelting workers: blood manganese-iron ratio as a novel tool for manganese exposure assessment*. Biomarkers, 2009. **14**(1): p. 3-16.
28. Gennart, J.P., et al., *Fertility of Male Workers Exposed to Cadmium, Lead, or Manganese*. American Journal of Epidemiology, 1992. **135**(11): p. 1208-1219.
29. Bader, M., et al., *Biomonitoring of manganese in blood, urine and axillary hair following low-dose exposure during the manufacture of dry cell batteries*. International Archives of Occupational and Environmental Health, 1999. **72**(8): p. 521-527.
30. Lauwerys, R., et al., *Fertility of Male Workers Exposed to Mercury-Vapor or to Manganese Dust - a Questionnaire Study*. American Journal of Industrial Medicine, 1985. **7**(2): p. 171-176.
31. Standridge, J.S., et al., *Effect of Chronic Low Level Manganese Exposure on Postural Balance: A Pilot Study of Residents in Southern Ohio*. Journal of Occupational and Environmental Medicine, 2008. **50**(12): p. 1421-1429.
32. Woolf, A., et al., *A child with chronic manganese exposure from drinking water*. Environmental Health Perspectives, 2002. **110**(6): p. 613-616.
33. Wasserman, G.A., et al., *Water manganese exposure and children's intellectual function in Araihasar, Bangladesh*. Environmental Health Perspectives, 2006. **114**: p. 124-129.

34. Ljung, K.S., et al., *Maternal and Early Life Exposure to Manganese in Rural Bangladesh*. Environmental Science & Technology, 2009. **43**(7): p. 2595-2601.
35. Bazzi, A., J.O. Nriagu, and A.M. Linder, *Determination of toxic and essential elements in children's blood with inductively coupled plasma-mass spectrometry*. Journal of Environmental Monitoring, 2008. **10**(10): p. 1226-1232.
36. Rollin, H.B., et al., *Examining the association between blood manganese and lead levels in schoolchildren in four selected regions of South Africa*. Environmental Research, 2007: p. 160-167.
37. Rollin, H., et al., *Blood manganese concentrations among first-grade schoolchildren in two South African cities*. Environmental Research, 2005. **97**(1): p. 93-99.
38. Aschner, M., *Manganese: Brain transport and emerging research needs*. Environmental Health Perspectives, 2000. **108**: p. 429-432.
39. Yokel, R.A., *Brain uptake, retention, and efflux of aluminum and manganese*. Environmental Health Perspectives, 2002. **110**: p. 699-704.
40. Davis, J.M., *Methylcyclopentadienyl manganese tricarbonyl: Health risk uncertainties and research directions*. Environmental Health Perspectives, 1998. **106**: p. 191-201.
41. Davis, J.M., et al., *The EPA health risk assessment of methylcyclopentadienyl manganese tricarbonyl (MMT)*. Risk Analysis, 1998. **18**(1): p. 57-70.
42. Roels, H., et al., *Relationship Between External and Internal Parameters of Exposure to Manganese in Workers From a Manganese Oxide and Salt Producing Plant*. American journal of industrial medicine, 1987. **11**(3): p. 297-305.
43. Jarvisalo, J., et al., *Urinary and blood manganese in occupationally nonexposed populations and in manual metal arc welders of mild-steel*. International archives of occupational and environmental health, 1992. **63**(7): p. 495-501.
44. Smyth, L., et al., *Clinical manganism and exposure to manganese in the production and processing of ferromanganese alloy*. Journal of occupational medicine, 1973. **15**(2): p. 101-9.
45. Klaassen, C., *Biliary-Excretion of Manganese in Rats, Rabbits, and Dogs*. Toxicology and applied pharmacology, 1974. **29**(3): p. 458-468.
46. Malecki, E., et al., *Biliary manganese excretion in conscious rats is affected by acute and chronic manganese intake but not by dietary fat*. The Journal of nutrition, 1996. **126**(2): p. 489-498.
47. Agency for Toxic Substances and Disease Registry (ATSDR), *Toxicological profile for manganese*. September 2012: Atlanta, GA.
48. Agency for Toxic Substances and Disease Registry (ATSDR). *Toxicological profile for Selenium.*, ATSDR, Editor. 2003, U.S. Department of Health and Human Services, Public Health Service: Atlanta, GA.
49. Goldhaber, S.B., *Trace element risk assessment: essentiality vs. toxicity*. Regulatory Toxicology and Pharmacology, 2003. **38**(2): p. 232-242.

50. Combs, G.F. and W.P. Gray, *Chemopreventive agents: Selenium*. Pharmacology & Therapeutics, 1998. **79**(3): p. 179-192.
51. Arthur, J.R., *The role of selenium in thyroid hormone metabolism*. Canadian Journal of Physiology and Pharmacology, 1991. **69**(11): p. 1648-1652.
52. Corvilain, B., et al., *Selenium and the thyroid: how the relationship was established*. American Journal of Clinical Nutrition, 1993. **57**(2): p. S244-S248.
53. Levander, O.A., *Nutrition and newly emerging viral diseases: An overview*. Journal of Nutrition, 1997. **127**: p. S948-S950.
54. McKenzie, R.C., T.S. Rafferty, and G.J. Beckett, *Selenium: An essential element for immune function*. Immunology Today, 1998. **19**(8): p. 342-345.
55. Ellis, D.R. and D.E. Salt, *Plants, selenium and human health*. Current Opinion in Plant Biology, 2003. **6**(3): p. 273-279.
56. Combs, G.E., *Food system-based approaches to improving micronutrient nutrition: The case for selenium*. Biofactors, 2000. **12**(1-4): p. 39-43.
57. Zimmermann, M.B. and J. Köhrle, *The Impact of Iron and Selenium Deficiencies on Iodine and Thyroid Metabolism: Biochemistry and Relevance to Public Health*. Thyroid, 2002. **12**(10): p. 867-878.
58. Beck, M.A., O.A. Levander, and J. Handy, *Selenium deficiency and viral infection*. Journal of Nutrition, 2003. **133**(5): p. 1463S-1467S.
59. Tanner, S.D., V.I. Baranov, and D.R. Bandura, *Reaction cells and collision cells for ICP-MS: a tutorial review*. Spectrochimica Acta Part B-Atomic Spectroscopy, 2002. **57**(9): p. 1361-1452.
60. Tanner, S.D. and V.I. Baranov, *Theory, design, and operation of a dynamic reaction cell for ICP-MS*. Atomic Spectroscopy, 1999. **20**(2): p. 45-52.
61. Lutz, T.M., P.M.V. Nirel, and B. Schmidt, *Whole Blood Analysis by ICP-MS*, in *Applications of Plasma Source Mass Spectrometry*. 1991, Royal Society of Chemistry. p. 96-100.
62. Burguera, J.L., et al., *Electrothermal atomic absorption spectrometry determination of molybdenum in whole blood*. Spectrochimica Acta - Part B Atomic Spectroscopy, 2002. **57**(3): p. 561-569.
63. Jarrett, J.M., et al., *Eliminating molybdenum oxide interference in urine cadmium biomonitoring using ICP-DRC-MS*. Journal of Analytical Atomic Spectrometry, 2008. **23**(7): p. 962-967.
64. Division of Laboratory Sciences, *Division of Laboratory Sciences Policies and Procedures Manual*. 2017, version 6.0, Centers for Disease Control and Prevention: Atlanta, GA.
65. *Occupational Safety and Health Administration, Occupational Safety and Health Standards, in 29 CFR part 1910, Subpart Z, Standard number 1910.1025, "Lead",. 1989.*
66. *Occupational Safety and Health Administration, Cadmium (OSHA 3136-06R 2004). 2004.*

67. *American Conference of Governmental Industrial Hygienists, Tlvs and Beis 2007: Based on the Documentation for Chemical Substances and Physical Agents & Biological Exposure Indices. 2007: American Conference of Governmental Industrial Hygienists.*
68. Baselt, R.C., *Disposition of Toxic Drugs and Chemicals in Man*. 2011, Seal Beach, CA: Biomedical Publications.
69. Nuttall, K.L., *Evaluating selenium poisoning*. *Annals of Clinical and Laboratory Science*, 2006. **36**(4): p. 409-420.
70. Centers for Disease Control and Prevention, *Fourth National Report on Human Exposure to Environmental Chemicals, Updated Tables, January 2017, Volume 1*. 2017, CDC: Atlanta, GA.