



Laboratory Procedure Manual

Analyte: **Folate Forms**

Matrix: **Serum**

Method: **Liquid Chromatography Tandem Mass Spectrometry**

Method No: **4013.06**

Revised: January 2023

as performed by: Nutritional Biomarkers Branch (NBB)
Division of Laboratory Sciences (DLS)
National Center for Environmental Health (NCEH)

contact: Zia Fazili-Qari, Ph.D.
Phone: 770-488-7581
Email: zxq0@cdc.gov

James L. Pirkle, M.D., Ph.D.
CLIA Laboratory Director
Centers for Disease Control and Prevention

Important Information for Users

CDC periodically refines these laboratory methods. It is the responsibility of the user to contact the person listed on the title page of each write-up before using the analytical method to find out whether any changes have been made and what revisions, if any, have been incorporated.

Images are included in this document as visual aids for certain topics. They are intended to be representative images only and should not be construed as absolute references. Discrepancies between the images in this document and the actual application design are not a cause for revisions to this document.

Public Release Data Set Information

This document details the Lab Protocol for testing the items listed in the following table for FOLFMS_L:

File Name	Variable Name	SAS Label
FOLFMS_L	LBDFOTSI	Serum total folate (nmol/L)
	LBDFOT	Serum total folate (ng/mL)
	LBXSF1SI	5-Methyl-tetrahydrofolate (nmol/L)
	LBDSF2SI	Folic acid (nmol/L)
	LBXSF3SI	5-Formyl-tetrahydrofolate (nmol/L)
	LBXSF4SI	Tetrahydrofolate (nmol/L)
	LBXSF5SI	5,10-Methenyl-tetrahydrofolate (nmol/L)
	LBXSF6SI	Mefox oxidation product (nmol/L)

1. Summary of Clinical Relevance and Principle

A. Clinical Relevance

Folate belongs to the group of water-soluble B vitamins that occur naturally in food. It is required in cellular one carbon metabolism and hematopoiesis [1]. Prolonged folate deficiency leads to megaloblastic anemia. Low folate status has been causally linked to an increased risk in women of childbearing age to have an offspring with neural tube defects. Low folate status also increases plasma homocysteine levels, a potential risk factor for chronic diseases such as cardiovascular disease or cognitive function. Potential roles of folate and other B vitamins in modulating the risk for diseases (e.g., heart disease, cancer, and cognitive impairment) are under investigation.

Serum folate is an indicator of short-term status. The measurement of folate forms circulating in serum may further elucidate the role of folate vitamers relative to various health outcomes. The measurement of total folate (tFOL), which is the sum of the individual folate forms, provides information on the folate status of the individual. The primary circulating folate vitamer in serum is 5-methyltetrahydrofolic acid (5-methylTHF), while the actual bioactive form of folate is tetrahydrofolic acid (THF). Pteroylglutamic acid (PGA, also called folic acid) is primarily derived from supplements and fortified foods. If the intake of PGA exceeds 200 µg per meal, unmetabolized PGA may appear in serum.

B. Test Principle

Five folate forms, 5-methylTHF, PGA, THF, 5-formyltetrahydrofolic acid (5-formylTHF), 5,10-methenyltetrahydrofolic acid (5,10-methenylTHF), and one oxidation product of 5-methylTHF called MeFox (pyrazino-s-triazine derivative of 4- α -hydroxy-5-methylTHF) are measured by isotope-dilution high performance liquid chromatography coupled to tandem mass spectrometry (LC-MS/MS) [2]. This is a modification of a previously published method [3, 4] to add the measurement of an oxidation product [5] and to scale down the amount of specimen needed and increase the sample throughput [2]. The assay is performed by combining specimen (150 µL serum) with ammonium formate buffer and an internal standard mixture. Sample extraction and clean-up is performed by automated 96-probe solid phase extraction (SPE) using 96-well phenyl SPE plates and takes ~1 h for a 96-well plate. Folate forms are separated within 4 min using isocratic mobile phase conditions and measured by LC-MS/MS (7 min to next injection). Quantitation is based on peak area ratios interpolated against a five-point aqueous linear calibration curve using 1/x² weighting. The following analytes are quantified (Table I):

Table I: Analytes and abbreviations

Compound	Abbreviation	
	Scientific literature	Database analyte code ¹
5-Methyltetrahydrofolic acid	5-methylTHF	MET
5-Formyltetrahydrofolic acid	5-formylTHF	FOT
Tetrahydrofolic acid	THF	THF
5,10-Methenyltetrahydrofolic acid	5,10-methenylTHF	MYT
Pteroylglutamic acid (folic acid)	PGA	PGA
Pyrazino-s-triazine derivative of 4- α -hydroxy-5-methyltetrahydrofolate	MeFox	MFO
Total folate (sum of folate forms)	tFOL	FOL3

¹ MET, PGA, and MFO are the major analytes. FOT, THF, and MYT are the minor analytes.

2. Safety Precautions

Consider all serum specimens as potentially positive for infectious agents including HIV, hepatitis B and hepatitis C. We recommend the hepatitis B vaccination series for all analysts working with whole blood and /or serum. Observe universal precautions; wear protective gloves, lab coat, and safety glasses during all steps of this method. Discard any residual sample material by autoclaving after analysis is completed. Place all disposable plastic, glassware, and paper (pipette tips, autosampler vials, gloves etc.) that contact serum/blood in a biohazard autoclave bag and keep these bags in appropriate containers until sealed and autoclaved. Use disposable bench diapers during sample preparation and serum/ blood handling and discard after use. Also wipe down all contaminated work surface with 10% bleach solution or other appropriate disinfectant when work is finished.

Handle acids (e.g., formic and acetic acid) and bases (e.g., ammonium hydroxide; produces strong fumes; handle only in chemical fume hood) used in sample and reagent preparation with extreme care; they are caustic and toxic. Handle organic solvents (e.g., methanol, acetonitrile) only in a well-ventilated area or, as required, under a chemical fume hood.

Reagents and solvents used in this study include those listed in Section 6. Safety data sheets (SDSs) for all chemicals are readily available in the SDS section as hard copies in the laboratory. SDSs for other chemicals can be viewed at <http://www.ilpi.com/msds/index.html> or at <http://intranet.cdc.gov/ossam/workplace-safety/safety-practices/chemical-safety/index.html>.

Observe universal precautions (i.e., PPE) during operation of automated liquid handlers (e.g., Hamilton, Caliper-Zephyr); keep instrument doors always locked when it is running; wipe down contaminated decks (10% bleach solution or other appropriate disinfectant) after work is finished. Either turn off the instruments or keep in standby mode.

Additional information on hazard identification, risk evaluation, and risk mitigation for this method can be found in the method risk assessment form.

3. Computerization; Data System Management

- (A) During sample preparation and analysis, samples are identified by their sample ID. The sample ID is a number that is unique to each sample that links the laboratory information to demographic data recorded by those who collected the sample.
- (B) The raw data file and respective batch file from the tandem mass spectrometer are collected using the instrument software and stored on the instrument workstation. The data file and batch file are transferred to the CDC network where the data file is processed into a results file that is also saved on the CDC network. Results are typically generated by auto-integration but may require in some cases manual integration. The results file (including analyte and internal standard names, peak areas, retention times, sample dilution factor, data file name, acquisition time, etc.) is imported into a STARLIMS database for review of the patient data, statistical evaluation of the QC data, and approval of the results. For details, refer to **Appendix B_D: JA-4013-DR-01-Computerization and Data System Management**.
- (C) For NHANES, data is transmitted electronically on a regular basis (approximately weekly for 3-week turnaround analytes). Abnormal values are confirmed by the analyst, and codes for missing data are entered by the analyst and are transmitted as part of the data file. NCHS makes arrangements for the abnormal report notifications to the NCHS Survey Physician.

- (D) The batch and the raw data file from the instrument workstation are typically backed up to the CDC network after a run is completed. This is the responsibility of the analyst under the guidance of the project lead. Files stored on the CDC network are automatically backed up nightly by ITSO support staff.

4. Specimen Collection, Storage, and Handling Procedures; Criteria for Specimen Rejection

- (A) We recommend that specimen donors fast prior to specimen collection but fasting is not required.
- (B) Serum folate assays are performed on fresh or frozen serum. Ascorbic acid (0.5%) is sometimes added to serum prior to storage to improve folate stability.
- (C) A 500- μ L serum specimen is required to allow for repeat analysis and for automated sample pipetting; a volume of 150 μ L is required for analysis.
- (D) Serum specimens may be collected with regular red-top Vacutainers. The appropriate amount of serum is dispensed into a Nalgene cryovial or plastic screw-capped vial labeled with the participant's ID.
- (E) Specimens collected in the field should be kept cold and protected from light. After processing, specimens should be frozen and shipped on dry ice by overnight mail. Once received, samples should be kept frozen during "in-processing", which is typically completed within less than 4 hours and then stored at $\leq -50^{\circ}\text{C}$ for up to 15 business day until samples are transferred to the testing laboratory for longer storage at deep frozen conditions typically around -70°C (-50°C to -90°C). Folates are stable for only a few weeks if the specimen is kept frozen (-10°C to -30°C). Up to three short (2 h) freeze-thaw cycles cause only minor ($<10\%$) folate degradation [6].
- (F) Specimens should generally arrive frozen. Specimens received at ambient temperature ($+15^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) are rejected. Refrigerated samples may be used provided they are brought promptly from the site where the blood was collected. Results from hemolyzed serum specimens should be interpreted with caution because they may have falsely elevated values. Specimens exposed to light for longer than 8 h may undergo 10–20% folate degradation [7]. Therefore, specimens intended for folate analysis should be processed and stored frozen promptly if analysis is not to be performed within 8 h of collection.
- (G) If the specimen vial is received either damaged or leaking, the specimen is rejected due to suspected contamination. If insufficient sample volume for at least one analysis ($<200 \mu\text{L}$ for manual pipetting) is received, the specimen is rejected. Serum specimens with a volume either below or above of what is typically expected ($<550 \mu\text{L}$ or $>770 \mu\text{L}$) are suspect if ascorbic acid has been added to the specimen during processing. However, these specimens are analyzed and assigned an appropriate description and comment code during data review ("no reportable", code 98).
- (H) Specimen handling conditions are outlined in the DLS Policies and Procedures Manual. The protocol discusses in general collection and transport of specimens and the special equipment required. If there is more than one test of interest in the specimen and it needs to be divided, the appropriate amount of serum or plasma should be transferred into a sterile Nalgene cryovial labeled with the participant's ID; avoid cross-contamination.

- (I) A series of standard comment codes are available in the STARLIMS database to identify any issues related to sample quality (e.g., not enough specimen for repeat analysis (code 22; set no reportable); lab error, spills, contamination etc. (code 23; set no reportable); or instrument error/failure (code 24; repeat analysis). These codes can be used, along with test descriptions, to document why a result was not reported (specimen rejection) or that a result should be interpreted with caution based on the sample quality.

5. Procedures for Microscopic Examinations; Criteria for Rejection of Inadequately Prepared Slides

Not applicable for this procedure.

6. Preparation of Reagents, Calibration (Standards), Controls, and All Other Materials; Equipment and Instrumentation

Chemicals and other materials used in the preparation of reagents, calibrators, and quality control materials are tracked as indicated in **Appendix B_B: JA-4013-R&S-02-Reagent Tracking**. To facilitate tracking of solvents and other chemicals transferred from the original manufacturer containment into a secondary container (e.g., solvent bottle), include the expiration date provided by the manufacturer or the lot number on the secondary container.

A. Reagent Preparation

Prepare all reagents with 0.22 µm filtered (cellulose nitrate filters) deionized water with a resistivity of at least 18 MΩ·cm. Use Class A volumetric glassware where a volumetric flask is specified. Though each reagent preparation specifies a total volume of reagent prepared, these directions may be scaled up or down to prepare larger or smaller quantities if desired.

1) 1% Ammonium formate buffer, pH 3.2

Into a 1-L reagent bottle, add 980 mL deionized water (measured by graduated glass cylinder) and 10 mL concentrated formic acid (measured by graduated glass pipette). This mixture is titrated with ammonium hydroxide (30%) to adjust the pH to 3.2. Transfer the buffer into a 1-L graduated glass cylinder and make up to final volume with deionized water. Transfer back into the 1-L reagent bottle and degas under vacuum for 3–5 min. This buffer can be stored at ambient temperature (+15°C to +30°C) for ~1 month.

Note: 10x concentrated ammonium formate buffer can be prepared as above (use 100 mL of concentrated formic acid instead of 10 mL measured by graduated glass cylinder). This buffer can be stored at ambient temperature (+15°C to +30°C) for 6 months. At the time of use, dilute 100 mL of the 10x buffer to 1 L using deionized water and degas under vacuum.

a) Conditioning solvent for 96-well SPE plates

1% ammonium formate buffer (pH 3.2) described above is used as is to condition the sorbent of the 95-well SPE plate (1.3 mL/well). Typically, 200 mL are used daily. Other solvents used to condition the sorbent are 0.5 mL/well each of acetonitrile and methanol, respectively.

- b) **Solvent #1:** 1% ammonium formate buffer (pH 3.2) with 0.5% ascorbic acid

Add 0.5% ascorbic acid (0.5 g/100 mL) to this buffer at the time of use (prepared daily for single use).

- c) Buffer for folate calibrator Mix A and ISTD Mix B: 1% ammonium formate buffer (pH 3.2) with 0.1% ascorbic acid

Add 0.1% ascorbic acid (0.1 g/100 mL) to this buffer at the time of use (prepared daily for single use).

- d) Wash buffer (**Solvent #2**): 0.05% ammonium formate (pH 3.4)

Dilute 50 mL of 1% ammonium formate buffer, pH 3.2 to 1 L with deionized water in a graduated glass cylinder. Transfer into a 1-L reagent bottle, check pH, and degas under vacuum for 3–5 min. This buffer can be stored at ambient temperature (+15°C to +30°C) for 2 weeks. Add ascorbic acid powder to a final concentration of 0.1% (0.1 g/100 mL) prior to use (typically 100 mL prepared daily for single use).

- 2) Sample elution buffer (**Solvent #3**): 49% deionized water, 40% methanol, 10% acetonitrile, 1% acetic acid, 0.5% ascorbic acid

Using graduated glass cylinders (500 mL and 100 mL), add 490 mL of deionized water, 400 mL of methanol, and 100 mL of acetonitrile to a 1-L reagent bottle and degas under vacuum for 3–5 min. At the time of use add 1% (final concentration) acetic acid (1 mL/100 mL) and 0.5% (final concentration) ascorbic acid powder (0.5 g/100 mL). Typically, 100 mL prepared daily for single use.

- 3) HPLC mobile phase: 49.5% deionized water, 40% methanol, 10% acetonitrile, 0.5% acetic acid

Using graduated glass cylinders (500 mL and 100 mL), add 495 mL of deionized water, 400 mL of methanol, and 100 mL of acetonitrile to a 1-L reagent bottle and degas under vacuum for 3–5 min. At the time of use add 5 mL of acetic acid (0.5% final concentration). Typically, 300 mL prepared daily for single use.

- 4) HPLC needle wash: 50% methanol and 50% water

Using a 500-mL graduated glass cylinder, add 250 mL methanol and 250 mL deionized water to a 1-L reagent bottle. Solution is used by the HPLC after each injection. Prepare this solution as needed.

- 5) HPLC rinse and shutdown solvent: 90% methanol and 10% water

Using a 500-mL graduated glass cylinder, add 450 mL methanol and 50 mL deionized water to a 1-L reagent bottle. Solution is used by the HPLC at the end of an analytical run. Prepare this solution as needed.

Note: HPLC mobile phase reagent bottle is rinsed once a week thoroughly with deionized water; washed once a month with warm soapy water and rinsed thoroughly with deionized water. HPLC reagent bottles for needle wash, rinse, and shutdown solutions are washed thoroughly with warm soapy water once a month and bottles are rinsed thoroughly with deionized water.

6) L-Ascorbic acid solution (1% or 1 g/dL, pH 2.7)

To a 50-mL falcon tube, add 0.5 g ascorbic acid and 40 mL deionized water and mix well to dissolve. Add deionized water to the 50-mL mark. The solution should be made fresh before use. If the ascorbic acid is used to dilute folate stock solutions, degas under a stream of nitrogen for a few minutes and filter using 0.45 μm (mixed cellulose esters) 10-mL sterile syringe filter just prior to use.

7) Phosphate buffer (0.1 M)

Dissolve 1.4 g potassium phosphate monobasic and 2.2 g dibasic in 95 mL of deionized water (measured by graduated glass cylinder), measure the pH and adjust to 7.2 with KOH. Make up to final volume (100 mL) with deionized water. Filter and degas under a stream of nitrogen for a few minutes. To obtain a 20 mM concentration, dilute 20 mL buffer to 100 mL with degassed deionized water (using graduated glass cylinder and pipette). Prepare fresh at the time of use.

8) PPG dilution solvent for mass spectrometer calibration

Dissolve 15.4 mg of ammonium acetate completely in 49.9 mL of water first. Mix 49.9 mL of methanol with 0.1 mL of formic acid and 0.1 mL of acetonitrile. Then mix the above two solutions together to make the final PPG dilution solvent. Use a 1/50 dilution (400 μL of PPG standard 2000 + 19.6 mL of PPG dilution solvent) for positive ion calibration of the SCIEX tandem mass spectrometer. This solution is stable when refrigerated (+2°C to +8°C) for 6 months. For negative ion calibration, PPG standard 3000 and diluent is supplied in the kit by the manufacturer (SCIEX).

B. Standards Preparation

1) Primary and intermediate stock solutions

The concentrations of primary stock solutions (for each individual folate) are calculated using molar absorptivity. Information on absorption maxima, absorption coefficients, and formulas to calculate the concentration are provided in **Appendix C. Table II** summarizes diluent information for primary and intermediate stock solutions. Information on the preparation of individual stock solutions is maintained on the shared network drive.

Table II: Diluents used for primary and intermediate stock solutions

Compound	Diluent for primary stock solution I	Diluent for intermediate stock solution II (100 µg/mL)	Diluent for intermediate stock solution III (20 µmol/L)
5-MethylTHF	20 mM phosphate buffer (pH 7.2) [D&F]	1% ascorbic acid [D&F]	0.5% ascorbic acid [D&F]
5-FormylTHF	20 mM phosphate buffer (pH 7.2) [D&F]	1% ascorbic acid [D&F]	0.5% ascorbic acid [D&F]
THF	20 mM phosphate buffer (pH 7.2) [D&F]	1% ascorbic acid [D&F]	1% ascorbic acid [D&F]
Compound	Diluent for primary stock solution I	Diluent for intermediate stock solution II (20 µmol/L)	
5,10-MethenylTHF	1 M HCl [D]	0.5 M HCl with 1.0% ascorbic acid [D&F]	--
PGA	20 mM phosphate buffer (pH 7.2) D&F]	Deionized water [D&F]	
MeFox	0.1 N NaOH [D&F]	Deionized water [D&F]	

[D] degassed with nitrogen; [F] filtered using 0.45 µm sterile filters (mixed cellulose esters membrane)

a) Preparation for 5-MethylTHF, 5-FormylTHF and THF stocks:

5-MethylTHF, 5-FormylTHF and THF reduced folate forms are treated the same way. The stock solutions of both unlabeled and ¹³C₅-labeled compounds (used as internal standards) are prepared in the same way as described below.

1. Prepare a **primary stock solution I (≤200 µg/mL)** in a volumetric flask by dissolving an accurately known mass (±0.2 mg) of the pure solid compound in degassed 20 mM phosphate buffer (pH 7.2, containing 0.1% cysteine), targeting a final concentration of ≤200 µg/mL (e.g., ≤5 mg in 25 mL). Vortex briefly to help dissolve the contents and make up to final volume.
2. Remove a small aliquot (~1 mL) of the primary stock solution I to a microcentrifuge vial to determine its concentration by UV spectrophotometry. From this aliquot, prepare two dilutions (e.g., 1/10 and 1/20), each in duplicate, measure the UV absorbance at the peak maximum (5-methylTHF 290 nm and 245 nm; 5-formylTHF 285 nm; THF 298 nm) using scan analysis against phosphate buffer as a blank, and calculate the concentration of primary stock solution I (**Appendix C**). For 5-methylTHF, the ratio of absorbance at 290/245 nm can be monitored (simple reads analysis at each wavelength) to ensure that no oxidation takes place. This ratio should exceed 3.3.
3. To the remaining primary stock solution I, add ascorbic acid powder to a final concentration of 1%. Vortex to help dissolve the ascorbic acid. Aliquot the primary stock solution I with 1% ascorbic acid into labeled cryovials (typically 1 mL/vial) and store deep frozen (-50°C to -90°C). The primary stock solution I is stable for at least 2 years.

4. Prepare an **intermediate stock solution II (100 µg/mL)** by diluting the primary stock solution I in a 25-mL volumetric flask using the diluent specified in the **Table II**. Aliquots of the intermediate stock solution II are stored deep frozen (-50°C to -90°C) in labeled cryovials (typically 1.2 mL/vial) and used approximately every two months to generate a fresh intermediate stock solution III (20 µmol/L). The intermediate stock solution II is stable at least 2 years.
5. Prepare an **intermediate stock solution III (20 µmol/L)** by diluting a portion of the intermediate stock solution II (100 µg/mL) in the diluent specified in **Table II** (typically 10 mL volume). Aliquots of the intermediate stock solution III (typically 0.2 mL/vial) are stored deep frozen (-50°C to -90°C) in labeled microcentrifuge vials and used to generate daily working solutions.

b) **Preparation for 5,10-MethenylTHF stocks:**

Although 5,10-methenylTHF is also a reduced folate, it is treated differently because it is only stable at acidic pH. At neutral pH it is in equilibrium with 5-formylTHF and at alkaline pH it converts to 10-formylTHF. The ¹³C₅-labeled compound is used as internal standard and is prepared in the same way as described below.

1. Prepare a **primary stock solution I (~100 µg/mL)** in a volumetric flask by dissolving an accurately known mass (±0.1 mg) of the pure solid compound in 1 M HCl, targeting a final concentration of ~100 µg/mL (e.g., 5 mg in 50 mL). Vortex briefly and keep the flask for ~10 min in a beaker with warm water (~70°C) to help dissolve the contents. Vortex a few times in between. After contents are completely dissolved, keep the volumetric flask at ambient temperature (15°C to 30°C) for ~30 min to let it cool down. Make up to final volume.
2. Remove a small aliquot (~1 mL) of the primary stock solution I to a microcentrifuge vial to determine its concentration by UV spectrophotometry. From this aliquot, prepare two dilutions (e.g., 1/10 and 1/20), each in duplicate, measure the UV absorbance at the peak maximum using scan analysis against 1 M HCl as a blank, and calculate the concentration of primary stock solution I (**Appendix C**).
3. To the remaining primary stock solution I, add ascorbic acid powder to a final concentration of 1%. Vortex to help dissolve the ascorbic acid. Aliquots of the primary stock solution I with 1% ascorbic acid are stored deep frozen (-50°C to -90°C) in labeled cryovials (typically 1 mL/vial) and used approximately every two months to generate a fresh intermediate stock solution II (20 µmol/L). The primary stock solution I is stable for at least 2 years.
4. Prepare an **intermediate stock solution II (20 µmol/L)** by diluting the primary stock solution I in a 10-mL volumetric flask using the diluent specified in the **Table II**. Aliquots of the intermediate stock solution II are stored deep frozen (-50°C to -90°C) in labeled microcentrifuge vials (typically 0.2 mL/vial) and used to generate daily working solutions.

c) **Preparation for PGA stocks:**

The solubility of PGA decreases as the pH decreases from alkaline to acidic [8]. To maintain optimum solubility, keep higher concentration stock solutions ($\mu\text{mol/L}$) at neutral (or alkaline) pH and ensure buffering when working at slightly acidic pH at much lower PGA concentrations (nmol/L). The $^{13}\text{C}_5$ -labeled compound is used as internal standard and is prepared in the same way as described below.

1. Prepare a **primary stock solution I (~50 $\mu\text{g/mL}$)** in a volumetric flask by dissolving an accurately known mass (± 0.1 mg) of the pure solid compound in degassed 20 mM phosphate buffer (pH 7.2), targeting a final concentration of ~ 50 $\mu\text{g/mL}$ (e.g., 5 mg in 100 mL). Vortex briefly to help dissolve the contents. Keep the volumetric flask in a luke-warm water bath (~ 30 °C) for ~ 20 min to ensure complete solubility. Keep the volumetric flask at ambient temperature (15°C to 30°C) for ~ 20 min to let it cool down. Make up to final volume.
2. Remove a small aliquot (~ 1 mL) of the primary stock solution I to a microcentrifuge vial to determine its concentration by UV spectrophotometry. From this aliquot, prepare two dilutions (e.g., 1/10 and 1/20), each in duplicate, measure the UV absorbance at the peak maximum using scan analysis against phosphate buffer as a blank, and calculate the primary stock solution I concentration (**Appendix C**).
3. Aliquot the remainder of the primary stock solution I into labeled cryovials (typically 1 mL/vial) that are stored deep frozen (-50 °C to -90 °C). The primary stock solution I is stable for at least 2 years and used approximately every two months to generate a fresh intermediate stock solution II (20 $\mu\text{mol/L}$).
4. Prepare an **intermediate stock solution II (20 $\mu\text{mol/L}$)** by diluting the primary stock solution I in a 10-mL volumetric flask using the diluent specified in the **Table II**. Aliquots of the intermediate stock solution II are stored deep frozen (-50 °C to -90 °C) in labeled microcentrifuge vials (typically 0.2 mL/vial) and used to generate daily working solutions.

d) **Preparation for MeFox stocks:**

MeFox is an oxidation product of 5-methylTHF. The $^{13}\text{C}_5$ -labeled compound is used as internal standard and is prepared in the same way as described below.

1. Prepare a **primary stock solution I (~100 $\mu\text{g/mL}$)** in a volumetric flask by dissolving an accurately known mass (± 0.1 mg) of the pure solid compound in degassed 0.1 N NaOH (pH ~ 12.4), targeting a final concentration of ~ 100 $\mu\text{g/mL}$ (e.g., 5 mg in 50 mL). Vortex briefly to help dissolve the contents and make up to final volume.
2. Remove a small aliquot (~ 1 mL) of the primary stock solution I to a microcentrifuge vial to determine its concentration by UV spectrophotometry. From this aliquot, prepare two dilutions (e.g., 1/10 and 1/20), each in duplicate, measure the UV absorbance at the peak maximum using scan analysis against 0.1 N NaOH as a blank, and calculate the primary stock solution concentration (**Appendix C**).

3. Aliquot the remainder of the primary stock solution I into labeled cryovials (typically 1 mL/vial) that are stored deep frozen (-50°C to -90°C). The primary stock solution I is stable for at least 2 years and used approximately every two months to generate a fresh intermediate stock solution II (20 µmol/L).
4. Prepare an **intermediate stock solution II (20 µmol/L)** by diluting the primary stock solution I in a 10-mL volumetric flask using the diluent specified in the **Table II**. Aliquots of the intermediate stock solution II are stored deep frozen (-50°C to -90°C) in labeled microcentrifuge vials (typically 0.2 mL/vial) and used to generate daily working solutions.

Note: Fresh individual primary stock solutions I for all folate forms and their respective ¹³C₅-labeled compounds are prepared approximately every 2 years. Individual intermediate stock solution II (100 µg/mL) for 5-methylTHF, 5-formylTHF and THF are prepared at the same time with primary stock solution preparation. Individual intermediate stock solutions II or III (20 µmol/L) are prepared approximately every 2 months. Analyte and internal standard stock solutions are always prepared together. In-house long-term storage stability data showed that 5-methylTHF and ¹³C₅-5-methylTHF stock solutions (100 µg/mL in 1% ascorbic acid and 10 µg/mL in 0.1% ascorbic acid) were stable for at least 9 years when stored deep frozen (-50°C to -90°C) [5]. Buffers & diluents are degassed with nitrogen & filtered using 0.45 µm sterile filters (mixed cellulose esters membrane) before use.

2) Mixed calibrator solutions

At the beginning of each run, prepare a mixed calibrator (Mix A) and a mixed internal standard solution (Mix B) using respective intermediate stock solutions II and III (20 µmol/L).

a) Calibrator mix (Mix A)

Contains a mixture of each standard prepared in 1% ammonium formate buffer, pH 3.2, with 0.1% ascorbic acid, as shown in **Table III**. The total volume of Mix A is 1.0 mL.

Table III: Information for calibrator Mix A

Mix A	5-MethylTHF	PGA	5-FormylTHF	MeFox	THF	5,10-MethenylTHF
Intermediate stock solution II 20 µmol/L (µL)	--	50	--	50	--	50
Intermediate stock solution III 20 µmol/L (µL)	100	--	50	--	50	--
1% Ammonium formate buffer with 0.1% ascorbic acid (µL)	650					
Concentration in Mix B (nmol/L)	2.0	1.0	1.0	1.0	1.0	1.0

b) Internal standard mix (Mix B)

The internal standard Mix B (total volume of 16 mL) contains a mixture of each internal standard prepared in 1% ammonium formate buffer, pH 3.2, with 0.1% ascorbic acid, as shown in **Table IV**.

Table IV: Information for internal standard Mix B

Mix A	¹³ C ₅ - 5-MethylTHF	¹³ C ₅ - PGA	¹³ C ₅ - 5-FormylTHF	¹³ C ₅ - MeFox	¹³ C ₅ - THF	¹³ C ₅ -5,10- MethenylTHF
Intermediate stock solution II 20 μmol/L (μL)	--	20	--	20	--	20
Intermediate stock solution III 20 μmol/L (μL)	80	--	20	--	20	--
1% Ammonium formate buffer with 0.1% ascorbic acid (mL)	15.82					
Concentration in Mix B (nmol/L)	100	25	25	25	25	25

3) Calibration standards

Prepare mixed calibrators S1 to S5 for the calibration curve in 1.5-mL microcentrifuge vials: add 50 μL mix A to 950 μL ammonium formate buffer (Solvent #1) to prepare the highest calibrator S5 (100 nmol/L 5-methylTHF, and 50 nmol/L each of PGA, MeFox, THF, 5-formylTHF, and 5,10-methenylTHF). Prepare calibrators S4 to S1 from calibrator S5 by using the amounts specified in **Table V**.

Table V: Information for mixed calibrators S1 to S5

Calibrator level	Calibrator mix	Solvent #1	Concentration (nmol/L): 5-methylTHF/PGA/5-formylTHF/MeFox/THF/5,10-methenylTHF
S5	50 μL Mix A	950 μL	100/50/50/50/50
S4	200 μL S5	800 μL	20/10/10/10/10
S3	40 μL S5	960 μL	4/2/2/2/2
S2	20 μL S5	980 μL	2/1/1/1/1
S1	10 μL S5	990 μL	1/0.5/0.5/0.5/0.5

C. Preparation of Quality Control Materials

Quality control materials for this assay are prepared in-house from blood products acquired from blood banks or from other volunteer blood donors. Approximate QC target values for serum tFOL are 7–10 (low), 25 (medium), and 50 nmol/L (high). The low QC aims to be close to the deficiency cutoff value of 7 nmol/L. The high QC aims to be in the top third of the population distribution (75th percentile is ~50 nmol/L). Because the United States has mandatory food fortification with folic acid, it is difficult to find donors with low serum folate levels and the prevalence of clinical folate deficiency is <1% in the US population.

In addition to tFOL concentrations of individual folate forms are considered. If specimens don't contain the approximate target values for the individual folate forms as shown in **Appendix D**, manipulation through spiking with standard compounds or dilution with physiologic sodium chloride solution or albumin solution (4%) is considered. Not all folate forms have to be low in the "low" pool, medium in the "medium" pool, and high in the "high" pool, just as long as there is sufficient distinction between the pools.

The serum is pooled, and pools are filtered through gauze before being dispensed to remove fibrinogen. Serum (usually 800 µL) is aliquoted into 2.0-mL Nalgene cryovials, capped, and frozen. The QC pools are aliquoted for single use and stored deep frozen (-50°C to -90°C) and are stable for at least 3 years. Ascorbic acid (0.5%) is added to the serum pools to ensure even better long-term folate stability.

Characterization limits are established by analyzing duplicates of each pool for at least 20 consecutive runs.

Information on the preparation of quality control pools is maintained on the shared network drive. For more detailed information on the preparation of QC materials, homogeneity testing, and characterization refer to **SOP NBB-OC-LABOP.01.01 QC Materials**.

D. Other Materials

With some exceptions, a material listed herein may be substituted with an equivalent product from a different manufacturer provided if it meets or exceeds the specifications of the product listed. In the case of standards, internal standards, chemicals, and reagents, the chemical and/or isotopic purity of the substitute must meet or exceed that of the listed product. In the case of the HPLC column and guard cartridge, equivalent performance must be demonstrated experimentally in accordance with DLS Policies and Procedures.

- (1) Consumables for automated sample dilution on Hamilton liquid handler
 - a) 96-well 2 mL sample/collection plate (Whatman)
 - b) 96-well collection plate seals (Agilent Technologies)
 - c) 1 mL and 0.3 mL plastic pipette tips (Hamilton)
 - d) 300 mL plastic reusable reagent and water troughs (Hamilton)
- (2) Consumables for automated solid phase extraction (SPE) on 96-probe SPE instrument (Caliper-Zephyr)
 - a) 96-well Bond Elute SPE plate (50 mg phenyl sorbent, Agilent Technologies)
 - b) Captiva 96-well filter plate (0.45 µm PVDF embedded well, Agilent Technologies) for efficient automated filtration with vacuum manifold (IST VacMaster-VCU)
- (3) General consumables
 - a) C-8(2) analytical HPLC column, 150 x 3.2 mm, 5 µm (Phenomenex)
 - b) 0.5 µm stainless frits A-102X (Upchurch Scientific)
 - c) PEEK tubing 0.005 and 0.007 ID (Upchurch Scientific)
 - d) HPLC Solvent glass inlet filters, purge frits, gold seal and outlet caps (Agilent Technologies)
 - e) Eppendorf pipette tips (0.5-10 µL, 10-100 µL, 100-1000 µL, and 10 mL)
 - f) Eppendorf Combitip advanced pipette tips (0.5 mL and 1 mL) for repeater pipette

- g) Gilson Pipetman positive displacement pipette tips (50 μ L, 100 μ L, and 1000 μ L)
 - h) Nunc 30 mm, 1-mL 96-well HPLC plate for 96-well autosampler (Fisher Scientific)
 - i) Fisherbrand silicone plate seal for 1 mL 96-well HPLC plate (Fisher Scientific)
 - j) HPLC solvent filter degasser, model FG-256 (Lazar Research Laboratories)
 - k) 0.22 μ m cellulose nitrate nonpyrogenic sterile water filter system 500 mL capacity (Corning)
 - l) 0.45 μ m mixed cellulose esters (MCE) membrane syringe filter (Millipore)
 - m) 2.0 mL polypropylene cryovials (Nalgene)
 - n) 1.0 mL disposable syringes (Hamilton)
 - o) Various glass beakers, volumetric flasks (class A), graduated glass cylinders, and bottles
- (4) Folate Standards
- a) PGA (Pteroylglutamic acid, free acid or Na₂-salt), 5-methylTHF ([6S]-5CH₃-H₄PteGlu, Ca-or Na₂-salt), 5-formylTHF ([6S]-5CHO-H₄PteGlu, Ca-or Na₂-salt), MeFox ([6S](pyrazino-s-triazine derivative), THF ([6S]-H₄PteGlu, free acid, Ca-or-Na₂-salt) and 5,10-methenylTHF ([6S]-5,10-CH=H₄PteGlu-Cl x HCl salt) (Merck & Cie, <https://www.emdgroup.com>)
 - b) ¹³C₅-PGA, ¹³C₅-5-methylTHF, ¹³C₅-5-formylTHF, ¹³C₅-MeFox, ¹³C₅-THF and ¹³C₅-5,10-methenylTHF (Merck & Cie.)
- (5) Chemicals and Solvents
- a) Ammonium hydroxide (28–30% as NH₃, Mallinckrodt Chemicals)
 - b) L-Cysteine (Sigma Aldrich)
 - c) Potassium phosphate dibasic and monobasic salts (Fisher Scientific)
 - d) Formic acid (>95%) reagent grade (Sigma Aldrich)
 - e) Acetic acid (99%) reagent grade (Fisher Scientific)
 - f) L-ascorbic acid (vitamin C min 99% purity, Spectrum Chemicals, Fisher Scientific)
 - g) Hydrochloric acid (36.5-38%, JT Baker)
 - h) Water, 18 M Ω /cm, HPLC grade (Aqua Solutions)
 - i) Methanol, acetonitrile HPLC certified solvent (Honeywell/Burdick & Jackson Laboratories)
 - j) Nitrogen ultrapure (>99.99 % purity, Air Products)

E. Instrumentation

In the case of simple laboratory instrumentation (e.g., pipettes, vortex mixer, analytical balance, etc.) a product listed herein may be substituted with an equivalent product from a different manufacturer provided if it meets or exceeds the specifications of the product listed. In the case of analysis instrumentation (e.g., HPLC components, tandem quadrupole mass spectrometer) equivalent performance must be demonstrated experimentally in accordance with DLS Policies and Procedures if a product substitution is made. For details, refer to **Appendix B_C: JA-4013-I-01-Instrument Comparison and System Verification**.

To provide adequate throughput for this method as well as backup instrumentation during times of repair and maintenance, we utilize multiple HPLC-MS/MS systems of the SCIEX type. Equivalent performance must be demonstrated in accordance with DLS policies and procedures when multiple analysis systems are used in parallel, even if they are of the exact same type.

- (1) Agilent 1260 Infinity HPLC system (Agilent Technologies) equipped with control pilot, degasser, binary pump, autosampler with thermostat, and column compartment
- (2) SCIEX 6500 QTrap triple quadrupole mass spectrometer with Turboionspray (TIS) as ion source in ESI mode equipped with Analyst software (SCIEX)
- (3) Peak Scientific Infinity 1031 nitrogen generator (Peak Scientific Instruments Ltd)
- (4) Caliper-Zephyr automated 96-probe liquid handler equipped with SPE (Perkin Elmer Inc.); can also be used for sample transfer for filtration on Captiva 96 well-plate filters
- (5) Microlab Starlet automated liquid handler (Hamilton)
- (6) Syringe pump (Harvard Apparatus)
- (7) Eppendorf repeater pipette (volume range from 1 μ L to 10 mL, Eppendorf)
- (8) Eppendorf single channel pipettes (10 μ L, 100 μ L, 200 μ L, 1000 μ L, and 10 mL; Eppendorf)
- (9) Gilson MICROMAN positive air displacement pipettes (20 μ L, 100 μ L, and 1000 μ L; Gilson)
- (10) Digiflex CX (ICN Biomedicals, Inc. Diagnostics Division)
- (11) Vortex Genie 2 mixer (VWR)
- (12) Magnetic stirrer (Baxter Scientific Products)
- (13) pH meter (Corning Pinnacle 530 or Accumet XL150 pH/MV; Fisher Scientific)
- (14) Analytical Balance Model AG104 (Mettler Instrument Corp.)
- (15) Cary 300 Bio UV-visible spectrophotometer (Agilent Technologies)

7. Calibration and Calibration Verification Procedures

A. Method Calibration

In-house studies showed that aqueous calibration provides equivalent results to calibration in serum and slopes for the two calibration curves (serum vs. water) were less than $\pm 5\%$ different for all folate forms [2]. Aqueous calibrators (at the beginning of each run) are carried through the entire sample processing procedure. Calculation of folate concentrations in QCs and unknown patient samples is based on results obtained from a daily 5-point calibration curve (S1 to S5). A blank sample (S0, containing internal standard mix) and a double blank (containing reagents only) are also included in each run. Area ratios of analyte to internal standard from single analysis of each calibrator are calculated and a linear regression equation ($1/x^2$ weighting) is generated. At the end of each run, the calibration curve is re-injected to assess potential calibrator drift. The measured concentrations of the calibrators at the beginning of each run should generally agree within 15% of their set values, although $>15\%$ agreement will be observed at concentrations approaching the LOD.

This method is linear from 0-220 nmol/L for each folate form [2, 3]. The calibration range is from 1–100 nmol/L for 5-methylTHF and from 0.5–50 nmol/L for all other folate forms. The removal of a calibration point(s) is not suggested for this assay. This assay uses a freshly prepared mixed calibration curve of all six analytes. In our experience excluding a calibrator point(s) for a particular analyte(s) didn't affect the accuracy of major analyte results, but the impact was noticeable on minor analyte results. We recommend repeating the sample analysis for the analyte(s) where a calibrator(s) point error is noticed. Samples with concentrations that exceed the calibration range are diluted with 0.1% ascorbic acid and re-analyzed.

Since calibrators are included in every run, there is no additional calibration verification required. However, for good laboratory practice, calibration verification is conducted at least twice a year using international reference materials to monitor the assay performance. For details, refer to **Appendix B_A: JA-4013-G-01-Calibration and Calibration Verification**.

In 2005, the National Institute of Standards and Technology (NIST) released a new three-level standard reference material (SRM) for homocysteine and folate in human serum, SRM 1955. This material was characterized by various mass-spectrometry-based methods used at NIST and the CDC. Because of the good agreement between the NIST and CDC methods for 5-methylTHF and PGA, NIST used the CDC results as part of the value assignment process. Good agreement was also found for tFOL between the CDC LC-MS/MS and microbiologic assay (level 1: 6.0 vs. 5.6; level 2: 13 vs. 14; level 3: 41 vs. 44).

In 2011, NIST released a one-level standard reference material for metabolites in human plasma, SRM 1950. This material was characterized by various mass-spectrometry-based methods used at NIST and the CDC. The CDC LC-MS/MS procedure matched the certified value [uncertainty] for 5-methylTHF (26.7 vs. 26.9 [0.70] nmol/L) and the reference value for PGA (4.03 vs. 3.42 [1.02] nmol/L). The SRM is stored deep frozen at -70°C (-50°C to -90°C). Prior to use SRM is allowed to thaw at ambient temperature ($+15^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) for at least 30 min protected from light.

In 2018, NIST released a new three-level standard reference material (SRM 3949) for folate vitamers in human serum. This material was characterized by various mass-spectrometry-based methods used at NIST and the CDC. Because of good agreement between the NIST and CDC methods for 5-methylTHF, PGA, 5-formylTHF, MeFox and THF; NIST used for most of the analytes the CDC results as part of the value assignment process. A good agreement was also found for tFOL between the NIST and CDC method results.

NIST SRMs are shipped on dry ice. Once received, the vials are stored in the dark under deep frozen conditions (-50°C to -90°C). Prior to use the material is thawed at ambient temperature (+15°C to +30°C), contents mixed gently, and used immediately. For economical use we aliquot the material into multiple cryovials and store them deep frozen until analysis (1 freeze/thaw cycle).

In 2006, the National Institute for Biological Standards and Control (NIBSC) issued the first WHO certified reference material (lyophilized, one level) for folate in human serum (03/178, established in 2004). The folate concentration in this material has been certified by HPLC-MS/MS through measurements provided by NIST and CDC. NIBSC material is shipped on dry ice. It's recommended to store unopened ampoules at -20°C. Accelerated degradation studies indicated that lyophilized material is adequately stable at -20°C for B12 and folate content. Once reconstituted, we store multiple aliquots under deep frozen conditions (-50°C to -90°C). Per WHO guidelines expiry date is not assigned for this material.

Details about our proficiency testing (PT) activities can be found in the proficiency testing form. We participate in external proficiency testing programs from the College of American Pathologists (CAP) K Ligand survey (3 times per year) and LN5 Calibration Verification Ligand survey (2 times per year). When possible, we also participate in other external proficiency testing program such as the UK NEQAS Haematinics survey (once or twice a year). General information on the handling, analysis, review, and reporting of proficiency testing materials is maintained on the shared network drive.

Information on method characteristics has been published [2–5]. Results from a series of in-house ruggedness testing experiments designed to assess how much method accuracy changes when certain experimental parameters are varied are presented in **Appendix E**.

We have performed in-house comparisons of the HPLC-MS/MS assay and the microbiologic assay; the latter still considered an accurate “reference point” for total folate (tFOL). For serum samples, there is excellent correlation and agreement between the two assays [2, 9, 10]. The microbiologic assay produces results that are within ±10% of the HPLC-MS/MS results.

B. Instrument Calibration

1) Tandem mass spectrometer

The calibration of the mass spectrometer is scheduled on a semi-annual basis as part of a preventive maintenance program and is performed by service engineer from SCIEX.

The tuning and mass calibration of the first (Q1) and third (Q3) quadrupoles of the mass spectrometer is performed using a solution of polypropylene glycol (PPG) by infusion and running the instrument in either Manual Tuning mode or using Automatic Mass Calibration. Resolution for Q1 is always set at “Unit” while the resolution for Q3 is sample matrix dependent and can be set either “Unit or High”. The method for the analysis of serum folate forms operates in positive ion mode with “Unit” resolution for both Q1 & Q3 quadrupoles.

The PPG calibration solution kit contains a PPG standard 2000 vial for positive ion mode calibration and a PPG standard 3000 vial for negative ion mode calibration. Preparation of the PPG standard 2000: SCIEX 6500 uses a 1:50 dilution (for calibration of positive ion mode add 0.4 mL of PPG standard to 19.6 mL of PPG dilution solvent). PPG 3000 is used directly without any dilution for calibration of negative ion mode.

After preventive maintenance (PM) or repair service is completed, the engineer performs tuning and mass calibration of the instrument and saves PPG calibration profiles and data in respective PM folder on the instrument computer. The analyst (and/or team lead) reviews the PPG data and peak responses (intensities and peak widths); verifies the instrument mass calibration table with the calibration data saved (pdf file) in PM folder; compares the data and PPG responses with previously conducted services (PM or repairs) and updates. Newly conducted calibration data, PPG peak widths and responses should either exceed or compare well with the previously conducted service data outcomes otherwise new instrument service (deep cleaning, repairs etc.) and updated mass calibration is requested. After service is completed, the analyst conducts a test run (including folate forms calibration curve, and 2 sets of QCs) to validate method performance. If the method performs is within acceptable responses for all analytes, the service call is closed, otherwise service is continued, and troubleshooting that includes additional instrument cleaning, repairs, and recalibration is conducted.

Rough pump oil (1L) in the mass spectrometer is replaced during PM as part of service and waste oil is discarded appropriately.

2) Hamilton Microlab Starlet liquid sample handler

Twice a year a Hamilton service engineer performs preventative maintenance including volume verification at 10 μ L and 1000 μ L. Additionally, pipetting accuracy may be checked using a gravimetric volume verification kit.

Daily and weekly maintenance of the system is executed through the instrument operated software that checks for deck cleaning, tip waste, and 8 channel tightness. For details, refer to **Appendix B_C: JA-4013-I-02-Maintenance and Verification of Hamilton Liquid-Handler.**

3) Caliper-Zephyr SPE instrument maintenance and function checks

A 96-probe automated SPE instrument is used for sample extraction and cleanup. Once per year a Perkin-Elmer service engineer performs preventative maintenance including arm adjustment and volume verification (10 μ L to 1.0 mL), and replacement of parts if needed. As a routine maintenance and function checks the analyst performs daily pressure check for vacuum pump (maintained at 15 \pm 5 psi during SPE); rinsing of solvent troughs with deionized water after each run (troughs replaced as needed); checking appropriate deck position for SPE supplies (tip boxes, solvent troughs, SPE plate and collection plate); checking of nitrogen gas cylinder regulator for pressure (maintained at 80 \pm 5 psi) and instrument regulator pressure (maintained at 60 \pm 5 psi). Analyst performs weekly instrument cleanup (dry wipe) and silicone lubrication of the 96-probe head. This instrument can also be used for non-volume-critical solvent transfers.

4) Pipettes (air displacement and positive displacement)

Pipette calibration verification is performed biannually. Typically, one calibration verification is done by a certified company and one by the analyst gravimetrically using a calibrated analytical balance.

5) Varian UV-Vis spectrophotometer

Calibration verification is performed three times per year by participation in the CAP Instrumentation survey. Additionally, every time the instrument is turned on there are internal diagnostics that are run. Calibration verification using certified filters is performed at least annually, generally twice per year. Calibration verification of the certified filters is performed externally every other year.

6) Balances

On site calibration is performed annually by a certified company. Calibration verification is performed by the analyst as needed using certified weights.

7) pH meter

Calibration verification is performed by the analyst prior to use as needed using commercial calibration solutions.

8. Procedure Operating Instructions; Calculations; Interpretation of Results

A typical run consists of the following sequence of samples: reagent blank (double blank), blank (contains internal standard mix), 5 calibrators, first set of QCs, 83 patient samples, and second set of QCs, for a total of 96 samples (96-well plate format). Three levels of serum QCs are analyzed in duplicate in each run as bench QC materials. In preparation for SPE, samples are typically pipetted by a liquid sample handler, but they can also be pipetted manually. SPE is carried out using an automated 96-probe instrument.

A. Preliminaries

- (1) Thaw frozen serum specimens (QCs and unknown patient samples), folate intermediate stock solutions (calibrator and internal standard); it takes about 40 min for the samples to reach ambient temperature.
- (2) Prepare buffers and mobile phase (can be prepared ahead of time).
- (3) Add 0.5% acetic acid to the pre-made mobile phase prior to use.
- (4) Prepare fresh sample solvent #1 (with 0.5% ascorbic acid), sample solvent # 2 (with 0.1% ascorbic acid) and sample solvent #3 (with 0.5% ascorbic acid and 1% acetic acid).
- (5) Prepare fresh 1% ammonium formate buffer (with 0.1 % ascorbic acid) for calibrator Mix A and internal standard Mix B.
- (6) Mark the 96 well-plate rows for the number of samples to be analyzed.
- (7) Prepare calibrator Mix A and internal standard Mix B as described in section 6.B.(2).

- (8) Prepare calibration standards S1 to S5 as described in section 6.B.(3). Vortex thoroughly for adequate mixing.
- (9) Vortex all thawed specimens thoroughly prior to pipetting and visually check for any unusual sample volume, specimen color, or debris/precipitate.

B. Automated Sample Pipetting using a Liquid Handler to Prepare for SPE (routine runs)

Daily maintenance and function checks are performed prior to sample pipetting. For details, refer to **Appendix B_C: JA-4013-I-02-Maintenance and Verification of Hamilton Liquid-Handler**.

- (1) Hamilton Microlab Starlet is used for automated pipetting of calibrators, IS Mix, water, buffer, and specimens from respective cryovials or troughs into a 2-mL 96-well plate that is then subjected to automated SPE sample extraction and clean-up.
- (2) Check and restock tip racks.
- (3) Fill two reagent troughs with solvent #1 and deionized water respectively and load the troughs onto the reagent carrier.
- (4) Put calibrators, internal standard mixture, QC samples, and unknown patient samples onto respective sample carriers.
- (5) Put a 2-mL 96-well collection plate (sample destination plate) in place.
- (6) Using Microlab Star Run Method Software the pipetting program is executed and dispenses according to the scheme in **Table VI** for a final volume of 550 μ L.
- (7) After the pipetting is completed, the sample plate is removed, covered with the 96-well plastic seal, and incubated at ambient temperature (15°C to 30°C) for 20 min prior to SPE for the internal standards to equilibrate with the endogenous folates.

Table VI: Pipetting Scheme

Well #	Sample Type	Internal Standard Mix B	Calibrator Mix	Sample Solvent #1	Water	QC or Patient Specimen	Final Volume
1	Reagent blank (Double blank)	--	--	400 µL	150 µL	--	550 µL
2	Calibrator S0 (Blank)	60 µL	--	340 µL	150 µL	--	550 µL
3	Calibrator S1	60 µL	150 µL	190 µL	150 µL	--	550 µL
4	Calibrator S2	60 µL	150 µL	190 µL	150 µL	--	550 µL
5	Calibrator S3	60 µL	150 µL	190 µL	150 µL	--	550 µL
6	Calibrator S4	60 µL	150 µL	190 µL	150 µL	--	550 µL
7	Calibrator S5	60 µL	150 µL	190 µL	150 µL	--	550 µL
8	Low QC – Set 1	60 µL	--	340 µL		150 µL	550 µL
9	Medium QC – Set 1	60 µL	--	340 µL		150 µL	550 µL
10	High QC – Set 1	60 µL	--	340 µL		150 µL	550 µL
11-92	Patient Samples	60 µL	--	340 µL		150 µL	550 µL
93	Low QC – Set 2	60 µL	--	340 µL		150 µL	550 µL
94	Medium QC – Set 2	60 µL	--	340 µL		150 µL	550 µL
95	High QC – Set 2	60 µL	--	340 µL		150 µL	550 µL

Notes: Use positive displacement pipettes wherever possible for the preparation of the calibration standards and Internal standards. Use positive displacement pipettes or an automated positive displacement pipettor for transferring the preliminary assay volume of standards, QCs, and patient samples into the sample plates.

Folate forms are light sensitive. Perform all sample preparation under yellow subdued light.

C. Manual Sample Pipetting to Prepare for SPE (occasional run)

- (1) Mark 96-well plate with date, study ID and mark the number of wells that need to be used to prepare samples for SPE.
- (2) To construct a 5-point calibration curve, follow the pipetting scheme shown in **Table VI**. A reagent blank that contains only reagents and a blank that contains reagents and only the internal standard mix (S0) is included in each run.
- (3) QC and patient samples (vial 8 and beyond): add 340 µL solvent #1, 60 µL of internal standard mix (Mix B) and 150 µL specimen for a final volume of 550 µL.
- (4) Seal the 96-well plate with silicone plate seal and mix the contents gently on plate shaker (typically 1 min.).
- (5) Incubate the sample plate at ambient temperature (15°C to 30°C) for 20 min for the internal standards to equilibrate with the endogenous folates prior to SPE.

D. Automated Solid Phase Extraction Method

- (1) A 96-probe SPE instrument (Caliper-Zephyr) is used for automated SPE. All SPE steps, such as SPE plate conditioning, sample loading, SPE plate washing, and sample elution are performed automatically. The instrument processes one 96-well plate extraction in about an hour which includes blanks, calibrators, 2 sets of QCs, and 82 unknown patient specimens.
- (2) Open the “Maestro Software” on the desktop and select the Folate SPE method.
- (3) The SPE sample plate, collection plate, conditioning solvents, wash, and elution buffers are placed on their respective deck positions on the SPE instrument.
- (4) Start the Folate SPE method within the software.
- (5) The SPE method performs the following steps:
 - a) **Condition:** SPE plate is conditioned with acetonitrile and methanol (0.5 mL each) in 2 steps each, followed by conditioning with 1.3 mL solvent #1 in 3 steps (0.5 mL x 2; and 0.3 mL x 1)
 - b) **Load:** 500 µL sample is loaded in 4 steps (125 µL x 4)
 - c) **Wash:** SPE plate is washed in 3 steps with 1.3 mL of solvent #2 (0.5 mL x 2; and 0.3 mL x 1)
 - d) **Elution:** Sample elution is carried out in 2 steps (0.3 mL and 0.2 mL) with solvent #3
- (6) The chemical waste (~2L/week) from this procedure consists of water, methanol, acetonitrile, ammonium formate, ascorbic acid, and serum (~4.0%) and is discarded appropriately.

E. Automated Sample Filtration for HPLC-MS/MS Analysis

- (1) After SPE is completed, samples are filtered using Captiva filter plates (96-well, 0.45 µm PVDF).
- (2) Caliper-Zephyr liquid handler is used to transfer 300 µL of the extracted sample into a fresh Captiva filter plate.
- (3) Captiva filter plate is placed on top of a fresh HPLC collection plate (30 mm Nunc plate).
- (4) The samples in the Captiva filter plate are filtered into the HPLC collection plate using a vacuum manifold (IST) at ≤5 psi pressure within <5 min.
- (5) The HPLC collection plates are sealed with the pre-slit seals and arranged on the HPLC autosampler for analysis.

Note: Once per year a Perkin-Elmer service engineer performs preventive maintenance of Caliper-Zephyr including arm adjustment and volume verification (10 μ L to 1.0 mL), and replacement of parts if needed.

As a routine maintenance the analyst performs a weekly instrument cleanup (dry wipe); silicone lubrication of the 96-probe head; vacuum pump pressure checks; rinsing of solvent troughs (used for solvents #1, # 2 and # 3) thoroughly with deionized water after each run, and air dried (replaced as needed); checks nitrogen gas cylinder. This instrument can also be used for non-volume-critical solvent transfers.

F. HPLC-MS/MS Instrument Preparation

- (1) The Agilent HPLC system coupled to the SCIEX MS/MS system is used to quantitate folate vitamers in extracted serum. Typical MS/MS method parameters for each folate vitamer and the corresponding internal standards are listed in **Appendix F**.
- (2) Prior to sample analysis, HPLC instrument function checks are performed:
 - a) Switch all HPLC components to "ON" mode.
 - b) Purge each solvent line with the respective solvent bottle.
 - c) Prime the HPLC column with a series of solvents as described in the next steps.
- (3) Methanol and water (90:10 v/v) are used for line purging for ~5 min at a flow rate of 5 mL/min. The column is primed with this solvent for ~20 min at a flow rate of 500 μ L/min.
- (4) The lines are purged with mobile phase for 5 min at a flow rate of 5 mL/min and column is primed \geq 15 min at a flow rate of 250 μ L/min. Pressure is recorded. HPLC system is ready for analysis; samples are analyzed (20 μ L injection volume) at a flow rate of 250 μ L/min for 7 min that includes column equilibration time.
- (5) Typically, once per week (or as needed) acetonitrile and water (65:35 v/v) is used for line purging (~5 min) at a flow rate of 5 mL/min. The column is primed with this solvent ~20 min at a flow rate of 500 μ L/min.
- (6) Typically, to increase sample throughput on column and/or improve peak resolution (occasionally for issues with a brand-new column) the column is reversed once it has run in one direction; primed as described below prior to its use for routine sample analysis.
- (7) Brand new column is primed with acetonitrile and water (65:35) at 800 μ L/min for \geq 60 min; with methanol and water at 800 μ L/min for ~30 min; equilibrated with mobile phase at 250 μ L/min for ~60 min prior to sample run. Typically, a new column is also primed with matrix by injecting a set of processed serum samples (15-20 samples) prior to its use for routine analysis.
- (8) Prior to sample analysis, MS/MS instrument maintenance and function checks are performed:
 - a) Put the instrument in standby mode and remove the curtain plate from source. Clean curtain plate first with water, then methanol, and wipe dried with lint-free Kimwipe.
 - b) Wipe orifice plate with methanol damped lint-free Kimwipe.

- c) Check ion spray needle for any blockage and clean if necessary.
 - d) Check instrument vacuum (1.0 to 2.0 e⁻⁵Torr) and gas generator for curtain (~60 psi), source (90-110 psi), and exhaust (~70 psi) status.
 - e) Load the appropriate instrument method and create an acquisition batch containing the sample sequence.
- (9) Prior to sample analysis, the HPLC-MS/MS instrument is equilibrated:
- a) At least 2 sample solvent injections are performed prior to sample analysis to verify that the system is performing properly.
 - b) Either a blank (S0, mixture of internal standards) or the low QC sample from the previous run can be re-injected prior to the analysis of the current run to check the instrument performance.
- (10) The sample plate is loaded on the autosampler, and the sample acquisition batch created in the earlier step is submitted for analysis.
- (11) The HPLC rinse method is loaded at the end of the sample batch. It runs isocratic (90:10 v/v of methanol and water) for 20 min in Q1 scan mode at the end of the batch to clean the HPLC column and MS/MS system. The data is recorded in an acquisition rinse batch file so that it can be reviewed later. If necessary, multiple batches can be submitted for analysis, each separated by the HPLC rinse method.
- (12) The HPLC shutdown method is loaded after the HPLC rinse method. It runs isocratic (90:10 v/v of methanol and water) for 5 min in MRM mode after the rinse method. Finally, the instrument goes to standby mode until the next run and the sample plate in the autosampler is maintained at ≤10°C during standby mode.
- (13) The chemical waste (~2 L/week) from this procedure consists of water, methanol, acetonitrile, acetic acid and is discarded appropriately.

G. Processing and Reporting a Run

- (1) The SCIEX Analyst software is used to review and process a run. For details, refer to **Appendix B_D: JA-4013-DR-02-Processing and Reporting a Run.**
- (2) A STARLIMS database is used for additional levels of data review by the analyst, project lead, QA officer, and supervisor and for data reporting.
- (3) Import the results file into the STARLIMS database for further data review, see **Appendix B_D: JA-4013-DR-01-Computerization and Data System Management.**
- (4) Calculate the results using Excel (exception, for R&D runs or troubleshooting):
 - a) The final integrated results can either be directly imported to the STARLIMS database (typical) or alternatively can be processed by importing into an Excel template sheet for final calculations and interpretation.
 - b) Transfer the peak areas for the analyte and the internal standard for each sample into the appropriate fields in the Excel sheet.

- c) The calibration curve with slope, intercept, and R2 is automatically generated (linear, not forced through zero, no weighting) based on area ratios.
- d) The concentrations of QCs and unknowns are automatically calculated using the slope and intercept information.
- e) The file is saved and maintained for documentation.
- f) For studies where data is not imported into the database, we directly transfer the results from the “**Analyst Result**” table into the Excel sheet for final summary and results interpretation since our results are from a weighted calibration curve ($1/x^2$ weighted).

H. Exporting a Run

The procedure to export a run to a STARLIMS database is described in section 3 and in **Appendix B_D: JA-4013-DR-01-Computerization and Data System Management**.

I. Calculations

Serum results for each folate species are directly reported as nmol/L. To obtain a serum tFOL result, the individual results from each folate form (generally excluding MeFox; see special method notes) have to be added up. If a vitamers result is less than the LOD, a fill value of LOD divided by the square root of the 2 [LOD/SQRT (2)] is used for summation. If one of the folate forms is missing, the serum tFOL result is also missing.

J. System Maintenance

The system maintenance consists mainly of the different prime, purge, and wash cycles described in Section 8 of this document and performed before and after each run. Column connections are checked for leaks daily and are wiped with a water-moistened tissue if any residues have built up. Solvent bottles are filled as needed and cleaned typically once in a month (or as needed). Filters in the solvent bottles are replaced as needed (typically every 6-8 months). The pre-column filters are replaced after ~500 injections. The curtain plate is cleaned on daily basis first with water, then wiped with lint free Kimwipes dabbed in methanol. The orifice plate is also wiped daily with methanol dabbed lint-free Kimwipes. Preventative maintenance is performed by service engineers on all major equipment (MS/MS, HPLC, Caliper-Zephyr, Hamilton, and Cary UV spectrophotometer) at least once a year.

Note: In preparation for an occasional power shutdown or repairs, the HPLC-MS/MS instruments are powered off following the steps:

- (1) Complete or stop any ongoing scans and put all instruments on standby
- (2) Close all software, shut-down the computer, and unplug the computer power cable
- (3) Close the air supply valve and vent the MS/MS (press vent button for ~3 sec)
- (4) Wait for ~15 minutes for the turbo pump to spin down
- (5) Power off all HPLC modules and MS/MS system

K. Special Method Notes

Serum tFOL may or may not include MeFox, depending on the study and investigator request. It is still of scientific debate whether MeFox is only generated *in vitro* or may already be present *in vivo* [11]. Including MeFox into the tFOL may slightly overestimate folate status, while excluding it may slightly underestimate status, however the difference between the two approaches is rather small (~5%) [11]. The 3 minor analytes THF, 5-formylTHF, and 5,10-methenylTHF, particularly the latter two, are usually <LOD in serum specimens, however they still have to be measured to accurately capture serum tFOL and they can reach concentrations >1 nmol/L in some samples [11].

9. Reportable Range of Results (AMR – Analytical Measurement Range)

This method is linear from the LOD (see section 7) to 100 nmol/L for 5-methylTHF, and 50 nmol/L for 5-formylTHF, PGA, MeFox, THF, and 5,10-methenylTHF. Samples with 5-methylTHF results <7 nmo/L (3 ng/mL) are re-analyzed for confirmation before results are released. Samples with folate concentrations exceeding the highest calibrator are re-analyzed after appropriate dilution with 0.1% ascorbic acid. There is no known maximum acceptable dilution. When possible, avoid small volume pipetting and minimize use of serial dilutions when generating diluted samples.

10. Quality Control (QC) Procedures

A. Blind Quality Controls

Blind QC specimens can be inserted into the mix of patient specimens. These QC specimens are generally prepared at two levels that would be encountered in patient samples; the labels used are identical to those used for patient samples. Alternatively, open label blind QC specimens can be used where the analyst knows that the sample is a blind QC, but they do not know what pool the sample is from. Open label blind QCs are only used if one can choose from at least 6 different pools and the analyte concentrations are similar to those found in patient samples. The frequency of blind QC specimens in a run is typically 1 in every 20 specimens analyzed. After a run is completed, used blind QC are removed from the run, marked with a black dot on the cap to indicate that the vial has been thawed, and returned to the blind QC box. This helps to identify which vials have been used. If a run needs to be repeated, the same blind QC can be inserted as in the initial run.

The use of blind QCs is optional but encouraged. Blind QCs are used in this method as a supplementary tool to assist in monitoring accuracy, precision, and aid in detecting errors; these are not used as part of the primary control procedures to determine if a run is out of control.

B. Bench Quality Controls

Bench QC specimens are prepared from a minimum of 2 pools that represent low and high levels of 6 different folate forms. This assay typically uses 3 serum pools, which represent low, medium, and high levels of 5-methylTHF, 5-formylTHF, PGA, MeFox, THF and 5,10-methenylTHF. These QC samples are prepared in the same manner as patient samples and analyzed in duplicate as part of each run (placed at the beginning and end of each run).

The QC results are checked after each run using of a multi-rule quality control system [12] based their characterization data, namely: the pool mean; the pooled within-run standard deviation associated with individual QC results measured in the same run (S_w); the standard deviation associated with individual QC results (S_i); and the standard deviation associated with run mean QC results (S_m). QC rules have been designed to accommodate the use of 1–3 different QC pools during a run, the use of 1–2 measurements of each pool per run, and as many instruments as needed. These QC rules are described in the DLS Policies and Procedures Manual and a relevant selection applicable to this assay is shown below. The system is declared “in control” if all individual QC results are within 2S limits; the run is accepted. If not, then the rules shown below are applied and the run is rejected if any condition is met; the run is declared “out of control”:

(1) One QC result per pool (due to accidental loss of duplicate)

a) If one of the three QC run **individual** results is outside a 2 S_i limit, reject run if:

1. 1 3S Rule—Run result is outside a 3 S_i limit or
2. 2 2S Rule—Two or more of the three run results are outside the same 2 S_i limit
3. 10 X-bar Rule—Current and previous nine run results are on the same side of the characterization mean or
4. R 4S Rule—Two consecutive standardized run results differ by more than 4 S_i

(2) Two QC results per pool (all three pools have duplicate results)

a) If one of the three QC run **means** is outside a 2 S_m limit, reject run if:

1. 1 3S Rule—Run mean is outside a 3 S_m limit or
2. 2 2S Rule—Two or more of the three run means are outside the same 2 S_m limit or
3. 10 X-bar Rule—Current and previous nine run means are on the same side of the characterization mean

b) If one of the six QC **individual** results is outside a 2 S_i limit, reject run if:

1. Outlier—One individual result is beyond the characterization mean $\pm 4 S_i$ or
2. R 4S Rule—Two or more of the within-run ranges in the same run exceed 4 S_w (i.e. 95 percent range limit)

Abbreviations:

S_i = Standard deviation of individual results (the limits are not shown on the chart unless run results are actually single measurements)

S_m = Standard deviation of the run means (the limits are shown on the chart)

S_w = Within-run standard deviation (the limits are not shown on the chart)

A QC program written in SAS is available from the DLS Quality Assurance Officer and should be used to apply these rules to QC data and generate Shewhart QC charts. No results for a given analyte are to be reported from an analytical run that has been declared “out of control” for that analyte as assessed by internal (bench) QC.

The initial limits are established by analyzing pool material in 20 consecutive runs. The SAS QC program is used to monitor the QC performance over time for potential shifts, trending, or changes in assay precision. For assays performed routinely, quarterly statistics (mean, SD, CV) are calculated for each pool and compared to the characterization target values. For assays performed infrequently, statistics are calculated at least annually. As more QC data become available (covering multiple lots of reagents, multiple analysts, etc.), the initial QC limits can be reevaluated and updated. QC limits can also be reevaluated and updated as a result of a non-conforming event when the assay shows a higher than expected out of control rate and the root cause investigation does not reveal a correctable course of action to bring the assay back into control. This needs to be documented by a CAPA in STARLIMS.

While a study is in progress, QC results are stored in a STARLIMS database. For runs that are not imported into the database (i.e., R&D, troubleshooting, research-type runs), QC results are stored electronically in the analyte-specific folder on the DLS network. At the conclusion of studies, complete QC records are prepared and submitted as a study QC report in STARLIMS for review by the laboratory chief, branch chief, and a DLS statistician.

C. Sample QC Criteria

Sample QC is set of criteria used to evaluate the quality of individual test result within run, and to evaluate the quality of the calibrators associated with the run. In addition to the sample QC criteria set forth in the DLS Policies and Procedures Manual that pertain to the reportable range of concentration results and calibration curves, sample QC criteria are also established for method-specific concentration and non-concentration data associated with an individual result.

The method-specific concentration and non-concentration parameters identified for sample QC evaluation, along with their associated thresholds and flagging protocols (‘Pass’, ‘Check’, ‘Warn’, ‘Fail’) are maintained and updated in the STARLIMS database, and sample QC assessment is performed and documented as part of run review process. A sample QC result flagged as ‘fail’ should not be reported. A sample QC result flagged ‘Warn’ or ‘Check’ should be reviewed both by the analyst and supervisor to determine if the quality of the result is suitable for reporting. Results that are flagged during sample QC evaluation may also be assigned one of a series of standard comment codes available in the STARLIMS database to identify the nature of the sample QC flag

The following parameters are subject to sample QC evaluation in this method. For details on how the sample QC criteria are used, refer to **Appendix B_D: JA-4013-DR-03-STARLIMS Data Review and Criteria**.

(1) Calibration curve R^2 should be:

- a) ≥ 0.98 for MET, MFO, FOT and MYT (pass; code 0)
- b) ≥ 0.95 for PGA and THF (pass; code 0)

- (2) Calibrator difference to the target and calibration drift (drift between front and back calibrator injections) should be:
 - a) Calibrator S1 within 30% for all analytes (pass; code 0)
 - b) Calibrators S2-S5 within 15% for MET, PGA and MFO (pass; code 0)
 - c) Calibrator S2 within 20% for FOT, THF and MYT (pass; code 0)
 - d) Calibrators S3-S5 within 15% for all analytes (pass; code 0)
- (3) Relative retention time (retention time quantitation ion/retention time ISTD) <minimum or >maximum threshold should be:
 - a) ≥ 0.98 to 1.02 min for MET and PGA (pass; code 0)
 - b) ≥ 0.985 to 1.02 min for MFO (pass; code 0)
 - c) ≥ 0.985 to 1.018 min for FOT (pass; code 0)
 - d) ≥ 0.963 to 1.06 min for THF (pass; code 0)
 - e) ≥ 0.982 to 1.02 min for MYT (pass; code 0)
- (4) ISTD peak area minimum threshold should be:
 - a) For major analytes: ≥ 10000 for MET, ≥ 7000 for PGA, and ≥ 9000 for MFO; (pass; code 0)
 - b) For minor analytes: ≥ 9000 for FOT, ≥ 6000 for THF, and ≥ 8000 for MYT; (pass; code 0)
 - c) If ISTD peak area is less than minimum threshold (check; enter code 23), then verify instrument result, troubleshoot, and address appropriately (perform instrument maintenance, repair, etc.). Finally, perform a test run to check instrument response and repeat sample analysis.
- (5) Patient sample result:
 - a) >High calibrator (fail; code 26); dilute/repeat (pass; code 97)
 - b) Result <LOD or <LOQ for MET (check or incomplete; code 37 or 33); Repeat/confirm (warn; code 37 and 33)
 - c) Result <LOD for PGA and MFO (check; code 37); Repeat/confirm (warn; code 37)
 - d) MET result \geq LOD and <7 nmol/L (incomplete; code 33); Repeat/confirm (warn; code 33)
 - e) No peak for major analyte results: MET (fail; code 26); PGA & MFO (warn; code 37); Repeat/confirm (warn; code 37)
 - f) Delta difference between repeat results $\leq 15\%$; otherwise, repeat/confirm
 - g) Signal/noise ratio <3 & raw result >LOD (warn; code 0); check dilution factor
 - h) Specimen volume less than expected for analysis (code 21); set no reportable (code 98)

- (6) Manual integration (check; no code assigned); verify instrument peak integration is correct; enter code '0' if result is correct.
- (7) Check and update results for appropriate comment codes, e.g., not enough specimen for repeat analysis (code 22; set no reportable [code 98]); lab error, spills, contamination etc. (code 23; set no reportable [code 98]); or instrument error/failure (code 24; repeat analysis).

11. Remedial Action if Calibration or QC Systems Fail to Meet Acceptable Criteria

For initial steps to investigate QC failures, see **Appendix B_D: JA-4013-DR-04-Out-of-Control Corrective Action**. Additional steps are provided as a general guideline for identifying possible problems resulting in "out of control" values for QC materials. The troubleshooting process should be done in consultation with the supervisor and may involve additional experiments beyond what is indicated below.

- (A) Look for sample preparation errors, i.e., if the analyst forgot to add internal standard, specimen, right volume of buffer etc.
- (B) Check the calibrations of the pipettes.
- (C) Check to make sure that the Instrument hardware is functioning properly. Check the autosampler for proper sample injections. Check the proper gas flow for curtain, exhaust, and source from the nitrogen generator. Make sure the mass spec calibrations are ok. Run PPG's in Q1 and Q3 scan mode to verify the instrument calibration.
- (D) Run folate standards in Q1 scan mode to see if molecular ion is detected.
- (E) If the steps outlined above do not result in correction of the "out of control" values for QC materials, consult the supervisor for other appropriate corrective actions.
- (F) Do not report analytical results for runs not in statistical control.

12. Limitations of Method; Interfering Substances and Conditions

- (A) The most common causes of error are intermittently inaccurate micro-pipettors or pipetting errors.
- (B) Calibrators, internal standards, and specimens should be vortex-mixed thoroughly before pipetting.
- (C) Handling calibrators and internal standards in stepwise sequential manner will minimize the chances of cross-contamination.
- (D) Working bench should be cleaned and small bench top waste bags should be emptied daily into the metal waste bins. The blue pads should be replaced weekly to keep the work area clean and free of contamination.
- (E) Change of gloves after the preparation of stock and working standards and internal standards are recommended to avoid any contamination.
- (F) All solvents should be degassed before use.
- (G) Buffers should be made fresh as needed and pH checked when prepared.

- (H) Ascorbic acid powder should be added to sample solvents #1 (0.5% w/v), #2 (0.1% w/v), and #3 (0.5% w/v) only prior to use.
- (I) Acetic acid should be added to sample solvent #3 (1% v/v) and to the mobile phase (0.5% v/v) only prior to use.
- (J) Hemolyzed serum samples may give falsely elevated values.
- (K) Samples (including calibrators) should be prepared in yellow subdued light.
- (L) Multiple freeze/thaw cycles of specimens for extended time at ambient temperature (+15°C to +30°C) will cause degradation of folates and should be avoided.
- (M) Nitrogen gas cylinder for 96-probe SPE system (Caliper-Zephyr) should be periodically monitored for gas. The pressure for out flow should be 40 *psi*, and the gas pressure to the instrument at the regulator should be adjusted to 5 *psi*. Change the cylinder before the gauge reads 500 *psi*.
- (N) HPLC system (lines and column) should be purged and primed properly.
- (O) Interference testing was performed on this method as part of its method validation and is documented in **Appendix A**. This method has also undergone a series of in-house ruggedness testing experiments designed to assess how much method accuracy changes when certain experimental parameters are varied. A total of five parameters judged to most likely affect the accuracy of the method have been identified and tested. Testing generally consisted of performing replicate measurements on a test specimen with the selected parameter set at a value substantially lower and higher than that specified in this method while holding all other experimental variables constant. **The ruggedness testing findings for this method are presented in Appendix E**. Please refer to the *DLS Policies and Procedures Manual* for further information on ruggedness testing.

13. Reference Ranges (Normal Values)

Clinical reference ranges reported for serum folate are 11–36 nmol/L with the microbiologic assay and 7–36 with chemiluminescence assays [13].

Post-fortification reference ranges for the U.S. population generated with the microbiologic assay for NHANES 2006–2010 are shown below [14].

Serum folate: 12.7–104 nmol/L (2.5th–97.5th percentile; *n* = 23,528)

In our hands, the LC/MS/MS method and the microbiologic assay give relatively comparable results [2, 9, 10]. Pfeiffer et al. also reported microbiologic assay-equivalent reference ranges for pre-fortification (NHANES 1988–1994) and early post-fortification (NHANES 1999–2004) periods, as well as reference ranges by population subgroups for all three time periods [14].

Reference ranges for individual folate vitamers based on LC-MS/MS were determined using data from NHANES 2011–2012 [11] and updated with data from NHANES 2011–2016 [15]. While these values are of interest to researchers, they cannot be interpreted clinically. Serum tFOL levels <7 nmol/L (3 ng/mL) represent a negative folate balance and are usually indicative of inadequate folate intake [16].

14. Critical Call Results (“Panic Values”)

Any NHANES samples with serum tFOL levels <7 nmol/L (3 ng/mL) are considered to require follow-up. Since survey data are transmitted approximately weekly to Westat, abnormal reports are automatically forwarded to the NCHS survey physician for follow-up. For smaller, non-NHANES studies, abnormal values are identified to the study principal investigator. Emails sent concerning abnormal results are maintained by the supervisor for the duration of the study. Most of these studies are epidemiological in nature.

15. Specimen Storage and Handling during Testing

Specimens should be brought and maintained at ambient temperature (+15°C to +30°C) during preparation and testing and then returned to deep frozen storage (-50°C to -90°C) as soon as possible.

16. Alternate Methods for Performing Test of Storing Specimens if Test System Fails

If only tFOL is of interest, the microbiologic assay could be performed instead of the LC-MS/MS method under some circumstances.

17. Test Result Reporting System; Protocol for Reporting Critical Calls (If Applicable)

Test results are reported to the collaborating agency at a frequency and by a method determined by the study coordinator. Generally, data from this analysis are compiled with results from other analyses and sent to the responsible person at the collaborating agency as an Excel file, generally through electronic mail or via ftp site. Data are transmitted via the CLIA Director after review by the Lab Supervisor, Branch Chief, and a CDC Statistician.

For NHANES 1999+, all data are reported electronically approximately weekly to WESTAT who then transfer the results to NCHS. For some smaller studies, electronic copies of a data report are sent and upon request hard copies can be sent as well.

18. Transfer or Referral of Specimens; Procedures for Specimen Accountability and Tracking

This protocol does not involve referral of specimens for testing the analytes of this method at another laboratory.

A STARLIMS database is used to keep records and track specimens for NHANES 1999+. If plasma or serum folate analyses are used for smaller, non-NHANES studies, records may be kept in Excel files on the network.

We recommend that records, including related QA/QC data, be maintained for 10 years after completion of the NHANES study. Only numerical identifiers should be used (e.g., case ID numbers). All personal identifiers should be available only to the medical supervisor or project coordinator. Residual serum from these analyses for non-NHANES studies are retained for at least 1 year after results have been reported and may be then returned or discarded at the request of the principal investigator. Very little residual material will be available after NHANES analyses are completed, however residual serum is retained for at least 2 years after results have been publicly released; at that point, samples with sufficient volume (>0.2 mL) are returned to NHANES and samples with insufficient volume may be autoclaved.


The exact procedure used to track specimens varies with each study and is specified in the study protocol or the interagency agreement for the study. Copies of these documents are kept by the supervisor. In general, when specimens are received, the specimen ID number is entered into a database and the specimens stored deep frozen (-50°C to -90°C). The specimen ID is read off the vial by a barcode reader used to prepare the electronic specimen table for the analytical system. When the analyses are completed, the result file is loaded into the database. The analyst is responsible for keeping records of specimens prepared incorrectly, those with labeling problems, and those with abnormal results, together with information about these discrepancies. In general, these are documented using codes in the STARLIMS database.

19. Summary Statistics and QC Graph

Please see following page.


**August 2021 – August 2023 Summary Statistics and QC Chart
 LBXSF1 (5-Methyl-tetrahydrofolate (nmol/L))**

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBF1H10	131	28SEP21	11SEP23	41.6729	0.9010	2.2
4758_LBF1L10	131	28SEP21	11SEP23	8.1931	0.1830	2.2
4759_LBF1M10	131	28SEP21	11SEP23	23.7073	0.5128	2.2




August 2021 – August 2023 Summary Statistics and QC Chart LBXSF2 (Folic acid, serum (nmol/L))

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBCFSH9	130	28SEP21	11SEP23	9.6411	0.3566	3.7
4758_LBCFSL9	130	28SEP21	11SEP23	0.6072	0.0403	6.6
4759_LBFISM10	130	28SEP21	11SEP23	5.3218	0.1904	3.6




**August 2021 – August 2023 Summary Statistics and QC Chart
 LBXSF3 (5-Formyl-tetrahydrofolic acid (nmol/L))**

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBC3FH7	122	28SEP21	11SEP23	5.5598	0.1809	3.3
4759_LBC3FM8	122	28SEP21	11SEP23	2.7966	0.0987	3.5




**August 2021 – August 2023 Summary Statistics and QC Chart
 LBXSF4 (Tetrahydrofolic acid, serum (nmol/L))**

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBCF4H7	122	28SEP21	28AUG23	4.9600	0.1800	3.6
4759_LBCF4M8	122	28SEP21	28AUG23	2.7137	0.1129	4.2




**August 2021 – August 2023 Summary Statistics and QC Chart
 LBXSF5 (5,10--Methenyl-tetthrofolc acid (nmol/L))**

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBCF5H7	123	28SEP21	28AUG23	4.9372	0.1827	3.7
4759_LBCF5M8	123	28SEP21	28AUG23	2.4983	0.0824	3.3



**August 2021 – August 2023 Summary Statistics and QC Chart
 LBXSF6 (Mefox oxidation product, serum (nmol/L))**

Lot	n	Start Date	End Date	mean	Standard Deviation	Coefficient of Variation
4760_LBCF6H7	123	28SEP21	28AUG23	2.0211	0.0922	4.6
4758_LBCF6L7	123	28SEP21	28AUG23	4.5276	0.1470	3.2
4759_LBCF6M8	123	28SEP21	28AUG23	1.3110	0.0668	5.1



20. Method Performance Documentation

Method performance documentation for this method including accuracy, precision, sensitivity, specificity, and stability is provided in **Appendix A**. The approval of this procedure by the Branch Chief and CLIA Director denote that the method performance is fit for the intended use of the method.

References

- (1) Bailey LB, Stover PJ, McNulty H, Fenech MF, Gregory III JF, Mills JL, Pfeiffer CM, Fazili Z, Zhang M, et al. Biomarkers of nutrition for development – folate review. *J Nutr.* 2015;145:1636S–80S.
- (2) Fazili Z, Whitehead Jr. RD, Paladugula N, Pfeiffer CM. A high-throughput LC-MS/MS method suitable for population biomonitoring measures five serum folate vitamers and one oxidation product. *Anal Bioanal Chem.* 2013;405:4549–60.
- (3) Pfeiffer CM, Fazili Z, McCoy L, Zhang M, and Gunter EW. Determination of folate vitamers in human serum by stable-isotope dilution tandem mass spectrometry and comparison to radioassay and microbiologic assay. *Clin Chem.* 2004;50:423–32.
- (4) Fazili Z, Pfeiffer CM. Measurement of folates in serum and conventionally prepared whole blood lysates: Application of an automated 96-well plate isotope-dilution tandem mass spectrometry method. *Clin Chem.* 2004;50:2378–81.
- (5) Fazili Z, Pfeiffer CM. Accounting for an isobaric interference allows correct determination of folate vitamers in serum by isotope dilution-liquid chromatography-tandem MS. *J Nutr.* 2013;143:108–113.
- (6) Fazili A, Sternberg MR, Paladugula N, Whitehead RD, Chen H, Pfeiffer CM. The loss of 5-methyl-tetrahydrofolate in human serum under suboptimal preanalytical conditions can only partially be recovered by an oxidation product. *J Nutr.* 2014;144:1873–9.
- (7) Mastropaolo W, Wilson MA. Effect of light on serum B12 and folate stability. *Clin Chem.* 1993;39(5):913.
- (8) Biamonte AR, Schneller GH. Study of folic acid stability in solutions of the B complex vitamins. *J Am Pharmaceut Assoc.* 1951;XL (7):313–20.
- (9) Fazili Z, Pfeiffer CM, Zhang M. Comparison of serum folate species analyzed by LC-MS/MS with total folate measured by microbiologic assay and Bio-Rad radioassay. *Clin Chem.* 2007;53:781–4.
- (10) Zhang M, Sternberg MR, Pfeiffer CM. Harmonizing the calibrator and microorganism used in the folate microbiologic assay leads to improvements in serum and whole blood folate comparability as shown in a CDC Round Robin study. *J Nutr.* 2018;148:807–17.
- (11) Pfeiffer CM, Sternberg MR, Fazili Z, Lacher DA, Zhang, M, Johnson CL, Hamner HC, Bailey RL, Radar JI, Yamini S, et al. Folate status and concentrations of serum folate forms in the US population: National Health and Nutrition Examination Survey 2011–2. *Br J Nutr.* 2015;113:1965–77.
- (12) Caudill SP, Schleicher RL, Pirkle JL. Multi-rule quality control for the age-related eye disease study. *Stat Med.* 2008;27:4094–4106.
- (13) Tietz NW. *Clinical Guide to Laboratory Tests.* Third Edition. 1995. WB Saunders Company.
- (14) Pfeiffer CM, Hughes JP, Lacher DA, Bailey RL, Berry RJ, Zhang M, Yetley EA, Rader JI, Sempos CT, Johnson CL. Estimation of trends in serum and RBC folate in the U.S. population from pre-

to postfortification using assay-adjusted data from the NHANES 1988-2010. *J Nutr.* 2012;142:886–93.

- (15) Fazili Z, Sternberg MR, Potischman N, Wang C-Y, Storaandt RJ, Yeung L, Yamini S, Gahche JJ, Juan WY, Qi YP, Paladugula N, Gabey G, Pfeiffer CM. Demographic, physiologic, and lifestyle characteristics observed with serum total folate differ among folate forms: cross-sectional data from fasting samples in the NHANES 2011-2016. *J Nutr.* 2020;150:851–60.
- (16) Institute of Medicine, Food and Nutrition Board. Dietary reference intakes: Thiamin, riboflavin, niacin, vitamin B6, folate, vitamin B12, pantothenic acid, biotin, and choline. Washington, D.C.: National Academy Press; 1998.

Acknowledgements

We gratefully acknowledge the contributions of Zia Fazili-Qari, Ph.D., Phuong Ngac, M.S., and Christine M. Pfeiffer, Ph.D. who assisted in developing the methodology and preparing the manuscript for this chapter.

Appendix A: Method Performance Documentation

A. Accuracy

Accuracy compared to Reference Material - fill in yellow shaded cells												
Mean concentration should be within $\pm 15\%$ of the nominal value except at 3^*LOD , where it should be within $\pm 20\%$												
Method name:	Folate Forms by LC-MS/MS											
Method #:	4013.04											
Matrix:	Serum											
Units:	nmol/L											
Reference material:	NIST & NIBSC materials											
Analyte:	Folate Forms											
5-MethylTHF			Measured concentration									Difference from nominal value (%)
Reference material	Replicate	Nominal value	Day 1	Day 2	Day 3	Day 4	Day 5	Mean	SD	CV (%)		
NIST SRM 1955 Level 1	1	4.26	4.36	4.75	4.78	4.97	4.66	4.71	0.21	4.46	10.6	
	2		4.37	4.75	4.78	4.99	4.72					
NIST SRM 1955 Level 2	1	9.73	10.00	10.39	10.50	11.00	10.58	10.5	0.32	3.02	7.7	
	2		10.02	10.33	10.52	10.84	10.62					
NIST SRM 1955 Level 3	1	37.1	37.30	39.34	37.71	39.01	38.14	38.2	0.80	2.08	2.9	
	2		36.98	38.79	37.96	38.89	37.62					
NIST SRM 1950	1	26.93	27.59	28.90	27.83	29.65	28.77	28.7	0.80	2.77	6.6	
	2		28.16	28.93	28.46	30.22	28.47					
NIBSC 03-178	1	9.75	10.65	10.00	10.40	10.05	10.06	10.2	0.25	2.42	4.2	
	2		10.40	10.05	10.16	9.87	9.95					

Accuracy using Spike Recovery - fill in yellow shaded cells											
Recovery = (final concentration – initial concentration)/added concentration											
Recovery should be 85-115% except at 3^*LOD where can be 80-120%											
Method name:	Folate Forms by LC-MS/MS										
Method #:	4013.04										
Matrix:	Serum										
Units:	nmol/L										
Analyte:	5-MethylTHF										
Replicate	Spike concentration	Sample 1				Sample 2				Mean recovery (%)	SD (%)
		Measured concentration			Recovery (%)	Measured concentration			Recovery (%)		
		Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			Day 1 (12/07/17)	Day 2 (12/08/17)		Mean	
Sample	1	0	12.5	11.8	12.40	0	17.5	17.1	17.45	96	2.3
	2		12.8	12.3			17.9	17.4			
	3		12.8	12.2			17.6	17.2			
Sample + Spike 1	1	5	17.5	16.8	17.22	5	22.1	21.5	22.08	93	
	2		17.5	17.0			22.6	21.8			
	3		17.6	16.9			22.3	22.2			
Sample + Spike 2	1	10	22.1	21.8	22.15	10	26.8	27.1	27.10	97	
	2		22.5	21.7			27.4	26.7			
	3		23.0	21.8			27.2	27.4			
Sample + Spike 3	1	20	32.0	31.4	32.15	20	36.0	34.9	36.20	94	
	2		33.2	31.5			37.0	35.5			
	3		33.2	31.6			37.2	36.6			

Folate Forms - Total & Individual - Serum
NHANES August 2021 - August 2023

Analyte:		MeFox											
Replicate		Sample 1					Sample 2					Mean recovery (%)	SD (%)
		Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration			Recovery (%)		
			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			
Sample	1	0	4.37	4.05	4.2		0	5.52	5.56	5.64		85	1.8
	2		4.68	3.95				5.67	5.49				
	3		4.23	3.86				5.48	6.10				
Sample + Spike 1	1	5	8.86	7.89	8.5	86	5	9.84	9.74	9.8	83		
	2		8.39	8.68				9.33	9.79				
	3		8.53	8.66				10.30	9.77				
Sample + Spike 2	1	10	11.8	12.6	12.7	85	10	14.0	14.5	14.1	84		
	2		12.5	13.4				13.8	13.3				
	3		12.7	13.2				14.3	14.5				
Sample + Spike 3	1	20	22.8	21.6	21.9	88	20	21.7	22.4	22.8	86		
	2		21.3	22.1				22.1	24.1				
	3		21.6	21.9				23.1	23.1				

Analyte:		Folic acid											
Replicate		Sample 1					Sample 2					Mean recovery (%)	SD (%)
		Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration			Recovery (%)		
			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			
Sample	1	0	0.33	0.43	0.41		0	0.55	0.62	0.60		107	2.8
	2		0.38	0.44				0.56	0.64				
	3		0.32	0.54				0.57	0.65				
Sample + Spike 1	1	5	5.60	5.34	5.6	103	5	5.85	5.81	5.8	105		
	2		5.35	5.74				6.04	5.64				
	3		5.50	5.78				6.01	5.63				
Sample + Spike 2	1	10	11.1	11.3	11.2	108	10	11.0	11.2	11.4	108		
	2		10.9	11.6				11.5	11.5				
	3		10.9	11.6				11.4	12.0				
Sample + Spike 3	1	20	22.4	22.5	22.4	110	20	20.8	21.8	21.5	105		
	2		23.1	22.3				21.6	21.3				
	3		21.8	22.5				20.9	22.8				

Analyte:		Tetrahydrofolate											
Replicate		Sample 1					Sample 2					Mean recovery (%)	SD (%)
		Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration			Recovery (%)		
			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			Day 1 (12/07/17)	Day 2 (12/08/17)	Mean			
Sample	1	0	0.39	0.18	0.27		0	0.62	0.37	0.45		97	1.4
	2		0.25	0.26				0.44	0.39				
	3		0.33	0.18				0.55	0.31				
Sample + Spike 1	1	5	4.92	5.04	5.11	97	5	5.25	5.48	5.16	94		
	2		4.83	5.40				4.94	5.08				
	3		5.19	5.27				4.69	5.54				
Sample + Spike 2	1	10	9.76	9.83	9.95	97	10	9.58	10.30	10.33	99		
	2		9.72	10.10				9.84	11.20				
	3		9.79	10.50				9.95	11.10				
Sample + Spike 3	1	20	19.40	19.90	19.73	97	20	19.30	20.70	19.85	97		
	2		18.80	20.70				19.20	20.50				
	3		19.40	20.20				18.80	20.60				

Folate Forms - Total & Individual - Serum
NHANES August 2021 - August 2023

Analyte:		5-FormylTHF											
Replicate	Spike concentration	Sample 1					Sample 2					Mean recovery (%)	SD (%)
		Day 1 (12/07/17)	Day 2 (12/08/17)	Mean	Recovery (%)	Day 1 (12/07/17)	Day 2 (12/08/17)	Mean	Recovery (%)				
Sample	1	0	0.14	0.14	0.14	0	0.14	0.14	0.14				
	2	0	0.14	0.14	0.14	0	0.14	0.14	0.14				
	3	0	0.14	0.14	0.14	0	0.14	0.14	0.14			97	2.4
Sample + Spike 1	1	5	5.14	4.72	5.00	97	5	4.8	4.81	4.85	94		
	2	5	4.67	5.42	5.00	97	5	4.98	4.79	4.85	94		
	3	5	4.88	5.14	5.00	97	5	4.76	4.94	4.85	94		
Sample + Spike 2	1	10	9.92	10.5	10.09	100	10	8.96	10.3	9.58	94		
	2	10	9.64	10.2	10.09	100	10	8.96	9.7	9.58	94		
	3	10	10.1	10.2	10.09	100	10	9.34	10.2	9.58	94		
Sample + Spike 3	1	20	19.6	20.6	20.1	100	20	18.7	19.5	19.5	97		
	2	20	19.8	21.0	20.1	100	20	18.7	20.1	19.5	97		
	3	20	19.6	20.0	20.1	100	20	18.8	20.9	19.5	97		

Analyte:		5, 10 -MethenylTHF											
Replicate	Spike concentration	Sample 1					Sample 2					Mean recovery (%)	SD (%)
		Day 1 (12/07/17)	Day 2 (12/08/17)	Mean	Recovery (%)	Day 1 (12/07/17)	Day 2 (12/08/17)	Mean	Recovery (%)				
Sample	1	0	0.14	0.14	0.14		0	0.14	0.14	0.14			
	2	0	0.14	0.14	0.14		0	0.14	0.14	0.14			
	3	0	0.14	0.14	0.14		0	0.14	0.14	0.14			
Sample + Spike 1	1	5	6.07	5.17	5.52	108	5	5.24	5.13	5.48	107		
	2	5	5.79	5.45	5.52	108	5	5.81	5.11	5.48	107		
	3	5	5.42	5.21	5.52	108	5	5.75	5.84	5.48	107		
Sample + Spike 2	1	10	10.50	10.80	10.92	108	10	11.40	10.50	10.77	106		
	2	10	11.00	11.50	10.92	108	10	11.30	9.80	10.77	106		
	3	10	11.00	10.70	10.92	108	10	10.50	11.10	10.77	106		
Sample + Spike 3	1	20	23.9	21.5	22.6	112	20	23.9	21.0	22.0	109		
	2	20	22.0	21.8	22.6	112	20	21.4	21.6	22.0	109		
	3	20	24.3	21.9	22.6	112	20	22.1	21.7	22.0	109		

B. Precision

Precision - fill in yellow shaded cells													
Total relative standard deviation should be $\leq 15\%$ (CV $\leq 15\%$)													
Method name:	Folate Forms by LC-MS/MS												
Method #:	4013.04												
Matrix:	Serum												
Units:	nmol/L												
Run dates:	2017-03-09 to 2017-03-21												
Analyte:	5-MethylTHF												
Quality material 1 (MS16431)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	24.5	23.6	24.05	0.202	0.202	1157	1	42.8	41.7	42.25	0.302	0.302	3570
2	25.2	24.3	24.75	0.202	0.202	1225	2	44.4	43.2	43.80	0.360	0.360	3837
3	23.9	24.4	24.15	0.063	0.063	1166	3	42.6	42.1	42.35	0.063	0.063	3587
4	24.9	25.0	24.95	0.003	0.003	1245	4	44.0	43.2	43.60	0.160	0.160	3802
5	24.4	24.2	24.30	0.010	0.010	1181	5	43.1	42.8	42.95	0.022	0.023	3689
6	24.1	24.0	24.05	0.003	0.003	1157	6	42.9	41.9	42.40	0.250	0.250	3596
7	25.2	25.0	25.10	0.010	0.010	1260	7	43.6	44.1	43.85	0.063	0.063	3846
8	24.2	24.2	24.20	0.000	0.000	1171	8	42.9	42.4	42.65	0.063	0.063	3638
9	24.7	25.2	24.95	0.063	0.063	1245	9	43.9	44.1	44.00	0.010	0.010	3872
10	24.7	24.7	24.70	0.000	0.000	1220	10	43.4	43.6	43.50	0.010	0.010	3785
Grand sum	490.4	Grand mean	24.52				Grand sum	862.7	Grand mean	43.14			
		Std Dev	0.41						Std Dev	0.69			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	1.1100	0.1110	0.3332	1.36			Within Run	2.6050	0.2605	0.5104	1.18		
Between Run	3.0420	0.3380	0.3369	1.37			Between Run	8.5205	0.9467	0.5858	1.36		
Total	4.15		0.4738	1.93			Total	11.13		0.7769	1.80		

Analyte:	MeFox												
Quality material 1 (LS16430)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	4.77	4.76	4.77	0.000	0.000	45.4105	1	2.02	2.09	2.06	0.001	0.001	8.4461
2	5.03	5.05	5.04	0.000	0.000	50.8032	2	2.34	2.15	2.25	0.009	0.009	10.0801
3	4.85	4.68	4.77	0.007	0.007	45.4105	3	2.09	2.02	2.06	0.001	0.001	8.4461
4	4.71	4.73	4.72	0.000	0.000	44.5568	4	2.00	2.12	2.06	0.004	0.004	8.4872
5	4.82	4.50	4.66	0.026	0.026	43.4312	5	2.10	2.13	2.12	0.000	0.000	8.9465
6	5.14	5.08	5.11	0.001	0.001	52.2242	6	2.33	2.01	2.17	0.026	0.026	9.4178
7	4.72	4.69	4.71	0.000	0.000	44.2741	7	2.06	2.02	2.04	0.000	0.000	8.3232
8	4.60	4.56	4.58	0.000	0.000	41.9528	8	2.05	2.08	2.07	0.000	0.000	8.5285
9	5.14	4.94	5.04	0.010	0.010	50.8032	9	2.17	2.18	2.18	0.000	0.000	9.4613
10	4.88	4.74	4.81	0.005	0.005	46.2722	10	2.11	2.14	2.13	0.000	0.000	9.0313
Grand sum	96.39	Grand mean	4.82				Grand sum	42.21	Grand mean	2.11			
		Std Dev	0.18						Std Dev	0.07			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.0989	0.0099	0.0995	2.06			Within Run	0.0836	0.0084	0.0914	4.33		
Between Run	0.5869	0.0652	0.1663	3.45			Between Run	0.0835	0.0093	0.0215	1.02		
Total	0.69		0.1938	4.02			Total	0.17		0.0939	4.45		

Folate Forms - Total & Individual - Serum
NHANES August 2021 - August 2023

Analyte:		Folic acid											
Quality material 1 (LS16430)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	0.62	0.62	0.62	0.0000	0.0000	0.769	1	9.79	9.49	9.64	0.022	0.023	185.86
2	0.63	0.67	0.65	0.0004	0.0004	0.845	2	9.83	10	9.92	0.007	0.007	196.61
3	0.63	0.7	0.67	0.0012	0.0012	0.884	3	10.3	9.54	9.92	0.144	0.144	196.81
4	0.62	0.63	0.63	0.0000	0.0000	0.781	4	10.2	9.91	10.06	0.021	0.021	202.21
5	0.64	0.61	0.63	0.0002	0.0002	0.781	5	10.1	9.92	10.01	0.008	0.008	200.40
6	0.64	0.61	0.63	0.0002	0.0002	0.781	6	9.63	9.79	9.71	0.006	0.006	188.57
7	0.61	0.65	0.63	0.0004	0.0004	0.794	7	9.83	10.4	10.12	0.081	0.081	204.63
8	0.61	0.61	0.61	0.0000	0.0000	0.744	8	9.66	9.64	9.65	0.000	0.000	186.25
9	0.60	0.64	0.62	0.0004	0.0004	0.769	9	10.2	10.5	10.35	0.023	0.023	214.25
10	0.61	0.57	0.59	0.0004	0.0004	0.696	10	9.70	9.95	9.83	0.016	0.016	193.06
Grand sum	12.52	Grand mean	0.63				Grand sum	198.38	Grand mean	9.92			
		Std Dev	0.02						Std Dev	0.22			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.0066	0.0007	0.0257	4.10			Within Run	0.6582	0.0658	0.2566	2.59		
Between Run	0.0075	0.0008	0.0092	1.48			Between Run	0.9074	0.1008	0.1323	1.33		
Total	0.01		0.0273	4.36			Total	1.57		0.2887	2.91		

Analyte:		Tetrahydrofolate											
Quality material 1 (MS16431)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	2.64	2.21	2.43	0.046	0.046	11.761	1	4.59	4.25	4.42	0.029	0.029	39
2	2.53	2.52	2.53	0.000	0.000	12.751	2	4.85	4.54	4.70	0.024	0.024	44
3	2.09	2.55	2.32	0.053	0.053	10.765	3	4.32	4.46	4.39	0.005	0.005	39
4	2.47	2.53	2.50	0.001	0.001	12.500	4	4.89	4.65	4.77	0.014	0.014	46
5	2.73	2.66	2.70	0.001	0.001	14.526	5	4.73	4.97	4.85	0.014	0.014	47
6	2.66	2.70	2.68	0.000	0.000	14.365	6	5.02	4.78	4.90	0.014	0.014	48
7	2.49	2.62	2.56	0.004	0.004	13.056	7	4.96	4.88	4.92	0.002	0.002	48
8	2.57	2.40	2.49	0.007	0.007	12.350	8	4.74	4.68	4.71	0.001	0.001	44
9	2.75	2.89	2.82	0.005	0.005	15.905	9	4.78	4.87	4.83	0.002	0.002	47
10	2.78	2.47	2.63	0.024	0.024	13.781	10	4.90	4.74	4.82	0.006	0.006	46
Grand sum	51.26	Grand mean	2.563				Grand sum	94.6	Grand mean	4.73			
		Std Dev	0.15						Std Dev	0.19			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.2841	0.0284	0.1686	6.58			Within Run	0.2239	0.0224	0.1496	3.16		
Between Run	0.3813	0.0424	0.0835	3.26			Between Run	0.6229	0.0692	0.1530	3.23		
Total	0.67		0.1881	7.34			Total	0.85		0.2140	4.52		

Folate Forms - Total & Individual - Serum
NHANES August 2021 - August 2023

Analyte:		5-FormylTHF											
Quality material 1 (MS16431)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	2.59	2.6	2.60	0.0000	0.0000	13.4681	1	5.18	5.33	5.26	0.006	0.006	55.230
2	2.69	2.75	2.72	0.0009	0.0009	14.7968	2	5.47	5.48	5.48	0.000	0.000	59.951
3	2.71	2.64	2.68	0.0012	0.0012	14.3113	3	5.33	5.35	5.34	0.000	0.000	57.031
4	2.67	2.80	2.74	0.0042	0.0042	14.9605	4	5.71	5.42	5.57	0.021	0.021	61.938
5	2.80	2.64	2.72	0.0064	0.0064	14.7968	5	5.52	5.39	5.46	0.004	0.004	59.514
6	2.78	2.67	2.73	0.0030	0.0030	14.8513	6	5.63	5.45	5.54	0.008	0.008	61.383
7	2.73	2.66	2.70	0.0012	0.0012	14.5261	7	5.34	5.39	5.37	0.001	0.001	57.566
8	2.70	2.69	2.70	0.0000	0.0000	14.5261	8	5.25	5.30	5.28	0.001	0.001	55.651
9	2.94	2.85	2.90	0.0020	0.0020	16.7621	9	5.48	5.63	5.56	0.006	0.006	61.716
10	2.77	2.86	2.82	0.0020	0.0020	15.8485	10	5.41	5.16	5.29	0.016	0.016	55.862
Grand sum	54.54	Grand mean	2.73				Grand sum	108.2	Grand mean	5.41			
		Std Dev	0.08						Std Dev	0.12			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.0422	0.0042	0.0650	2.38			Within Run	0.1232	0.0123	0.1110	2.05		
Between Run	0.1166	0.0130	0.0661	2.42			Between Run	0.2660	0.0296	0.0928	1.72		
Total	0.16		0.0927	3.40			Total	0.39		0.1447	2.67		

Analyte:		5,10-MethenylTHF											
Quality material 1 (MS16431)							Quality material 2 (HS16432)						
Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2	Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	2.54	2.67	2.61	0.0042	0.0042	13.5721	1	5.17	4.88	5.03	0.021	0.021	50.501
2	2.74	2.64	2.69	0.0025	0.0025	14.4722	2	5.20	5.55	5.38	0.031	0.031	57.781
3	2.37	2.38	2.38	0.0000	0.0000	11.2813	3	4.61	4.47	4.54	0.005	0.005	41.223
4	2.93	2.20	2.57	0.1332	0.1332	13.1585	4	5.54	4.89	5.22	0.106	0.106	54.392
5	2.63	2.70	2.67	0.0012	0.0012	14.2045	5	5.26	5.00	5.13	0.017	0.017	52.634
6	2.73	2.60	2.67	0.0042	0.0042	14.2045	6	5.15	5.09	5.12	0.001	0.001	52.429
7	2.74	2.76	2.75	0.0001	0.0001	15.1250	7	5.24	5.48	5.36	0.014	0.014	57.459
8	2.71	2.69	2.70	0.0001	0.0001	14.5800	8	5.36	5.42	5.39	0.001	0.001	58.104
9	3.00	2.67	2.84	0.0272	0.0272	16.0745	9	5.45	5.40	5.43	0.001	0.001	58.861
10	2.74	2.79	2.77	0.0006	0.0006	15.2905	10	5.14	5.27	5.21	0.004	0.004	54.184
Grand sum	53.23	Grand mean	2.66				Grand sum	103.6	Grand mean	5.18			
		Std Dev	0.13						Std Dev	0.26			
	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)				Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)		
Within Run	0.3470	0.0347	0.1863	7.00			Within Run	0.4003	0.0400	0.2001	3.86		
Between Run	0.2911	0.0323	0.0000	0.00			Between Run	1.2322	0.1369	0.2201	4.25		
Total	0.64		0.1863	7.00			Total	1.63		0.2974	5.74		

C. Stability

Stability - fill in yellow shaded cells							
The initial measurement can be from the same day for all stability experiments.							
Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions							
Describe condition: QC material thawed 3 times (2 hrs at room temperature) and re-frozen at -70°C (3 freeze-thaw cycles)							
Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typicall at room temperature)							
Describe condition: QC material (not yet processed) stored at room temperature for 3 hrs, then processed for analysis							
Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler							
Describe condition: (1): Processed samples were kept in the autosampler overnight prior to analysis next day							
(2): Processed samples were kept at -20°C for 3 days (over the weekend) prior to analysis							
Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis							
Describe condition: QC material stored at -70°C for 1.5 years							
All stability sample results should be within $\pm 15\%$ of nominal concentration							
Method name:	Folate Forms by LC-MS/MS						
Method #:	4013.04						
Matrix:	Serum						
Units:	nmol/L						
Run date:	2017-08-31 (HZ:6500)						

Analyte:	5-MethylTHF						
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	25.10	25.00	25.30	25.50	25.00	22.98	25.10
Replicate 2	25.10	25.20	25.70	25.40	24.90	23.21	25.10
Replicate 3	25.50	25.10	24.80	25.00	24.90	23.89	25.50
Mean	25.23	25.10	25.27	25.30	24.93	23.36	25.23
% difference from initial measurement	--	-0.53	0.13	0.26	-1.19	--	8.02
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	44.50	43.70	44.20	44.00	44.20	41.11	44.50
Replicate 2	44.10	43.80	43.40	43.60	44.30	41.58	44.10
Replicate 3	44.70	44.50	44.10	43.40	43.20	40.98	44.70
Mean	44.43	44.00	43.90	43.67	43.90	41.22	44.43
% difference from initial measurement	--	-0.98	-1.20	-1.73	-1.20	--	7.79

Analyte:		MeFox					
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)*	Processed sample stability (2)**	Initial measurement	Long-term stability
Replicate 1	1.36	1.35	1.46	1.39	1.46	1.21	1.36
Replicate 2	1.34	1.32	1.41	1.42	1.37	1.23	1.34
Replicate 3	1.33	1.35	1.31	1.34	1.38	1.25	1.33
Mean	1.34	1.34	1.39	1.38	1.40	1.23	1.34
% difference from initial measurement	--	-0.25	3.72	2.98	4.47	--	9.4
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	2.12	2.07	2.09	2.05	2.10	1.88	2.12
Replicate 2	1.98	2.01	1.99	2.24	2.17	1.90	1.98
Replicate 3	2.15	2.02	2.17	1.94	2.20	1.87	2.15
Mean	2.08	2.03	2.08	2.08	2.16	1.88	2.08
% difference from initial measurement	--	-2.40	0.00	-0.32	3.52	--	10.6
Analyte:		Folic acid					
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)*	Processed sample stability (2)**	Initial measurement	Long-term stability
Replicate 1	5.75	5.62	5.68	5.73	5.74	5.26	5.75
Replicate 2	5.58	5.67	5.76	5.67	5.68	5.24	5.58
Replicate 3	5.69	5.58	5.6	5.77	5.82	5.31	5.69
Mean	5.67	5.62	5.68	5.72	5.75	5.27	5.67
% difference from initial measurement	--	-0.88	0.12	0.88	1.29	--	7.6
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	10.30	10.10	10.30	10.40	10.20	9.23	10.30
Replicate 2	10.50	10.10	10.00	10.20	10.30	9.23	10.50
Replicate 3	10.70	10.20	10.60	10.40	10.30	9.41	10.70
Mean	10.50	10.13	10.30	10.33	10.27	9.29	10.50
% difference from initial measurement	--	-3.49	-1.90	-1.59	-2.22	--	13.0

Analyte:		Tetrahydrofolate					
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	2.69	2.41	2.70	2.93	2.49	2.60	2.69
Replicate 2	2.40	2.52	2.54	2.83	2.48	2.46	2.40
Replicate 3	2.61	2.37	2.29	2.60	2.55	2.71	2.61
Mean	2.57	2.43	2.51	2.79	2.51	2.59	2.57
% difference from initial measurement	--	-5.19	-2.21	8.57	-2.34	--	-0.9
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	4.65	4.74	4.77	4.88	4.81	4.56	4.65
Replicate 2	5.01	4.64	4.45	4.75	4.62	4.54	5.01
Replicate 3	4.83	4.47	4.74	4.89	4.86	4.63	4.83
Mean	4.83	4.62	4.65	4.84	4.76	4.58	4.83
% difference from initial measurement	--	-4.42	-3.66	0.21	-1.38	--	5.5
Analyte:		5-FormylTHF					
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	2.83	2.79	2.86	2.68	2.82	2.54	2.83
Replicate 2	2.73	2.81	2.84	2.77	2.76	2.65	2.73
Replicate 3	2.72	2.67	2.8	2.79	2.81	2.54	2.72
Mean	2.76	2.76	2.83	2.75	2.80	2.58	2.76
% difference from initial measurement	--	-0.12	2.66	-0.48	1.33	--	7.2
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	5.40	5.44	5.76	5.38	5.57	5.13	5.40
Replicate 2	5.39	5.55	5.40	5.61	5.58	5.20	5.39
Replicate 3	5.56	5.55	5.54	5.54	5.55	5.18	5.56
Mean	5.45	5.51	5.57	5.51	5.57	5.17	5.45
% difference from initial measurement	--	1.16	2.14	1.10	2.14	--	5.4

Analyte: 5,10-MethenylTHF							
Quality material 1 (MS16431)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)*	Processed sample stability (2)**	Initial measurement	Long-term stability
Replicate 1	2.54	2.55	2.68	2.88	2.83	2.38	2.54
Replicate 2	2.77	2.46	2.79	2.75	2.70	2.38	2.77
Replicate 3	2.82	2.60	2.67	2.86	2.93	2.32	2.82
Mean	2.71	2.54	2.71	2.83	2.82	2.36	2.71
% difference from initial measurement	--	-6.40	0.12	4.43	4.06	--	14.7
Quality material 2 (HS16432)							
	Initial measurement	Three freeze-thaw cycles	Bench-top stability	Processed sample stability (1)	Processed sample stability (2)	Initial measurement	Long-term stability
Replicate 1	5.38	5.27	5.37	5.36	5.55	4.63	5.38
Replicate 2	5.28	5.21	4.92	5.31	5.66	4.72	5.28
Replicate 3	5.29	5.42	5.21	5.49	5.42	4.61	5.29
Mean	5.32	5.30	5.17	5.39	5.54	4.65	5.32
% difference from initial measurement	--	-0.31	-2.82	1.32	4.26	--	14.2

D. LOD, Specificity, and Fit for Intended Use

LOD, specificity and fit for intended use - fill in yellow shaded cells				
Method name:	Folate Forms by LC-MS/MS			
Method #:	4013.04			
Matrix:	Serum			
Units:	nmol/L			
Analytes	Mass Spectrometer Instrument (2017-01-18)	Limit of Detection (LOD)	Interferences successfully checked in at least 50 human samples	Accuracy, precision, LOD, specificity and stability meet performance specifications for intended use
5-MethylTHF	Hazel (6500)	0.13	yes	yes
Folic acid	Hazel (6500)	0.14	yes	yes
5-FormylTHF	Hazel (6500)	0.20	yes	yes
Tetrahydrofolate	Hazel (6500)	0.25	yes	yes
5, 10-MethenylTHF	Hazel (6500)	0.20	yes	yes
MeFox	Hazel (6500)	0.10	yes	yes

Appendix B: Job Aids

A. General:

1) JA-4013-G-01-Calibration and Calibration Verification

a) Calibration

This assay is calibrated daily by a 5-point calibration curve (linear and 1/x² weighted). Calibration is performed with the following concentrations:

- 5-MethylTHF at 1, 2, 4, 20, 100 nmol/L
- PGA, 5-formylTHF, MeFox, THF, and 5,10-methenylTHF at 0.5, 1, 2, 10, 50 nmol/L

The reportable range for this method is LOD – high calibrator. Because the daily calibration range starts at a concentration slightly higher than the LOD, the expanded calibration range LOD – high calibrator is verified periodically. Previous versions of this method (#4001 and #4012) were calibrated from 0 to high calibrator and generated comparable results to method #4013 throughout the calibration range. Furthermore, for the main analyte 5-methylTHF, all results obtained so far have far exceeded the LOD. LODs for serum folate forms analyzed by SCIEX 6500 are:

Analytes	LOD (nmol/L)	LOQ (nmol/L)
5-MethylTHF (MET)	0.13	0.43
Folic acid (PGA)	0.14	0.46
5-FormylTHF (FOT)	0.20	0.66
Tetrahydrofolate (THF)	0.25	0.83
5,10-MethenylTHF (MYT)	0.20	0.66
MeFox (MFO)	0.10	0.33

Serum MET results <7 nmol/L should be reanalyzed for confirmation (critical limit). Serum results >high calibrator (>100 nmol/L for MET and >50 nmol/L for all other folate forms) should be diluted with 0.1% ascorbic acid and reanalyzed.

b) Calibration verification

According to the updated CLIA regulations from 2003 (see also statement of the Joint Commission on Accreditation of Healthcare Organizations), the requirement for calibration verification is met if the test system's calibration procedure includes three or more levels of calibration materials, and includes a low, mid, and high value, and is performed at least once every six months. All these conditions are met with the calibration procedure of this assay, and therefore no additional calibration verification is required by CLIA.

Perform calibration verification whenever any of the following occur:

- All of the reagents used for a test procedure are changed to new lot numbers, unless the laboratory can demonstrate that changing reagent lot numbers does not affect the range used to report patient test results, and control values are not adversely affected by reagent lot number changes.
- There is major preventative maintenance or replacement of critical parts that may influence the test's performance. This includes when the laboratory sends a

test system to the manufacturer for repairs. The laboratory must check the calibration of a repaired test system before resuming patient testing and reporting results.

- Control materials reflect an unusual trend or shift, or are outside of the laboratory’s acceptable limits, and other means of assessing and correcting unacceptable control values fail to identify and correct the problem.
- The laboratory has determined that the test system’s reportable range for patient test results should be checked more frequently.

Since 2008, we participate in the **College of American Pathologists (CAP) LNS Ligand calibration verification/linearity PT challenge** twice a year. However, since the LC-MS/MS assay doesn’t have a peer group, this challenge can only be evaluated for linearity and not for calibration verification.

c) Periodic calibration verification activities

Analysis of international standard reference material (SRM) multiple times per year can be used to satisfy calibration verification requirements if the SRM material covers the reportable range of the assay and has at least 3 levels. Reference materials are stored deep frozen (-50°C to -90°C).

- (1) **NIST SRM 1955 Homocysteine and Folate in Frozen Human Serum** (NIST SRM 1955 Serum folate_Hcy.pdf) has been available since fall 2005, but is no longer available from NIST. This is a 3-level material that approximately covers the reportable range of the assay for serum folate:

Assay (Lab)	Analyte – Type of value	NIST SRM 1955 – Serum Folate (nmol/L)		
		Level 1	Level 2	Level 3
LC-MS/MS (CDC & NIST)	5-MethylTHF – Certified (target ± expanded uncertainty)	4.26 ± 0.25	9.73 ± 0.24	37.1 ± 1.4
LC-MS/MS (CDC & NIST)	PGA – Reference (target ± expanded uncertainty)	0.49 ± 0.17	1.05 ± 0.16	1.07 ± 0.24
LC-MS/MS (CDC)	5-Formylthf – Information (target ± 2 SD)	1.2 ± 0.4	2.3 ± 0.8	3.6 ± 1.3
LC-MS/MS (CDC)	Total folate – Information (target ± 2 SD)	6.0 ± 0.4	13 ± 1	41 ± 2
Microbiologic assay (CDC)	Total folate – Information (target ± 2 SD)	5.6 ± 1.2	14 ± 3	44 ± 11

Analyze each level as unknown against working standards. The yearly mean measured concentration for 5-methylTHF should be within ±15% of the target values. The yearly mean measured concentration for PGA should be within ±30% of the target values (concentrations ≤1 nmol/L).

- (2) **NIST SRM 3949 Folate Vitamers in Frozen Human Serum** ([NIST SRM 3949 Serum folate_Hcy.pdf](#)) has been available since fall 2018. This is a 3-level material that approximately covers the reportable range of the assay for serum folate forms:

Assay (Lab)	Analyte – Type of value	NIST SRM 3949 – Serum Folate		
		Level 1	Level 2	Level 3
LC-MS/MS (CDC & NIST)	5-MethylTHF – Certified (target ± expanded uncertainty)	14.69 ± 2.18	45.71 ± 4.07	29.3 ± 3.77
LC-MS/MS (CDC & NIST)	PGA – Certified (target ± expanded uncertainty)	1.00 ± 0.32	6.75 ± 1.54	4.67 ± 1.01
LC-MS/MS (CDC)	THF – Non-certified (target ± expanded uncertainty)	1.14 ± 0.18	1.53 ± 1.12	1.39 ± 0.99
LC-MS/MS (CDC)	MeFox – Non-certified (target ± expanded uncertainty)	1.58 ± 0.78	1.96 ± 0.40	2.22 ± 0.55
LC-MS/MS (CDC)	5-FormylTHF – Non-certified (target ± expanded uncertainty)	n/a	n/a	7.33 ± 4.01
LC-MS/MS (CDC)	Total folate – Non-certified (target ± expanded uncertainty)	17.0 ± 0.4	56.0 ± 0.8	41.8 ± 0.5

Analyze each level as unknown against working standards. The yearly mean measured concentration for 5-methylTHF, MeFox (levels 1, 2, 3), PGA (levels 2, & 3), and 5-formylTHF (level 3) should be within ±15% of the target values. The yearly mean measured for THF and for concentrations ≤ 1.0 nmol/L should be within ±30% of the target values.

(3) **NIST SRM 1950 Metabolites in Frozen Human Plasma** [NIST SRM 1950 Metabolites in human plasma](#), a 1-level serum-based reference material, has been available since 2011 August.

- Certified concentration (uncertainty) for 5-methylTHF is 26.93 (0.70) nmol/L.
- Reference concentration (uncertainty) for PGA is 3.42 (1.02) nmol/L.
- Information value concentration for total folate by microbiologic assay is 30.6 nmol/L.
- Analyze as unknown against working standards. The yearly mean measured concentration for 5-methylTHF, and PGA should be within 15% of the target value.

(4) **NIBSC Code 03/178 1st International Standard For Vitamin B12 and Serum Folate** ([NIBSC 03-178 Serum folate and vitamin B12.pdf](#)) is another 1-level serum-based reference material.

- Assigned content of 12.1 nmol/L total folate (when reconstituted with 1.0 mL distilled/de-ionized water), made up of 9.75 nmol/L 5-methylTHF (CV 5.5%), 1.59 nmol/L 5-formylTHF (CV 4.2%), and 0.74 nmol/L PGA (CV 31.6%) as measured by LC-MS/MS; 12.1 nmol/L is equivalent to 5.33 ng/mL.
- Analyze reconstituted material as unknown against working standards. The yearly mean measured concentration for 5-methylTHF should be within 15% of the target value. The yearly mean measured concentration for PGA should be within 30% of the target value (concentration ≤1 nmol/L).

d) Verification of expanded calibration range to cover LOD to lowest calibrator


Additional aqueous calibration materials between the LOD and the lowest calibrator from the routine calibration curve (5-methylTHF: 1.0–100 nmol/L; other folate forms: 0.5–50 nmol/L) are tested periodically to verify the reportable range. Triplicate calibrators are prepared for the entire calibration range and carried through the entire sample preparation procedure. QC samples are analyzed together with the additional calibrators:

Analytes	Concentration 1 (nmol/L)	Concentration 2 (nmol/L)	Concentration 3 (nmol/L)
MET	0.125	0.250	0.50
PGA	0.063	0.125	0.25
MFO	0.063	0.125	0.25
FOT	0.063	0.125	0.25
MYT	0.063	0.125	0.25
THF	0.063	0.125	0.25

Two calibration curves, the routine, and the expanded curve (5-methylTHF: 0.125–100 nmol/L; other folate forms: 0.063–50 nmol/L) are evaluated for slope, intercept, R^2 . Slopes are expected to be interchangeable (overlapping 95% CI) and R^2 are expected to be >0.99 . Calculated results for the QC samples using the two calibration curves are expected to agree within $\pm 5\%$ for concentrations >1 nmol/L, within $\pm 15\%$ for concentrations between the LOQ and ≤ 1 nmol/L and within $\pm 30\%$ for concentrations $< LOQ$.

B. Reagents & Standards

1) JA-4013-R&S-01-Sample Preparation Flowchart



2) JA-4013-R&S-02-Reagent Tracking

For details on labeling requirements, see the Division of Laboratory Sciences Safety and Quality Joint Labeling Requirement.

a) Reagent Preparation

Date of reagent preparation is logged on reagent bottle or tags attached to reagent bottles. The reagent bottle also has a label indicating for how long the reagent can be used (e.g., Expires 1 month after preparation). Chemicals used to prepare the reagent are documented on reagent tracking sheets in the laboratory. The following information is included on the reagent tracking sheet:

- (1) Name of chemical, supplier, catalog number, and lot number
- (2) Method and reagent(s) in which it is used
- (3) CHaTS barcode (if available)
- (4) Date received, started, and stopped use
- (5) Any important notes regarding these chemicals

b) Frequently Prepared Reagents:

- (1) 1% Ammonium formate buffer (pH 3.2)
- (2) 0.05% Ammonium formate buffer (pH 3.4)
- (3) Elution solvent (49% water, 40% methanol, 10% acetonitrile, 1% acetic acid with 0.5% ascorbic acid)
- (4) HPLC mobile phase (49.5% water, 40% methanol, 10% acetonitrile, and 0.5% acetic acid)
- (5) HPLC rinse & shutdown solvent (90% methanol and 10% water)
- (6) 0.1% and 0.5% Ascorbic acid (used as powder)

c) Infrequently Prepared Reagents:

- (1) 0.1 M Phosphate buffer pH 7.2 (using potassium phosphate monobasic and dibasic salts)
- (2) 1 M Hydrochloric acid (HCl)
- (3) 0.5 M Hydrochloric acid (HCl)
- (4) 0.1N Sodium hydroxide (NaOH)
- (5) Cysteine (~0.1%)

d) Standards Preparation

- (1) Stock Solutions:

All calibrators and internal standards (ISTDs) specified in the APM originate from single-analyte reagents. Information available on chemicals used to prepare stock solutions is tracked in individual written or electronic records (network shared drive) generated when the stock solution is prepared.

(2) Intermediate Standard Solutions:

Intermediate standard and ISTD solutions are prepared by diluting primary stocks as specified in the APM. Aliquots of working standards and ISTD solutions for single use are kept deep frozen (-50°C to -90°C). Information available on the stock solutions used to prepare the working standard solutions is tracked in individual electronic records (network shared drive) generated when the working standard is prepared.

e) Quality Control Materials

Chemicals and other materials (e.g., human serums, antioxidants) used to prepare quality control materials are tracked in individual written or electronic records (network shared drive) generated when the material is prepared. Quality control materials are aliquoted for single use and kept deep frozen (-50°C to -90°C).

Chemicals may be used to amend (i.e., spike) analyte into a quality control material. In these cases, available information on the chemicals used is tracked in the record (either written or electronic).

f) Other Materials

To facilitate tracking of solvents and other chemicals transferred from the original manufacturer containment into a secondary container (e.g., solvent bottle), expiration provided by the manufacturer, or the lot number is included on the secondary container.

C. Instruments

1) JA-4013-I-01-Instrument Comparison and System Verification

a) Instrument Comparison

When a method is analyzed on multiple instruments, an initial instrument comparison must be conducted to establish analytical comparability. All calibrators, quality control materials, and at least 30 samples that span the measurement range should be analyzed on each system. The same preparation should be analyzed on all systems on the same day. However, this may not be feasible. In this case, it is acceptable to prepare different aliquots of the samples and perform the analysis as close to the same time as possible. The results from each instrument should be plotted against the results from the original instrument. The parameters assessed are correlation (Pearson $r > 0.95$), regression fit (r^2), and slope (m). The same data analysis and documentation procedure as described below under System Verification can be used.

b) System Verification

Following the initial instrument comparison, semi-annual system verification is performed to ensure that the systems are maintaining comparability.

Requirement: According to DLS Policies and Procedures (section 12.2), if a DLS method is run on multiple instruments or at multiple sites, a set of at least five samples spanning the reportable range of the analytes must be run at least once every six months. The Pearson correlation coefficient of the compared results should be greater than 0.95, and if not, appropriate corrective action should be taken. In special situations, the Division Director may give written approval that the methods are sufficiently similar for the intended use of the data.

Procedure: No separate sample preparation is needed. A set of samples spanning approximately the reportable range that were analyzed on the primary instrument, are re-analyzed shortly on the secondary instrument. The time delay should be within the processed sample stability parameters determined during method validation (e.g., sample kept in autosampler or refrigerated/frozen for a certain time).

Data analysis:

- (1) Identify a subset of results ($n \geq 5$) from the two analyses and describe any inclusion or exclusion criteria applied (e.g., only include samples with analyte results $\geq 3x$ the LOD and \leq highest calibrator).
- (2) Determine the Pearson correlation coefficient. Note: Pearson correlation is a parametric test that requires normally distributed data. Most nutritional biomarkers show right-skewed analyte concentration distributions benefiting from a log-transformation to yield data that approximates a normal distribution. While Pearson analysis verifies correlation and not concordance, high concordance is expected for instrument comparisons because most critical variables are the same (measurement technique, sample preparation, operator, calibration, etc.) and only 1 variable changed (instrument).

- (3) Optional: Assess Lin's rho coefficient for concordance. Perform regression and Bland-Altman bias analysis. Assess whether a similar proportion of samples is <LOD on both instruments. Note: Most nutritional biomarkers show non-constant variance (constant CV with increasing concentration); thus, weighted Deming regression and relative (%) Bland-Altman analysis are generally preferred.

Documentation, review, and approval:

- (1) Summarize the results in a spreadsheet that contains the raw data, the data analyses, and the summary information and request review by the supervisor.
- (2) General supervisor reviews the data and, if acceptable, approves the data. Convert summary information tab to a PDF and electronically sign the PDF in the designated field. Save the signed PDF to electronic QA Manual.

2) JA-4013-I-02-Maintenance and Verification of Hamilton Liquid-Handler

a) Daily Maintenance

- (1) Load instrument software method to execute the daily maintenance for deck and waste – Click **[Yes]** & continue. Empty tip waste; click **[OK]** to continue; (avoid clicking on 'Cancel' - aborts the daily maintenance)
- (2) Execute the 1000 µL channel tightness check; ensure 8 teaching needles are installed. Click **[Yes]** & continue. Execute the 1000 µL channel cLLD check; ensure 8 teaching needles are installed. Click **[Yes]** and continue until maintenance method is completed.

b) Weekly Maintenance

- (1) Load instrument software to execute the daily maintenance for deck & waste – Click **[Yes]**
- (2) Using auto-load all racks are moved & put back (can empty deck manually) – Click **[OK]** & continue. Open Front cover, clean the deck – Click **[OK]** & continue
- (3) Open front cover, empty & clean the tip waste – Click **[OK]** & continue
- (4) Check the laser scanner window (clean using lint free cloth or Q tips lightly soaked in 70% ethanol) – Click **[OK]**
- (5) Upon prompting: Execute the 1000 µL channel tightness check; ensure 8 teaching needles are installed – Click **[Yes]** & continue. Execute the 1000 µL channel cLLD check; ensure 8 teaching needles are installed – Click **[Yes]** and continue until maintenance method is completed.

c) Instrument Maintenance (2 times/year) using Volume Field Verification (VFV) Kit

- (1) Load instrument software to execute the daily maintenance for deck and waste – Click **[Yes]**
- (2) Open method & at deck layout find tip carrier and weighing unit locations
- (3) Double click Microlab STAR Method Editor icon → **[File]** → **[Open]** → Liquid Class Developer → Use ***“Ultimate Liquid Class Validation. Method”***
- (4) On worktable check the tip racks for respective tips (300 µL, 1000 µL etc.) are at appropriate deck locations (Carrier-Tip_CAR_480_A00,), restock if needed.
- (5) Following manufacturer’s instructions assemble the weighing unit & its components; plug in the corresponding cables from the weighing unit to the balance control; from the balance control to the computer and power source.
- (6) Place appropriate test liquid in the given vial; put it on the weighing unit. Go to the **“Microlab Star Run”** method icon on desktop; double click and open a workstation window.

- (7) Open **[File]**, locate the target file (same as step 1) → double click to open method → Click **"Start"** → **"Initialize"** will occur automatically → Prompt window appears **"Test Type Choice"** → **Validation** → Select **"test & tip type"** → In tip count window → Look at **"Labware"** for tip size → Ensure that **"First"** and the **"Last"** standard volume tip counts are correct → Click **"OK"** → Run starts at position 1 → Select **"liquid type"** → Select **"liquid class to test"** → Check **"liquid transfer information"**. Enter density volume (μL), liquid density (g/mL); program calculates total times transferred/channel → Select **"default pipetting settings"** Click **"Yes"** → For Anti-Droplet Control (ADC) Select **"No"** → For tip re-use Select **"No"** → Enter aspiration/dispense submerge depths → Select **"OK"** → Select aspiration/dispense "mix volume & number of mix cycles" → Enter appropriate values for these selections & click **"OK"** → Tips are picked up automatically & weights of the set liquid volume measured with each tip during respective aspirating/ dispensing steps.
- (8) Auto-generated data are reviewed on instrument computer; results imported to network (summarized & saved for reference).

D. Data Review

1) JA-4013-DR-01-Computerization and Data System Management

a) Sample Identification

During sample preparation and analysis, samples are identified by their sample ID. The sample ID is a number that is unique to each sample that links the laboratory information to demographic data recorded by sample collectors.

b) Data Collection

Instrument raw data files and respective batch files collected using Analyst software (SCIEX 6500) are stored on mass spectrometer workstation.

c) Data Back-up

Raw batch files and data files from instrument computers are generally transferred via ISLE, KVM-switch, or encrypted thumb drive to NBB shared network drive.

d) Peak Review and Integration

The raw data (calibration curve data, patient data; QC data) are processed into results using the Quantitation Wizard software (Analyst); results are saved and reviewed by team lead for peak integrations, chromatograms, slope/intercept and R^2 . Finally, the processed and reviewed result files are saved in the import data file folder on the shared network drive.

e) Data Import

The result files containing patient data as well as QC data stored on the shared network drive are further exported into the STARLIMS Database for QC and statistical evaluation.

f) STARLIMS Data Review

Level I – Analyst:

- Double click the STARLIMS icon on desktop
- Under 'Run-based Tasks', select 'Pending Runs Assigned to My Labs'
- Choose 'Show Pending Tests' and select '4013 (Folate Vitamers (serum) - Caliper)'
- Click on 'Add' and select the Instrument
- Run# and Equipment ID will be populated
- [0] Run Instrument Macro – select the excel result file to run macro for STARLIMS import
- [1] Upload Instrument File – import the post-macro result file to STARLIMS
- [2] Mark Null Results – click this button which replaces the null with “|” (a pipe tab)
- [3] Evaluate Sample QC – check the sample QC flags according to the defined criteria
- [4] Evaluate Run QC – evaluate bench QC via the DLS SAS Multi-Rule System QC program to determine QC pass/fail

- [5] Set Run QC Statuses – set analytes pass/fail based on SAS out-of-control assessment
- [6] Attach SAS QC file – upload both the SAS input file (.csv) and output file (.pdf)
- Enter run bench QC (SAS) information in Run Comments column
- Click on 'Manage Attachments' and upload the data review checklist for the run
- Click 'Finish Results' located under the test workflow steps and notify Project Lead

Level II – Project Lead:

- Double click the STARLIMS icon on desktop
- Under 'Run-based Tasks', select 'Run Approval'
- Choose 'Show Pending Tests' and select '4013 (Folate Vitamers (serum) - Caliper)'
- Review analyst run sheet, data review checklist, peak integration, and Sample QC evaluation
- Review and confirm Run QC evaluation
- Assess blind QC results - click on 'Blind QC Results Only' tab, 'Assess Blind QC', 'Final Result' and 'Use Default Characterization Sets', 'Proceed to Next Step', 'OK'
- Print blind QC report - click on 'BQC Reports', 'All data displayed in the datagrid', 'A paper-based report from template', 'NBB Blind QC Report v2', 'OK', 'Proceed to the SSRS Report', Save PDF
- Enter bench QC (SAS) and blind QC evaluation status in the Run Comments column
- Set results final - in 'All Results (S)' tab, click on 'Set Final' Wizard, select 'Process all samples displayed in the datagrid' and 'Run the Set Final Wizard' and click 'Proceed'
- Choose Set final criteria - check 'Required Sample QC Passed' and 'Required Run QC Passed'; check 'Pass' and 'Warn' for 'Allowable Results Statuses for Set Final'; choose date range to cover runs that may include the previous analysis of these samples; click 'Proceed'
- Resolve samples with retest results and set final
- Submit sample IDs and repeat instructions to the analyst to schedule the repeats
- Click on 'Manage Attachments' and upload the blind QC report and final data review checklist for the run
- In Run Approval tab, click Release Run and notify QA Officer (level III review)

General Supervisor (Lab Chief):

- Conduct random “spot checks” to verify proper handling of lab results
- Discuss with Team Lead or QA Officer course of action on difficult questions
- Results set reportable, released, and reported by QA officer

2) JA-4013-DR-02-Processing and Reporting a Run

a) Reviewing the run:

- (1) After a batch or run is acquired, the data are reviewed in Analyst software. Chromatograms for all folate forms are checked for retention time, peak shape, peak separation, intensity, and any potential interferences.
- (2) Data are typically processed on a network computer. Batch and result files are copied from the instrument computer onto the network drive for data processing.

b) Quantitation and integration of the data file in Analyst software:

- (1) Under the left bar menu, select “Quantitate” and open “Quantitation Wizard”.
- (2) In the quantitation wizard menu:
 - Under “Available Data Files” select the data file for processing, under “Available Samples” select the number of individual samples for processing, use the arrow key (=>) to transfer the samples under “Selected Samples”, and click “Next”
 - Open the dropdown list next to “Setting to Use” and select “Summary Table”, select “None” under the “Default Query” setting, and click “Next”
 - Select “Choose Existing Method”, open the dropdown list next to “Method” and select “FOL_SeQuant.qmf”, and click “Finish” to process data file.
- (3) An auto integrated results table is generated.
- (4) In the top bar menu, select “File” and “Save As”. Save the quantitation result file (.rbd) using the following format: “FolateSem m ddyy”, where “**Se**” defines “**Serum**” and m ddyy (month/date/year) is the sample preparation date.
- (5) Review the integration of each sample (blanks, calibrators, QCs, and unknowns).
- (6) Check calibration curve accuracy:
 - Accuracy from target value for S1 should be within 30% for all analytes
 - Accuracy from target value for S2-S5 should be within 15% for major analytes
 - Accuracy from target value for S2-S5 should be within 20% for minor analytes
- (7) Review integration parameters (minimum peak height & peak width, RT window, noise percent, smoothing width, peak splitting factor, etc.) and make necessary integration corrections using either manual or auto integration options. **Note:** auto integration is preferred over manual integration.
- (8) Integration parameter changes can be applied to all samples by right click on the integrated peak and select “Update Method”.
- (9) Check QC results for limits and accuracy.
- (10) Save any changes to the results table.

- (11) When data review and integration is completed, the saved Analyst “rdb” file on the network is imported to STARLIMS database:
 - In Analyst results table, right click and select “Table Settings” then select “Serum Results for Import”.
 - In the top bar menu, select “File” then select “Export” and save the results in a designated folder on the network drive.

3) JA-4013-DR-03-STARLIMS Data Review and Criteria

Processed instrument results are imported into STARLIMS database for data review and sample QC evaluation. The quality of individual test result within run is subjected to the sample QC criteria and action is taken based on the flagged result.

- a) For flagged results on calibration curve see details on **Table 1**.
- b) For flagged results on relative retention time see details on **Table 2**.
- c) Evaluate bench QC and blind QC:
 - (1) Bench QC evaluation is performed by Level I (Analyst)
 - If pass → set Run QC “Pass”
 - If fail → set Run QC “Fail”; set code 61; repeat the entire run
 - (2) Blind QC evaluation is performed by Level II (Project Lead)
 - If pass → no action; continue with reporting
 - If fail (of 4 if 1 or 2 QC are [$<$ or $>3SD$]) → repeat 10% samples from the run (including blind QC) to confirm results
- d) For flagged results on patient sample see details on **Table 3**.
 - (1) Pass (codes: 0, 33, 37, or 97) – result passed all criteria and can be reported
 - (2) Warn, Check, and Fail – result needs to be reviewed or repeated to confirm
- e) Check sample reproducibility:

When sample is repeated, multiple results will be shown in “Set Final” Wizard window. The “Retest Delta” is difference between the results and should be within $\pm 15\%$. Set final on the result that meets sample QC criteria.

Table 1. Calibration Curve and Blank Sample QC Criteria				
Analyte	Flagged Result	Result Status	DLS Comment Code	Action
MET, PGA, MFO, FOT, THF & MYT	Numeric Blank Value in Reagent Blank (RB) and Blank (S0)	Warn	None	Check instrument result <ul style="list-style-type: none"> For all analytes & repeat the run if RB contaminated For area in MET, PGA & MFO from 'S0' Area in FOT, THF, & MYT expected from 'S0'/override, & enter code "0"
	Calibrator (S1) diff. to target: >30%	Check	0	Check instrument result; repeat the run
	Calibrator drift (within-run): >15%	Warn	0	Check instrument result; check % diff. from target value
MET, PGA & MFO	Calibrator (S2-S5) diff. to target: >15%	Check	0	Check instrument result; repeat the run
	Calibrator drift (within-run): >15%	Warn	0	Check instrument result; check % diff. from target value
FOT, THF & MYT	Calibrator (S2) diff. to target: >20%	Check	0	Check instrument result; repeat the run
	Calibrator drift (within-run): >20%	Warn	0	Check instrument result; check % diff. from target value
	Calibrator (S3-S5) diff. to target: >15%	Check	0	Check instrument result; repeat the run
	Calibrator drift (within-run): >15%	Warn	0	Check instrument result; check % diff. from target value
MET, MFO, FOT & MYT	$R^2 < 0.98$	Fail	26	Check instrument result; repeat the run
	$R^2 0.98 \geq$ to < 0.99	Check	None	Check calibration result; enter DLS code "0"
PGA, & THF	$R^2 < 0.95$	Fail	26	Check instrument result; repeat the run
	$R^2 0.95 \geq$ to < 0.98	Check	None	Check calibration result; enter DLS code "0"

Table 2. Relative Retention Time Sample QC Criteria				
Analyte	Flagged Result	Result Status	DLS Comment Code	Action
MET & PGA	Result \geq LOQ & Rel. RT: ≤ 0.980 or ≥ 1.02	Fail	26	Check instrument result/repeat the run
	Results $<$ LOQ & Rel. RT: ≤ 0.980 or ≥ 1.02	Check	None	Check instrument result/resolve or repeat
MFO	Results \geq LOQ & Rel. RT: ≤ 0.985 or ≥ 1.02	Fail	26	Check instrument result/repeat the run
	Result $<$ LOQ & Rel. RT: ≤ 0.985 or ≥ 1.02	Check	None	Check instrument result/resolve or repeat
FOT	Results \geq LOQ & Rel. RT: ≤ 0.985 or ≥ 1.018	Fail	26	Check instrument result/repeat the run
	Results $<$ LOQ & Rel. RT: ≤ 0.985 or ≥ 1.018	Check	None	Check instrument result/resolve or repeat
THF	Results \geq LOQ & Rel. RT: ≤ 0.963 or ≥ 1.06	Fail	26	Check instrument result/repeat the run
	Results $<$ LOQ & Rel. RT: ≤ 0.963 or ≥ 1.06	Check	None	Check instrument result/resolve or repeat
MYT	Results \geq LOQ & Rel. RT: ≤ 0.982 or ≥ 1.02	Fail	26	Check instrument result/repeat the run
	Results $<$ LOQ & Rel. RT: ≤ 0.982 or ≥ 1.02	Check	None	Check instrument result/resolve or repeat

Table 3. Patient Sample QC Criteria


Analytes	Flagged Result	Result Status	DLS Comment Code	Action
MET	Null, 0 or No Peak	Fail	26	Check instrument result, repeat/confirm
	<LOD	Check	37	Check instrument result, repeat/confirm
	<LOQ	Incomplete (Result 1) / Warn (Result 2)	33	Check instrument result, repeat/confirm
	≥LOD <7 nmol/L	Incomplete (Result 1) / Warn (Result 2)	33	Check instrument result, repeat/confirm
	≥7 nmol/L & ≤100 nmol/L	Pass	0	None
	>100 nmol/L	Fail	26	Check instrument result, dilute/repeat
	≥7 nmol/L & ≤100 nmol/L; Retest # ≥1 & dilution >1.07	Pass	97	None
PGA & MFO	Null, 0 or No Peak	Check	37	Check instrument result, repeat/confirm
	<LOD	Check	37	Check instrument result, repeat/confirm
PGA, MFO, FOT, THF, & MYT	≥LOD & ≤ 50 nmol/L	Pass	0	None
	>50 nmol/L	Fail	26	Check instrument result, dilute/repeat
	≥LOD & ≤ 50 nmol/L Retest # ≥1 & dilution >1.07	Pass	97	None
FOT, THF & MYT	Null, 0 or No Peak	Pass	37	None
	<LOD	Pass	37	None
	Signal/Noise (S/N) <3, & Raw result >LOD	Warn	0	Check instrument result & dilution factor (if >1.07); report undiluted result
MET, PGA, MFO, FOT, THF, & MYT	Retest Delta: ≥ 0.85 2 nd result within ±15% of 1 st result	Pass	0/33/37	None
	2 nd result not within ±15% of 1 st result	Pass	0/33/37	Check instrument result; repeat 3rd time
MET, PGA, MFO, FOT, THF, & MYT	Manual Integration	Check	None	Check instrument result integration; update code "0"

4) JA-4013-DR-04-Out-of-Control Corrective Action

- a) QC performance is evaluated by SAS - run comments updated appropriately (pass/fail)
- b) SAS QC failed analytes are investigated, and appropriate corrective action measures as indicated are applied.
- c) Verify calibrator response, area ratios, slope/Intercept, R^2 (>0.989) of the analyte in question and other analytes.
- d) Verify QC status of other analytes in the run. Verify if there were problems with sample preparation and extraction; check chromatograms, peak shape, integration, relative retention time, MS instrument response, or interferences; apply appropriate code.
- e) Verify if the issue is caused due to an error in one of the calibrators and eliminating such calibrator from the curve can correct the issue; discuss and seek approval from supervisor.
- f) Verify if QC in question passed for other analytes; if appropriate, eliminate failed QC and re-run SAS; update comments (pass/fail); discuss and seek approval from supervisor.
- g) If the QC failure is true (potential statistical issue); apply code 61 and repeat samples in next run.
- h) Enter appropriate comments in STARLIMS database (User fields 1 & 2).

5) JA-4013-DR-05-STARLIMS Data Review Flowchart

The following instructions reflect general steps for data review, approve, sample reconciliation, and report in STARLIMS.



Appendix C: Information on Absorption Maxima, Absorption Coefficients, and Formulas to Calculate the Folate Concentration

Compound	Diluent for primary stock solution I	Absorption maximum (nm)	Molar extinction coefficient (L/mol*cm)	Reference	Molecular weight (g/mol)	
					Unlabeled	¹³ C5-Labeled
5-MethylTHF	20 mM phosphate buffer with 0.1% cysteine (pH 7.2)	290	31,700	1	459.46	464.46
5-FormylTHF	20 mM phosphate buffer with 0.1% cysteine (pH 7.2)	285	37,200	2	473.44	478.44
THF	20 mM phosphate buffer with 0.1% cysteine (pH 7.2)	298	25,000	3	445.43	450.43
5,10-MethenylTHF	1 M HCl (pH 0)	288 & 348	13,500 & 26,500	3–5	455.45	460.45
Folic acid (PGA)	20 mM phosphate buffer (pH 7.2)	282 & 346	27,600 & 7,200	3	441.40	446.40
MeFox	0.1 N NaOH (pH 13)	280	19,365	6	473.40	482.40

a) References for UV-spectrophotometric determination of folate concentration:

- (1) Gupta VS, Huennekens M. Arch. Biochem. Biophys. 1967;120:712.
- (2) Uyeda K, Rabinowitz JC. J. Biol. Chem. 1965;240:1701.
- (3) Rabinowitz JC. In: Boyer PD, Lardy H, and Myrbaeck K (eds.). The Enzymes, vol. 2. 2nd ed., Academic Press, New York 1960. p. 185.
- (4) Rabinowitz JC. In: Colowick SP and Kaplan NO (eds.). Methods in Enzymology, vol. 6, Academic Press, New York 1963. p. 814.
- (5) Huennekens FM, Ho PPK, Scrimgeour KG. In: Colowick SP and Kaplan NO (eds.). Methods in Enzymology, vol. 6, Academic Press, New York 1963, p. 806.
- (6) Personal communication with Jean-Pierre Knapp at Merck Cie, March 2012.
- (7) Blakley R.L. The biochemistry of folic acid and related pteridines. In: Neuberger A, Tatum EL, eds. Frontiers of biology. Amsterdam: North Holland Publishing Company, 1969:91–5.

b) Formulas to calculate the concentration of folate stock solution based on molar absorptivity:

$$\text{Conc. (ppm or } \mu\text{g / mL)} = [\text{Absorbance (cm}^{-1}) \times \text{dilution} \times 1000 \times \text{MW (g mol}^{-1})] / \epsilon \text{ max (L mol}^{-1} \text{ cm}^{-1})$$

$$\text{Conc. (}\mu\text{mol / L)} = [\text{Absorbance (cm}^{-1}) \times \text{dilution} \times 1000 \times 1000] / \epsilon \text{ max (L mol}^{-1} \text{ cm}^{-1})$$

with $\epsilon \text{ max}$ = molar extinction coefficient and conc = concentration

Example calculation for folic acid (PGA):

Absorbance	Dilution	$\epsilon \text{ max}$	MW
0.661	10	27600	441.4

$$\text{Concentration (ppm or } \mu\text{g/mL)} = 0.661 \times 10 \times 1000 \times 441.4 / 27600 = 105.7$$

$$\text{Concentration (}\mu\text{mol/L)} = 0.661 \times 10 \times 1000 \times 1000 / 27600 = 239.5$$

c) Conversion factors from conventional (ng/mL) to SI units (nmol/L) for different folate forms:

	5-MethylTHF	5-FormylTHF	PGA	THF	5,10-MethenyTHF	MeFox
Unlabeled form	2.176	2.112	2.266	2.245	2.196	2.112
¹³ C5-labeled form	2.153	2.09	2.24	2.22	2.172	2.09

Appendix D: Approximate QC Pool Target Concentrations (nmol/L) for the Various Folate Vitamers and the 5-methylTHF Oxidation Product (MeFox)

QC level	5-MethylTHF	5-FormylTHF	PGA	THF	5,10-MethenylTHF	MeFox	TFOL
Low	10*	1	1	1	1	1	15*
Medium	20	2.5	5	2.5	2.5	2.5	30
High	50	5	10	5	5	5	70

*As low as possible

Appendix E: Ruggedness Testing

A previous version of this method (method 4001) has undergone a series of in-house ruggedness testing experiments designed to assess how much method accuracy changes when certain experimental parameters are varied. A total of six parameters judged to most likely affect the accuracy of the method have been identified and tested. Testing generally consisted of performing replicate measurements on a test specimen with the selected parameter set at a value substantially lower and higher than that specified in this method while holding all other experimental variables constant. Please refer to Chapter 20 of the 2017 DLS Policies and Procedures Manual for further information on ruggedness testing. Because the method basics did not change by going to method 4013, there is no need to repeat the ruggedness testing with this method.

Folate is an important nutrient involved in one carbon cellular metabolism. Serum folates are measured to determine folate status. We use an isotope-dilution tandem mass spectrometric method in multiple reaction-monitoring mode (MRM) coupled with liquid chromatography (LC-MS/MS) for quantitative measurements of 5 different folate species.

A. Variations in sample preparation

- (1) Principle: The buffers used for sample preparation, solid-phase extraction and analyte elution use formic acid, ammonium hydroxide (to adjust pH), ascorbic acid (as antioxidant), and acetic acid (acid modifier). The changes in buffer pH, concentration of formic acid, ascorbic acid, and acetic acid are critical for analyte and/or sample matrix recovery during sample preparation and solid-phase extraction and will affect analyte sensitivity and potentially affect the results.
- (2) Proposal: To vary and test the sample preparation and solid-phase extraction conditions.
 - pH of ammonium formate buffer (solvent #1)
 - Formic acid concentration in ammonium formate buffer (solvent #1)
 - Ascorbic acid concentration in ammonium formate buffer (solvent #1)
 - Ammonium formate concentration in intermediate wash step during which matrix compounds are eluted but analytes are retained (SPE wash buffer)
 - Ascorbic acid concentration in SPE elution buffer (solvent # 3)
 - Acetic acid concentration in SPE elution buffer (solvent # 3)
- (3) Findings:
 - Varying the pH of the ammonium formate sample preparation buffer does not appear to affect folate species results in serum samples.
 - Varying the formic acid concentration in the ammonium formate sample preparation buffer does not appear to affect folate species results in serum samples.
 - Varying the ascorbic acid concentration in the ammonium formate sample preparation buffer does not appear to affect folate species results in serum samples.
 - Varying the ammonium formate concentration in the SPE wash buffer does not appear to affect the folate species results in serum samples.
 - Varying the ascorbic acid concentration in the SPE elution buffer does not appear to affect folate species results in serum samples.

- Varying the acetic acid concentration in the SPE elution buffer does not appear to affect folate species results in serum samples.

B. Ruggedness testing results for serum folate vitamers by LC-MS/MS

Factor	Method specifies	Results ^a (nmol/L)	Lower level	Results ^a (nmol/L)	Higher level	Results ^a (nmol/L)
1. pH of ammonium formate buffer (solvent #1)	3.2	MET: 22.3 FOT: 1.8 PGA: 1.8	3.0	MET: 22.5 FOT: 2.1 PGA: 1.8	3.4	MET: 22.2 FOT: 1.8 PGA: 1.9
2. Formic acid concentration in ammonium formate buffer (solvent #1)	1%	MET: 22.7 FOT: 2.1 PGA: 2.0	0.8%	MET: 21.6 FOT: 2.1 PGA: 1.9	1.2%	MET: 21.5 FOT: 2.1 PGA: 1.9
3. Ascorbic acid in ammonium formate buffer (solvent #1)	0.5%	MET: 22.4 FOT: 1.6 PGA: 1.9	0.3%	MET: 22.3 FOT: 1.6 PGA: 2.1	0.7%	MET: 23.1 FOT: 1.6 PGA: 2.2
4. Ammonium formate concentration in SPE wash buffer	0.05%	MET: 22.2 FOT: 2.1 PGA: 2.0	0.04%	MET: 22.0 FOT: 2.3 PGA: 2.0	0.06%	MET: 21.7 FOT: 2.0 PGA: 2.0
5. Ascorbic acid concentration in SPE elution buffer	0.5%	MET: 22.3 FOT: 1.3 PGA: 1.8	0.3%	MET: 22.5 FOT: 1.4 PGA: 1.9	0.7%	MET: 21.7 FOT: 1.5 PGA: 1.7
6. Acetic acid concentration in SPE elution buffer	1%	MET: 22.0 FOT: 2.0 PGA: 1.9	0.8%	MET: 21.5 FOT: 2.0 PGA: 1.9	1.2%	MET: 21.6 FOT: 1.8 PGA: 2.0

^a Results are shown for the medium QC sample.

Abbreviations for folate vitamers in serum: MET (5-methylTHF); PGA (Folic acid); FOT (5-formylTHF), THF (tetrahydrofolate); MYT (5,10-methenylTHF).

Appendix F: Typical MRM Method Parameters (Analysis in Positive Ion Mode)

Analyte (Transition)	Tr (min)	*DP (V)	CE (V)	CXP (V)	EP (V)
5-MethylTHF (m/z 460.2 \rightarrow m/z 313.2)	2.37	90	28	13	10
¹³ C ₅ -5-MethylTHF (m/z 465.2 \rightarrow m/z 313.2)	2.37	90	28	13	10
PGA (m/z 442.2 \rightarrow m/z 295.2)	3.17	70	18	13	10
¹³ C ₅ -PGA (m/z 447.2 \rightarrow m/z 295.2)	3.17	70	19	13	10
MeFox (m/z 474.4 \rightarrow m/z 284.2)	3.22	90	48	13	10
¹³ C ₅ -MeFox (m/z 479.4 \rightarrow m/z 284.4)	3.22	90	48	13	10
5-FormylTHF (m/z 474.4 \rightarrow m/z 299.2)	3.16	80	43	13	10
¹³ C ₅ -5-FormylTHF (m/z 479 \rightarrow m/z 299.2)	3.16	80	43	13	10
THF (m/z 446.2 \rightarrow m/z 299.2)	2.33	85	25	13	10
¹³ C ₅ THF (m/z 451.1 \rightarrow m/z 299.1)	2.33	85	25	13	10
5,10-MethenylTHF (m/z 456.1 \rightarrow m/z 412.2)	2.30	155	41	13	10
¹³ C ₅ -5,10-MethenylTHF (m/z 461.1 \rightarrow m/z 416.2)	2.30	155	42	13	10

Tr (min) retention time; DP, declustering potential; CE, collision energy; CXP, collision cell exit potential; EP, Entrance potential; The general instrument parameters used for LC/MS/MS detection and quantitation of all four analytes in multiple reaction mode (MRM) were as follows: resolution Q1 and Q3: unit; dwell time: 110 msec; ion spray voltage: 5500 V; source temperature: 450 °C; curtain gas: 35 psi; gas 1: 55 psi; gas 2: 60 psi; CAD gas: 8 psi (medium)

*DP voltages are subject to change with instrument sensitivity. These voltages can be optimized in Q1 & Q3 based on analyte response.