



Laboratory Procedure Manual

Analytes: **Aromatic Amines**

Matrix: **Urine**

Method: **GC Tandem Mass Spectrometry**

Method No: **2020**

Revised:

As performed by:

**Tobacco and Volatile Branch
Division of Laboratory Sciences
National Center for Environmental Health**

Contact:

**Dr. Tiffany Seyler
Phone: 770-488-4527
Fax: 770-488-0181
Email: tvh2@cdc.gov**

**Dr. James Pirkle, Director
Division of Laboratory Sciences**

Important Information for Users

The Centers for Disease Control and Prevention (CDC) periodically refines these laboratory methods. It is the responsibility of the user to contact the person listed on the title page of each write-up before using the analytical method to find out whether any changes have been made and what revisions, if any, have been incorporate

Division of Laboratory Sciences

Laboratory Protocol



Analytes: *o*-toluidine, 2,6-dimethylamine, *o*-anisidine, 1-aminonaphthalene, 2-aminonaphthalene, 4-aminobiphenyl

Matrix: Urine

Method: Gas Chromatography Tandem Mass Spectrometry (GC-MS/MS)

Method code: 2020.02

Branch: Tobacco and Volatile Branch

Prepared By: Tiffany Seyler, PhD, MPH _____
Author's name Signature Date

Supervisor: Lanqing Wang, PhD _____
Supervisor's name Signature Date

Branch Chief: Benjamin Blount, PhD _____
Branch Chief's name Signature Date

Date current version of method first used in lab: _____
Date

Director's Signature Block:

Reviewed: _____
Signature Date

Procedure Change Log

Procedure: Aromatic amines in Urine
DLS Method Code: 2020.02

Date	Changes Made	By	Rev'd By (Initials)	Date Rev'd
08/13/2016	Section 2.h: Disposal of wastes-10% bleach solution was replaced by Lysol® I.C Quaternary Disinfectant Cleaner solution Section 6.a: Handling information for solvents and reagents- Acetic Anhydride is no longer used.	sm		08/13/2017
06/20/2016	Section 8.d.1: Analyzing and storing GC/MS-MS data- Indigo ASCENT™ was used to process chromatographic data and quantification in place of Agilent MassHunter™ software.	sm		06/20/2016
08/22/2016	Section 8.d.1: Analyzing and storing GC/MS-MS data- Repeat manager database was added	ra		08/22/2017
08/31/2016	Appendix A.A.1: Native standard • Table A1: Stock B was added (OTOL only). • Table A2: Standard 13-17 were added (OTOL only) to extend OTOL dynamic range.	ts		08/31/2016
04/06/2017	Section 6.g: Instrumentation • GC column 1 changed from DB-17MS to DB-FFAP. • Post column is added for backflush (listed as column 1 in Instrument Operation Parameter). • GC temperature gradient- final temperature changed from 280°C to 240°C due to new GC column 1. • GC temperature gradient- post run hold for 5 minutes at 240°C during backflush.	sm		04/06/2017
12/16/2016	Appendix A.B.1: Quality control materials • Table B1: QCLow_2016 and QCHigh_2016 were added	ts		12/16/2016
05/18/2017	Method cross validation for column DB-17MS and DB-FFAP and all changed parameters that were made to GC method (for sample in-reinjections).	sm		05/18/2017
11/11/2017	Completed method re-validation Appendix B: all method validation data were update	ts		11/16/2017
11/29/2017	Appendix C: Method performance documentations are added.	ts		12/29/2017
02/16/2018	Appendix B: updated LoDs	ts		02/20/2018
04/18/2018	Appendix B: updated LoDs, table B4	ts		06/20/2018

Public Release Data Set Information

This document details the Lab Protocol for testing the items listed in the following table:

DATA FILE NAME	VARIABLE NAME	ANALYTE DESCRIPTION
AA_H & AAS_H	URX1NP	1-Aminonaphthalene, urine (pg/mL)
	URX2NP	2-Aminonaphthalene, urine (pg/mL)
	URX4BP	4-Aminobiphenyl, urine (pg/mL)
	URXANS	o-Anisidine, urine (pg/mL)
	URXDMN	2,6-Dimethylaniline, urine (pg/mL)
	URXOTD	o-Toluidine, urine (pg/mL)

1. CLINICAL RELEVANCE AND SUMMARY OF TEST PRINCIPLE

A. *Analytes*

o-toluidine (OTOL), C₇H₉N; Mol Wt 107.15
2,6-dimethylaniline (26DM), C₈H₁₁N, Mol Wt 121.18
o-anisidine (OANS), C₇H₉NO, Mol Wt 123.15
1-aminonaphthalene (1-AMN), C₁₀H₉N, Mol Wt 143.19
2-aminonaphthalene (2-AMN), C₁₀H₉N; Mol Wt 143.19
4-Aminobiphenyl (4-ABP), C₁₂H₁₁N; Mol Wt 169.22

B. *Clinical Relevance*

o-toluidine (OTOL), o-anisidine (OANS), 2-6-dimethylaniline (26DM), 2-aminonaphthalene (2-AMN), and 4-aminobiphenyl (4-ABP) are classified as carcinogens or possible carcinogens (1-6). They are present in mainstream and sidestream tobacco smoke, with the latter containing up to thirty times as much 4-ABP as mainstream smoke (7-8). Aromatic amines are metabolized mainly in the liver where they are N-glucuronidated, N-acetylated, or form N-hydroxyarylamines via oxidation. Aromatic amines are believed to exert their carcinogenic effect by reaction of the n-hydroxylamine metabolite with DNA in the target organ (e.g. bladder) (7-8). Free or conjugated forms (N-acetylated or glucuronidated) are excreted in urine directly from the bladder. Consequently, urinary concentrations of aromatic amines are effective surrogate measures of the carcinogenic metabolite at the target tissues (bladder, liver, kidney, pancreas, spleen, thyroid, etc.).

C. *Assay Principle*

OTOL, OANS, 26DM, 1-AMN, 2-AMN, and 4-ABP are quantified by an isotope-dilution gas chromatographic, tandem mass spectrometric method (ID GC-MS/MS). Urine samples are collected and stored at approximately -70±10°C. ¹³C and ²H internal standards are added, and the samples are hydrolyzed, cleaned up, and extracted on support liquid extraction (SLE) cartridges. The analytes are then derivatized to form pentafluoropropionamides, and analyzed by GC/MS/MS, using multiple reaction monitoring (MRM). The analyte concentrations are derived from the ratio of the integrated peaks of native to labeled ions by comparison to a standard curve.

D. *Special Precaution*

Because of the sensitive nature of these assays and the off-gassing of the target analytes by smokers, all analysts performing this method must be nonsmokers, and measurements must be performed in a smoke-free environment.

2. SAFETY PRECAUTIONS

A. *Reagent Toxicity or Carcinogenicity*

Many aromatic amines (AAs) are carcinogenic. Care should be taken to avoid inhalation or dermal exposure. Use a chemical fume hood when working with AAs. Appropriate use of personal protection including lab coats, gloves, and safety glasses are required when preparing or handling neat materials, standard solutions, extraction solutions, or collected samples.

B. *Radioactive Hazards*

None.

C. *Microbiological Hazards*

This assay involves human urine samples. Universal precautions must be followed. Analysts working directly with the specimens must use proper technique and avoid any direct contact with the samples. Lab coats, gloves, and safety glasses (as required) should be worn while handling the specimens.

D. *Mechanical Hazards*

There are no unusual mechanical hazards associated with this method. Analysts should know and follow the manufacturer's recommendations concerning the safe handling of instruments and other equipment. High voltages are found within certain areas of the mass spectrometer and care must be taken when working in those areas.

E. *Protective Equipment*

Standard chemical laboratory personal safety equipment is required including lab coats, gloves, and safety glasses.

F. *Training*

Training for sample preparation, sample handling, and equipment operation is required.

G. *Personal Hygiene*

Follow standard precaution and comply with all established laboratory safety practices. Care should be taken when handling chemicals to avoid inhalation or dermal exposure. Lab coats, gloves, and safety glasses should be worn at all times handling standards or samples.

H. ***Disposal of Wastes***

Dispose all waste materials in compliance with laboratory, federal, state, and local regulations. Solvents and reagents should always be disposed of in an appropriate container that has been clearly marked for waste products and temporarily stored in a chemical fume hood. Place disposable laboratory supplies such as vials, pipette tips, syringe, etc. that directly contact AAs or samples in a biohazard autoclave bag or similar approved storage container. Unshielded needles, glass pipets and disposable syringes with attached needles must be placed in a sharps container and autoclaved when the container is full. Wipe down all surfaces potentially exposed to biological samples with Lysol® I.C Quaternary Disinfectant Cleaner after each sample preparation. Non-disposable glassware or other equipment that comes into contact with biological samples must be rinsed with bleach before reuse.

3. **COMPUTERIZATION; DATA-SYSTEM MANAGEMENT**

A. ***Software and knowledge requirements***

This method has been validated using an automated sample preparation system - followed by gas chromatography tandem mass spectrometry. The Agilent GC Triple Quad 7000C is controlled by MassHunter™ software. Indigo ASCENT™ is utilized in chromatographic and MS quantitation analysis. Proficiency is required for the analytical software package of automation system, GC, and mass spectrometer used in the analysis. Further, statistical analysis of results requires proficiency in a standard statistical analysis software package. The Statistical Analysis System (SAS Institute, Cary, NC) is one such.

B. ***Sample information***

Typically samples are analyzed in runs of 32-64 samples -including one water blank, one quality control (QC) low, one QC high and unknowns. Each run is identified as "AYYMMDDxxxx" (YearYearMonthMonthDateDateAnalyst'userID, e.g. A160211leo9). Each sequence file contains such information as Run ID, sample ID, sample file name, date of analysis, analyst, and sample volume. The GC/MS/MS relative response data are transferred electronically into the database for each sample and associated calibrators, QCs, and blanks.

C. ***Data maintenance***

Check data entered into the database for transcription or transmission errors. Routinely back-up the database on a weekly basis or as needed (software update, etc.).

D. ***Information security***

Information security is provided at multiple levels. The data systems used in this work are accessed via computers that require individual login and passwords and that default to locked conditions during extended periods of nonuse. Sensitive portions of custom software are protected with additional password requirements. In addition, on the Chamblee campus of CDC, all systems and equipment have restricted access with security personnel approving all entry. Furthermore, the individual laboratory building has multiple levels of controlled access including the requirement for key cards to access the building itself, and also the individual floors where the equipment is located. Confidentiality of the results is protected by use of blind coded ID numbers (no clinical specimen are ever labeled with personal identifiers).

4. **COLLECTION, STORAGE, AND HANDLING PROCEDURES;
CRITERIA FOR SPECIMEN REJECTION**

A. ***Special requirements***

There are no special requirements such as fasting or adherence to special diets for this assay.

B. ***Sample collection***

The specimen for these analyses is human urine. Based on the relatively short physiological half-lives of these analytes, urine samples integrating over longer time periods are preferred over spot urine samples.

C. ***Sample handling***

Specimens should be frozen prior to shipment, must be sent and received frozen where they will be stored at $-70\pm 10^{\circ}\text{C}$ until analysis. All samples are vortex thoroughly prior to preparation (more details in sample preparation section).

D. ***Sample quantity***

The sample size is 2.00 ml of urine. This sample volume is required to quantify the analyte concentrations listed for the Limits of Detection in Section 9a. However, smokers have much higher concentrations of these analytes, so detectable concentrations can likely be measured in a reduced volume sample collected from a smoker (1 ml).

E. ***Sample rejection criteria***

Criteria for defining a sample as unacceptable include (1) use of improper collection materials or techniques leading to possible background

contamination; and (2) sample volume below 1 ml.

F. Long-term stability

The OTOL, OANS, 26DM, 1-AMN, 2-AMN, and 4-ABP are stable in urine samples stored in glass or cryogenic polypropylene tubes at low temperatures, $-70\pm 10^{\circ}\text{C}$ for at least 2-10 years.

5. PROCEDURES FOR MICROSCOPIC EXAMINATIONS; CRITERIA FOR REJECTION OF INADEQUATELY PREPARED SLIDES

Not applicable for this procedure.

6. PREPARATION OF REAGENTS, CALIBRATORS (STANDARDS), CONTROLS, AND ALL OTHER MATERIALS; EQUIPMENT AND INSTRUMENTATION

Note: Class A glassware such as pipets and volumetric flasks are used unless otherwise stated.

A. Handling information for solvents and reagents

- (1) o-toluidine, 2,6-dimethylaniline, o-anisidine, 1-aminonaphthalene, 2-aminonaphthalene, 4-aminobiphenyl. These chemicals are known or suspected carcinogens, and suitable protective clothing, gloves and eye/face protection must be used. It can be dangerous if inhaled, swallowed or absorbed through the skin, and should only be used in a chemical safety hood. If contact occurs, flush area immediately with copious amounts of water. If inhaled, remove to fresh air or consult a physician.
- (2) Sodium hydroxide. This chemical is a caustic base that is corrosive to all tissues. It generates considerable heat when mixed with water or an acid. It is nonflammable but would be harmful if inhaled or swallowed. Protective clothing and safety glasses must be worn while working with this reagent.
- (3) Hexane. This is a flammable solvent. The vapor or mist is irritating to the skin, eyes, mucous membranes and upper respiratory tract. Protective clothing and safety glasses should be used. Use chemical fume hood when working with this solvent.
- (4) Toluene. This is a flammable liquid and may form explosive vapors. Toluene vapor is heavier than air and may travel some distance to an ignition source. Toluene forms an irritating vapor. As a liquid it is a skin

irritant, and may be absorbed through the skin. Large volumes of toluene should be handled with gloves in a chemical fume hood.

- (5) 1.07 M Trimethylamine. Aldrich 98% #T7, 276-1 FW 95.6, stored in a desiccator. Weigh 1.0 gm trimethylamine-HCl at 98% (0.98 gm TMA-HCl) and dissolve in 2.0 ml B&J water, neutralized with 10 M NaOH, approx. 1-3 μ l, and extract with 5.0 ml hexane (Burdick & Jackson #216-4 high purity solvent). Transfer 3.0 ml of the hexane extract to a clean 10 ml vial and add 3.0 ml additional hexane. Final concentration is 1.0 M, reagent is stored in refrigerator up to 6 months.
- (6) Methylene chloride. This solvent is chemically stable and relatively unreactive. It is not flammable, but the vapor can be irritating to the eyes, nose and throat. Skin or eye contact with the liquid should be avoided. Flush exposed tissue copiously with water if any contact should occur. Evaporation of significant volumes of this solvent must be performed in the Savant evaporator, or in a chemical fume hood.
- (7) Ethanol. This is a flammable solvent whose vapors can form ignitable and explosive mixtures with air at normal room temperatures. Exposure to liquid, vapors, fumes, or mist may cause irritation, redness, pain, and possible corneal damage. Prolonged contact may lead to defatting, redness, pain, itching, inflammation, and possible secondary infections. Use appropriate gloves, lab coat, safety glass, and recommended PPE (personal protective equipment) when handling this reagent.

B. **Stock reagent preparation**

Native standards are weighed and diluted in ethanol to make Stock A. Chemical purity, lot and stock concentrations would be different each time a new stock is made. An example of Stock A is listed below.

Date 1mg/mL Stock Made	Analyte	Vendor	Chemical Purity (%)	Purity Test	Lot	Analyte Weight (mg)	Stock A Conc.* (mg/mL)
03/31/15	OTOL	Fluka	99.5	GC	1421704V	10.720	1.0666
07/28/14	26DM	Aldrich	98.5	GC	STBB5452V	7.840	0.7722
07/28/14	OANS	Fluka	99.7	GC	SZBB343XV	9.630	0.9601
03/31/15	1AMN	Fluka	99.9	GC	SZB8268XV	10.100	1.0090
03/31/15	2AMN	Sigma	98.0	TLC	SLBD0851V	11.090	1.0868
07/28/14	4ABP	Aldrich	99.0	TLC	128K1386V	10.830	1.0722

*Made in 10-mL volumetric flask, Stock A in ethanol, further dilution in hexane

Date 1mg/mL Made	Analyte	Vendor	Chemical Purity (%)	Isotopic Purity (%)	Purity Test	Lot	Actual Weight (mg)	Stock DA Final Con.* 10 ml (mg/mL)
01/29/15	OTOL- ¹³ C ₆	MI	98.0	98.0	Mass Spec	713	2.360	0.2267
07/28/14	PTOL-D ₉	MI	98.6	99.0	GC	18	9.230	0.9010
07/28/14	MTOL-D ₉	MI	98.4	99.0	GC	310	9.070	0.8836
07/28/14	26DM-D ₆	TRC	98.0	95.3	Mass Spec	3-MIC-149-1	9.000	0.8405
07/28/14	OANS-D ₇	MI	99.8	99.4	HPLC & NMR	533	10.180	1.0099
01/29/15	2ABP-D ₉	MI	99.0	99.6	HPLC	932	10.430	1.0284
07/28/14	1AMN-D ₉	MI	98.1	99.1	HPLC & NMR	145	10.790	1.0490
12/21/15	2AMN-D ₇	CDN	98.0	95.0	NMR & Mass Spec	Z-330	3.170	0.3170
07/28/14	3ABP-D ₉	TRC	98.0	98.7	Mass Spec	12-SDJ-187-1	8.700	0.8415
07/28/14	4ABP-D ₉	CI	99.0	99.6	HPLC & NMR	11-9733	9.570	0.9436

*Made in 10-mL volumetric flask, Stock A in ethanol, further dilution in hexane, concentration is adjusted by purity as described in Appendix A.

Toluene. Burdick & Jackson, ACS/HPLC Grade GC 99.96% purity.

Hexane. Burdick & Jackson, High purity GC 99.63% purity.

Ethanol. Pharmco, ACS Grade GC 99.9% purity.

C. **Calibrators**

Prepare one complete set of calibration standards (50-100 mL) at one time for use over a period of several years. Prepare labeled ¹³C or ²H (D) internal standards at the same time or as needed and then aliquot, seal, and store these solutions at approximately -20°C±4°C. The internal standard spiking (ITSD) solution for unknown samples was prepared in ethanol, dispensed in 4.0 ml aliquots; and stored at approximately -20°C±4°C. Derivative standards and ISTD spiking solutions were prepared as described below. Analyze at least three calibration curves to confirm acceptable linearity (R² > 0.98). A standard curve is run with every analytical run/batch. A total of 14 standards were prepared ranging from 0 to 200 pg/μl (details in Appendix A) for all analytes except for OTOL, which has 17 standards and a higher range (1,300 pg/μl)

Derivatization of Hexane Standards. 5 mL of each standard is aliquoted via 5-mL volumetric flask. A known amount of internal standard of labelled-analytes is added, followed by TMA and PFAA (pentafluoropropionic acid anhydride, Pierce #65193M MW 310.0). The standards were capped, vortexed, held at room temperature for 30-40 minutes, transferred to silanized tubes, and dried completely in a Savant Speedvac System. Reconstitution to 5 mL toluene is done by 2x2 mL rinse of the silanized tubes, transferred to 5-mL volumetric flasks, and bring to final volume mark. Final standard solutions are aliquoted to 5x1 mL into high recovery amber vials, capped, labelled, and stored at -20°C±4°C. Larger volume, up to 25 ml, can be prepared with volumetric flask.

Details and actual concentrations for each analyte of the currently standard curve are in Appendix A.

Concentrated Internal Standard Spiking Solution:

Internal standard (ISTD) spiking solution for unknown sample preparation was made by diluting the DA and DC stock solutions in ethanol. The ITSD spiking solution was dispensed into 4.0 mL aliquots and stored frozen at $-20\pm 4^{\circ}\text{C}$. Appendix A provided the details of how the current ISTD spiking solution was prepared.

D. Controls

- (1) **Quality control (QC) materials.** There are two quality control pools for the urinary aromatic amines assay: low and high. The QC pools were prepared in-house from urine collected from non-smokers, filtered with 0.2 μm filters and spiked with standard stock solution containing 6 analytes: *o*-toluidine, 2,6-dimethylaniline, *o*-anisidine, 1-aminonaphthalene, 2-aminonaphthalene, and 4-aminobiphenyl. The pool concentrations were made at approx. 100-125 pg/mL and 400-500 pg/mL, respectively. Details and exact concentrations for each analyte can be found Appendix A. Each pool was mixed well, dispensed (with constant stirring) in 2.4 mL aliquots into 5mL cryovials with screw cap, and stored frozen at $-70\pm 10^{\circ}\text{C}$. One box of each QC is kept in a convenient freezer for daily analysis.
- (2) **Proficiency testing (PT) materials.** At this time, there are no external PT programs or certified reference materials for aromatic amines in human urine. Therefore, we developed an in-house PT program administered by a QC officer (more details are in Section 10c).

E. Other Material and Supplies

Materials, supplies and their sources that used during the development and validation are listed below. Materials and supplies for use with this method should be equivalent to those listed if obtained from other sources.

- Pipettes and disposable tips capable of accurately dispensing the following volumes: 50 μl to 200 μl , 1 ml, 5 ml (Hamilton)
- Gas-tight syringes capable of accurately dispensing the following ranges of volumes: 1-100 μL , 5-500 μL , 0.1-2.5 mL (Hamilton)
- High recovery 3.7 mL glass vials (ChemGlass)
- PTFE coated cap (Wheaton)
- Disposable silanized test tubes, 13x100 mm (Fisher Scientific)
- 300 μl insert 12x32 mm amber vials (Wheaton)

- Assorted glassware

F. **Equipment.**

- Commercial Hamilton STAR automation system
- Commercial tandem mass spectrometer such as the Agilent triple quad G7000C (or comparable)
- Commercial gas chromatography system (GC) such as the Agilent 7968 GC system (or comparable).
- Commercial autosampler system (AS) such as the Agilent G4514A system, GC Sampler 80 (or comparable)
- Commercial Thermo Savant SpeedVac SPD 2010 (vacuum evaporator).
- Digital block heater, VWR block heater

G. **Instrumentation**

- (1) Gas chromatograph

Instrument Operating Parameters

Instrument Operation Parameter (Agilent 7890A, 7890B)	Setting
Injection port	Pulsed splitless
Injection port temperature	250°C
Injection volume	1 µL
Injection Pulse Pressure	30 psi until 0.35 min
Septum Purge Flow	3 mL/min
Purge flow to split vent	60 mL/min at 0.35 min
Column 1 (Agilent J&W, DB-FFAP, 122-3232)	Constant flow 1.68 ± 0.05 mL/min (analytical run); Constant flow -2.20 ± 0.05 mL/min (post-run)
Column 2 (Agilent J&W, Inert fused silica, 160-7265-5)	Constant flow 1.85 ± 0.05 mL/min (analytical run); Constant flow 10.5 ± 0.05 mL/min (post-run)
Liner (Agilent 5183-4693)	4mm ID, 900 µL splitless, single taper, glass wool, deactivated
MSD Transfer Line	280°C
Oven program (analytical run)	80°C (hold 2 min); 30°C/min ramp to 180°C; 15°C/min ramp to 240°C
Oven program (post-run)	240°C (hold 5 min)

Mass Spectrometer Source/Gas Parameters

Instrument Operation Parameter (Agilent 7000C)	Setting
Mass spectrometer mode	Positive ion detection
Ionization source	Electron impact (EI)
Source temperature	280°C
Quad temperature	150°C
Collision gas	Nitrogen at 1.5 mL/min
Collision pressure	35 psi
Quench gas	Helium at 2.25 mL/min
Detector setting	Gain at 20 (x 10 ⁵)

Multiple Reaction Monitoring Parameters

The ion pairs and compound dependent parameters are listed below.

Analyte	Precursor ion	Product ion	Collision energy	Dwell time
	(m/z)	(m/z)	(V)	(ms)
OTOL	253	116	30	55
	253	106	25	17
	259	122	30	35
	259	112	25	35
26DM	267	148	10	65
	267	120	25	65
	273	154	15	65
	273	126	30	65
OANS	269	135	30	35
	269	190	15	55
	276	139	35	35
	276	157	15	35
1AMN	289	144	20	45
	289	127	30	45
	295	150	20	45
	296	151	20	45
2AMN	289	144	20	45
	289	127	30	45
	296	151	20	45
	295	150	20	45
4ABP	315.1	170	30	65
	315.1	141	40	65
	324	179	20	65
	324	150	45	65

*Ion units in m/z. **DT = Dwell time units in milliseconds. CE = collision energy (in Volts).

The Agilent Triple Quad specific operational variables CE and DT refer to collision energy and dwell time.

(2) GC-MS/MS Instrument Control Program

An instrument control program for the Agilent GC/MS Triple Quad created using the MassHunter™ software that incorporates the above parameters is used for data acquisition.

7. CALIBRATION AND CALIBRATION VERIFICATION PROCEDURES

A. *Calibration Curve*

(1) Data Collection

A calibration curve is constructed at the beginning of each study using response factors (i.e., peak area ratio of analyte to labeled internal standard) versus calibration standard concentration of 14-17 calibrators (pg/μl). All 14-17 standard solutions are injected 3 times and subsequently analyzed on the GC/MS-MS to evaluate the linearity response for each analyte. A calibration curve is obtained with each analytical run and after a major service on the instrument (such as oil change, filament replacement, source cleaning, etc.).

(2) Calculation and Evaluation of Curve Statistics

The slope and intercept of the calibration curve are generated using linear regression with 1/x weighting. This analysis can be performed using the instrument's data analysis software or other suitable data analysis software (such as Indigo ASCENT™). The resulting plot should be examined for linearity over the entire calibration range (R-square \geq 0.98). Determine the slope and intercept of the calibration curve by linear least squares fit. Any deviations from this procedure (e.g., using a quadratic fit) must have a valid scientific justification and be approved by a supervisor.

B. *Usage of Curve*

The calibration range was chosen based on AA levels from previously measured urine samples from smokers and nonsmokers. Quantification can only be reported for values that fall within the calibration range (between highest and lowest calibrator levels).

For sample results that are higher than the highest calibrator, the analysis can

be repeated with a lower volume of sample to bring the result within the calibration range.

C. Calibration Verification

The accuracy of the calibration curve is verified by using testing calibrators. The testing calibrators were prepared using native standards purchased from a different vendor. If a different vendor is not available, a different lot of native standards from the same vendor can be used. Three levels of testing calibrators (low, medium, and high) were prepared and tested to verify calibration accuracy. This accuracy test is performed each time a new standard calibration set, or new internal standard stocks or sources are prepared and used for analyte quantitation.

8. PROCEDURE OPERATING INSTRUCTIONS; CALCULATIONS; INTERPRETATION OF RESULTS

A. Hamilton- volume verification

Volume verification is done monthly at: 45 μL , 1000 μL , and 2000 μL . The temperature is measured and used to determine the density of the liquid used. The Ultimate Liquid Class Validation_Venus2.med program and an integrated balance on the Hamilton is used to determine the volume dispensed by the Hamilton. If the %error between the dispensed and set values is less than 5%, then the volume delivered is considered accurate. If the difference is more than 5%, then a service call will be placed with Hamilton. No aliquot can be performed until a service engineer from Hamilton services certifies the volume verification.

B. Sample Preparation

An analytical run consists of 1 water blank, 1 QC low, 1 QC high, and unknown urine samples.

- (1) Thaw samples at room temperature if they are frozen.
- (2) De-cap samples and load onto Hamilton sample carriers. The following sequence is utilized: water blank, QC low, unknown samples 1-xx, and QC high. This will result in one run/batch with a total of 1 water blank, 1 QC low, 1 QC high, and maximum 29 unknown samples.
- (3) Place 1-32 barcoded high recovery vials in one-four sample delivery racks on the Hamilton deck carriers. The barcode is AYYMMDDxx, with xx being 01 to 32.
- (4) Place eight 0.5-mL vial of ISTD spiking solution in position 1-8 of

Hamilton carrier.

- (5) Scan in the barcodes from the original sample vials and receiving barcoded high recovery vials to generate a Hamilton output file that is used to make a sequence file and run sheet.
- (6) Aliquot 2.0 mL of each sample to a corresponding barcode high recovery vial.
- (7) Add 45 μL of ISTD and 50 μL of 10 M NaOH to each barcoded vial on the receiving rack.
- (8) Cap, vortex, and transfer the 1-32 barcoded vials to the dry heat block, pre-heat at 90°C. Hydrolyze overnight for approximately 15 hours.
- (9) Remove the samples from the heating block. Let cool to room temperature. Vortex samples.
- (10) Load the samples back onto the Hamilton carriers as before. The following sequence is utilized: water blank, QC low, unknown samples 1-29, and QC high.
- (11) Load each sample to a corresponding SLE cartridge. Wait approximately 10 minutes to let the sample absorb onto the cartridge gravitationally.
- (12) Deliver approximately 10 mL DCM to each cartridge, collecting the eluent into 13x100 silanized tube placing directly below each SLE cartridges.
- (13) Transfer the siliconized tubes to a Savant, dry down to approximately 250 μL .
- (14) Add 3 μL of TMA and 3 μL of PFPA to each sample for derivatization. Cap and vortex, and leave at room temperature for 30 minutes.
- (15) Transfer the derivatized sample to a 12x32 mm amber vial with 300 μL insert. Dry down to completion in a Savant.
- (16) Reconstitute with 10 μL toluene.
- (17) Cap vials with aluminum crimp caps.
- (18) Samples are immediately ran on GC-QQQ or store at $-20\pm 4^\circ\text{C}$.

C. **Sample Analysis**

This assay uses a GC coupled with a triple quadrupole mass spectrometer to quantitate AAs in human urine. The Hamilton output file that was generated during the sample preparation process is converted to a MassHunter sequence file. This sequence file is used to run the prepared batch on an Agilent GC triple quad 7000C.

The analytes are first resolved from other potential interferences on an Agilent J&W DB-17MS column (or comparable column such as Restek, Thermo can be used). Afterwards, further selectivity is accomplished using a triple quadrupole mass spectrometer operated under positive electron impact ionization and multiple reaction monitoring (MRM) mode. Comparison of the area ratio (native analyte area/isotope labeled analyte area) with previously generated calibration curve yields individual analyte concentrations.

Before the run:

- (1) Analyze a toluene solvent blank to check for system contamination in the GC/MS-MS system.
- (2) A typical run order is as follow: a hexane standard curve, toluene solvent blank, analytical samples, and a toluene solvent blank.

D. **Processing of Data**

Process all the raw data files using instrument's quantification software (or comparable software package).

- (1) Analyzing and Storing the GC/MS-MS Data
 - a. Upon completion of an analytical run, quantitation will be done via Indigo ASCENT™ quantitation software: selection and integration of quantification and confirmation peaks for all native analytes and the internal standards; sample QC verification such as retention time, internal standard counts, carry over, etc.
 - b. Review the automated integrations of peaks to ensure correct integration. Manually re-integrate if the integration was chosen incorrectly.
 - c. Verify and certify the quantitation results.
 - d. Download Indigo result files to the corresponding run sheet file folder.

- e. Upload the Indigo result files to Repeat Manager database per protocol.
- f. All Hamilton output files and report data files are stored in the TEB share drive folder:
[\\cdc\project\CCEHIP_NCEH_DLS_ERATB_TEBL\AA](#), organized according to study and run sheet name.

(2) Evaluation of Calibration Curves

The y-intercept of each calibration curve should not be significantly different from zero ($p \geq 0.05$); if it is, the source of bias should be identified. An R^2 of ≥ 0.98 is acceptable. Through visual inspection, check to see if any single standard is an outlier. If removal of a point changes the slope or intercept by more than 10% it should be considered an outlier. If either the highest or lowest standard is removed, the reporting limits must be adjusted to reflect the new reporting range.

(3) Evaluation of Quality Control Material

After the completion of a run, the calculated results from the analysis of quality control samples are compared to the established quality control limits to determine if the run is “in control”. Quality control procedures implemented in this method are defined by the Division’s Policies and Procedures Manual (for more information see: Caudill SP, Schleicher RL, and Pirkle JL (2008) Multi-rule quality control for the age-related eye disease study, Stat Med, 27: 4094-4106.). QC samples are subjected to the complete analytical process. The data from these materials are then used to estimate method precision and to assess the magnitude of any time-associated trends. The concentrations of these materials should cover the expected concentration range of the analytes for the method.

- g. If both the low and the high QC results are within the 2σ limits, then accept the run.
- h. If one of two QC results is outside the 2σ limits, then apply the rules below and reject the run if any condition is met.
 - i. Extreme Outlier – Run result is beyond the characterization mean $\pm 4 S_i$
 - ii. 3S Rule - Run result is outside a $3S_i$ limit
 - iii. 2S Rule - Both run results are outside the same $2S_i$ limit
 - iv. 10 X-bar Rule – Current and previous 9 run results are on same side of the characterization mean
 - v. R 4S Rule – Two consecutive standardized run results differ by more than $4S_i$. Note: Since runs have a single result per pool for

2 pools, comparison of results for the R 4S rule will be with the previous result within run or the last result of the previous run. Standardized results are used because different pools have different means.

(4) Calculation

Calibration curve of each analyte provides the concentration in pg/total volume. The sample result must be reported as ng/L, therefore, the ng/total volume must be corrected as follow: Final concentration (ng/L) = [(calculated concentration from calibration curve)*(volume of ITSD spiking solution in each sample)*(concentration of ITSD spiking solution)]/[(concentration of ISTD in standard)*(sample volume)]

9. REPORTABLE RANGE OF RESULTS

People are exposed to aromatic amines from a variety of sources. A broad range of urinary aromatic amine levels can be expected, extending from <LOD to greater than 10 µg/L. As described below, we designed our method to quantify aromatic amines at the concentrations typically found in smokers and some non-smokers. If an unusually high value is observed that is greater than the highest standard on the curve, that sample is flagged as outside the calibration range and a more dilute sample is re-analyzed if sufficient sample exists (see criteria given in section 10).

A. **Limit of Detection**

The detection limits for AAs in human samples are determined according to the guideline for determination of limits of detection by the Clinical and Laboratory Standard Institute (CLSI. Protocols for Determination of Limits of Detection and Limits of Quantitation: Approved Guideline. CLSI document EP17-A. Wayne, PA: Clinical and Laboratory Standards Institute, 2004).

B. **Accuracy**

Neat reference and internal standards are obtained from commercial sources. Stock solution concentrations are based on stated purity using gravimetric analysis.

Accuracy was determined by spiking known amounts of AA standard solution into hexane (accuracy in solution) and urine (accuracy in matrix). The accuracy was calculated by the following formula.

$$\%bias = 100 * (observed AA level - expected AA level) / expected AA level$$

Criteria for accuracy passing is the same as PT (Section 10c below).

C. **Precision**

The relative %RSD values calculated over 4-5 runs in 5 days include both within-day and between day error. Acceptable RSD values should consistently fall below 20% for all the analytes. If higher RSD values are obtained, the origin should be investigated and corrective action discussed with supervisor.

D. **Analytical Specificity**

A high degree of analytical specificity is achieved with this approach. Correct retention times, correct ion mass-to-charge ratios, and correct precursor/product ion transitions help ensure a very high degree of specificity and minimize the influence from any potential interference.

An established range of ratios of the response of quantitation ion to that of confirmation ion of QC samples is used to determine if an unknown sample test positive for a given analyte.

E. **Recovery**

Sample matrix effects for each analyte are evaluated. Spiking same amounts of isotope labeled AA internal standards in urine samples. Prepare these spiked samples according to sample preparation procedures. % recovery is calculated as the ratio of the responses of labeled internal standards (area count) in the urine samples to the responses of internal standard in the standards. The average recovery for all AAs is ranging from 30% to 50%.

F. **Linearity Limits**

The AA calibration curves established are linear over the concentration ranges from the low and high standard with R² values greater than or equal to 0.98. The lower reportable limit is either the LoD or the lowest standard concentration, whichever is higher. The upper reportable limit is the highest standard concentration. A residual plot of the calibrators is checked to confirm linearity.

G. **Ruggedness test**

Ruggedness testing was performed to assess the potential of important analytical variables to affect results. Each of these variables was systematically varied to examine their influence, if any, on the analytical results and was optimized to achieve sensitivity and high throughput.

10. **QUALITY ASSESSMENT AND PROFICIENCY TESTING**

A. Quality Assessment

Quality assessment procedures follow standard practices.

Examine the blank to check for possible contamination in the system or extraction solution or reagents. Next, evaluate standards obtained under the same instrument conditions as that for samples to check the performance of the GC/MS-MS system. If the retention times and peak intensities of the analytes are within the acceptable ranges described in section 8.b., then evaluate sample blanks, QCs and unknown samples in the run sequence.

Compare the QC results obtained from the run with the acceptance criteria to assure the proper operation of the analysis. If a QC result is “out of control”, the cause of the failure should be determined. No results from the associated batch may be reported.

B. Establishing QC Limits

As per division policy, acceptable QC concentration limits must be calculated from the concentration results observed in at least 20 characterization runs. During the 20 characterization runs, previously characterized QCs or pools with target values assigned by outside laboratories should be included to evaluate the analysis. The process of limits calculation is performed using the laboratory database and the SAS division QC characterization program (for more information see: Caudill SP, Schleicher RL, and Pirkle JL (2008) Multi-rule quality control for the age-related eye disease study, *Stat Med*, 27: 4094-4106.).

C. Proficiency Testing

- (1) PT is conducted twice a year. At this time there are no external PT programs or certified reference materials for aromatic amines in human urine. Therefore, we developed an in-house PT program administered by a quality control officer. PT samples are spiked urine samples at three different levels cover the calibration range (more details in Appendix A).
- (2) When a PT test is conducted, five PT samples are blind-coded by a QC officer. The five blind-coded PT samples are prepared as unknown samples by an analyst or a team as described in Section 8a. A correct determination, within $\pm 20\%$ of the known concentration, on at least 80% of the samples must be achieved to be considered proficient.
- (3) Performance in the PT program along with documentation of remedial action taken for unacceptable performance is to be documented in a QC

Manual in the laboratory that is available for review.

11. REMEDIAL ACTION IF CALIBRATION OR QC SYSTEMS FAIL TO MEET ACCEPTABLE CRITERIA

If the calibration or QC fails, all operations are suspended until the source or cause of failure is identified and corrected. Analytical results are not reported. After calibration and/or quality control have been reestablished, analytical runs may be resumed.

A. *Internal Standard Response*

If the area counts of the internal standards of the check standard fall below the established optimal absolute counts, or if the chromatograph peak start to have tailing, this indicates that the instrumental sensitivity has fallen below acceptable limits. The following steps should be taken, and the instrument sensitivity is checked after each step is performed. Once sensitivity has been reestablished, further steps are not necessary.

- (1) Replace injection liner & o-ring, and septum
- (2) Clean the injection port and needle guide
- (3) Trim the GC column
- (4) Replace the syringe wash vials and rinse or replace the syringe
- (5) Bake out the ion source and quad
- (6) Bake out the MMI
- (7) Tune the mass spectrometer
- (8) Clean the mass spectrometer ion source
- (9) Clean the mass spectrometer quads (this must be done by a service engineer)
- (10) If the sensitivity is lowered due to band broadening, inspect all
- (11) GC connections, leak test, and consider changing the analytical column.

B. Calibration Regression

If the linearity of the calibration curve criterion 0.98 is not met, check if the standards are prepared correctly or if an instrument malfunction has occurred. If no error is found in standard preparation, check if the detector is saturated. Also check if GC delivery pressure is deviated-. Other instrument specific factors that could cause calibrations problems, such as leak, should be checked and corrective action is taken as needed.

C. Analyte in Standards or QC Materials

If an unexpectedly large amount of analyte is measured in one of the calibration standards or QC materials, but is not seen in the remainder of the samples, this indicates a contamination of this particular sample. The source of this contamination should be investigated to prevent repeat occurrences, but no further action is required.

D. Analyte in All Samples

If an unexpectedly large amount of analyte is present in all measurements for a particular day, it is likely that one or more of the solvents or reagents used are contaminated. If necessary, prepare new solvents and/or reagents.

E. QC Sample Outside of Control Limits

Verify the integrity of the QC material if the result of QC sample falls outside the control limits. Check if the proper amount of internal standard was added to that sample. Also confirm that the integration was performed correctly. No analytical results can be reported for runs that QC is outside of the control limits.

11. LIMITATIONS OF METHOD; INTERFERING SUBSTANCES AND CONDITIONS

Some plastic labware, solvents, air and water may contain trace amounts of AAs which could contaminate urine samples. We guard against biased data by screening reagents and laboratory materials that come in contact with samples. In addition, the specificity of detection by GC/EI-MS/MS helps to avoid background chemical interferences. It is highly unlikely that another substance would have the same mass transitions, retention times, and relative abundances of different MRMs as any of our analytes.

12. REFERENCE RANGES (NORMAL VALUES)

The study population typically includes both smokers and non-smokers, therefore, a large range of urinary aromatic amine levels are expected. Current literature reported mean levels for nonsmokers and smokers are listed in the

table below. (9-10, 12). We plan to apply the method described in this document to characterize population-based ranges for U.S. smokers and non-smokers.

Analyte	Nonsmoker (ng/L)	Smoker (ng/L)
o-TOL	55	117
2-AMN	5	12
4-ABP	2	9

13. CRITICAL CALL RESULTS (“PANIC VALUES”)

Not applicable for this procedure.

14. SPECIMEN STORAGE AND HANDLING DURING TESTING

Samples are received refrigerated or frozen and stored frozen at $-70\pm 10^{\circ}\text{C}$. If the entire sample volume is not used in the assay, then residual urine is refrozen and stored at $-70\pm 10^{\circ}\text{C}$. All analytes were tested for stability for at least 5 thawing-re-freezing cycles as specified in section h in the Appendix B.

15. ALTERNATE METHODS FOR PERFORMING TEST OR STORING SPECIMENS IF TEST SYSTEM FAILS

If a problem with the method exists, samples are held in the freezer until it can be resolved. If necessary, extracted and derivatized samples ready for analysis can be sealed and stored at $-20^{\circ}\text{C}\pm 4^{\circ}\text{C}$ up to one year before they are analyzed.

16. TEST RESULT REPORTING SYSTEM; PROTOCOL FOR REPORTING CRITICAL CALLS (IF APPLICABLE)

Not applicable.

17. TRANSFER OR REFERRAL OF SPECIMENS; PROCEDURES FOR SPECIMEN ACCOUNTABILITY AND TRACKING


Any residual urine is stored at $-70\pm 10^{\circ}\text{C}$.

18. SUMMARY STATISTICS AND QC GRAPHS

See following pages.


2013-2014 Summary Statistics and QC Chart for 1-Aminonaphthalene, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	46	06APR17	16AUG17	445.000	28.758	6.5
4566	46	06APR17	16AUG17	106.591	8.452	7.9
6323	85	13SEP17	25JAN19	295.424	15.974	5.4
6322	85	13SEP17	25JAN19	100.893	5.706	5.7




2013-2014 Summary Statistics and QC Chart for 2,6-Dimethylaniline, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	40	07APR17	16AUG17	503.775	29.528	5.9
4566	40	07APR17	16AUG17	121.273	11.706	9.7
6323	123	13SEP17	25JAN19	425.374	23.127	5.4
6322	123	13SEP17	25JAN19	145.772	8.411	5.8




2013-2014 Summary Statistics and QC Chart for 2-Aminonaphthalene, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	47	06APR17	16AUG17	467.596	29.027	6.2
4566	47	06APR17	16AUG17	108.972	7.177	6.6
6323	122	13SEP17	25JAN19	281.139	17.297	6.2
6322	122	13SEP17	25JAN19	91.437	6.281	6.9




2013-2014 Summary Statistics and QC Chart for 4-Aminobiphenyl, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	47	06APR17	16AUG17	431.340	19.378	4.5
4566	47	06APR17	16AUG17	106.747	4.019	3.8
6323	122	13SEP17	25JAN19	316.557	18.826	5.9
6322	122	13SEP17	25JAN19	104.941	7.159	6.8




2013-2014 Summary Statistics and QC Chart for o-Anisidine, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	41	13JUN17	16AUG17	494.659	26.163	5.3
4566	41	13JUN17	16AUG17	153.146	9.843	6.4
6323	122	13SEP17	15NOV18	381.885	17.088	4.5
6322	122	13SEP17	15NOV18	138.648	8.180	5.9



2013-2014 Summary Statistics and QC Chart for o-Toluidine, urine (pg/mL)

Lot	N	Start Date	End Date	Mean	Standard Deviation	Coefficient of Variation
4565	43	06APR17	16AUG17	479.186	43.314	9.0
4566	43	06APR17	16AUG17	158.953	26.781	16.8
6323	114	13SEP17	25JAN19	443.167	36.197	8.2
6322	114	13SEP17	25JAN19	182.211	29.761	16.3



19. REFERENCES

1. Hammond SK, Coghlin J, Gann PH, Paul M, Taghizadeh K, Skipper PL, Tannenbaum SR. Relationship between Environmental Tobacco Smoke Exposure and Carcinogen-Hemoglobin Adduct Levels in Nonsmokers. *J Natl Cancer Inst* 1993; 85:474-8.
2. Perera FP, Santella RM, Brenner D, Poirier MC, Munshi AA, Fischman HK, Ryzin JV. DNA Adducts, Protein Adducts, and Sister Chromatid Exchange in Cigarette Smokers and Nonsmokers. *J Natl Cancer Inst*. 1987; 79:449-56.
3. Skipper PL, Tannenbaum SR. Molecular Dosimetry of Aromatic Amines in Human Populations. *Env. Health Perspec*. 1994; 102:17-21.
4. Wishnok JS. Environmental Carcinogens: In Vivo Monitoring Using GC/MS. *Anal Chem* 1992; 64:1126-35.
5. Skipper PL, Peng X, Soohoo CK, Tannenbaum SR. Protein Adducts as Biomarkers of Human Carcinogen Exposure. *Drug Metab Rev* 1994; 26:111-24.
6. Maclure M, Katz RB, Bryant MS, Skipper PL, Tannenbaum SR. Elevated Blood Levels of Carcinogens in Passive Smokers. *Am J. Pub Health* 1989; 79:1381-4.
7. Bryant MS, Skipper PL, Tannenbaum SR, Maclure M. Hemoglobin Adducts of 4-Aminobiphenyl in Smokers and Nonsmokers. *Cancer Research* 1987; 47:602-8.
8. Bryant MS, Skipper PL, Wishnok JS, Stillwell WG, Glogowski JA, Tannenbaum SR. Determination of Haemoglobin Adducts of Aromatic Amines by Gas Chromatography-Mass Spectrometry. *IARC Public Sci* 1993; 109:281-92.
9. Riedel K, Scherer G, Engl J, Hagedorn HW and Tricker AR. Determination of three carcinogenic aromatic amines in urines of smokers and nonsmokers. 2006; *J. Anal. Tox.* 30: 187-195.
10. Seyler TH and JT Bernert, Analysis of Urinary 4-Aminobiphenyl Metabolites in Non-smokers and Smokers by LC and GC Tandem Mass Spectrometry, *Biomarkers*, 2011; 16(3).
11. Caudill SP, Schleicher RL, Pirkle JL, Multi-rule quality control for the age-related eye disease study, *Stat Med*. 2008; 10;27(20):4094-106.

12. Airoidi L. 4-Aminobiphenyl-Hemoglobin Adducts and Risk of Smoking-Related Disease in Never smokers and Former Smokers in the European Prospective Investigation into Cancer and Nutrition Prospective Study. *Cancer Epi Biomarkers Prev.* 2005;14:2118- 2124.
13. Taylor JK, Quality Assurance of Chemical Measurements, Lewis Publishers, Boca Raton, Florida, 1987.

Appendix A

Standard Materials

A. Stock reagent preparation

Each analyte was ordered separately in pure form (vendors and purity shown below).

(1) Native standard

Stock solutions of individual analyte and subsequent combined analyte dilutions are prepared in ethanol. The following table shows an example how the stock solutions (A) and combined analyte dilutions (B, C, & D) are made. B contained only OTOL to make standards with higher levels of OTOL.

Table A1:

Date 1mg/mL Stock Made	Analyte	Vendor	Chemical Purity (%)	Purity Test	Lot	Analyte Weight (mg)	Stock A Conc. (mg/mL)	Vol. of Stock A in Dilution B (μ L)	Conc. of Analyte in Dilution B (ng/ μ L)	Vol. of Stock A in Dilution C (μ L)	Conc. of Analyte in Dilution C (ng/ μ L)	Vol. of Stock A in Dilution D (μ L)	Conc. of Analyte in Dilution D (ng/ μ L)
03/31/15	OTOL	Fluka	99.5	GC	1421704V	10.720	1.0666	500.000	106.664	235.000	10.026	95.000	1.013
07/28/14	26DM	Aldrich	98.5	GC	STBB5452V	7.840	0.7722			325.000	10.039	130.000	1.004
07/28/14	OANS	Fluka	99.7	GC	SZBB343XV	9.630	0.9601			260.000	9.985	105.000	1.008
07/28/14	QNL	Aldrich	97.5	GC	MKBL8695V	9.330	0.9097			275.000	10.006	110.000	1.001
03/31/15	1AMN	Fluka	99.9	GC	SZB8268XV	10.100	1.0090			245.000	9.888	100.000	1.009
03/31/15	2AMN	Sigma	98.0	TLC	SLBD0851V	11.090	1.0868			230.000	9.999	90.000	0.978
07/28/14	4ABP	Aldrich	99.0	TLC	128K1386V	10.830	1.0722			230.000	9.864	90.000	0.965

Final volume of stock A was 10 mL in ethanol

Final volume of stock B was 5 mL in ethanol

Final volume of Dilution C was 25 mL in ethanol

Final volume of Dilution D was 100 mL in ethanol

Table A2:

Prepare one complete set of calibration standards in hexane (B&J, lot DH930) as described in detail below. A total of 18 standards were prepared ranging from 0 to 1,300 pg/ μ L.

The concentration of AA in standard 1-17 is listed in table A3. Standard 14-17 contain only OTOL to accommodate higher levels of OTOL that were observed in unknown samples.

Table A3:

Analyte	Std. No. 1 (pg/ μ L)	Std. No. 2 (pg/ μ L)	Std. No. 3 (pg/ μ L)	Std. No. 4 (pg/ μ L)	Std. No. 5 (pg/ μ L)	Std. No. 6 (pg/ μ L)	Std. No. 7 (pg/ μ L)	Std. No. 8 (pg/ μ L)	Std. No. 9 (pg/ μ L)	Std. No. 10 (pg/ μ L)	Std. No. 11 (pg/ μ L)	Std. No. 12 (pg/ μ L)	Std. No. 13 (pg/ μ L)	Std. No. 14 (pg/ μ L)	Std. No. 15 (pg/ μ L)	Std. No. 16 (pg/ μ L)	Std. No. 17 (pg/ μ L)
OTOL	0.5	2.5	5.1	7.1	10.1	20.1	25.1	30.1	40.1	50.1	75.2	100.3	213.3	426.6	853.3	1,066.60	1,279.90
26DM	0.5	2.5	5.0	7.0	10.0	20.1	25.1	30.1	40.2	50.2	75.3	100.4	200.8				
OANS	0.5	2.5	5.0	7.1	10.1	20.0	25.0	30.0	39.9	49.9	74.9	99.9	199.7				
1AMN	0.5	2.5	5.0	7.1	10.1	19.8	24.7	29.7	39.6	49.4	74.2	98.9	197.8				
2AMN	0.5	2.4	4.9	6.8	9.8	20.0	25.0	30.0	40.0	50.0	75.0	100.0	200.0				
4ABP	0.5	2.4	4.8	6.8	9.6	19.7	24.7	29.6	39.5	49.3	74.0	98.6	197.3				

All standards were made at a volume of 50 mL, except standard 1 and 2 which were made at a volume of 100 mL.

Standard Number	Standard Conc. (pg/ μ L)	Final Standard Volume (mL)	Dilution Used	Vol. of Dilution Added (μ L)
0	0.000	0.000	NA	50.000
1	0.500	100.000	D	50.000
2	2.500	100.000	D	250.000
3	5.000	50.000	D	250.000
4	7.000	50.000	D	350.000
5	10.000	50.000	D	500.000
6	20.000	50.000	C	100.000
7	25.000	50.000	C	125.000
8	30.000	50.000	C	150.000
9	40.000	50.000	C	200.000
10	50.000	50.000	C	250.000
11	75.000	50.000	C	375.000
12	100.000	50.000	C	500.000
13	200.000	50.000	C	1000.000
13	213.328	50.000	B	100.000
14	426.656	50.000	B	200.000
15	853.312	50.000	B	400.000
16	1066.640	50.000	B	500.000
17	1279.968	50.000	B	600.000

1) Accuracy testing solutions

To test the accuracy of our calibration curve, accuracy testing solutions were prepared using native stocks that were purchased from different vendors. If a different vendor was not commercially available, then a different lot from the same vendor was purchased and used.

Table A4: Accuracy testing stocks and solutions

Date stock made	Analyte	Vendor	Lot	Purity (%)	Actual Weight	Stock A conc. (mg/mL)	Volume stock A used in C2 (µL)	Dilution C2 conc. (ng/µL)	Volume stock A used in D2 (µL)	Dilution D2 conc. (ng/µL)
7/28/2014	oTOL	TRC	14-ANR-105	97.0	11.540	1.1194	220.000	9.851	90.000	1.007
7/28/2014	26DM	TRC	12-ABY-5-1	98.0	8.410	0.8242	245.000	8.077	120.000	0.989
7/28/2014	oANS	TCI	SDGQC	99.7	9.400	0.9372	265.000	9.934	105.000	0.984
7/28/2014	1AMN	Fluka	7080X	99.9	10.180	1.0170	210.000	8.543	100.000	1.017
7/28/2014	2AMN	TRC	6-GNM-1-2	95.0	14.620	1.3889	180.000	10.000	70.000	0.972
7/28/2014	4ABP	Sigma	056K1361	99.0	10.630	1.0524	240.000	10.103	95.000	1.000

Stock A was made in 10 mL ethanol.

Dilution C2 was made in 25 mL ethanol.

Dilution D2 was made in 100 mL ethanol.

Testing ID	Conc. (pg/µL)	Final Volume (mL)	Stock Used	Stock Volume (µL)
ACC-1	2.500	100.000	D	250.000
ACC-2	7.000	50.000	D	350.000
ACC-3	25.000	25.000	D	650.000
ACC-4	50.000	50.000	C	250.000
ACC-5	200.000	25.000	C	500.000

Table A5: Accuracy testing solution concentrations

Analyte	ACC-1	ACC-2	ACC-3	ACC-4	ACC-5
oTOL	2.519	7.052	26.193	49.253	197.011
26DM	2.510	7.027	25.714	40.385	161.539
oANS	2.520	7.057	25.585	49.671	198.682
1AMN	2.503	7.008	26.442	42.713	170.853
2AMN	2.423	6.783	25.278	50.000	200.002
4ABP	2.412	6.755	25.994	50.514	202.055

2) Isotopically Labeled Internal Standards

All concentrations are calculated with correction for percentage of chemical and isotopic purity.

Table A6:

Date 1mg/mL Made	Analyte	Vendor	Chemical Purity (%)	Isotopic Purity (%)	Purity Test	Lot	Actual Weight (mg)	Stock DA Final Con.* 10 ml (mg/mL)
01/29/15	OTOL- ¹³ C ₆	MI	98.0	98.0	Mass Spec	713	2.360	0.2267
07/28/14	PTOL-D ₉	MI	98.6	99.0	GC	18	9.230	0.9010
07/28/14	MTOL-D ₉	MI	98.4	99.0	GC	310	9.070	0.8836
07/28/14	26DM-D ₆	TRC	98.0	95.3	Mass Spec	3-MIC-149-1	9.000	0.8405
07/28/14	OANS-D ₇	MI	99.8	99.4	HPLC & NMR	533	10.180	1.0099
01/29/15	2ABP-D ₉	MI	99.0	99.6	HPLC	932	10.430	1.0284
07/28/14	1AMN-D ₉	MI	98.1	99.1	HPLC & NMR	145	10.790	1.0490
12/21/15	2AMN-D ₇	CDN	98.0	95.0	NMR & Mass Spec	Z-330	3.170	0.3170
07/28/14	3ABP-D ₉	TRC	98.0	98.7	Mass Spec	12-SDJ-187-1	8.700	0.8415
07/28/14	4ABP-D ₉	CI	99.0	99.6	HPLC & NMR	II-9733	9.570	0.9436

a. Internal standard spiking solution for standard solutions

Table A7:

Date made	Analyte	Stock A Final Con. (mg/mL)	Volume Stock A used (μ L)	¹ Final DC Conc. (mg/mL)	¹ Final DC Conc. (ng/ μ L)	² Final DD conc (ng/ μ L)	³ Final conc of ISTD in ea STD solution (ng/mL or pg/ μ L)
12/21/2015	OTOL- ¹³ C ₆	0.227	1100.000	0.010	9.973	0.997	39.891
12/21/2015	PTOL-D ₉	0.901	300.000	0.011	10.812	1.081	43.247
12/21/2015	MTOL-D ₉	0.884	300.000	0.011	10.603	1.060	42.411
12/21/2015	26DM-D ₆	0.841	300.000	0.010	10.087	1.009	40.346
12/21/2015	OANS-D ₇	1.010	300.000	0.012	12.118	1.212	48.474
12/21/2015	2ABP-D ₉	1.028	300.000	0.012	12.341	1.234	49.365
12/21/2015	1AMN-D ₉	1.049	300.000	0.013	12.588	1.259	50.351
12/21/2015	2AMN-D ₇	0.317	1000.000	0.013	12.680	1.268	50.720
12/21/2015	3ABP-D ₉	0.842	250.000	0.008	8.415	0.842	33.661
12/21/2015	4ABP-D ₉	0.944	300.000	0.011	11.324	1.132	45.295

¹DC was made in 25 mL ethanol

²DD was made by 1:10 dilution of DC in hexane

³200 μ l of DD was added to each 5 mL standard solution

Toluene. Burdick & Jackson Microsolve VLSI MS80863-4 GC 99.9% purity.

Hexane. Burdick & Jackson 216-4 High Purity

a. Derivatization of calibrators and accuracy testing solutions

Aliquot 5mL of each calibration standard, 0-13, and accuracy testing solutions,

ACC-1 to ACC-5. An aliquot of hexane solvent blank was also added.

Add 200µL of ITSD DD solution in each aliquot.

Add 25 µL of TMA and 25 µL of PFPA, wait 30 minutes.

Transfer to siliconized tubes, dry completely in a Savant, reconstituted to 5 ml toluene using volumetric flasks: rinsing each siliconized tube 2x2 ml with toluene and bring up to final 5 ml.

Aliquot to 5 sets of 1 mL in high recovery amber vial.

b. Internal standard spiking solution for unknown samples:

Internal standard (ISTD) spiking solution for unknown sample preparation was prepared from stock DC. Spiking solution was dispensed in 4.0 ml aliquots and stored frozen at $-20\pm 4^{\circ}\text{C}$. The exact concentration of each isotopically labeled analyte in the ITSD spiking solution is shown below.

Table A8:

Date made	Analyte	Final DC Conc. (ng/µL)	Volume of Stock DC used (µL)	Spiking solution conc. (ng/mL) or (pg/µL)	Total amount ISTD added in ea sa (pg)
12/21/2015	OTOL- ¹³ C ₆	9.973	1000.000	9.973	448.776
12/21/2015	PTOL-D ₉	10.812	1000.000	10.812	486.528
12/21/2015	MTOL-D ₉	10.603	1000.000	10.603	477.124
12/21/2015	26DM-D ₆	10.087	1000.000	10.087	453.895
12/21/2015	OANS-D ₇	12.118	1000.000	12.118	545.329
12/21/2015	2ABP-D ₉	12.341	1000.000	12.341	555.357
12/21/2015	1AMN-D ₉	12.588	1000.000	12.588	566.445
12/21/2015	2AMN-D ₇	12.680	1000.000	12.680	570.600
12/21/2015	3ABP-D ₉	8.415	1000.000	8.415	378.682
12/21/2015	4ABP-D ₉	11.324	1000.000	11.324	509.566

Dilution DC was made in 25 mL ethanol.

Final spiking solution was made in 1L ethanol.

45 µl is spiked in each unknown sample.

B. Controls

1. Quality control materials

There are two quality control pools for the urinary aromatic amines assay. Pools QC Low and “QC High contain relatively low and high levels of urinary aromatic amines, respectively. The QC pools were prepared in-house from urine collected from non-smokers, filtered with 0.2µm filters and spiked with combined analyte standard stock solution. Each pool was mixed well, dispensed in 2.4mL aliquots into 5mL cryovials, and stored frozen at -

70±10°C.

Table B1: Concentration of AA in QC pools

Pool ID	Stock Dilution used for Spiking	Amount Spiked (µL)	Final Volume of Urine (mL)	Final Pool Desired Conc. (pg/mL)
QC_Low2015	D*	150	1200	125
QCHigh_2015	C*	60	1200	500
QCLow_2016	D*	350	3,200	100
QCHigh_2016	C*	125	3,200	400

*Dilution D and C are the same Dilution stocks used to make standard solution (Table A1)

Date Img/mL Stock Made	Analyte	Vendor	Lot	Weight (mg)	Stock A Conc. (mg/mL)	Vol Sotck A used in Dil C (µL)	Dilution C Conc. (ng/µL)	Vol Sotck A used in Dil D (µL)	Dilution D Conc. (ng/µL)	QC_Low2015 (pg/mL)	QCHigh_2015 (pg/mL)	QCLow_2016 (pg/mL)	QCHigh_2016 (pg/mL)
08/31/15	OTOL	Fluka	1421704V	10.720	1.0666	235.000	10.026	95.000	1.013	126.664	501.321	110.831	391.657
07/28/14	26DM	Aldrich	STBB5452V	7.840	0.7722	325.000	10.039	130.000	1.004	125.489	501.956	109.803	392.153
07/28/14	OANS	Fluka	SZBB343XV	9.630	0.9601	260.000	9.985	105.000	1.008	126.015	499.258	110.263	390.045
08/31/15	IAMN	Fluka	SZB8268XV	10.100	1.0090	245.000	9.888	100.000	1.009	126.124	494.405	110.358	386.254
08/31/15	2AMN	Sigma	SLBD0851V	11.090	1.0868	230.000	9.999	90.000	0.978	122.267	499.937	106.984	390.576
07/28/14	4ABP	Aldrich	I28K1386V	10.830	1.0722	230.000	9.864	90.000	0.965	120.619	493.198	105.542	385.311

Stock A was made in 10 mL ethanol.

Dilution C was made in 25 mL ethanol.

Dilution D was made in 100 mL ethanol.

2. Limit of Detection Materials

There are three limit of detection pools for the urinary aromatic amines assay. Pools “LoD1” and “LoD2” and “LoD3” contain relatively low, medium and high levels of urinary aromatic amines, respectively. The LoD pools were prepared in-house from urine collected from non-smokers, filtered with 0.2µm filters and spiked with standard stock solution. Each pool was mixed well, dispensed in 2.4mL aliquots into 5mL cryovials with screw cap, and stored frozen at -60°C to -70°C.

Table B2: Concentration of LoD pools

Pool ID	Stock Dilution used for Spiking	Amount Spiked (µL)	Final Volume of Urine (mL)	Final Pool Desired Conc. (pg/mL)
LoD0		0	300	TBD, background matrix
LoD1	E	60	300	20
LoD2	E	120	300	40
LoD3	E	180	300	60

Date 1mg/mL Stock Made	Analyte	Vendor	Lot	Weight (mg)	Stock A Conc. (mg/mL)	Vol Sotck A used in Dil D (µL)	Dilution D Conc. (ng/µL)	Dilution E conc. (pg/µL)	LOD 1 (pg/mL)	LOD 2 (pg/mL)	LOD 3 (pg/mL)
08/31/15	OTOL	Fluka	1421704V	10.720	1.0666	95.000	1.013	101.331	20.266	40.532	60.798
07/28/14	26DM	Aldrich	STBB5452V	7.840	0.7722	130.000	1.004	100.391	20.078	40.156	60.235
07/28/14	OANS	Fluka	SZBB343XV	9.630	0.9601	105.000	1.008	100.812	20.162	40.325	60.487
08/31/15	1AMN	Fluka	SZB8268XV	10.100	1.0090	100.000	1.009	100.899	20.180	40.360	60.539
08/31/15	2AMN	Sigma	SLBD0851V	11.090	1.0868	90.000	0.978	97.814	19.563	39.126	58.688
07/28/14	4ABP	Aldrich	128K1386V	10.830	1.0722	90.000	0.965	96.495	19.299	38.598	57.897

Stock A was made in 10 mL ethanol.

Dilution C was made in 25 mL ethanol.

Dilution D was made in 100 mL ethanol.


Dilution E was made by diluting 1 ml of Dilution D to 10 mL of ethanol.




Appendix B


Method Validation Data

a. Linearity Limits


We have confirmed linear responses for all analytes across a broad range of analyte concentrations relevant to urinary levels of aromatic amines ($R^2 \geq 0.99$). As shown below, the linear range of the analytical method extends from the lowest to highest calibrator, 0.5 to 200 $\mu\text{g}/\mu\text{L}$ for all analytes except OTOL, which extend to 1,200 $\mu\text{g}/\mu\text{L}$. Show below are calibration curves for all analytes.












A weighting factor of 1/x is used for all analytes base on the residual plot for each analyte shown below.










b. Calibration curve- Matrix vs non-matrix based calibrator



To verify matrix equivalency, 14 standard solutions were prepared as specified in Appendix A in hexane (non-matrix) and in urine (matrix) from the same stock solutions at the same time the hexane standards were made. The urine standards were carried through sample preparation as described in Section 8a for urine samples.



Each calibrator set was run 3 times, and the average of three slopes was obtained. The difference between the two average slopes of hexane and matrix set was less than 5% for all 5 analytes. The results indicate that there is no significant difference observed between the two calibration curves prepared in hexane or in urine.


Table B3:

Analyte	Hexane and urine calibration range (pg/μL)	Run date	Hexane calibration curve, slope			Hexane, averaged slopes	Urine calibration curve, slope			Urine, averaged slopes	Difference in averaged slopes (%)
			Run 1	Run 2	Run 3		Run 1	Run 2	Run 3		
OTOL	0.507 - 1280	11/6 - 11/9/17	0.01811	0.01673	0.01753	0.01746	0.01768	0.01688	0.01598	0.01685	3.49
26DM	0.502 - 100	10/16 - 10/18/17	0.01613	0.01567	0.01563	0.01581	0.01643	0.01573	0.01611	0.01609	-1.76
OANS	0.504 - 99.9	10/16 - 10/18/17	0.01164	0.01144	0.01160	0.01156	0.01153	0.01099	0.01141	0.01131	2.19
1AMN	0.504 - 99.9	10/16 - 10/18/17	0.03636	0.03736	0.03735	0.03702	0.03553	0.03610	0.03574	0.03579	3.33
2AMN	0.489 - 100	10/16 - 10/18/17	0.02788	0.02850	0.02849	0.02829	0.02924	0.02948	0.02924	0.02932	-3.64
4ABP	0.482 - 98.6	10/16 - 10/18/17	0.00699	0.00698	0.00673	0.00690	0.00712	0.00716	0.00697	0.00708	-2.62









c. Limit of Detection

The method detection limits for AAs in human samples are determined according to the guideline for determination of limits of detection by the Clinical and Laboratory Standard Institute (CLSI. Protocols for Determination of Limits of Detection and Limits of Quantitation: Approved Guideline. CLSI document EP17-A. Wayne, PA: Clinical and Laboratory Standards Institute, 2004). Four urine pools: one non-smoker blank pool (LoD0), non-smoker blank spiked approximately at 20, 40, and 60 pg/mL (LoD1, LoD2, and LoD3 respectively) for 6 analytes: *o*-toluidine, 2,6-dimethylaniline, *o*-anisidine, 1-aminonaphthalene, 2-aminonaphthalene, and 4-aminobiphenyl. Preparing details are in Appendix A, Section B2. These four pools are being used to estimate the LoDs. LoDs were obtained from 50 independent runs, and the results are listed below.

Table B4:

Analytes	LOD (pg/mL)
OTOL	111.22
26DM	15.67
OANS	7.02
1AMN	1.29
2AMN	2.79
4ABP	1.75

c. Accuracy

All native compounds were purchased from Sigma and affiliates and used to prepare calibration curve. Native standards purchased from a different vendor, mostly TRC, were used to prepared accuracy testing solutions.

Accuracy was determined by spiking known amounts of AA standard solution into hexane (accuracy in solution) and urine (accuracy in matrix). The spiked urine standards were prepared the same procedure as for unknown samples.

$$\%bias = 100 * (observed\ AA\ level - expected\ AA\ level) / expected\ AA\ level$$

Accuracy in solution- All analytes used to prepare the calibration curves and testing calibrators were purchased from different manufactures as listed in table below. If a different vendor was not available, a different lot was purchased from the same vendor.

Accuracy testing solutions were run with a hexane calibration curve, and the results are shown below.

Table B5: Accuracy in solution (%bias)

Analyte	%bias				
	ACC-1	ACC-2	ACC-3	ACC-4	ACC-5
OTOL	-7.43	-22.71	-2.22	-1.63	-7.48
26DM	3.78	1.84	1.80	-1.03	0.13
OANS	-13.19	4.86	-7.44	-12.39	-16.53
1AMN	-4.83	-1.25	2.55	-12.24	-3.52
2AMN	5.28	6.54	2.55	5.18	0.12
4ABP	-9.06	0.15	11.84	-8.43	-14.83

Accuracy in matrix was tested with three separate runs in three days, at three levels, and in triplicate at each level. All %bias except one are less than 15%, and accuracy in matrix tests are acceptable.

Table B6: Accuracy in matrix (%bias)

Analyte	Expected Conc (pg/mL)	Avg Calc Conc (pg/mL)	%-bias
OTOL	506.66	478.08	5.64
	101.33	96.77	4.50
	50.67	56.35	-11.23
26DM	506.96	524.09	-3.38
	152.09	165.04	-8.52
	50.70	54.72	-7.94
OANS	509.06	545.71	-7.20
	152.72	159.90	-4.70
	50.91	57.11	-12.19
1AMN	110.76	108.13	2.37
	106.61	94.51	11.35
	458.10	457.46	0.14
2AMN	489.07	525.85	-7.52
	146.72	160.45	-9.36
	48.91	53.46	-9.30
4ABP	109.03	103.12	5.42
	109.17	104.17	4.58
	349.32	321.83	7.87

d. Precision

Several pools were used to determine precision of our aromatic amines assay. The QC pools were prepared in house from urine collected from non-smokers and spiked with native standard solution (Appendix A, Section B.2). The QCs samples were prepared and analyzed. The CV% values for intra-day (n=5) and inter-day (n=5) are calculated and listed in the table below. All CV% values are less than 10% and thus acceptable.

Table B7: Inter and intra-day precision (n = 5)

Analyte	Sample ID	Intra-day (N=5)		Inter-day (N=5)	
		Avg Measured Conc (pg/mL)	CV%	Avg Measured Conc (pg/mL)	CV%
OTOL	QC_Low2016	188.5	7.73	185.4	5.60
	QC_Low2015	164.5	5.40	159.8	6.26
	QCHigh_2016	448.8	4.74	461.3	5.18
	QCHigh_2015	503.2	5.55	527.3	5.75
26DM	QC_Low2016	137.4	2.82	140.8	3.60
	QC_Low2015	113.9	6.22	115.6	5.88
	QCHigh_2016	410.5	3.97	421.4	5.61
	QCHigh_2015	487.1	3.50	503.0	4.37
OANS	QC_Low2016	140.6	2.14	145.3	1.84
	QC_Low2015	149.5	3.61	157.9	2.62
	QCHigh_2016	385.6	3.94	410.5	1.38
	QCHigh_2015	506.5	2.78	526.1	3.44
1AMN	QC_Low2016	108.5	2.35	109.8	3.37
	QC_Low2015	120.0	2.19	118.2	4.48
	QCHigh_2016	309.1	3.14	310.9	2.89
	QCHigh_2015	491.9	2.19	492.8	3.34
2AMN	QC_Low2016	102.3	1.93	105.7	1.97
	QC_Low2015	116.7	1.38	121.0	4.31
	QCHigh_2016	314.3	2.11	331.4	5.45
	QCHigh_2015	505.6	2.44	524.4	3.84
4ABP	QC_Low2016	108.5	1.76	113.6	2.99
	QC_Low2015	112.7	2.51	116.5	5.12
	QCHigh_2016	345.4	2.96	366.5	3.22
	QCHigh_2015	468.3	1.13	479.1	3.31

Matrix blank contained oTOL, oANS, and 2AMN.

f. Carry over

Carryover was examined by comparing successive pairs of injections of the highest calibrator, 200 pg/μl (ng/ml) or high QC samples followed by solvent blank, toluene. No carry over was observed in the solvent blank after any injection of the highest calibrator or QC. As a precaution, a toluene blank injection is done followed a calibration curve, QC high samples, and an analytical run. Between each individual injection, the injection syringe is washed with 6x3 μl toluene (injection volume is 1 μl).

g. Thaw-refreeze and storage stability

Samples are shipped frozen and stored at -70°C±10°C until processed (thawed and aliquoted) or prepared. Residual samples are re-frozen and stored at -70°C±10°C. Samples might be thawed again if repeated analysis is needed.

Repeated thaw and re-freeze cycles were carried out, and the results indicate that all analytes were stable following at least 5 thaw-refreeze cycles (T/RF) (sample loss \leq 5-10%). Results for individual analytes are listed below.

Table B8: Freeze/thaw stability after 5 cycles

Analyte	Control Sample Conc (pg/mL)	%-bias after five T/RF cycles
OTOL	218.0	9.09
26DM	157.5	5.46
OANS	158.2	8.71
1AMN	109.5	16.39
2AMN	101.4	10.93
4ABP	111.6	5.75

Prepared urine samples are routinely stored at approximately $-20^{\circ}\text{C}\pm 4^{\circ}\text{C}$ until they are analyzed on a GC/MS-MS system. During GC/MS-MS analysis, samples are stored in a cooled sample drawer at approximately 10°C . Results from repeated injections of samples stored in the cool tray at 10°C indicated that all analytes levels are not be significantly impacted by this short term storage conditions: 24 hour (overnight) after the first injection and 72 hours after the first injection (over a weekend). These conditions were chosen to ensure the longest time (that a sample would be left in the cool tray before being taken out and stored again) were tested.

Table B9: Short-term storage stability at 10°C

Analyte	Day 1 Conc (pg/mL)	%-bias to Day 1 in 24 hrs	%-bias to Day 1 in 72 hrs
OTOL	623.0	5.53	9.66
26DM	459.6	7.11	8.98
OANS	427.1	1.60	5.15
1AMN	269.2	-6.87	-12.71
2AMN	318.8	-0.81	-1.54
4ABP	347.0	2.06	1.46

Table B10: Short-term storage stability at -20°C

Analyte	Day 1 Conc (pg/mL)	%-bias to Day 1 in one week	%-bias to Day 1 in 13 weeks
OTOL	57.0	-6.55	6.72
	67.6	11.61	5.86
	82.2	-6.38	-12.25
	168.0	-0.40	-7.79
	490.2	-1.75	-2.26
26DM	19.8	10.28	10.92
	43.8	1.95	13.38
	68.1	9.89	12.45
	149.6	1.57	-0.63
	543.4	11.58	8.16
OANS	44.9	-3.25	-20.39
	70.5	5.88	9.05
	92.8	10.51	6.87
	153.5	10.39	7.81
	544.4	2.28	5.14
1AMN	13.6	-14.37	-3.62
	34.1	1.30	6.84
	43.8	-21.19	-15.84
	94.1	-22.87	-21.25
	449.4	-5.08	-11.64
2AMN	14.3	-5.54	-6.70
	35.4	-7.76	-6.46
	53.4	-6.35	-6.27
	96.1	-4.72	-3.42
	528.8	2.93	3.71
4ABP	21.5	5.36	-0.32
	42.5	2.65	0.90
	62.4	6.43	2.17
	120.2	11.81	10.04
	476.0	3.15	0.00

Table B11: Long-term stability of un-processed QC samples stored at -70°C for 2 years

Analyte	Nominal Conc (pg/mL)	%-bias
OTOL	485.9	-1.48
26DM	502.6	2.01
OANS	505.8	1.44
1AMN	458.1	-8.36
2AMN	479.8	-6.62
4ABP	446.1	-4.80

Table B12: Short-term stability for length of time needed to handle un-processed samples.

Analyte	Nominal Conc (pg/mL)	%-bias
OTOL	485.9	-5.24
26DM	502.6	-3.27
OANS	505.8	3.37
1AMN	458.1	5.08
2AMN	479.8	0.63
4ABP	446.1	3.99

h. Ruggedness test

QC AND UNKNOWN SAMPLES WERE USED IN RUGGEDNESS TEST BY VARYING THE FOLLOWING PARAMETERS:

- i. Injection Pulse Pressure: 25 psi, 30 psi (final method condition), and 35 psi. This could happen occasionally if the carry gas, He, is not pressurized enough during injection time. Injection pressure plays an important role on how much analyte being pushed on to analytical column.
- ii. Injection Port Temperature: 240°C, 250 °C (final method condition), and 260 °C. The injection port temperature could be varied slightly due to the insulator or heating element around the injection port. The injection port temperature affect the evaporation of injected sample, thus overall assay sensitivity.
- iii. Source temperature: 270°C, 280 °C (final method condition), and 290 °C. The source temperature could be varied slightly due to the flow of the GC column. Source temperature could affect the ionization of the analytes, thus the overall sensitivity of the assay.
- iv. Amount of PFPA used during derivatization: 2 µl, 3 µl (final method condition), 4 µl. The amount of PFPA has a crucial role in the efficiency of the assay derivatization.

- v. GC analytical column lots (country it was made, year it was made, the work order/polymer lot or batch, the run of the work order, the column of the work order, internal notification): GC column performance plays a crucial roles in the overall data quality: chromatographic peak, sensitivity and thus quantification of each final result. This test was done differently than other parameters (no %different between tested and final method condition). The same samples were injected on different columns of different lots, and the % different was calculated and listed. There were 3 lots tested, including one special order. The other 2 lots were spanned approximately one year.

Table B13: Ruggedness test results (%different between tested and final method condition)

Analyte	Measured Conc (pg/mL) at 250°C inlet (control)	240°C:250°C	260°C:250°C
		%-bias	%-bias
OTOL	214.7	-2.02	5.27
26DM	150.7	1.02	-2.65
OANS	133.6	-4.33	-3.49
1AMN	97.9	-9.15	-0.71
2AMN	93.0	-9.37	1.89
4ABP	105.6	-13.20	-6.91

Analyte	Measured Conc (pg/mL) at 30 psi pulsed pressure (control)	25psi:30psi	35psi:30psi
		%-bias	%-bias
OTOL	214.7	15.61	-0.20
26DM	150.7	-3.67	-1.09
OANS	133.6	-5.81	-13.27
1AMN	97.9	-0.96	-2.87
2AMN	93.0	1.08	0.07
4ABP	105.6	-8.10	-12.17

Analyte	Measured Conc (pg/mL) with 3 µL PFPA (control)	2µL:3µL	4µL:3µL
		%-bias	%-bias
OTOL	448.1	-15.12	1.40
26DM	456.0	9.23	7.98
OANS	403.4	2.86	7.00
1AMN	314.9	6.21	1.02
2AMN	315.8	2.00	2.43
4ABP	337.0	0.32	0.78

Analyte	Measured Conc (pg/mL) with 15hr hydrolysis (control)	14hr:15hr	16hr:15hr
		%-bias	%-bias
OTOL	448.1	-2.73	-14.49
26DM	456.0	4.60	9.74
OANS	403.4	1.20	3.52
1AMN	314.9	-1.39	8.64
2AMN	315.8	0.95	1.44
4ABP	337.0	-2.26	0.38

	OTOL	26DM	OANS	1AMN	2AMN	4ABP
Urine samples tested	Lot #1:Lot #2, %-bias	Lot #1:Lot #2, %-bias	Lot #1:Lot #2, %-bias	Lot #1:Lot #2, %-bias	Lot #1:Lot #2, %-bias	Lot #1:Lot #2, %-bias
Sample 1	-12.59	10.68	7.36	1.05	-0.78	-2.22
Sample 2	4.54	-0.21	-2.16	1.09	-0.22	-1.02
Sample 3	-1.64	-11.54	1.72	7.80	3.16	2.56
Sample 4	11.94	-2.67	-3.41	3.79	0.51	1.73
Sample 5	3.07	-4.01	-0.67	13.74	6.16	4.79
Sample 6	1.16	-5.13	0.45	8.81	1.19	1.72
Sample 7	12.78	-12.72	-0.13	6.46	2.62	-1.49
Sample 8	-4.92	-9.90	3.89	2.52	1.78	2.17
Sample 9	-7.84	-0.75	-1.36	9.28	1.54	-0.24
Sample 10	-0.42	2.21	-3.49	3.45	0.66	-0.08
Sample 11	-6.17	-9.94	-2.83	7.23	-1.86	5.45
Sample 12	-2.91	-0.56	-2.87	0.28	3.82	0.42
Sample 13	4.08	-2.78	4.40	4.44	-2.44	0.13
Sample 14	-7.27	1.34	1.10	2.39	4.18	3.72
Sample 15	0.42	-15.33	4.49	0.78	3.71	-11.10
Sample 16	10.68	-0.76	2.29	6.02	-0.30	0.02

	OTOL	26DM	OANS	1AMN	2AMN	4ABP
Urine samples tested	Lot #1:Lot #3, %-bias	Lot #1:Lot #3, %-bias	Lot #1:Lot #3, %-bias	Lot #1:Lot #3, %-bias	Lot #1:Lot #3, %-bias	Lot #1:Lot #3, %-bias
Sample 17	12.61	-9.02	-4.67	6.98	-0.24	-6.86
Sample 18	16.36	-3.12	2.54	6.51	0.07	-5.98
Sample 19	16.33	5.74	8.13	7.09	-16.65	-16.30
Sample 20	-14.85	-0.03	0.68	7.01	-2.18	-6.37
Sample 21	-7.09	-8.38	-10.71	16.07	-7.54	-4.61
Sample 22	1.97	-5.85	4.52	12.22	0.62	-0.71
Sample 23	18.66	-3.37	6.80	-10.32	-7.83	0.63
Sample 24	-21.82	-2.87	1.17	8.51	-9.39	-6.92
Sample 25	8.38	-2.22	-1.44	20.45	-5.99	-7.07
Sample 26	1.31	-6.26	6.46	7.47	-1.91	-1.51
Sample 27	20.14	-1.64	-3.48	-1.06	-3.44	9.40
Sample 28	9.83	-4.40	2.24	5.27	5.31	-7.42
Sample 29	-0.49	-11.39	0.35	4.26	1.97	-3.96
Sample 30	-14.65	-10.20	-3.50	-3.01	-4.60	-16.75
Sample 31	-5.62	-6.70	-1.51	0.97	3.90	-5.75
Sample 32	-14.46	-4.36	-1.90	9.78	2.69	-0.14

i. Instrument cross validation

There are 2 instruments used to analyze prepared samples, both are Agilent QQQ 7000C GC/MS-MS named Gonzo (main) and Animal (back up) respectively. 5 smokers and 5 nonsmokers were used for instrument cross validation.

Analyte	Pearson Value
OTOL	0.9782
26DM	0.9716
OANS	0.9985
1AMN	0.9993
2AMN	0.9998
4ABP	0.9998

Appendix C

Method Performance Documentation

Method performance documentation for this method including accuracy, precision, specificity, and stability is provided in this appendix. **The signatures of the Branch Chief and Director of the Division of the Laboratory Sciences on the first page of this procedure denote that the method performance is fit for the intended use of the method.**

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **OTOL**

Replicate	Sample 1					Sample 2					
	Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration			Recovery (%)	
		Day 1	Day 2	Mean			Day 1	Day 2	Mean		
Sample	1	0	119.31	247.25	177.1	0	88.11	55.71	63.7		
	2		157.16	141.29				50.24		57.80	
	3		165.89	231.92				61.79		68.61	
Sample + Spike 1	1	506.655	646.83	685.59	655.19	94.36	116.42	169.07	176.07	166.12	87.96
	2		619.61	734.63				163.24	165.82		
	3		608.12	636.38				171.77	150.72		
Sample + Spike 2	1	101.331	237.84	370.17	273.91	95.50	223.88	255.12	281.61	277.61	95.54
	2		241.17	257.53				263.44	277.61		
	3		219.22	317.51				307.45	280.46		
Sample + Spike 3	1	50.6655	201.94	336.48	233.49	111.23	805.95	945.11	946.63	902.46	104.07
	2		219.36	268.93				930.08	819.24		
	3		179.09	195.16				841.79	931.88		

Mean recovery (%)	SD (%)
98.1	8.2

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **26DM**

Replicate	Sample 1					Sample 2					
	Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration			Recovery (%)	
		Day 1	Day 2	Mean			Day 1	Day 2	Mean		
Sample	1	0	7.68	7.57	8.7	0	3.18	3.91	4.2		
	2		9.01	7.04				3.07		4.37	
	3		10.41	10.27				3.75		6.99	
Sample + Spike 1	1	55.71	58.00	54.56	58.02	88.60	121.29	125.61	103.56	121.40	96.62
	2		60.70	58.24				135.19	115.45		
	3		59.91	56.71				138.18	110.43		
Sample + Spike 2	1	107.14	113.68	102.76	109.70	94.31	144.1	174.26	147.63	154.11	104.03
	2		114.42					153.02	147.44		
	3		116.12	101.53				147.67	154.67		
Sample + Spike 3	1	385.72	402.13	383.81	403.50	102.36	416.19	432.54	429.72	420.74	100.08
	2		408.27	394.23				435.00	425.08		
	3		432.34	400.25				402.57	399.56		

Mean recovery (%)	SD (%)
97.7	5.7

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **OANS**

Replicate	Sample 1					Sample 2					Mean recovery (%)	SD (%)
	Spike concentration	Day 1	Day 2	Mean	Recovery (%)	Spike concentration	Day 1	Day 2	Mean	Recovery (%)		
Sample	0	30.20 31.32 28.84	32.66 26.46 29.59	29.8		0	75.01 73.60 (128.591)	71.54 72.14 70.54	72.6		101.0	8.1
Sample + Spike 1	58.48	89.94 86.49 90.06	89.06 92.78 92.45	90.13	103.09	50.91	117.26 118.44 111.25	111.07 115.38 123.08	116.08	85.48		
Sample + Spike 2	112.46	145.02 155.43 153.45	151.76 132.27 131.02	144.82	102.24	101.81	179.09 178.11 172.73	168.54 178.72 176.89	175.68	101.28		
Sample + Spike 3	404.86	492.42 481.03 495.74	414.61 429.40 435.30	458.08	105.77	203.62	284.22 293.82 300.42	(498.286) 294.48 293.63	293.31	108.41		

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **1AMN**

Replicate	Sample 1					Sample 2					Mean recovery (%)	SD (%)
	Spike concentration	Day 1	Day 2	Mean	Recovery (%)	Spike concentration	Day 1	Day 2	Mean	Recovery (%)		
Sample	0	1.59 1.35 1.70	2.27 1.22 1.64	1.6		0	1.30 1.32 1.17	3.21 2.97 3.02	2.2		92.2	5.5
Sample + Spike 1	110.76	99.27 97.10 102.32	122.08 119.18 118.64	109.8	97.63	52.88	47.47 50.21 49.84	46.73 42.79 47.05	47.3	85.45		
Sample + Spike 2	106.61	90.88 100.82 95.71	99.04 84.24 106.11	96.1	88.65	101.7	94.07 94.26 98.75	96.97 94.86 84.78	93.9	90.25		
Sample + Spike 3	458.1	421.48 422.24 421.61	501.63 500.44 487.14	459	99.86	366.11	326.39 338.34 352.28	327.79 314.35 363.91	337	91.51		

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **ZAMN**

Sample	Replicate	Sample 1				Sample 2				Mean recovery (%)	SD (%)		
		Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration				Recovery (%)	
			Day 1	Day 2	Mean		Day 1	Day 2	Mean				
Sample	1	0	2.51	3.75	2.51		0	1.72	2.30	1.83		103.2	6.2
	2		2.54	1.60				1.61	1.74				
	3		2.16	2.51				1.86	1.75				
Sample + Spike 1	1	489.07	508.44	543.22	528.36	107.52	111.68	116.07	117.59	113.31	99.82		
	2		511.95	546.92				111.07	117.76				
	3		520.32	539.28				102.76	114.61				
Sample + Spike 2	1	146.721	161.90	164.00	162.96	109.36	99.14	99.39	98.16	98.19	97.19		
	2		153.25	168.88				98.24	99.21				
	3		157.48	172.25				99.66	94.48				
Sample + Spike 3	1	48.907	55.95	60.69	55.96	109.30	306.55	301.70	282.44	295.67	95.85		
	2		51.29	57.95				306.00	292.66				
	3		52.17	57.74				302.21	289.03				

Accuracy using Spike Recovery

Recovery = (final concentration – initial concentration)/added concentration
Recovery should be 85-115% except at 3*LOD where can be 80-120%

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: **Urine**
Units: pg/mL
Analyte: **4ABP**

Sample	Replicate	Sample 1				Sample 2				Mean recovery (%)	SD (%)		
		Spike concentration	Measured concentration			Recovery (%)	Spike concentration	Measured concentration				Recovery (%)	
			Day 1	Day 2	Mean		Day 1	Day 2	Mean				
Sample	1	0	5.84	5.83	5.28		0	4.39	1.93	2.82		91.9	3.4
	2		5.92	3.15				4.07	1.49				
	3		5.74	5.22				3.47	1.57				
Sample + Spike 1	1	109.03	108.25	107.91	108.40	94.58	54.72	53.11	48.60	50.00	86.22		
	2		104.73	114.45				48.89	50.80				
	3		101.71	113.38				51.63	46.99				
Sample + Spike 2	1	109.17	108.67	110.67	109.45	95.42	105.24	104.52	90.55	97.49	89.96		
	2		111.38	107.67				102.68	89.61				
	3		105.59	112.71				107.60	89.97				
Sample + Spike 3	1	349.32	334.07	330.11	327.13	92.13	378.85	384.06	327.53	355.81	93.17		
	2		340.64	321.86				362.76	327.06				
	3		328.93	307.14				395.87	337.57				

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: OTOL

Quality material 1	A170828		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	224.88	201.83	295.26	280.28	283.608	236.12	161.53	160.8
Replicate 2	262.30	187.39	200.69	202.34	258.932	253.09	161.53	178.8
Replicate 3	203.32	205.41	199.98	291.86	217.414	263.71	161.53	167.1
Mean	230.17	198.21	231.98	258.2	253.318	250.98	161.53	168.9
% difference from initial measurement	--	-13.9	--	11.3	--	-0.9	--	4.6

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	622.97	542.71	622.97	459.99	590.785	564.95	492.29	488.2
Replicate 2	426.25	538.93	426.25	413.58	613.187	591.12	492.29	485.5
Replicate 3	470.00	489.84	470.00	560.16	574.191	578.63	492.29	505.4
Mean	506.41	523.8	506.41	477.9	592.721	578.2	492.29	493.0
% difference from initial measurement	--	3.4	--	-5.6	--	-2.4	--	0.2

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: 26DM

Quality material 1	A170828		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	174.26	179.97	174.26	175.36	125.611	126.28	119.27	103.6
Replicate 2	153.02	146.80	153.02	147.34	124.631	118.98	119.27	115.4
Replicate 3	147.67	147.42	147.67	155.25	120.34	131.98	119.27	110.4
Mean	158.32	158.06	158.32	159.3	123.53	125.75	119.27	109.8
% difference from initial measurement	--	-0.2	--	0.6	--	1.8	--	-7.9

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	459.57	399.25	459.57	429.72	525.347	487.59	503.6	514.9
Replicate 2	454.49	418.56	454.49	425.08	531.809	510.95	503.6	477.4
Replicate 3	457.57	412.54	457.57	399.56	500.118	471.89	503.6	485.3
Mean	457.21	410.1	457.21	418.1	519.09	490.1	503.6	492.6
% difference from initial measurement	--	-10.3	--	-8.6	--	-5.6	--	-2.2

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: OANS

Quality material 1	A170828		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	151.21	140.70	151.21	148.89	155.149	161.82	155.14	143.3
Replicate 2	143.32	145.12	143.32	136.86	156.647	161.11	155.14	149.0
Replicate 3	151.22	147.58	151.22	142.93	160.037	162.48	155.14	151.8
Mean	148.59	144.47	148.59	142.9	157.28	161.80	155.14	148.0
% difference from initial measurement	--	-2.8	--	-3.8	--	2.9	--	-4.6

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	427.06	377.25	427.06	399.84	475.003	491.15	502.18	511.7
Replicate 2	402.21	388.93	402.21	375.67	501.076	483.17	502.18	496.5
Replicate 3	404.62	380.46	404.62	376.16	490.281	462.83	502.18	487.3
Mean	411.30	382.2	411.30	383.9	488.79	479.05	502.18	498.5
% difference from initial measurement	--	-7.1	--	-6.7	--	-2.0	--	-0.7

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: 1AMN

Quality material 1	A170828		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	90.88	88.94	90.88	89.68	107.104	112.18	110.76	122.1
Replicate 2	100.82	92.49	100.82	93.38	109.298	101.80	110.76	119.2
Replicate 3	95.71	90.51	95.71	99.67	101.758	100.82	110.76	118.6
Mean	95.80	90.65	95.80	95.9	106.05	104.93	110.76	120.0
% difference from initial measurement	--	-5.4	--	0.1	--	-1.1	--	8.3

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	269.21	261.05	269.21	288.82	444.881	429.71	454.24	501.6
Replicate 2	300.47	273.78	300.47	311.55	434.998	443.18	454.24	500.4
Replicate 3	329.36	285.60	329.36	296.55	424.66	419.29	454.24	487.1
Mean	299.68	273.5	299.68	299.0	434.85	430.7	454.24	496.4
% difference from initial measurement	--	-8.7	--	-0.2	--	-0.9	--	9.3

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: ZAMN

Quality material 1	A170828		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	99.39	99.27	99.39	96.55	116.067	113.28	111.63	117.6
Replicate 2	98.24	97.77	98.24	92.69	114.422	114.74	111.63	117.8
Replicate 3	99.66	90.82	99.66	99.27	121.123	116.57	111.63	114.6
Mean	99.10	95.95	99.10	96.2	117.20	114.87	111.63	116.7
% difference from initial measurement	--	-3.2	--	-3.0	--	-2.0	--	4.5

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	318.76	255.50	318.76	282.44	475.08	466.81	477.84	521.6
Replicate 2	316.33	264.26	316.33	292.66	473.418	489.88	477.84	503.4
Replicate 3	315.29	284.68	315.29	289.03	481.723	466.57	477.84	509.5
Mean	316.79	268.14	316.79	288.0	476.74	474.42	477.84	511.5
% difference from initial measurement	--	-15.4	--	-9.1	--	-0.5	--	7.0

Stability

Freeze and thaw stability = Assess for a minimum of 3 freeze-thaw cycles; conditions should mimic intended sample handling conditions

Describe condition: example: five times frozen at -70°C and then thawed (5 freeze-thaw cycles)

Bench-top stability = Assess short-term stability for length of time needed to handle study samples (typically at room temperature)

Describe condition: example: original samples (not yet prepared for instrument analysis) frozen at -70°C, thawed, left at room temperature for approx. 4 hours before prepared

Processed sample stability = Assess short-term stability of processed samples, including resident time in autosampler

Describe condition: example: processed samples (ready for instrument analysis) stored in the cool auto-sampler tray at 10°C for 24 hour before reinjection

Long-term stability = Assess long-term stability that equals or exceeds time between date of first sample collection and date of last sample analysis

Describe condition: example: QC samples made in May 2015, stored at -70±10°C for 2 years

All stability sample results should be within ±15% of nominal concentration

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: 4ABP

Quality material 1	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	108.67	98.89	99.04	105.81	108.252	110.65	110.29	107.9
Replicate 2	111.38	103.30	84.24	102.69	105.265	109.31	110.29	114.4
Replicate 3	105.59	102.62	106.11	106.50	111.863	108.81	110.29	113.4
Mean	108.55	101.60	96.46	105.0	108.46	109.59	110.29	111.9
% difference from initial measurement	--	-6.4	--	8.8	--	1.0	--	1.5

Quality material 2	A170808		A170808		A170828		ABC170501	
	Initial measurement	Three freeze-thaw cycles	Initial measurement	Bench-top stability	Initial measurement	Processed sample stability	Initial measurement	Long-term stability
Replicate 1	347.03	300.35	347.03	330.11	426.174	430.10	441.33	473.9
Replicate 2	336.31	314.18	336.31	321.86	423.264	430.90	441.33	463.3
Replicate 3	337.61	323.48	337.61	307.14	435.255	437.58	441.33	465.2
Mean	336.96	312.7	336.96	319.7	428.23	432.86	441.33	467.5
% difference from initial measurement	--	-7.2	--	-5.1	--	1.1	--	5.9

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: OTOL

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	196.4439198	169.26	182.85	184.7413737	184.7413737	66869.67847
2	191.0568	205.456	198.26	51.83424016	51.83424016	78611.20028
3	186.75	162.061	174.41	152.3866803	152.3866803	60834.55686
4	184.422	171.538	177.98	41.499364	41.499364	63353.7608
5	169.118	165.294	167.21	3.655744	3.655744	55915.69287
6	211.86	198.473	205.17	44.80294225	44.80294225	84186.58544
7	153.917	182.033	167.98	197.627364	197.627364	56431.20125
8	218.799	172.885	195.84	527.023849	527.023849	76708.17793
9	163.18	178.374	170.78	57.714409	57.714409	58329.56746
10	221.08	192.69	206.88	201.441249	201.441249	85601.15138

Grand sum 3694.68672 Grand mean 184.734336

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	2925.454431	292.5454431	17.10395987	9.26
Between Run	4306.074878	478.4527642	9.641247875	5.22
Total	7231.529309		19.63413109	10.63

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	414.03	438.82	426.42	153.6473889	153.6473889	363668.9566
2	533.46	412.06	472.76	3684.589324	3684.589324	447011.2555
3	460.03	423.46	441.75	334.4692323	334.4692323	390278.1735
4	483.59	386.14	434.87	2374.028176	2374.028176	378218.6154
5	396.21	455.81	426.01	888.0102003	888.0102002	362975.0044
6	418.56	420.17	419.37	0.64561225	0.64561225	351736.5226
7	428.21	433.45	430.83	6.8644	6.8644	371232.4244
8	349.81	450.52	400.17	2535.525316	2535.525316	320267.2558
9	419.93	451.52	435.72	249.450436	249.450436	379710.8084
10	427.87	408.14	418.01	97.318225	97.318225	349456.3601

Grand sum 8611.80972 Grand mean 430.590486

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	20649.09662	2064.909662	45.44127708	10.55
Between Run	6392.044413	710.227157	0	0.00
Total	27041.14103		45.44127708	10.55

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: **26DM**

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	115.13	117.66	116.40	1.600225	1.600225	27095.59205
2	119.99	120.1	120.05	0.003025	0.003025	28821.60405
3	139.04	154.76	146.90	61.7796	61.7796	43159.22
4	117.76	122.1	119.93	4.7089	4.7089	28766.4098
5	128.35	115.76	122.06	39.627025	39.627025	29794.84605
6	125.28	125.64	125.46	0.0324	0.0324	31480.4232
7	127.55	125.8	126.68	0.765625	0.765625	32093.11125
8	127.77	120.86	124.32	11.937025	11.937025	30908.43845
9	150.91	183.24	167.08	261.307225	261.307225	55828.11125
10	107.8	92.77	100.29	56.475225	56.475225	20114.16245
Grand sum	2538.27	Grand mean	126.9135			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	876.47255	87.647255	9.362011269	7.38
Between Run	5921.188905	657.9098783	16.88583168	13.30
Total	6797.661455		19.30747437	15.21

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	503.91	533.86	518.89	224.250625	224.250625	538483.2865
2	504.78	490.64	497.71	49.9849	49.9849	495430.4882
3	466.59	480.5	473.55	48.372025	48.372025	448489.7341
4	505.63	493.99	499.81	33.8724	33.8724	499620.0722
5	514.46	506.1	510.28	17.4724	17.4724	520771.3568
6	499.07	501.57	500.32	1.5625	1.5625	500640.2048
7	535.09	495.02	515.06	401.401225	401.401225	530563.3061
8	538.91	512.1	525.51	179.694025	179.694025	552311.0101
9	502.36	476.66	489.51	165.1225	165.1225	479240.0802
10	457.72	484.34	471.03	177.1561	177.1561	443738.5218
Grand sum	10003.3	Grand mean	500.165			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	2597.7774	259.77774	16.11762203	3.22
Between Run	5987.5161	665.2795667	14.23906294	2.85
Total	8585.2935		21.50647933	4.30

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: **OANS**

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	153.94	149.2	151.57	5.6169	5.6169	45946.9298
2	149.31	151.24	150.28	0.931225	0.931225	45165.15125
3	161.19	149.34	155.27	35.105625	35.105625	48214.44045
4	146.67	151.04	148.86	4.774225	4.774225	44315.62205
5	151.97	168.44	160.21	67.815225	67.815225	51331.28405
6	141.68	149.28	145.48	14.44	14.44	42328.8608
7	152.21	151.93	152.07	0.0196	0.0196	46250.5698
8	146.29	151.07	148.68	5.7121	5.7121	44211.4848
9	149.34	156.12	152.73	11.4921	11.4921	46652.9058
10	150.72	148.4	149.56	1.3456	1.3456	44736.3872
Grand sum	3029.38	Grand mean	151.469			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	294.5052	29.45052	5.426833331	3.58
Between Run	296.47678	32.94186444	1.321238897	0.87
Total	590.98198		5.585355156	3.69

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	516.13	476.84	496.49	385.926025	385.926025	492994.7105
2	492.49	478.06	485.28	52.056225	52.056225	470983.6513
3	537.93	531.95	534.94	8.9401	8.9401	572321.6072
4	486.47	452.17	469.32	294.1225	294.1225	440522.5248
5	526.86	513.84	520.35	42.3801	42.3801	541528.245
6	460.06	454.75	457.41	7.049025	7.049025	418438.6681
7	472.39	484.93	478.66	39.3129	39.3129	458230.7912
8	487.53	542.9	515.22	766.459225	766.459225	530892.9925
9	487.01	472.78	479.90	50.623225	50.623225	460598.4221
10	470.95	483.61	477.28	40.0689	40.0689	455592.3968
Grand sum	9829.65	Grand mean	491.4825			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	3373.87645	337.387645	18.3681149	3.74
Between Run	11003.05312	1222.561458	21.03774956	4.28
Total	14376.92957		27.92802449	5.68

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: **1AMN**

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	94.47	109.31	101.89	55.0564	55.0564	20763.1442
2	102.74	106.9	104.82	4.3264	4.3264	21974.4648
3	110.52	118.28	114.40	15.0544	15.0544	26174.72
4	102.53	99.92	101.23	1.703025	1.703025	20493.00125
5	104.92	119.3	112.11	51.6961	51.6961	25137.3042
6	105.5	106.63	106.07	0.319225	0.319225	22499.56845
7	112.27	115.53	113.90	2.6569	2.6569	25946.42
8	100.44	95.9	98.17	5.1529	5.1529	19274.6978
9	105.4	116.43	110.92	30.415225	30.415225	24604.27445
10	96.16	109.74	102.95	46.1041	46.1041	21197.405
Grand sum	2132.89	Grand mean	106.6445			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	424.96935	42.496935	6.518967326	6.11
Between Run	604.012545	67.112505	3.508245288	3.29
Total	1028.981895		7.403021005	6.94

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	426.58	437.99	432.29	32.547025	32.547025	373740.6425
2	409.8	437.26	423.53	188.5129	188.5129	358755.3218
3	470.74	472.25	471.50	0.570025	0.570025	444615.0701
4	459.56	476.12	467.84	68.5584	68.5584	437748.5312
5	449.14	426.81	437.98	124.657225	124.657225	383644.2013
6	453.96	455.61	454.79	0.680625	0.680625	413658.7925
7	434.8	424.29	429.55	27.615025	27.615025	369017.8141
8	423.27	447.88	435.58	151.413025	151.413025	379451.1613
9	403.41	439.24	421.33	320.947225	320.947225	355029.5113
10	421.06	458.47	439.77	349.877025	349.877025	386786.5105
Grand sum	8828.24	Grand mean	441.412			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	2530.757	253.0757	15.90835315	3.60
Between Run	5556.48132	617.3868133	13.49650165	3.06
Total	8087.23832		20.86219683	4.73

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: **2AMN**

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	103.38	97.3	100.34	9.2416	9.2416	20136.2312
2	106.58	106.85	106.72	0.018225	0.018225	22776.18245
3	115.18	116.22	115.70	0.2704	0.2704	26772.98
4	109.14	113.59	111.37	4.950625	4.950625	24804.32645
5	107.33	104.28	105.81	2.325625	2.325625	22389.39605
6	110.37	111	110.69	0.099225	0.099225	24502.33845
7	97.67	107.2	102.44	22.705225	22.705225	20985.85845
8	104.7	106.09	105.40	0.483025	0.483025	22216.21205
9	99.92	103.07	101.50	2.480625	2.480625	20602.47005
10	112.4	113.87	113.14	0.540225	0.540225	25599.05645
Grand sum	2146.14	Grand mean	107.307			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	86.2296	8.62296	2.936487698	2.74
Between Run	489.20662	54.35629111	4.781910241	4.46
Total	575.43622		5.611561775	5.23

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	428.74	427.92	428.33	0.1681	0.1681	366933.1778
2	443.6	450.4	447.00	11.56	11.56	399618
3	462.64	513.27	487.96	640.849225	640.849225	476200.1641
4	464.61	471.91	468.26	13.3225	13.3225	438534.8552
5	463.3	457.43	460.37	8.614225	8.614225	423871.8665
6	465.14	467.8	466.47	1.7689	1.7689	435188.5218
7	459.38	458.91	459.15	0.055225	0.055225	421628.2621
8	450.97	484.35	467.66	278.5561	278.5561	437411.7512
9	448.34	449.88	449.11	0.5929	0.5929	403399.5842
10	457.43	474.31	465.87	71.2336	71.2336	434069.7138
Grand sum	9200.33	Grand mean	460.0165			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	2053.44155	205.344155	14.32983444	3.12
Between Run	4552.291105	505.8101228	12.25695655	2.66
Total	6605.732655		18.85675314	4.10

Precision

Total relative standard deviation should be $\leq 15\%$ ($CV \leq 15\%$)

Method name: Aromatic amines in human urine.
Method #: 2020.02
Matrix: Urine
Units: pg/mL
Analyte: **4ABP**

Quality material 1

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	100.45	99.82	100.14	0.099225	0.099225	20054.03645
2	99.82	105.58	102.70	8.2944	8.2944	21094.58
3	113.99	111	112.50	2.235025	2.235025	25310.25005
4	108.12	111.26	109.69	2.4649	2.4649	24063.7922
5	108.33	106.14	107.24	1.199025	1.199025	22998.69045
6	109.08	110.48	109.78	0.49	0.49	24103.2968
7	106.7	100.53	103.62	9.517225	9.517225	21472.13645
8	104.92	100.99	102.96	3.861225	3.861225	21199.46405
9	105.39	108.3	106.85	2.117025	2.117025	22831.70805
10	103.96	105.95	104.96	0.990025	0.990025	22031.10405
Grand sum	2120.81	Grand mean	106.0405			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev (%)
Within Run	62.53615	6.253615	2.500722895	2.36
Between Run	267.305745	29.70063833	3.4239614	3.23
Total	329.841895		4.239944182	4.00

Quality material 2

Run	Result 1	Result 2	Mean	SS 1	SS 2	2*mean^2
1	415.54	412.6	414.07	2.1609	2.1609	342907.9298
2	413.47	405.16	409.32	17.264025	17.264025	335077.5385
3	475.49	461.01	468.25	52.4176	52.4176	438516.125
4	436.05	423.1	429.58	41.925625	41.925625	369069.3613
5	428.65	432.55	430.60	3.8025	3.8025	370832.72
6	420.59	417.79	419.19	1.96	1.96	351440.5122
7	435.7	433.08	434.39	1.7161	1.7161	377389.3442
8	431.05	431.24	431.15	0.009025	0.009025	371772.0221
9	426.53	413.52	420.03	42.315025	42.315025	352842.0013
10	429.12	419.1	424.11	25.1001	25.1001	359738.5842
Grand sum	8561.34	Grand mean	428.067			

	Sum squares	Mean Sq Error	Std Dev	Rel Std Dev
Within Run	377.3418	37.73418	6.142815315	1.44
Between Run	4759.00862	528.7787356	15.66915051	3.66
Total	5136.35042		16.83022453	3.93

LOD, specificity and fit for intended use

Method name: Aromatic amines in human urine.
 Method #: 2020.02
 Matrix: Urine
 Units: pg/mL

Analytes	Limit of Detection (LOD)	Interferences successfully checked in at least 50 human samples	Accuracy, precision, LOD, specificity and stability meet performance specifications for intended use
OTOL	111.22	yes ¹	yes
26DM	15.67	yes ²	yes
OANS	7.02	yes ³	yes
1AMN	1.29	yes	yes
2AMN	2.79	yes	yes
4ABP	1.75	yes	yes

N=50

¹ interference is observed in Qual channel, approximately 10% in general population

² interference is observed in Qual channel, approximately 6% in general population

³ interference is observed in both Quant and Qual channel, approximately 6% in general population