



National Health and Nutrition Examination Survey (NHANES)

Home Urine Collection



August 2013

TABLE OF CONTENTS

<u>Chapter</u>		<u>Page</u>
1	HOME URINE COLLECTION	1-1
	1.1 Background.....	1-1
2	EQUIPMENT AND SUPPLIES.....	2-1
	2.1 Description of Examination Room in the MEC.....	2-1
	2.2 Description of Equipment and Supplies	2-2
	2.2.1 Equipment	2-2
	2.2.2 Supplies.....	2-2
	2.3 Start of Stand Procedures.....	2-4
	2.3.1 Room Setup Procedures	2-4
	2.3.2 Inventory Management	2-4
	2.4 Inventory Management during the Stand.....	2-5
	2.5 End of Stand Procedures.....	2-5
	2.6 Urine Collection Kit Storage	2-6
3	MEC PROTOCOL.....	3-1
	3.1 Overview.....	3-1
	3.2 Eligibility Criteria.....	3-2
	3.3 Pre-Examination Procedures.....	3-2
	3.4 Recruiting the Participant to Accept the Home Collection.....	3-3
	3.5 Examination Procedures	3-4
	3.5.1 Procedures When the SP Agrees to Take a Kit.....	3-4
	3.5.2 Labeling the Kits.....	3-10
	3.5.3 Urine Collection and Shipping Instructions.....	3-13
	3.6 Packing the Urine Collection Kit for the SP to Take Home	3-18
	3.7 Final Steps	3-18
	3.8 Coding Not Done Exams	3-20
4	USPS LABEL PRINTING UTILITY.....	4-1
	4.1 Warehouse Responsibility	4-1
	4.2 USPS Label Utility	4-1

TABLE OF CONTENTS (Continued)

<u>Chapter</u>		<u>Page</u>
5	FIELD FOLLOW-UP MANAGEMENT SYSTEM PROTOCOL	5-1
	5.1 Overview.....	5-1
	5.2 Supplies and Equipment	5-2
	5.3 Log-on to the FFMS	5-3
	5.4 Access the HUC Module	5-5
	5.5 Appointment Manager	5-5
	5.6 Call History.....	5-8
	5.7 Reports.....	5-11
	5.7.1 Reminder Postcard	5-12
	5.7.2 HUC Response Rate Report.....	5-14
	5.7.3 Check Status Report.....	5-16
	5.8 Supervisor.....	5-17
	5.8.1 Appointment Monitoring Report.....	5-17
	5.8.2 HUC Module Security.....	5-21
6	CHECK PROTOCOL.....	6-1
	6.1 Overview.....	6-1
7	TRACKING AND DELIVERY STATUS	7-1
	7.1 Overview.....	7-1
	7.2 Email Sent Weekly to Analysis Lab.....	7-3

List of Appendixes

<u>Appendix</u>		
A	“What To Do” Urine Collection and Shipping Instructions English and Spanish.....	A-1
B	Kit Labeling Instructions	B-1
C	Talking Points English and Spanish	C-1
D	Inventory of Home Urine Collection Supplies	D-1
E	Reminder Postcards English and Spanish.....	E-1

TABLE OF CONTENTS (Continued)

List of Exhibits

Exhibits

2-1	HUC Exam Room.....	2-1
2-2	Confidentiality label	2-5
3-1	Home Urine Collection exam screen.....	3-5
3-2	HUC Screen—SP answers “Yes.”	3-6
3-3	USPS shipping label tracking number bar code	3-7
3-4	USPS tracking numbers scanned into ISIS.....	3-8
3-5	Hard-stop edit	3-9
3-6	Dymo label examples.....	3-10
3-7	Shipping label and bar code linkage	3-11
3-8	Image label placement	3-12
3-9	Reprint labels	3-12
3-10	Evening Label on Styrofoam Lid.....	3-15
3-11	Morning Label on Styrofoam Lid.....	3-17
3-12	Status Screen - Complete.....	3-19
3-13	Status Screen – Not Done comments.....	3-20
3-14	Comment codes.....	3-21
4-1	Shipment Tracking Label Maker blank	4-2
4-2	Shipment Tracking Label Maker completed.....	4-3
4-3	Sample shipping label.....	4-4
5-1	FFMS connect to nhprodweb window.....	5-4
5-2	FFMS Log-in screen	5-4

TABLE OF CONTENTS (Continued)

List of Exhibits

Exhibits

5-3	System default opening screen.....	5-5
5-4	HUC Appointment Manager button.....	5-5
5-5	Appointment Manager HUC screen.....	5-6
5-6	Day opening tab	5-6
5-7	Appointments by month screen	5-7
5-8	Appointment status legend.....	5-7
5-9	HUC Call History search screen.....	5-8
5-10	Call History screen.....	5-9
5-11	Example of contact attempts.....	5-10
5-12	Contact status codes.....	5-10
5-13	HUC Reports screen and links.....	5-12
5-14	HUC Reminder Postcard screen	5-13
5-15	Response Rate Report first screen	5-14
5-16	Response Rate Report multiple stands selected.....	5-14
5-17	HUC Response Rate Report example	5-15
5-18	Check Status Report search screen	5-16
5-19	Check Status Report example	5-16
5-20	Supervisor main screen	5-17
5-21	Appointment Monitoring Report – select criteria	5-18
5-22	Appointment Monitoring Report example.....	5-19
5-23	Appointment Status change options.....	5-19

TABLE OF CONTENTS (Continued)

List of Exhibits

Exhibits

5-24	Appointment Status reason options	5-20
5-25	Call History displays at bottom of SPs window	5-20
5-26	HUC Module security screen.....	5-21
5-27	HUC Module security - assign menus	5-22
5-28	HUC Module security - available menus.....	5-23
5-29	HUC Security Module – assign supervisor status.....	5-24
6-1	English letter – check enclosure	6-2
6-2	Spanish letter – check enclosure	6-3
7-1	Sample tracking email	7-1
7-2	Example Status Dispatched to Sort Facility.....	7-2
7-3	Example Status Depart USPS Sort Facility	7-2
7-4	Example Status Delivered.....	7-3
7-5	Example emailed file sent to analysis lab	7-4

1. HOME URINE COLLECTION

1.1 Background

Heart disease and stroke are the leading causes of death in the United States. Too much sodium can increase blood pressure as well as the risk of heart attack and stroke. The Dietary Guidelines for Americans recommends an intake of $\leq 2,300$ mg of sodium per day for adults in general. The recommendation for individuals in the following population groups however, is $< 1,500$ mg of sodium per day:

- Individuals 51 years of age or older;
- African Americans;
- Individuals with high blood pressure;
- Individuals with diabetes; and
- Individuals with chronic kidney disease.^{1,2}

Americans eat on average about 3,300 mg of sodium a day.² The 1,500 mg recommendation applies to about half of the U.S. population overall and the majority of adults. Eating less sodium can help prevent, or control, high blood pressure.¹ The Federal Government is improving data collection on sodium, including the amount of sodium people consume, and their knowledge, behaviors and health outcomes.²

Home Urine Collection allows modeling of 24-hour sodium intake based on equations developed from the 2011 NHANES Urine Calibration Study. This study collected 24-hour urine specimens where each void was collected separately. A specimen from the first morning void plus a specimen from the morning, afternoon, and evening was tested for sodium along with a specimen from the 24-hour collection. Preliminary results of modeling the results from this study indicate the combination of the first morning void, plus the prior evening's void, best predicts the 24-hour excretion of sodium. A home urine collection of the first morning void and the prior evening's void will allow an alternative and less burdensome data collection modality to estimate sodium intake if 24-hour urine collection fails to be feasible on the NHANES.

The 2013 Home Urine Collection component of NHANES is sponsored by the CDC National Center for Health Statistics (NCHS).

¹ <http://www.cdc.gov/salt/>.

² <http://www.cdc.gov/VitalSigns/pdf/2012-02-vitalsigns.pdf>.

2. EQUIPMENT AND SUPPLIES

2.1 Description of Examination Room in the MEC

The Home Urine Collection component examination room is located in trailer 3 of the MEC. The room is equipped with a desk, a shelf above the desk, storage bins on the wall, a tech chair, and an SP chair. The HUC room has a storage area behind the curtain where kits and extra supplies are stored. Exhibit 2-1 shows the exam room.

Exhibit 2-1. HUC Exam Room



2.2 Description of Equipment and Supplies

The following equipment and supplies have been provided for the Home Urine Collection component of NHANES.

2.2.1 Equipment

- Dymo label printer.

2.2.2 Supplies

Consumable supplies include:

- Home Urine Collection kit (pre-assembled in the NHANES warehouse):
 - White cardboard shipping box with Styrofoam insert;
 - Two, 2-ounce Refrigerant gel pack pouch (Cold Ice™);
 - One 500 milliliter plastic specimen cup with screw on lid;
 - One 7" x 8," blue, double zip plastic bag containing a biohazard sticker and contents label;
 - One 8"X 8," clear, zip closable plastic bag containing an absorbent pad;
 - Two strips of packing tape; and
 - Temperature monitoring strip (affixed to the inside lid of Styrofoam container).
- 10" x 13" zip closable plastic bag;
- Preprinted Avery USPS shipping labels;
- Metered USPS postage stamps;
- Hard-copy color "WHAT TO DO Urine Collection and Shipping Instructions" (Appendix A) in English and Spanish.
- Permanent markers, heavy point;
- Dymo labels—1" x 2 1/8,";
- Jumbo paper clips;

- Poly tote hand bag with handle 20" x 20" x 5";
- Confidentiality label English;
- Confidentiality label Spanish;
- Sani Cloth towelettes;
- 7" Twist Ties;
- Temperature monitoring strips (extras/replacements);
- Purell hand soap;
- SORTKWIK finger moistener; and
- 3/4" Scotch tape.

Non-Consumable supplies include:

- Metal footstool;
- E-Z Fold stepstool - 12" black;
- Scissors;
- Tape dispenser; and
- Temperature monitor strip storage box w /lid.

Non-inventoried items:

- Wall aid – “Label the Kit” summary (Appendix B);
- Wall aid - Laminated hard-copy “Talking Points” English and Spanish (Appendix C);
and
- Container of extra HUC kit supplies.

2.3 Start of Stand Procedures

2.3.1 Room Setup Procedures

At the beginning of the stand, these procedures are followed to prepare the room and supplies for use. The FES will connect the computer, keyboard, mouse, dymo printer, and bar code wand.

1. Place the tech and SP chair in the upright position.
2. Clean the area (i.e., wipe countertop).
3. Remove supplies from the storage shelves and bin, including the HUC display samples and place them on the counter top. Store the extra supplies on the shelves and/or storage bin.

2.3.2 Inventory Management

The HUC component is managed by the health technologists. The par level for the component's supplies is set at the levels shown in Appendix D. A mid-stand inventory request for consumable supplies may be necessary to complete the stand.

2.3.2.1 Start of Stand Procedures

At the start of a stand, you will receive cases of kits that have each been pre-labeled at the warehouse as moon kits and cases of kits that have been pre-labeled at the warehouse as sun kits. On each sun and moon kit, there is an informational label located on the outside of the lid of the Styrofoam insert inside the cardboard shipping box. The label is affixed to the lid at the warehouse. This label contains space for the SP to write the date and time of the specimen collection and information about their collection (Section 3.5.3).

At the start of a stand, and as necessary on a weekly basis, MEC examiners make up a set number (as determined by the health tech in charge of the component) of HUC packets and tote bags. To make a packet, attach a stamp to three preprinted Avery USPS shipping labels. Put the three stamped shipping labels and one "WHAT TO DO Urine Collection and Shipping Instructions" sheet in a 10" x 13" clear Ziploc bag. These packets will be stored in a wall bin for easy access during exams. Next, apply a

confidentially label (Exhibit 2-2) to the upper middle area of each tote bag. Fold each tote bag and stack them in a bin in groups of 25.

Exhibit 2-2. Confidentiality label

<p>In order to protect your confidentiality, please do not write your name anywhere on the urine collection kit.</p> <p>However, it is very important that you urinate into the correct specimen cup. Please make sure that the picture on this label matches the label on the specimen cup and on the shipping box.</p>	<p>Para proteger su confidencialidad, por favor no escriba su nombre en ninguna parte de los materiales para tomar la muestra de orina.</p> <p>Sin embargo, es muy importante que orine dentro del envase correcto para tomar la muestra. Por favor asegúrese de que el símbolo de esta etiqueta coincide con la etiqueta del envase para la muestra y el de la caja para el envío.</p>
--	---

2.4 Inventory Management during the Stand

During the stand, store the stamps in the safe at the coordinator's station. When it is time to make more packets each week, the tech should ask the MEC manager for stamps. Toward the end of the stand, track how many additional packets will be needed to complete the stand. Avoid putting postage on shipping labels that will not be used.

If a mid-stand inventory request for consumable supplies is required to complete the stand, request what you think is needed plus a little cushion. Try to avoid over ordering since stamps and shipping labels are printed as needed and often cannot be used at another stand. An emergency shipment can be sent during the last 2 weeks of a stand if necessary. Please allow 3 days to receive the shipment since stamps and shipping labels are printed on demand.

2.5 End of Stand Procedures

At the end of every stand, discard unused preprinted Avery shipping labels. Count any unused stamps and return them to the warehouse via the "transfer inventory to warehouse inventory

sheet." Count the remaining kits as part of the end of stand inventory process and leave them onboard the MEC to continue onto the next MEC site.

NOTE: Kits cannot be stored in the Audiometry room during transit.

2.6 Urine Collection Kit Storage

A limited number of urine collection kits are shipped from the warehouse with the start of stand supplies. Each kit contains a temperature sensing strip and therefore cannot be stored in the belly compartments. Store the kits in the Audiometry room (not in the booth) and at the designated off-site space. At mid-stand inventory, submit a request for the number of supplies needed for the rest of the stand.

3. MEC PROTOCOL

3.1 Overview

The purpose of the Home Urine Collection is to collect and store urine specimens for potential future testing to estimate the daily intake of sodium. Participants do not receive any results.

All primary sample persons (SPs) aged 20 - 69 years old who have consented to future research, are not pregnant, and have vessel 45 marked as filled in the urine processing module of the Laboratory Subsystem, are eligible for the home urine collection exam. The participant is asked to take kits home and collect the first morning void and a void from the prior evening (collected between 5:30 p.m. and bedtime). The participant is instructed to collect a full void both times and mail it to the contract lab as soon as possible - ideally within 14 days after the MEC exam. However, if they cannot collect it within 14 days, then they are asked to collect and ship it at their earliest convenience. The kits contain instructions on how to collect and mail the urine specimen to the testing laboratory. Each participant who agrees to participate will take home three kits even though only two samples are required. Since most people urinate more in the morning, the SP is provided with two morning kits. One of the morning kits is for overflow only. Most people only need to use one morning kit.

The contents of each kit are assembled at the warehouse and placed inside the shipping box. All other supplies are sent to the field individually and are assembled in the MEC as needed.

In the MEC, the examiner opens the kits, explains the contents, and provides verbal, written, and visual instructions to reinforce the procedures. A self-addressed, stamped shipping box containing a U.S. Postal Service (USPS) Priority Mail shipping label is provided. The tracking number on the shipping label is used to track the urine shipment to the contract laboratory. This process is described in Chapter 7.

Most MEC examination staff members, including the MEC manager, are trained to perform HUC. Exceptions are the lab staff, coordinator, and physician. The exam takes place in trailer #3.

3.2 Eligibility Criteria

All primary SPs aged 20 - 69 years old are eligible for the Home Urine Collection (HUC) component. There are no safety exclusion criteria for HUC. The coordinator is only able to assign an SP to the HUC component if the SP has consented to future research, is not pregnant, has submitted a urine specimen and (1) the urine specimen volume is sufficient to fulfill the SPs urine processing protocol; or (2) the urine specimen volume is not sufficient to fulfill the SPs protocol, but the lab has enough urine to process vessel no. 45. If the SP submitted a small quantity of urine, but vessel no. 45 has not been processed by the MEC laboratory staff, the coordinator will not be able to assign the SP to HUC.

When the examiner attempts to recruit an SP who does not speak English and the examiner does not speak the language of the SP, a translator who does speak the language of the SP assists the examiner. The interpreter stays with the examiner and the SP for the entire examination.

The successful recruitment of the SP is an important aspect of this examination. Recruitment requires careful patience on the part of the examiner. Encourage the SP to accept the kits and follow through with the collection and shipment of the home specimen. To assist the examiner, talking points (Appendix C) and labeling summary (Appendix B) documents have been developed.

3.3 Pre-Examination Procedures

The Home Urine Collection ISIS application allows the examiner to (1) record that the participant consented to take kits; (2) scan the USPS tracking bar code into the system; and (3) print the labels for the kits.

Using the standard ISIS format, open the HUC application and log the SP into the exam by scanning the bar code on the SP ID bracelet or manually typing the SP ID when the SP arrives in the examination room. The HUC application only consists of two screens: one exam screen and the section status screen. Once the exam screen appears, recruit the SP to participate in the component.

3.4 **Recruiting the Participant to Accept the Home Collection**

Invite the SP to participate in the home urine collection with the assistance of the following talking points:

English:

- I would like to invite you to participate in a Home Urine Collection exam.
- We want to collect urine samples and store them for future studies. Examples include the analysis of how much salt you eat and how much iodine you get in your diet.
- I will ask you to collect two urine samples from home and then mail them directly to the laboratory.
- I will provide you with the urine collection kits, pre-addressed postage-paid shipping boxes, and detailed instructions to make it a very easy process for you.
- You will need to take the samples to the post office to mail it. When you mail the samples, we will send you a check for \$50.00.
- Would you like to participate in this exam?

Spanish:

- Deseo invitarle a participar en un examen en el cual usted mismo recolectará unas muestras de su orina en su hogar.
- Queremos recolectar muestras de orina para guardarlas para futuros estudios. Algunos ejemplos incluyen el análisis sobre cuánta sal usted consume y cuánto yodo hay en su dieta.
- Le voy a pedir que recolecte dos muestras de su orina en su hogar y que las mande por correo directamente al laboratorio.
- Le daré los materiales para recolectar las muestras de su orina, unas cajas de envío con nuestra dirección y con franqueo pre-pagado y también las instrucciones detalladas para que el proceso sea muy fácil para usted.
- Usted tendrá que llevar las muestras a la Oficina de Correos para enviarlas. Cuando envíe las muestras, le mandaremos un cheque por \$50.00.
- ¿Le gustaría participar en este examen?

If an SP initially refuses the component, ask questions to determine the reason for the refusal and try to address any of their specific concerns. Provide reassurance and encouragement. Show the SP the urine collection kits and describe the procedures. If he or she still refuses, code the exam as an SP refusal (Section 3.8).

3.5 Examination Procedures


3.5.1 Procedures When the SP Agrees to Take a Kit

Conveying clear and concise specimen collection instructions and shipping procedures is essential to the response rate and to specimen integrity. Once the SP agrees to participate in the home urine component, you will need:

- One moon urine collection kit
- Two sun urine collection kits
- One pre-assembled packet that contains:
 - One “WHAT TO DO: Urine Collection and Shipping Instructions”
 - Three stamped USPS shipping labels
- Two tote bags that have the confidentiality label already applied

Exhibit 3-1 illustrates the Home Urine Collection exam screen. It includes one drop-down list and three text boxes. Access the ISIS application drop-down box adjacent to the question, “Did the SP take a kit home with them?” The drop-down box options are either “Yes” or “No.” If No is selected, the reason must be identified on the Section Status screen.

Exhibit 3-1. Home Urine Collection exam screen



If the SP agrees to take a kit, select “Yes.” The three text boxes will activate. This is where the tracking numbers will be entered for the three shipping boxes (Exhibit 3-2).

Exhibit 3-2. HUC Screen—SP answers “Yes.”

The screenshot shows a software window titled "Stand:615 Session:615810 11/15/2010 08:30 am - 12:30 pm". The main content area is titled "Home Urine Collection" and displays the following information: "SP ID: 784174 Name: GONZALEZ, SIENNA N Age: 30 years Gender: Female Date: 07/10/2013 Time: 12:05 PM". Below this, a question is asked: "Did the SP take a kit home with them?" with a dropdown menu currently set to "Yes". Underneath, there are three text input fields labeled "Scan evening airbill barcode here", "Scan morning airbill barcode here", and "Scan 2nd morning airbill barcode here". At the bottom of the window, there are navigation buttons: "End of Section", "Close Exam", and "Finish", along with a status bar showing "Ready", "MEC Layer: 4/8/2013", "Application: 10.1.0", "Not connected to Coordinator", and "12:05 PM".

Use the bar code wand to scan the USPS bar codes into the text box. Exhibit 3-3 illustrates the shipping label's tracking bar code.

Exhibit 3-3. USPS shipping label tracking bar code



You will have three USPS labels in the pre-assembled packet. Exhibit 3-4 shows the tracking numbers once the bar code has been scanned into ISIS.

Exhibit 3-4. USPS tracking numbers scanned into ISIS

Home Urine Collection

File View Utilities Reports Window Help

Home Urine Collection: Stand:615 Session:615841 11/18/2010 01:30 pm - 05:30 pm

SP ID: 720141 Name: CHAN, JACKIE S Age: 20 years Gender: Male Date: 12/27/2012 Time: 11:56 AM

Home Urine Collection

Did the SP take a kit home with them? Yes

Scan evening airbill barcode here 123456789123456789123456789123

Scan morning airbill barcode here 789456123789456123789456123789

Scan 2nd morning airbill barcode here 894563211234567897894561231238

1 of 2


End of Section Close Exam Finish

Ready | MEC Layer: 11/28/2012 | Application: 10.1.0 | Connected to Coordinator | SP assigned to room | 12:01 PM

Start NHANES - Lab Vessel Pri... Lab Subsystem Home Urine Collection 12:01 PM

All USPS tracking numbers contain 30 numerals. If the tracking number field does not contain the requisite number of characters, a hard-edit stop (Exhibit 3-5) will alert the examiner to review the field and the application will not advance to the next screen.

Exhibit 3-5. Hard-edit stop



When a hard edit occurs, delete the tracking number entry, and re-scan the bar code. If the bar code on the envelope cannot be scanned, set this label aside and retrieve another. Return the label to the warehouse with an explanation that the bar code is not able to be scanned.

Notice in Exhibit 3-4 that the text beside each data entry field indicates what kit the tracking number will be linked to. Before the kits and shipping labels are scanned, they are not assigned to anything. Once the shipping labels are scanned however, it is important to keep track of the order in which the labels are scanned. The first label scanned will be linked to the evening kit. The second label scanned will be linked to the first morning kit and the third label scanned will be linked to the 2nd morning/overflow kit.

Once the tracking number on the USPS label is scanned, it also becomes linked to the bar code labels that go on the specimen cups that will print out later. This is described in detail in Section 3.5.2. After scanning the shipping labels, apply the first shipping label scanned to the moon box, the second shipping label scanned to one of the two sun boxes (it doesn't matter which one), and the third shipping label scanned to the remaining sun box. It is still important to keep the boxes in order.

3.5.2 Labeling the Kits

Experience has demonstrated that SPs tend to write their name somewhere on the kit. To assist the SP, a specimen cup identification system using computer-generated images has been developed to assure that the SP is using the correct specimen cup.

When the USPS bar codes have been successfully scanned into the application, select the arrow in the bottom right corner to launch the section status screen and print the labels.



Do not close the section status screen. Leave this final screen open while you explain the collection and shipping process to the SP.

When the screen is advanced to the section status screen, the Dymo printer will automatically print 10 labels. Three of the labels contain the SP ID number and bar code. The other seven labels containing matching pictorial (picture stickers) images will print using the household data to ensure that no two SPs in the dwelling unit are given the same picture image. An example of these two label types is illustrated in Exhibit 3-6. It is important to note that under no circumstance should the specimen cup display any identifiable SP personal information, including the age, first name, last name, or initials.

Exhibit 3-6. Dymo label examples



When the labels print, it is important to remember that each bar code label is linked to a specific shipping label. You will need to put the bar code label on a specific cup according to how you applied the shipping labels to the kits (Exhibit 3-7).

Exhibit 3-7. Shipping label and bar code linkage



Apply the bar code sticker ending in 157 with the moon picture to the moon cup. Apply the bar code sticker ending in 158 with the sun picture to the sun cup that's in the box that has the USPS label that was scanned in the first morning data field. The bar code sticker ending in 159 with the sun picture should be applied to the sun cup that's in the box that has the USPS label that was scanned in the second morning data field.

NOTE: The labels are not waterproof. Please put a piece of scotch tape over each bar code label. The tape will preserve the bar code if the cup leaks during transit.

You do not need to explain this linkage to the SP. Only ask them to put the moon cup back in the moon box and the sun cup back in a sun box. There are back-end procedures in place to handle situations where the SP inadvertently puts the sun cup in the wrong sun box.

Apply one picture sticker to the large tote bag beside the confidentiality label. Apply a picture sticker to each of the specimen cups beside the SP ID label on the side of the specimen cup. (Do NOT place the picture sticker on the lid of the specimen cup.) Apply a picture sticker to the bottom of each shipping box. Exhibit 3-8 illustrates the label placement. There is a wall card in the HUC room that details how to label the kits (Appendix B).


Exhibit 3-8. Image label placement



To reprint a label set, use the back navigation arrow to move back to the first Home Urine Collection screen and then advance the screen forward again.

A pop-up box (Exhibit 3-9) displays asking if you want to print another label set. Select “Yes” to reprint the labels or “No” to move forward without reprinting the labels.

Exhibit 3-9. Reprint labels



3.5.3 Urine Collection and Shipping Instructions

Present the urine collection and shipping instructions at the same time the kit is being labeled. The “WHAT TO DO” instruction document (Appendix A) contains all of the information the SP needs to successfully understand the procedure and collect and ship the sample. Follow this guideline and talking points as you explain the procedures to the SP:

Show the SP the “WHAT TO DO, Urine Collection and Shipping Instructions.”

- Collection and shipping instructions are included in the kit you are taking home.
- If you have any questions, call the toll-free telephone number on the instructions.

General Overview:

- Collect two urine samples – **one** in the evening between 5:30 p.m. and bedtime and **one** the **following** morning.
- Collect your urine sample as soon as possible and ideally within the next 14 days. However, if you cannot collect it within 14 days, then please collect and ship it at your earliest convenience. [EXAMINER: ASK FEMALE SPS NOT TO COLLECT THE SAMPLES DURING THEIR PERIOD.]

[NOTE: Tell all female SPS, regardless of age, about collecting samples during their period. If an SP collects the evening sample and starts her period overnight, she should throw the sample away and call the toll-free number to request another kit.]

- Whenever possible, collect the evening sample when you will be able to mail it the next day along with the morning sample.
- When you get home, place the ice packs in the freezer overnight or for at least 6 hours. Freeze the ice packs flat.

Label the parts of the kit and show the SP his or her unique picture sticker.

- This is your special sticker. If other members in your home participate, they will have a different sticker.

[NOTE: As you open each kit to label the cup and explain the contents, remove the ice packs and leave them on the counter. The ice packs will be put in a bag together at the end of the exam as explained in Section 3.6.]

Show the SP the difference between the evening and morning cups.

Evening Collection:


- Urinate (pee) once into the cup with the sticker of the moon in the evening between 5:30 p.m. and bedtime.
- Urinate a full void (until your bladder is empty).
- Screw the lid on the specimen cup **tightly** and package it in the shipping box immediately.
- Place the cup in the small blue Ziploc bag. Squeeze the air out of the bag and **seal** it.
- Put the **sealed** Ziploc bag in the larger clear Ziploc bag. Squeeze the air out of the bag and **seal** it.
- Put two of the frozen ice packs in the box – one on the bottom and one on the side.
- Put the sample in the Styrofoam box with the moon symbol.
- Put the Styrofoam lid on the package.
- On the label on top of the Styrofoam container, write the date and time and answer the questions. [EXAMINER: GO OVER THE LABEL WITH THE SP.]

Exhibit 3-10. Evening Label on Styrofoam Lid

Please write the date and time you collected the urine specimen.
Por favor anote la fecha y la hora en que tomó la muestra de orina.

Date/Fecha _____

Time/Hora _____ AM
 PM



Evening Collection:
Muestra de orina de la tarde:

All the urine is collected in the cup. / El envase contiene toda la orina.

Only a few drops are missing from the cup. / Únicamente unas pocas gotas hacen falta en el envase.

More than a few drops are missing from the cup. / Más de unas pocas gotas hacen falta en el envase.

[NOTE: When answering the questions at the bottom of the label, “missing” covers all reasons why the entire void may not be in the cup (i.e., spills, too much for the cup to hold).]

- Close the box with the tape provided.
- Store the urine inside the Styrofoam container with the frozen ice pack overnight until you collect the specimen in the morning.

Morning Collection:


- The following morning, when you first wake up, urinate (pee) the first void into the cup with the sticker of the sun.
- Urinate in the cup until your bladder is empty. You are only collecting **one** sample that contains a full morning void. If that one sample fits in one cup, throw the other kit away. However, if you fill the first cup and you still need to go more, stop mid-stream and finish in the second morning cup.
- If you do not use the second cup with a sun on it, throw it away. Do not mail back an empty cup.
- Screw the lid on the specimen cup **tightly**.
- If you used two cups, package each cup separately.
- EXAMINER: REVIEW PACKAGING INSTRUCTIONS
- EXAMINER: GO OVER THE LABEL WITH THE SP. POINT OUT THAT THIS LABEL IS SLIGHTLY DIFFERENT FROM THE LABEL ON THE EVENING KIT.
 - Please mark whether or not you are returning one or two morning cups on the label. If you use two morning cups, mark “**2 cups**” on both boxes.

Exhibit 3-11. Morning Label on Styrofoam Lid

Please write the date and time you collected the urine specimen.
Por favor anote la fecha y la hora en que tomó la muestra de orina.

Date/Fecha _____

Time/Hora _____ AM
 PM



Morning collection - I am returning:
Muestra de orina de la mañana – Estoy enviando:

1 cup only / 1 envase únicamente

2 cups / 2 envases

All the urine is collected in the cup(s). / El envase contiene (los envases contienen) toda la orina.

Only a few drops are missing from the cup(s). / Únicamente unas pocas gotas hacen falta en el envase (los envases).

More than a few drops are missing from the cup(s). / Más de unas pocas gotas hacen falta en el envase (los envases).

Mail the samples:

- Take the packages inside the Post Office and hand it to a teller or, take the packages inside the Post Office and put it in the large priority mail-drop box. (The large priority drop box is usually located near the self-serve kiosk.
- If the urine specimen has not been received by the laboratory in 14 days, then a reminder postcard will be sent to you. If you have already sent your urine then ignore this postcard.
- Do not mail the large outer plastic bag with your name on it to the laboratory.
- We will send you a check for \$50.00 as soon as you ship the package.

3.6 Packing the Urine Collection Kit for the SP to Take Home

After explaining the collection procedures and answering any questions that the SP may have, pack the urine collection kits using the following procedures.

1. Place the “WHAT TO DO Urine Collection and Shipping Instructions” and ice packs into the 10” x 13” Ziploc bag and seal the plastic bag.
2. SPs carry their kits home in a tote bag that we provide them. Two bags will be required for the three kits per SP.
3. Write the SP’s first and last name at the top of the tote bag that has the SP’s special picture sticker.
4. Place the two morning kits in one bag together. Place the moon kit and the 10” x 13” Ziploc bag containing the “WHAT TO DO Urine Collection and Shipping Instructions” and ice packs in the other tote bag.
5. Using a twist tie, tie one handle from each bag together. Make sure that the bag with the SP’s name is facing out so the coordinator can easily identify who the bag belongs to.

3.7 Final Steps

After you package the kits, click “Finish” on the status screen to release the SP from the component. Exhibit 3-12 illustrates the HUC status screen.

Exhibit 3-12. Status Screen - Complete

Home Urine Collection: Stand:430 Session:430122 12/15/2004 05:30 pm - 09:30 pm

SP ID: 921039 Name: DELANEY, JACOB Age: 17 years Gender: Male Date: 06/25/2009 Time: 08:21 PM

Status Screen

Status

Complete
 Partial
 Not Done

Comments

Other text

2 of 2


End of Section Close Exam Finish

Restock the shelf above the computer with one moon and two sun kits, escort the SP to the next exam, and place the tote bags in the old respiratory health room in trailer 1.

3.8 Coding Not Done Exams

Choose and enter the appropriate comment code when the HUC section status is Not Done. Exhibit 3-13 illustrates the Not Done comment text box and Exhibit 3-14 lists all the possible comment codes and their intended use.

Exhibit 3-13. Status Screen – Not Done comments



The screenshot shows a form titled "Status" with three radio button options: "Complete", "Partial", and "Not Done". The "Not Done" option is selected. Below the status options, there is a "Comments" section with a drop-down menu. The menu is open, showing a list of comment codes: "SP refusal", "safety exclusion", "SP refusal", "no time", "physical limitation", "communication problem", and "equipment failure". The "SP refusal" option is highlighted in blue. To the right of the drop-down menu, there is a text box labeled "Other text" which is currently empty.

To record a comment in the Comment text box, click on the drop-down list, use the scroll bar to view all choices and select the desired choice.

Exhibit 3-14. Comment codes

Comment code	Use when:
Safety exclusion	Not applicable.
SP refusal	The SP refuses to take a kit home.
No time	Not applicable.
Physical limitation	Not applicable.
Communication problem	SP is unable to understand and follow the instructions for the component due to cognitive impairment or other problem, and is unable to complete the test.
Equipment failure	Not applicable.
Language Barrier	SP is unable to understand and follow the instructions for the component due to language and interpreting difficulties, and is unable to complete the test.
SP ill/emergency	The SP became ill or an emergency occurred and the test was not performed on the SP.
Interrupted	An exam is interrupted, usually for a MEC-wide emergency, and cannot be completed by the SP.
Other, specify	If the above reason for a status code of Not Done is not explained by one of the comment codes, the examiner must choose Other, specify and record a comment in the text field.

When finished, select the bright blue arrow in the bottom right corner or select [Enter].

4. USPS LABEL PRINTING UTILITY

4.1 Warehouse Responsibility

Warehouse staff members prepare the HUC kits. This includes printing the USPS shipping labels. The shipping label routes the kit from the SP to the contract laboratory. The label contains a tracking number which is used to track the shipment from the SP's home to the analysis laboratory. The kits, stamps, and shipping labels are sent to the field individually and are assembled in the MEC as needed.

4.2 USPS Label Utility

The purpose of the USPS Label Printing Utility is to provide the warehouse staff with an easy way to print USPS mailing labels when preparing the HUC kits.

To print the USPS shipping labels, click on the **USPS Label Generation** icon on the desktop. The application screen is illustrated in Exhibit 4-1.

Select a Return Address (this is the Westat Mailing Address) and Destination Address (this is the address of the Contract Lab) by selecting a choice from the two “saved” drop-down lists.

Exhibit 4-1. Shipment Tracking Label Maker blank

The screenshot shows a software window titled "Shipment Tracking Label Maker". It features two main sections for address entry: "Return Address" and "Destination Address". Each section includes a "Saved" dropdown menu and text input fields for Name, Company, Address 1, Address 2, City, State, Zip 5, and Zip 4. At the bottom of the window, there are input fields for "Weight" (with the value "12") and "Quantity" (with the value "3"), and two buttons labeled "Print" and "Clear".

The window displays with the correct data in the various fields as illustrated in Exhibit 4-2.

Exhibit 4-2. Shipment Tracking Label Maker completed

The screenshot shows a window titled "Shipment Tracking Label Maker" with a blue header bar. The window is divided into two main sections: "Return Address" and "Destination Address".

Return Address:

- Saved: NHANES Warehouse (dropdown menu)
- Name: NHANES Warehouse
- Company: (empty)
- Address 1: 15835 Gaither Road
- Address 2: GA-W10
- City: Gaithersburg
- State: MD
- Zip 5: 20877
- Zip 4: (empty)

Destination Address:

- Saved: Department of Pediatrics (dropdown menu)
- Name: Microalbumin Laboratory
- Company: Department of Pediatrics
- Address 1: 420 Delaware St. S.E.
- Address 2: MMC 491
- City: Minneapolis
- State: MN
- Zip 5: 55455
- Zip 4: (empty)

Weight: 12

Quantity: 3 (spin button)

Buttons: Print, Clear

Enter the quantity of labels and select the print button.

Once the print button has been selected, the system will connect to the USPS web site and generate and print the requested number of shipping labels.

A sample shipping label is illustrated in Exhibit 4-3.

Exhibit 4-3. Sample shipping label

P

USPS PRIORITY MAIL®

WESTAT
WESTAT
1450 RESEARCH BOULEVARD
ROCKVILLE MD 20850

SHIP TO:
R SCHMOE
DUKE ST
EPHRATA PA 17522-2031

ZIP - e/ USPS DELIVERY CONFIRM

420 17522 9101 8052 1390 7842 7057 94

Electronic Rate Approved #805213907

Cut on dotted line

Instructions

1. Please use a laser or laser-quality printer.
2. Adhere shipping label to package with tape or other secure fasteners. **DO NOT TAPE OVER BARCODE.** Be sure all edges are secure. Self-adhesive label is recommended.
3. Place label so that it does not wrap around the edge of the package.
4. Affix Priority Mail service postage. There is no extra fee for Delivery Confirmation service.
5. Stamped packages weighing 16 ounces or more may not be placed in Postal Service collection boxes. For information on pickup options, go to the Pickup page on www.usps.com.
6. Each shipping label may be printed on one page. Be sure to use the correct paper size and weight. Be sure to use the correct postage meter. For more information, go to the www.usps.com website.
7. Please use this shipping label on the "Ship date" selected when you requested the label.

Delivery Confirmation™ Service Number
9101 8052 1390 7842

Priority Mail® with electronic confirmation service
Print Date: 07/24/08
Electronic Confirmation Service Fee: 0.00

From: WESTAT
WESTAT
1450 RESEARCH BOULEVARD
ROCKVILLE MD 20850

To: R SCHMOE
217 DUKE ST
EPHRATA PA 17522-2031

*Digital Priority Mail Service postage only. There is no fee for electronic postage. Delivery Confirmation is not available by phone or the website, only.

UNITED STATES
POSTAL SERVICE

Thank you for shipping with the United States Postal Service!
Check the status of your shipment on the [Track & Confirm](http://www.usps.com) page at www.usps.com

5. FIELD FOLLOW-UP MANAGEMENT SYSTEM PROTOCOL

5.1 Overview

The HUC module in the FFMS (Field Follow-up Management System) supports the tracking and management of the HUC component. A HUC appointment is created daily at the end of the MEC sessions for each eligible HUC participant.

- If the SP does not agree to take a kit or does not receive the component, the appointment status is set to Final Non-Response (FNR) with a comment of “Primary Appt. Not Complete.”
- If the SP agrees to take a kit, the appointment status is set to “Scheduled, Tentative.”
 - A nightly procedure creates a PF appointment record for 14 days following the MEC Exam appointment. The PF Appointment Status is set to “Scheduled, Tentative.”
- Once the mailperson scans the kit into the USPS system showing that the kit was picked up, the PF appointment is set to “Scheduled, Firm.”
 - “Scheduled, Firm” records get checks processed.
- If 14 days have elapsed since the SP was initially given the kit and there is no tracking activity, then the PF appointment is “Not Scheduled, Broken.”
 - A nightly procedure looks to see if there are appointments that are “Scheduled, Tentative” and 14 days past their appointment date and then marks the record as broken which triggers the reminder postcard.
- Once the USPS system marks the air bill as “receipted,” meaning they delivered it, the PF appointment status is “Complete.”
 - This triggers the creation of efiles that are emailed to the laboratory, as well as the expectations of data to be received back.
- If after 30 days the air bill has not been marked as “delivered,” then the PF appointment status is marked as FNR (Final Non-Response), SP refusal.
 - A nightly procedure looks for all records where the status is “Scheduled, Tentative” and >30 days past MEC appt. and marks them as FNR.

If the participant has not sent the sample to the lab within 14 days from the day he or she was recruited in the MEC, then the FFMS generates a reminder postcard list and updates the FFMS

contact history indicating that the reminder postcard was sent. Support staff print the labels and mail the postcards.

Once the air bill is marked by the USPS as “Processed through Sort Facility” (meaning they scanned the kit into their system), a remuneration check for \$50 is mailed to the participant. This process is described in Chapter 6.

Use the **Call History** option to document contact(s) with an SP and/or other activities like sending a replacement kit.

Use the **Reports** option to (1) print labels for the reminder postcards; (2) calculate and display the response rate; and (3) access information about the check status for an SP.

Use the **Supervisor** option to (1) view a report that lists SPs by appointment status and stand and (2) assign employees to the HUC component within the FFMS. The Interview Scheduler link is not applicable to the HUC component.

5.2 Supplies and Equipment

Store enough kits and shipping supplies to respond to requests for replacement kits. The warehouse staff provides all these supplies, except for the reminder postcards and Dymo labels.

Equipment

- Dymo label printer.

Supplies

- Home Urine Collection kit (pre-assembled in the NHANES warehouse):
 - White cardboard shipping box with Styrofoam insert;
 - Two, 2 ounce refrigerant gel pack pouch (Cold Ice™);
 - One 500 milliliter plastic specimen cup with screw-on lid;
 - One 7” x 8”, blue, double zip plastic bag containing a biohazard sticker and contents label;
 - One 8”X 8”, clear, zip closable plastic bag containing an absorbent pad;

- Two strips of packing tape; and
- Temperature monitoring strip (affixed to the inside lid of Styrofoam container).
- 10" x 13" zip closable plastic bag;
- Preprinted Avery USPS shipping labels;
- Metered USPS postage stamps;
- Hard-copy color "WHAT TO DO Urine Collection and Shipping Instructions" (Appendix A) in English and Spanish.
- Permanent markers, heavy point;
- Dymo labels—1" x 2¹/₈";
- Poly tote hand bag with handle 20" x 20" x 5";
- Confidentiality label, English;
- Confidentiality label, Spanish;
- 12" x 12" x 24" shipping box; and
- Reminder postcard (Appendix E) in English and Spanish.

Reminder postcards are printed once per year, so notify the component person when it is time to reorder; allow 6-weeks lead time. The Dymo labels are printed by ISIS as needed. The tracking numbers on the Avery USPS shipping labels expire over time and are printed by the warehouse as needed. The cost of postage varies by location and is printed by the warehouse as needed.

5.3 Log-on to the FFMS


To launch the FFMS, click on the **FFMS** icon on the desktop. The first screen is the **Connect to nhprodweb** screen, shown below. Enter the **User name** and **Password** and then select the **OK** button to move forward to the login screen or the **Cancel** button to exit. Exhibit 5-1 illustrates the Log-on window.

Exhibit 5-1. FFMS connect to nhprodweb window



The FFMS **Log-in** screen displays and is shown below in Exhibit 5-2. Enter the **User ID** and **Password** and press either **Enter** or click on the **Log-in** button to launch the application.

Exhibit 5-2. FFMS Log-in screen




The session times out after 60 minutes.

5.4 Access the HUC Module

After logging in, the system defaults to the **Dietary Appointment Manager** screen. Select the **HUC** header (far right) to access to the HUC module. Exhibit 5-3 illustrates the default opening screen.

Exhibit 5-3. System default opening screen




5.5 Appointment Manager

The **Appointment Manager** module is designed to track scheduled appointments. However, this functionality is not necessary for the HUC component. It is useful to look at these data by day, week, or month to get a sense of the number of SPs by appointment type.


Select the HUC **Appointment Manager** button as shown below in Exhibit 5-4.

Exhibit 5-4. HUC Appointment Manager button




The **HUC Appointment Manager** screen displays and it is shown below in Exhibit 5-5.

Exhibit 5-5. Appointment Manager HUC screen



The screen will open to the **Day** tab. **Day** and **Week** tabs list all SPs who are scheduled and their appointment status. If there are no appointments, then the screen will include only the statement, “There are no appointments for ‘today’s date,’ for the given criteria.” Exhibit 5-6 illustrates the Day opening tab screen.

Exhibit 5-6. Day opening tab



The appointment list includes a default time of 9:00 AM, SP ID, SP name, SP age, appointment status, and Rem:.

Select the **Month** tab to view a monthly calendar of all appointments by appointment type. Exhibit 5-7 is an example of the appointments by month screen. Use the links on the right and left side of the teal top border to go to the previous or next month.

Exhibit 5-7. Appointments by month screen

Field Follow-Up Management System - Version 9.2.0 You signed in as Barbara Lindstrom [Log Out](#)

Dietary Hepatitis C STD FFQ PAM PSA Allergy ER ROF HUC

Contact Manager Appointment Manager Call History Reports Supervisor

HUC Appointment Manager

Filter Criteria (Optional)

Interviewer:
 Appt. Status:

Day Week **Month** New Appt. Reschedule Appt. Monitor Not Sched.

Select A Month to View Appts, Day to View Details

Jun	July 2009						Aug
Sun	Mon	Tue	Wed	Thu	Fri	Sat	
			1 9:00 AM Mary Smith 1 9:00 AM Mary Smith 2	2 9:00 AM Mary Smith 3	3 9:00 AM Mary Smith 5	4 9:00 AM Mary Smith 6	

Names are displayed in color and the color represents their appointment status. The appointment status legend (Exhibit 5-8) is listed below.

Exhibit 5-8. Appointment status legend

Green	Scheduled, Firm Appt. Mailperson scanned the kit into their system showing that it was picked up. Can also include cases where the call history has been updated. Processed through Sort Facility (USPS data)
Black	Scheduled, Tentative Appt. SP said "Yes" in the MEC exam. Can also include cases where the call history has been updated. Electronic Shipping Info Received (USPS data)
Red	Not Scheduled, Broken Appt. 14 days have elapsed since the SP was initially given the kit and there has been no tracking activity. Can also include cases where the call history has been updated.
Blue	Appt. Completed USPS system marks the air bill as "received" meaning they delivered the kit. Delivered (USPS data)
Brown	Appt. Final-Nonresponse 30 days have elapsed since the SP was initially given the kit and there has been no tracking activity.

5.6 Call History

Document all contacts with an SP. To access the call history, click on his or her underlined name on the monthly appointment schedule or select the **Call History** button on the top of the module. If the link on the monthly calendar in Appointment Manager is selected, then the Call History screen for that SP displays.

If the **Call History** button on the top of the module is selected, then only blank text fields for SP ID, Last Name, and First name display. Enter the relevant information and select the **Search** button to access the call history for that SP. Exhibit 5-9 is an example of the HUC Call History search screen.

Exhibit 5-9. HUC Call History search screen

The screenshot displays the 'Field Follow-Up Management System - Version 9.2.0' interface. At the top, it indicates the user is signed in as 'Barbara Lindstrom' with a 'Log Out' link. A navigation bar contains tabs for 'Dietary', 'Hepatitis C', 'STD', 'FFQ', 'PAM', 'PSA', 'Allergy', 'ER', 'ROF', and 'HUC'. Below this, a secondary navigation bar features buttons for 'Contact Manager', 'Appointment Manager', 'Call History', 'Reports', and 'Supervisor'. The main heading is 'HUC Call History'. A search form titled 'Enter Search Criteria' includes input fields for 'SP ID:', 'Last Name:', 'First Name:', 'Suffix:', and 'Stand:', along with 'Search' and 'Clear' buttons.

This module tracks the date and time of each call and the contact code if contact was unsuccessful.

Once the search criteria have been successfully entered and the search executed, the Call History screen displays. The screen displays the history of the calls for the SP identified on the **Call Information** screen. An example of the Call History screen is illustrated in Exhibit 5-10.

Exhibit 5-10. Call History screen

HUC Call History

Enter Search Criteria

SP ID: 687317 Last Name: PEREZ First Name: 95F Suffix: Stand: 615

Current Appointment: Sunday - 2/24/2013 9:00 AM Initial Exam: Tuesday - 11/16/2010

Current Appt Status: Appt. Final-Non-Response, Updated on 7/23/2013

Age: 50 Gender: Female Language: English Interpreter Required: No Send Materials in: ENGLISH

Phone No: (242)342-3423 Location: Home Comment: SP Will Call In:

Attempt No.	Date/Time	Contact Made	Contact Status	Other, Reason	Rescheduled
1	7/23/2013 9:41 AM	No	Reminder Postcard Sent		

Record Call History for 95F PEREZ (6150182)

Appointment Date/Time: 2/24/2013 9:00 AM

Contact Made with SP: Yes No

Reschedule Appointment: Yes No

Contact Status:

Other, Reason:

The features of this screen include:

- **Search Criteria bar.** Use this feature to access information about another SP. Enter the SP ID or the first and last name in the search criteria boxes at the top of the screen.
- **Contact Information box.** This top box provides contact information specific to the SP. It includes current appointment date and time, initial exam date, current appointment status and last update date, phone number, location, comments, age, gender, language, and appointment status for all components (C = complete, S = scheduled [either tentative or firm], B = broken, FNR = Final non-response).
- **Record Call History table.** This table lists the history of call attempts for the SP. Record all contacts by marking **Contact Made with SP**, Yes or No, and **Reschedule Appointment**, Yes or No, and recording **Contact Status** and **Other, Reasons**.
- **Update Call Info button.** This button updates the contact status information, after it has been entered into the FFMS (described below).

Documenting these status codes allows all interviewers to see what has transpired with previous call attempts to a particular SP. Assign a contact status for every attempt to contact the SP. Two contacts were made to the particular SP shown in Exhibit 5-11.

Exhibit 5-11. Example of contact attempts

Attempt No.	Date/Time	Contact Made	Contact Status	Other, Reason	Rescheduled
1	6/22/2009 6:12 PM	No	Reminder Postcard Sent		
2	7/7/2009 1:02 PM	No	Other	USPS envelope sent	

The contact status will automatically update with the date and time the reminder postcard was sent once the SP's postcard has been marked as sent.

Select a **Contact Status** from the list of possible choices. Exhibit 5-12 illustrates the various codes.

Exhibit 5-12. Contact status codes

Appointment Date/Time: 6/18/2009 9:00 AM

Contact Made with SP: Yes No

Reschedule Appointment Yes No

Contact Status:


Other, Reason:

- Reminder Postcard Sent
- Replacement Kit Sent
- Other
- Appt. Rescheduled (DO NOT USE)

The contact codes are described below.

- **Reminder Postcard Sent** – The contact status will automatically update with the date and time the reminder postcard was sent once the SP’s postcard has been marked as sent.
- **Replacement Kit Sent** – Select this option when a replacement kit has been sent.
- **Other** – If none of the above contact codes applies, select “Other” and type in the reason.
- **Appt. Rescheduled (DO NOT USE)** – This choice is not applicable to the HUC.

After entering the information, select the **Update Call Info** button before exiting the module.




The button is highlighted in pink and is located above the updated information.

5.7 Reports

Use the **Reports** functionality to print labels for the reminder postcards, run the HUC Response Rate Report, and investigate the status of a particular payment check. Exhibit 5-13 illustrates the Reports screen and links.

Exhibit 5-13. HUC Reports screen and links



To access one of the three options, click on the desired link.

5.7.1 Reminder Postcard

At least daily, access the Reminder Postcard list, print labels, and mail the postcards. Select the [Reminder Postcard](#) link to access this functionality. Exhibit 5-14 illustrates the HUC Reminder postcard screen.

Exhibit 5-14. HUC Reminder Postcard screen

HUC Reports

Reminder Postcard

Print Label(s)
Mark Postcards as Mailed
 Print All

SP ID	Name	Age	Stand	Initial Exam	Language	Generate Label
184869	ERIC JONES	41	410		English	<input type="checkbox"/>
243260	JASON BALK	60	410		English	<input type="checkbox"/>
334755	JOSEPH OT	10	410		English	<input type="checkbox"/>
360722	CORY DEWICK	49	410		English	<input type="checkbox"/>
381736	JENNIFER DELACH	59	410		English	<input type="checkbox"/>
386651	HOPE CARROLL	14	410		English	<input type="checkbox"/>
392089	MICHAEL COSTANZO	42	410		English	<input type="checkbox"/>
411725	STEPHEN BREYER	33	410		English	<input type="checkbox"/>
437875	REUBEN OT	10	410		English	<input type="checkbox"/>
534770	BRIAN CATANZERITE	8	410		English	<input type="checkbox"/>
603252	JENNY BLOCK	42	410		English	<input type="checkbox"/>

The **Reminder Postcard** screen includes buttons for **Print Label(s)**, and **Mark Postcards as Mailed**, a check box for **Print All**, and columns for SP ID, Name (first, last), Age, Stand number, Initial Exam date, Language, and a checkbox for Generate Label.


To individually mark each SP, check each **Generate Label** checkbox or, to print labels for all SPs in the list, select the **Print All** checkbox.

Use the language displayed for each SP to select the correct English or Spanish postcard. Affix the label to the postcard and mail the postcard through Westat's routine mail system.

5.7.2 HUC Response Rate Report

Select the **HUC Response Rate Report** link to access this functionality. The first screen is illustrated in Exhibit 5-15.

Exhibit 5-15. Response Rate Report first screen



HUC Reports

Home Urine Collection Response Rate Report


Select Stand(s):

- 245 (10/15/2008 to 11/22/2008)
- 246 (11/06/2008 to 12/15/2008)
- 247 (11/30/2008 to 01/26/2009)
- 248 (01/11/2009 to 02/15/2009)

Filter

The Response Rate Report includes a **Select Stand(s)** filter. Select or highlight the stand or stands. To select more than one stand, hold down the **Ctrl** key and click on the desired stands. Exhibit 5-16 illustrates an example where multiple stands were selected.

Exhibit 5-16. Response Rate Report multiple stands selected



Home Urine Collection Response Rate Report

Select Stand(s):

- 256 (08/02/2009 to 09/06/2009)
- 257 (08/22/2009 to 09/27/2009)
- 258 (09/20/2009 to 10/19/2009)
- 259 (10/10/2009 to 11/09/2009)

Filter

Once the stands have been selected, click on the **Filter** button to generate the report.

Exhibit 5-17 is an example of a HUC Response Rate Report.

Exhibit 5-17. HUC Response Rate Report example

	All Ages		6 - 11 yrs		12 - 19 yrs		20 - 49 yrs		50 - 69 yrs		70 - 74 yrs		75+ yrs	
	#	%	#	%	#	%	#	%	#	%	#	%	#	%
SPs Eligible for MEC Urine Collection	261		54		43		113		47		4		4	
SPs Eligible - no urine collected	3	1.1	0	0	1	2.3	0	0	2	4.3	0	0	0	0
SPs Eligible - no vessel 45 collected	0	0	0	0	0	0	0	0	0	0	0	0	0	0
SPs Eligible for Home Urine Collection	258	98.9	54	100	42	97.7	113	100	45	95.7	4	100	4	100
SPs with Not Done HUC Status	49	19	11	20.4	7	16.7	18	15.9	11	24.4	2	50	2	50
SPs Provided Home Urine Collection Kit	209	81	43	79.6	35	83.3	95	84.1	34	75.6	2	50	2	50
Samples Outstanding <= 10 Days	39	18.7	8	18.6	8	22.9	20	21.1	3	8.8	0	0	0	0
Samples Outstanding > 10 Days	5	2.4	0	0	2	5.7	3	3.2	0	0	0	0	0	0
Broken Appointments with Reminder Postcards Sent	0	0	0	0	0	0	0	0	0	0	0	0	0	0
Samples in Transit	14	6.7	4	9.3	2	5.7	8	8.4	0	0	0	0	0	0
Samples Received by Contract Lab	164	78.5	34	79.1	25	71.4	72	75.8	31	91.2	2	100	2	100
Number of SPs that returned the evening void	14	6.7	4	9.3	2	5.7	8	8.4	0	0	0	0	0	0
Number of SPs that returned 1 container for the morning void	14	6.7	4	9.3	2	5.7	8	8.4	0	0	0	0	0	0
Number of SPs that returned 2 container for the morning void	164	78.5	34	79.1	25	71.4	72	75.8	31	91.2	2	100	2	100
Samples Marked as Final Non-Response	1	0.5	1	2.3	0	0	0	0	0	0	0	0	0	0

The report breaks down the response rate by number and percent for the following age groups: All ages, 6-11, 12-19, 20-49, 50-69, 70-74, and 75+. It includes the following data:

- SP eligible for MEC urine collection by:
 - SPs eligible – no urine collected
 - SPs eligible – no vessel 45 collected
- SPs eligible for home urine collection (SP eligible for MEC Urine Collection minus SPs eligible – no urine collected and SPs eligible – no vessel 45 collected)
- SPs with Not Done HUC status
- SPs provided home urine collection kit
 - Samples outstanding ≤ to 10 days
 - Samples outstanding > 10 days
 - – Broken appointments with reminder postcard sent
 - Samples in transit
 - Samples received by contract lab
 - The number of SPs that return the evening void.
 - The number of SPs that return 1 container for their morning void.
 - The number of SPs that return 2 containers for their morning void.
 - Samples marked as Final Non-Response

5.7.3 Check Status Report

Use the **Check Status Report** functionality to check the status of HUC payments.

Select the **Check Status Report** link to access this functionality. This search screen is illustrated in Exhibit 5-18.

Exhibit 5-18. Check Status Report search screen

Enter the SP ID or the first and last name in the search criteria boxes at the top of the screen. Select the **Search** button to initiate the search. Select the **Clear** button to clear the fields.

The Report displays after the search button has been selected. An example is illustrated in Exhibit 5-19.

Exhibit 5-19. Check Status Report example


SP ID	Name	Component	Appt Status	Check Date	Batch No	Check Amount
111111	MARY SMITH1	HUC	Appt. Completed	06/15/09	100058	\$40.00
111111	MARY SMITH1	Dietary	Appt. Completed	07/09/09	101173	\$30.00
111111	MARY SMITH1	ER	Not Scheduled			
111111	MARY SMITH1	FCBS	Appt. Completed	07/09/09	10757	\$15.00

The screen includes the SP ID, Name, Component, Appt Status, Check Date, Batch No, and Check Amount.

5.8 Supervisor

Use the **Supervisor** functionality to (1) view a report that lists SPs by appointment status and stand and (2) assign employees to the HUC component within the FFMS. There are three links: **Appointment Monitoring Report**, **Interviewer Scheduler**, and **HUC Module Security**. The Interview Scheduler link is not applicable to the HUC component. The main screen is illustrated in Exhibit 5-20.

Exhibit 5-20. Supervisor main screen



To access one of the two applicable options, click on the desired link.

5.8.1 Appointment Monitoring Report

Use the **Appointment Monitoring Report** to view a report that lists SPs by appointment status and stand. To access the module, select the link. The search criteria screen displays and it is illustrated in Exhibit 5-21.

Exhibit 5-21. Appointment Monitoring Report – select criteria

The screenshot shows a web interface for a HUC Supervisor. At the top, there is a purple header with the text "HUC Supervisor". Below this is a white header with the text "Appointment Monitoring Report". The main content area is a form titled "Select Search Criteria" with a blue header. The form has three columns: "SP ID:", "APPOINTMENT STATUS:", and "STAND(S):". The "SP ID:" column contains a text input field. The "APPOINTMENT STATUS:" column contains five radio button options: "Scheduled", "All Broken Appts.", "Complete Appts.", "All Final, Non-Response Appts.", and "Not Scheduled". The "STAND(S):" column contains a list box with four items: "Year 10 Stands", "Year 11 Stands", "Year 12 Stands", and "Year 13 Stands". At the top right of the form, there are "Search" and "Clear" buttons.

Select the **APPOINTMENT STATUS** by clicking on the desired radio button and then select or highlight the stand or stands. To select more than one stand, hold down the **Ctrl** key and click on the desired stands. When finished, select the **Filter** button to run the report.

The appointment status for all SPs who meet the criteria displays with each SP contained in a separate box. Exhibit 5-22 is an example of an appointment monitoring report.

Exhibit 5-22. Appointment Monitoring Report example

Field Follow-Up Management System - Version 10.4.0 You signed in as Pat Soto [Log Out](#)

Dietary Hepatitis C FFQ PAM PSA FCBS HUC 24 UC

Contact Manager Appointment Manager Call History Reports Supervisor

HUC Supervisor

Appointment Monitoring Report for Stand(s): 615

Select Search Criteria Search Clear

SP ID:	APPOINTMENT STATUS:	STAND(S):
<input type="text"/>	<input type="radio"/> Scheduled <input type="radio"/> All Broken Appts. <input type="radio"/> Complete Appts. <input checked="" type="radio"/> All Final, Non-Response Appts. <input type="radio"/> Not Scheduled	Year 10 Stands Year 11 Stands Year 12 Stands Year 13 Stands

Update Changes

SP ID	Name	Participant No.	Age	Initial Exam	Current Appt.	Stand ID
687317	95F PEREZ	615-20-15-1-1	50	11/16/2010	2/24/2013 9:00:00 AM	615

Appt. Status: Reason: Other Text:

Attempt No.	Date/Time	Contact Made	Contact Status/Description	Contact By	Rescheduled
1	7/23/2013 9:41 AM	No	Reminder Postcard Sent	Stephen Bernas	

The criteria box remains on the top of the screen. The teal header includes the SP ID, Name, age, initial exam date, current appointment date, and stand number. The lower section includes text boxes for Appt Status, Reason, and Other Text.

Change the appointment status, reason, and/or enter free text. Exhibit 5-23 shows the appointment status change options.

Exhibit 5-23. Appointment Status change options

Appt. Status:

- Scheduled, Firm Appt.
- Scheduled, Tentative Appt.
- Not Scheduled
- Not Scheduled, Broken Appt.
- Appt. Completed
- Appt. Final-Non-Response

To change the appointment status, select or highlight the correct choice from the drop-down list.

The default reason for HUC exams coded as final non-response is **SP Refusal**. Do not change this to any of the other options. Exhibit 5-24 illustrates the various reason options.

Exhibit 5-24. Appointment Status reason options

Reason: SP Refusal

- SP Refusal
- Language Problem
- Hearing Impaired
- Unable to contact SP
- End of Feasibility
- Other

Any **Call History** attempts display under the lower section as illustrated by Exhibit 5-25.

Exhibit 5-25. Call History displays at bottom of SPs window

Field Follow-Up Management System - Version 10.4.0 You signed in as Pat Soto [Log Out](#)

Dietary Hepatitis C FFQ PAM PSA FCBS **HUC** 24 UC

[Contact Manager](#) [Appointment Manager](#) [Call History](#) [Reports](#) [Supervisor](#)

HUC Supervisor

Appointment Monitoring Report for Stand(s): 615

Select Search Criteria Search Clear

SP ID:	APPOINTMENT STATUS: <input type="radio"/> Scheduled <input type="radio"/> All Broken Appts. <input type="radio"/> Complete Appts. <input checked="" type="radio"/> All Final, Non-Response Appts. <input type="radio"/> Not Scheduled	STAND(S): Year 10 Stands Year 11 Stands Year 12 Stands Year 13 Stands
--------	---	--

Update Changes

SP ID	Name	Participant No.	Age	Initial Exam	Current Appt.	Stand ID
687317	95F PEREZ	615-20-15-1-1	50	11/16/2010	2/24/2013 9:00:00 AM	615
Appt. Status:		Appt. Final-Non-Response		Reason: SP Refusal		Other Text:
Attempt No.	Date/Time	Contact Made	Contact Status/Description	Contact By	Rescheduled	
1	7/23/2013 9:41 AM	No	Reminder Postcard Sent	Stephen Bernas		

Select the **Update Changes** button to save the changes to the database.

Update Changes

5.8.2 HUC Module Security

Use the **HUC Module Security** link to assign employees to the HUC component within the FFMS and to give them permission to access the various modules within the HUC component. To access the module, select the link. Exhibit 5-26 illustrates the HUC Module security screen.


Exhibit 5-26. HUC Module security screen

The screenshot displays the 'HUC Module Security' interface. At the top, the system title is 'Field Follow-Up Management System - Version 9.2.0' and the user is signed in as 'Barbara Lindstrom'. The navigation bar includes modules like 'Dietary', 'Hepatitis C', 'STD', 'FFQ', 'PAM', 'PSA', 'Allergy', 'ER', 'ROF', and 'HUC'. Below this, there are buttons for 'Contact Manager', 'Appointment Manager', 'Call History', 'Reports', and 'Supervisor'. The main content area is titled 'HUC Module Security' and contains two columns of employee names. The left column, 'Assign Emp to HUC', lists: Apodaca, Kay; Arevalo, Maria; Ayala, Yanira; Benson, Robert; Binzer, Greg; Caceres, Becky; and Carrillo, Washington. The right column, 'Employees Assigned to HUC', lists: Ayala, Sandra; Bernas, Stephen; Bulaong, Julia; Glascoe, Terri; Hall, Janice; Hillard, Debbie; and Lash, Sharon. A red arrow points from the left to the right, and another points from the right to the left. Below the left list is a checkbox for 'Add as Supervisor'. Below the right list is a note: '* Employees with HUC as Default'. There are also buttons for 'Set As Default Module' and 'Assign Menu's'.

Two employee text boxes display. The left text box includes the names of employees who are eligible to be assigned to the HUC module. The right text box includes the names of employees who have been assigned to the HUC FFMS module. Use the bright red arrows to move the employees between the two boxes.

Assign specific menus to specific individuals who have been assigned to the module. This screen is illustrated in Exhibit 5-27.

Exhibit 5-27. HUC Module security - assign menus



To assign specific menus to specific individuals select or highlight the employee from the right box and then select the **Assign Menu's** button. The menu options text boxes display. These additional text boxes are illustrated in Exhibit 5-28.

Exhibit 5-28. HUC Module security - available menus

HUC Module Security

Employees:

Assign Emp to HUC:

- Apodaca, Kay
- Arevalo, Maria
- Ayala, Yanira
- Benson, Robert
- Binzer, Greg
- Caceres, Becky
- Carrillo, Washington

Add as Supervisor

Employees Assigned to HUC:

- Ayala, Sandra
- Bernas, Stephen
- Bulaong, Julia
- Glascoc, Terri
- Hall, Janice
- Hillard, Debbie
- Lash, Sharon

* Employees with HUC as Default

Buttons: Set As Default Module, Assign Menu's

HUC Menus for Lindstrom, Barbara:

Available Menus:

Menus Assigned to Lindstrom, Barbara:

- Appointment Manager *
- Call History
- Contact Manager
- Reports
- Supervisor

To move one or more of the menu options from one box to the other, click on the bright red arrow.

To designate someone as a supervisor, highlight the name and add a checkmark to the **Add as Supervisor** checkbox. This option is illustrated in Exhibit 5-29.

Exhibit 5-29. HUC Security Module – assign supervisor status

Employees:

Assign Emp to HUC:

- Apodaca, Kay
- Arevalo, Maria
- Ayala, Yanira
- Benson, Robert
- Binzer, Greg
- Caceres, Becky
- Carrillo, Washington

Employees Assigned to HUC:

- Ayala, Sandra
- Bernas, Stephen
- Bulaong, Julia
- Glascoe, Terri
- Hall, Janice
- Hillard, Debbie
- Lash, Sharon

Add as Supervisor

* Employees with HUC as Default

Buttons: Set As Default Module, Assign Menu's

Use the **Set as Default Module** button to set the module that is displayed when the application opens. However, if the application is opened in the same browser window that it was opened in previously, meaning logged out and back in without closing the window, then the application remembers the module that was displayed when logging out regardless of the default.

6. CHECK PROTOCOL

6.1 Overview

Remuneration checks for \$50 are generated once the USPS scans one of the three kits associated with an SP into the USPS system. The PF appointment is “Scheduled, Firm.” The ISIS check application is used to send check information to Westat accounting for the generation of checks to respondents who sent in their urine sample.

The check application generates a unique batch ID each time it generates a group of checks. The check extraction application first determines the list of participants who need checks. Then the check application retrieves the name of the participant from the nh_person table along with the address. The check is mailed from the same address sources used by the HUC kit generation page. This information is placed in the Access database specified by Accounting and emailed to them. Accounting then cuts the necessary checks. The SP ID and check batch number appear in the memo field of each check for tracking purposes.

The email to Accounting is generated and checks are mailed twice a week, on Tuesday and Friday. The home office dietary phone staff are responsible for mailing the checks. They place the check in an envelope with a letter (Exhibits 6-1 and 6-2) and mail the checks.

Exhibit 6-1. English letter – check enclosure

Dear Participant:

Thank you for participating in the National Health and Nutrition Examination Survey (NHANES) by sending us your urine collection kits.

Enclosed is a check for \$50 as a token of appreciation. If additional family members also received urine kits, they will receive their check once we receive their kits.

Thank you once again for taking part in this important study. If you have any questions regarding this payment, please call **1-888-458-4762** Monday through Friday between the hours of 8:00 AM and 4:30 PM Eastern Time.

Sincerely,

Sarah C. Zweig
Study Manager

Enclosure

Exhibit 6-2. Spanish letter – check enclosure

Estimado(a) participante:

Gracias por participar en la Encuesta Nacional de Examen de Salud y Nutrición (NHANES) al mandarnos los materiales de su muestra de orina.

Adjunto hay un cheque por \$50 dólares como muestra de agradecimiento. Si otros familiares también recibieron materiales para una muestra de orina, ellos recibirán su cheque una vez que nosotros recibamos sus materiales.

Gracias nuevamente por participar en este importante estudio. Si desea hacer alguna pregunta relacionada con este pago, por favor llame al **1-888-458-4762** de lunes a viernes entre las 8:00 de la mañana y las 4:30 de la tarde, hora del Este.

Atentamente,

Sarah C. Zweig
Administradora del estudio

Se adjunta cheque

7. TRACKING AND DELIVERY STATUS

7.1 Overview


Tracking emails are sent daily to home office staff (Exhibit 7-1). HUC shipments are sent by USPS Priority Mail. If one of the kits belonging to an SP is scanned by the USPS, all three HUC tracking numbers linked to that SP show up in the tracking report. The “Delivery Status” column indicates the latest status of the package. The home office component person investigates all shipments where the “No (number) of Days Pending” is greater than 8 days. He or she emails the analysis lab and includes the SP ID and number of days pending. If the analysis lab indicates that they received the sample, then the home office system analysis is notified and the database is updated accordingly.

Exhibit 7-1. Sample tracking email

From: FedexTracking@westat.com						Sent: Mon 7/22/2013 12:33
To: Pat Soto; Kimberly Sun; Stephen Bernas; Lakshmi Venkataraman; Debbie Hillard; Renee Storandt; Sarah Zweig; Ken Giffin						
Cc:						
Subject: NHANES Airbill tracking						
Lab airbills in transit:						
Lab						
Airbill No	Lab ID	Shipped DT	Container/SP ID	No of Days Since Last Status	e-mail address	Delivery Status
420554559405501699320003343091	16	7/18/2013	123456159 / 123456	4	microalb@umn.edu	Processed - Dispatched to Sort Facility
420554559405501699320003343114	16	7/18/2013	123456158 / 123456	4	microalb@umn.edu	Processed - Dispatched to Sort Facility
420554559405501699320003343121	16	7/18/2013	123456157 / 123456	4	microalb@umn.edu	Processed - Dispatched to Sort Facility
420554559405501699320003349581	16	7/22/2013	987654157 / 987654	0	microalb@umn.edu	
420554559405501699320003349611	16	7/22/2013	987654158 / 987654	0	microalb@umn.edu	
420554559405501699320003349628	16	7/22/2013	987654159 / 987654	0	microalb@umn.edu	

Selecting any tracking number will link you to the USPS web site where track and confirm information is available. Exhibits 7-2 through 7-4 illustrate three examples.

Exhibit 7-2. Example Status Dispatched to Sort Facility



Search USPS.com or Track Packages

Quick Tools | Ship a Package | Send Mail | Manage Your Mail | Shop | Business Solutions

Track & Confirm

GET EMAIL UPDATES | PRINT DETAILS



YOUR LABEL NUMBER	SERVICE	STATUS OF YOUR ITEM	DATE & TIME	LOCATION	FEATURES
420554559405501699320003343121 	Priority Mail®	Dispatched to Sort Facility	July 18, 2013, 4:50 pm	SARATOGA SPRINGS, NY 12866	Expected Delivery By: July 20, 2013 USPS Tracking / Delivery Confirmation™
		Acceptance	July 18, 2013, 12:16 pm	SARATOGA SPRINGS, NY 12866	

Exhibit 7-3. Example Status Depart USPS Sort Facility



Search USPS.com or Track Packages

Quick Tools | Ship a Package | Send Mail | Manage Your Mail | Shop | Business Solutions

Track & Confirm

GET EMAIL UPDATES | PRINT DETAILS


YOUR LABEL NUMBER	SERVICE	STATUS OF YOUR ITEM	DATE & TIME	LOCATION	FEATURES
420554559405501699320005006376 	Priority Mail®	Depart USPS Sort Facility	July 22, 2013	KANSAS CITY, MO 64121	Expected Delivery By: July 24, 2013 USPS Tracking / Delivery Confirmation™
		Processed through USPS Sort Facility	July 22, 2013, 5:36 pm	KANSAS CITY, MO 64121	
		Acceptance	July 22, 2013, 10:28 am	OLATHE, KS 66061	
		Electronic Shipping Info Received	June 26, 2013		

Exhibit 7-4. Example Status Delivered


The screenshot shows the USPS.COM website interface. At the top left is the USPS logo. To the right is a search bar with the text "Search USPS.com or Track Packages". Below the logo is a navigation bar with buttons for "Quick Tools", "Ship a Package", "Send Mail", "Manage Your Mail", "Shop", and "Business Solutions". The main heading is "Track & Confirm". Below this are two buttons: "GET EMAIL UPDATES" and "PRINT DETAILS". The tracking history is presented in a table with the following columns: YOUR LABEL NUMBER, SERVICE, STATUS OF YOUR ITEM, DATE & TIME, LOCATION, and FEATURES.

YOUR LABEL NUMBER	SERVICE	STATUS OF YOUR ITEM	DATE & TIME	LOCATION	FEATURES
420554559405501699320003346627 ⓘ	Priority Mail®	Delivered	July 19, 2013, 8:27 am	MINNEAPOLIS, MN 55455	Expected Delivery By: July 19, 2013 USPS Tracking / Delivery Confirmation™
		Out for Delivery	July 19, 2013, 8:18 am	MINNEAPOLIS, MN 55414	
		Sorting Complete	July 19, 2013, 8:08 am	MINNEAPOLIS, MN 55414	
		Arrival at Post Office	July 19, 2013, 6:41 am	MINNEAPOLIS, MN 55414	
		Depart USPS Sort Facility	July 18, 2013	SPRINGFIELD, MA 01152	
		Processed at USPS Origin Sort Facility	July 18, 2013, 12:43 am	SPRINGFIELD, MA 01152	
		Dispatched to Sort Facility	July 17, 2013, 5:11 pm	GANSEVOORT, NY 12831	
		Acceptance	July 17, 2013, 11:49 am	GANSEVOORT, NY 12831	

7.2 Email Sent Weekly to Analysis Lab

The analysis lab emails a file weekly. This file contains a list of every SP ID marked as received in the database in the past week. The contact lab staff use this file to reconcile the samples they received with the list of samples that were marked by the USPS as delivered. Exhibit 7-5 is an example of this email.

Exhibit 7-5. Example emailed file sent to analysis lab



No action is required.

Appendix A

“What To Do”

**Urine Collection and Shipping Instructions
English and Spanish**

“WHAT TO DO”

Urine Collection and Shipping Instructions

Collect one evening and one morning urine sample and ship it to the laboratory using these instructions.

Collect your urine samples as soon as possible - ideally within the next 14 days. However, if you cannot collect it within 14 days then please collect and ship it at your earliest convenience.

We will send you a check for \$50.00 as soon as you ship the packages.

The bag you received in the mobile examination center includes:

- Large tote bag with a confidentiality statement, your name written on the bag, and your personal picture sticker
- “WHAT TO DO” Urine Collection and Shipping Instructions sheet
- The home urine collection kits

Your Personal Picture Sticker:

Your picture sticker is unique to you and it is found on:

- The large tote bag with your name written on it
- Your urine collection cups
- The bottom of your shipping boxes

The picture sticker identifies your kits from kits assigned to other people in your household. When collecting and shipping your urine, make sure you use the collection cups and shipping boxes marked with your unique picture sticker.



Home Urine Collection Kit Contents:

- A.** Cardboard shipping box with postage-paid label
- B.** Styrofoam container inside the shipping box
- C.** Temperature monitoring strip
- D.** Ice packs
- E.** Plastic urine collection cup with your picture sticker
- F.** Small blue Ziploc bag
- G.** Larger clear Ziploc bag



Kit Preparation:

1. Remove the ice packs from your large tote bag and place them in the freezer compartment of your home freezer for at least 6 hours before you collect your samples. Freeze the ice packs flat.
2. When you are ready to collect your sample, remove the correct kit(s) from the large tote bag that has your name and picture symbol on it.
 - For the evening sample, remove the kit with a picture of a moon on the outside of the box.
 - For the morning sample, remove the two kits with a sun on the outside of the box.
3. Open the white cardboard box(es).
4. Remove the Styrofoam lid(s).
5. Remove the collection cup from the box.



Collect the Evening Urine Sample:

Evening Sample: Collect the sample between 5:30 p.m. and bedtime. Pick an evening when you will be able to collect the second sample on the following morning.

1. Before you collect your urine, verify that the picture sticker on the cup matches the picture sticker on the tote bag.
2. When collecting the evening sample, make sure you are using the evening cup with a picture of the moon on it.
3. Urinate (pee) into the collection cup until you have emptied your bladder.
4. Screw the cap on the urine cup tightly and wipe off the outside of the cup.
5. On the label located on the top of the Styrofoam container, write the date and time you collected the urine and answer the questions.



Collect the Morning Urine Sample:

Morning Sample: Collect one sample when you first wake up for the day. Collect it the morning after you collect the evening sample.

Whenever possible, mail all samples (morning and evening samples) on the same day you collect the morning sample.

1. Before you collect your urine, verify that the picture sticker on the cup matches the picture sticker on the tote bag that has your name written on it.
2. You will need the morning cups that have a picture of a sun on it.
3. Urinate (pee) into the collection cup until you have emptied your bladder. If you think you will have more urine than the cup can hold, stop midstream and finish the collection in the second cup morning cup.
4. Only collect one morning sample. If all of your urine fits in one cup, throw the second cup away.
5. Screw the cap on the urine cup tightly and wipe off the outside of the cup.
6. On the label located on the top of the lid of the container, write the date and



time you collected the urine and answer the questions. (If you needed to use two cups, mark "2 cups" on both labels.)

Pack the Urine Samples:

Pack each sample in its shipping container with two ice packs immediately after collection, even if you do not mail the package until later in the day.

Discard these instructions and the large bag with your name on it. Do not add additional packaging.

1. Before packing the samples, double-check the sun and moon symbols on the cup, Styrofoam container lid, and bottom of the shipping box to make sure they match.
2. If you used two morning cups, each cup will be packaged separately.
3. Place the urine cup into the small blue Ziploc bag.
4. Push the air out of the Ziploc bag and seal it.
5. Place the sealed Ziploc bag in the larger clear Ziploc bag. Push the air out of the bag and seal it.
6. Place the sample into the Styrofoam shipping box. Remember to match the moon or sun symbol on your cup with the moon or sun symbol on your box.
7. Place two of the frozen ice packs from your freezer into the box – one on the bottom and one on the side.

A white label with a sun symbol and text for recording collection date, time, and quantity.

Please write the date and time you collected the urine specimen.
Por favor anote la fecha y la hora en que tomo la muestra de orina

Date/Fecha 1/14/13

Time/Hora 7:30 AM PM

Morning collection - I am returning:
 1 cup only
 2 cups

All the urine is collected in the cup(s).
 Only a few drops are missing from the cup(s).
 More than a few drops are missing from the cup(s).

8. Place the Styrofoam lid onto the top of the container inside the box.
9. Verify that the picture sticker on the bottom of the shipping box matches the picture sticker on the tote bag that has your name written on it.

10. Tape the box closed with the tape provided.



Ship the Urine Sample:

Use the US Postal Service (USPS) to send the package.

Take the packages inside the Post Office and hand it to a teller.



Or, take the packages inside the Post Office and put it in the large priority mail drop box. (The large priority drop box is usually located near the self-service kiosk.)



If you have any questions call 1-888-458-4762.

“QUÉ HACER”

Instrucciones para recolección y envío de la muestra de orina

Recolecte una muestra de orina por la noche y una por la mañana y use estas instrucciones para enviarlas al laboratorio.

Tome sus muestras de orina tan pronto como sea posible – lo ideal es dentro de los próximos 14 días. Sin embargo, si no puede tomar las muestras dentro de los 14 días siguientes, entonces tómelas y envíelas tan pronto como pueda.

Le mandaremos un cheque por \$50.00 dólares tan pronto como envíe los paquetes.

La bolsa que recibió en el centro móvil de examen incluye:

- Una bolsa grande con una declaración de confidencialidad, su nombre escrito en la bolsa y una calcomanía impresa con un símbolo elegido únicamente para usted.
- Hoja “QUÉ HACER” Instrucciones para recolección y envío de las muestras de orina
- Materiales para tomar las muestras de orina en el hogar

Calcomanía con su símbolo personal:

El símbolo impreso en la calcomanía es únicamente para usted y se encuentra en:

- La bolsa grande que lleva su nombre escrito
- Sus envases para las muestras de orina
- La caja de envío para devolver los materiales

El símbolo distingue sus materiales de los materiales asignados a otras personas de su hogar. Cuando tome y envíe su muestra de orina, asegúrese de usar los envases para la muestra y las cajas marcadas con el símbolo asignado únicamente para usted.



Los materiales para tomar la muestra de orina en el hogar contienen:

- A.** Caja de cartón para envío, con etiqueta de franqueo prepagado
- B.** Caja de poliestireno dentro de la caja de envío
- C.** Cinta para monitorear la temperatura
- D.** Bolsas de hielo-gel refrigerante
- E.** Envase de plástico para poner la muestra, marcado con su símbolo.
- F.** Bolsa pequeña con cierre azul (Ziploc)
- G.** Bolsa grande transparente con cierre (Ziploc)



Preparación de los materiales:

1. Saque las bolsas de hielo-gel refrigerante de su bolsa grande y póngalas en el congelador de su casa por al menos 6 horas antes de recolectar la muestra. Las bolsas de hielo-gel refrigerante deben ponerse planas y horizontalmente en el congelador.
2. Cuando esté listo para recolectar su muestra, saque de la bolsa grande los materiales que llevan su nombre y su símbolo.
 - Para recolectar la muestra de orina en la noche, saque el envase que tiene la imagen de la luna afuera de la caja.
 - Para recolectar la muestra en la mañana, saque los dos envases que tienen la imagen del sol afuera de la caja.
3. Abra la(s) caja(s) de cartón blanco.
4. Saque la(s) tapa(s) de poliestireno.
5. Saque de la caja el envase para la muestra de orina.



Recolección de la muestra de orina en la noche:

Muestra de la noche: Recolecte la muestra entre las 5.30 de la tarde y la hora de dormir. Escoja una noche en la que pueda recolectar la segunda muestra a la mañana siguiente.

1. Antes de tomar la muestra de orina, verifique que el envase lleva el mismo símbolo de la bolsa.
2. Cuando recolecte la muestra de la noche, asegúrese de usar el envase para la noche el cual tiene una imagen de la luna.
3. Orine dentro del envase de plástico hasta que haya vaciado la vejiga.
4. Ajuste bien y apriete la tapa del envase de plástico y limpie la parte de afuera del envase.
5. En la etiqueta que está en la tapa de la caja de poliestireno, anote la fecha y la hora en que se tomó la muestra y conteste las preguntas.



Recolección de la muestra de orina en la mañana:

Muestra de la mañana: Recolecte una muestra cuando recién se levante. Recolecte la orina la mañana después de haber recolectado la muestra de la noche.

Cuando le sea posible, envíe sus muestras (la de la mañana y la de la noche) el mismo día que recolectó la muestra de la mañana.

1. Antes de recolectar la muestra, verifique que el símbolo del envase es igual al símbolo que aparece en la bolsa con su nombre.
2. Necesitará los envases de la mañana que tienen la figura del sol.
3. Orine dentro del envase de plástico hasta que haya vaciado la vejiga. Si tiene mucha orina y cree que el envase se podría desbordar, detenga el flujo de orina y use el segundo envase de la mañana para terminar de recolectar la muestra.
4. Recolecte únicamente una muestra. Si usa únicamente un envase para recolectar toda la orina, bote el otro envase.
5. Ajuste bien y apriete la tapa del envase de plástico y limpie la parte de afuera del envase.
6. En la etiqueta que está en la tapa del envase, anote la fecha y la hora en que recolectó la muestra y conteste las preguntas. (Si necesita usar los dos envases, escriba "2 cups" en ambas etiquetas.)



Empaque las muestras de orina:

Empaque cada muestra en su caja de envío con dos bolsas de hielo-gel refrigerante inmediatamente después de tomar la muestra, incluso si no las va a enviar hasta más tarde en el día.

Puede botar estas instrucciones y la bolsa grande con su nombre. No agregue material adicional de empaque.

1. Antes de empacar las muestras, vuelva a verificar que las imágenes del sol y la luna del envase son iguales a la de la tapa de la caja de poliestireno y a la de la caja de envío.
2. Si usó dos envases en la mañana, tendrá que empacar cada envase por separado.
3. Ponga el envase de plástico con la orina en la bolsa pequeña con cierre azul (Ziploc).
4. Saque el aire de la bolsa con cierre (Ziploc) y ciérrela bien.
5. Ponga la bolsa (Ziploc) cerrada dentro de la bolsa transparente Ziploc más grande. Saque el aire de la bolsa y ciérrela bien.
6. Ponga la muestra en la caja de envío de poliestireno. Recuerde verificar que las figuras del sol y la luna de los envases son iguales a las figuras del sol y la luna de su caja.
7. Ponga dos bolsas de hielo-gel de su congelador dentro de la caja de poliestireno – una al fondo de la caja y otra al lado.



Please write the date and time you collected the urine specimen.
Por favor anote la fecha y la hora en que tomó la muestra de orina.

Date/Fecha 1/14/13
Time/Hora 7:30 AM PM

Morning collection - I am returning:
 1 cup only
 2 cups

All the urine is collected in the cup(s).
 Only a few drops are missing from the cup(s).
 More than a few drops are missing from the cup(s).

8. Para tapar la caja de poliestireno que se encuentra dentro de la caja de envío, use la tapa de la caja de poliestireno.
9. Verifique que el símbolo que aparece en el fondo de la caja de envío es igual al símbolo de la bolsa grande que tiene su nombre escrito.

10. Use la cinta adhesiva que le dimos para cerrar la caja de envío.



Envíe las muestras de orina:

Use el Servicio Postal de Estados Unidos (USPS) para mandar el paquete.

Lleve los paquetes a la Oficina de Correos y entréguelos a un empleado.



O, lleve los paquetes dentro de la Oficina de Correos y póngalos en el buzón para paquetes grandes y envíos prioritarios. (El buzón para paquetes grandes y envíos prioritarios está habitualmente ubicado cerca del quiosco de autoservicio).




Si desea hacer alguna pregunta llame al 1-888-458-4762.

Appendix B

Kit Labeling Instructions

**Participant Unique Symbol:
(Located on cup, bottom of shipping box, tote bag)**



**Evening Barcode:
(Located on evening cup)**



**Morning Barcode:
(Located on morning cups)**



**Evening Kit Moon Sticker:
(Located on cup, top of Styrofoam lid,
bottom of shipping box)**



**Morning Kit Sun Sticker:
(Located on cup, top of Styrofoam lid,
bottom of shipping box)**



LABELING THE KIT

The kits will be prelabeled with sun and moon stickers. Before beginning the exam, you will need one evening kit and two morning kits. The printer will print 3 barcode stickers and 7 SP specific picture stickers.

- Open each urine collection kit. Confirm that the temperature monitoring strip is affixed to the inside lid of the Styrofoam shipping container.
- On each cup, apply an SP picture sticker and a Dymo barcode label to the center of each specimen cup. Place the stickers side-by-side. [REMEMBER TO MATCH THE SUN AND MOON STICKERS.]
- On the bottom of each shipping box, apply an SP picture sticker.
- On the large plastic tote bag:
 - Apply the appropriate language confidentiality label to the large plastic tote bag if necessary.
 - Write the SP's name with permanent black marker.
NOTE: This is the only location on the kit where the SP name is identified.
 - Apply one SP picture sticker beside the confidentiality label.

Appendix C

Talking Points
English and Spanish

Talking Points

Show the SP the “WHAT TO DO, Urine Collection and Shipping Instructions.”

- Collection and shipping instructions are included in the kit you are taking home.
- If you have any questions, call the toll-free telephone number on the instructions.

General Overview:

- Collect two urine samples – **one** in the evening between 5:30 p.m. and bedtime and **one** the **following** morning.
- Collect your urine sample as soon as possible and ideally within the next 14 days. However, if you cannot collect it within 14 days then please collect and ship it at your earliest convenience. [EXAMINER: ASK FEMALE SPS NOT TO COLLECT THE SAMPLES DURING HER PERIOD.]
- Whenever possible, collect the evening sample when you will be able to mail it the next day along with the morning sample.
- When you get home, place the ice packs in the freezer overnight or for at least 6 hours. Freeze the ice packs flat.

Label the parts of the kit and show the SP their unique picture sticker.

- This is your special sticker. If other members in your home participate, they will have a different sticker.

Show the SP the difference between the evening and morning cups.

Evening Collection:

- Urinate (pee) once into the cup with the sticker of the moon in the evening between 5:30 and bedtime.
- Urinate a full void (until your bladder is empty).
- Screw the lid on the specimen cup **tightly** and package it in the shipping box immediately.
- Place the cup in the small blue Ziploc bag. Squeeze the air out of the bag and **seal** it.
- Put the **sealed** Ziploc bag in the larger clear Ziploc bag. Squeeze the air out of the bag and **seal** it.
- Put two of the frozen ice packs in the box – one on the bottom and one on the side.
- Put the sample in the Styrofoam box with the moon symbol.
- Put the Styrofoam lid on the package.
- On the label on top of the Styrofoam container, write the date and time and answer the questions. [EXAMINER: GO OVER THE LABEL WITH THE SP.]
- Close the box with the tape provided.
- Store the urine inside the Styrofoam container with the frozen ice pack overnight until you collect the specimen in the morning.

Talking Points

Morning Collection:

- The following morning, when you first wake up, urinate (pee) into the cup with the sticker of the sun.
- Urinate in the cup until your bladder is empty. You are only collecting **one** sample that contains a full morning void. If that one sample fits in one cup, throw the other kit away. However, if you fill the first cup and you still need to go more, stop mid-stream and finish in the second morning cup.
- If you do not use the second cup with a sun on it, throw it away. Do not mail back an empty cup.
- Screw the lid on the specimen cup **tightly**.
- If you used two cups, package each cup separately.
- EXAMINER: REVIEW PACKAGING INSTRUCTIONS
- EXAMINER: GO OVER THE LABEL WITH THE SP.
 - If you use two morning cups, mark “**2 cups**” on both boxes.

Mail the samples:

Take the packages inside the Post Office and hand it to a teller or, take the packages inside the Post Office and put it in the large priority mail drop box. (The large priority drop box is usually located near the self-serve kiosk.)

If the urine specimen has not been received by the laboratory in 14 days, then a reminder postcard will be sent to you. If you have already sent your urine then ignore this postcard.

Do **not** mail the large outer plastic bag with your name on it to the laboratory.

We will send you a check for \$50.00 as soon as you ship the package.

Talking Points

Muéstrele a la persona a quien le da la muestra, las páginas “QUÉ HACER, Instrucciones para recolección y envío de muestras de orina”.

- Las instrucciones para la recolección y el envío de las muestras están incluidas en el paquete de materiales que usted se llevará a su hogar.
- Si tiene alguna pregunta, llame a la línea gratuita que se encuentra en las instrucciones.

Información general:

- Recolecte dos muestras de orina – **una** en la noche entre las 5:30 de la tarde y la hora de dormir, y la **otra** a la mañana **siguiente**.
- Recolecte la muestra de orina tan pronto como sea posible; lo ideal es dentro de los próximos 14 días. Sin embargo, si no puede tomar las muestras dentro de los siguientes 14 días, entonces tómelas y envíelas tan pronto como pueda. [EXAMINADOR: PIDA A LAS MUJERES PARTICIPANTES NO RECOLECTAR LA MUESTRA DURANTE SU PERIODO.]
- Cuando sea posible, recolecte la muestra de la noche únicamente cuando pueda enviarla al día siguiente junto con la muestra de la mañana.
- Cuando llegue a casa, ponga las bolsas de hielo gel refrigerante en el congelador de su casa por la noche o al menos 6 horas antes de recolectar la muestra. Congele las bolsas de hielo gel refrigerante planas y horizontalmente.

Marque las partes del paquete de materiales para tomar las muestras y muéstrele a la persona a quien da la muestra (SP) la etiqueta con su imagen única.

- Este es su símbolo especial. Si hay otros miembros de su hogar que también están participando, tendrán otro símbolo único para cada uno de ellos.

Muéstrele a la persona a quien da la muestra las diferencias entre los envases de recolección de la mañana y los de la noche.

Muestra de la noche:

- Orine una vez en el envase con el símbolo de la luna en la noche entre las 5:30 de la tarde y la hora de dormir
- Orine hasta que haya vaciado la vejiga.
- Cierre y apriete **bien** la tapa del envase y empáquela inmediatamente en la caja para el envío.
- Ponga el envase en la bolsa pequeña con cierre azul (Ziploc). Saque el aire de la bolsa y **ciérrela bien**.
- Ponga la bolsa **cerrada** (Ziploc) dentro de la bolsa Ziploc más grande. Saque el aire de la bolsa y **ciérrela** bien.
- Ponga dos bolsas de hielo gel refrigerante en la caja – una al fondo y otra al lado.
- Ponga la muestra dentro la caja de envío de poliestireno con el símbolo de la luna.
- Tape la caja de poliestireno.
- Apunte la fecha, hora y conteste las preguntas de la etiqueta que está encima de la caja de poliestireno. [EXAMINADOR: EXPLÍQUELE A LA PERSONA A QUIEN DA LA MUESTRA LA INFORMACIÓN QUE TIENE QUE PONER EN LA ETIQUETA]
- Use la cinta adhesiva que le dimos para cerrar la caja.
- Guarde la muestra de orina junto con la bolsa de hielo gel refrigerante dentro de la caja de poliestireno durante la noche, hasta que recolecte la muestra de orina de la mañana.

Talking Points

Muestra de la mañana:

- A la mañana siguiente, cuando recién se levante, orine en el envase con el símbolo del sol.
- Orine hasta que haya vaciado la vejiga. Debe recolectar únicamente **una** muestra que contenga todo el contenido de orina de su vejiga. Si la muestra cabe en un solo envase, bote el otro envase. Sin embargo, si llena el primer envase, pero todavía tiene más orina, detenga el flujo de orina y use el segundo envase de la mañana para terminar de recolectar la orina.
- Si no usa el segundo envase con el símbolo del sol, bótelo. No envíe ningún envase vacío.
- Cierre y apriete **bien** la tapa del envase de muestra.
- Si usa dos envases, empáquelos por separado.
- EXAMINADOR: REVISE LAS INSTRUCCIONES DE EMPAQUE
- EXAMINADOR: EXPLÍQUELE A LA PERSONA A QUIEN DA LA MUESTRA, LA INFORMACIÓN QUE TIENE QUE PONER EN LA ETIQUETA
- Si usa dos envases para la mañana, escriba “**2 cups**” en ambas cajas.

Envíe las muestras:

Lleve los paquetes a de la Oficina de Correos y entréguelos al cartero o, póngalos en el buzón para paquetes grandes y envíos prioritarios. (El buzón para paquetes grandes y envíos prioritarios está habitualmente ubicado cerca del quiosco de autoservicio).

Si el laboratorio no recibe las muestras en 14 días, se le enviará una tarjeta para recordárselo. Si usted ya ha enviado su muestra de orina, ignore esa tarjeta.

No envíe al laboratorio la bolsa plástica exterior con su nombre.

Le mandaremos un cheque por \$50.00 dólares tan pronto como envíe el paquete.

Appendix D

Inventory of Home Urine Collection Supplies

Home Urine Collection

Consumable

Part #	Description	Par	Unit	Count
EXAM SUPPLIES				
ConLab-Eng	Confidential Label - English - White Handle Bag Label	1	box	<input type="text"/>
ConLab-Span	Confidential Label - Spanish - White Handle Bag Label	1	box	<input type="text"/>
GOJ-9639 12CT	Handsoap-Waterless (Purell)	1	bottle	<input type="text"/>
HUK -MOON	Home Urine Kit - MOON	120	each	<input type="text"/>
HUK - SUN	Home Urine Kit - SUN	240	each	<input type="text"/>
HUC-IS- Eng/Span	Instruction Sheet - Eng/Span Combo - Stapled Packet	180	each	<input type="text"/>
S-596	Mini-Grip Bag 10" x 13" - Zip-Lock Bag, Clear	4	bundle	<input type="text"/>
AVE 5126	Online Shipping Label - Half Sheet - 5 1/2" x 8 1/2"	360	each	<input type="text"/>
12-23W	Poly Tote Hand Bag - White - 20 x 20 x 5	2	case	<input type="text"/>
TM0103-C	Temperature Monitoring Strip - Range C: 54-65C	20	each	<input type="text"/>
06-66997	Towelettes, Sani-Cloth (Cleaner)	1	bottle	<input type="text"/>
26-32	Twist Tie - 7" White (Trashbag Tie)	1	box	<input type="text"/>
Postage Stamp	USPS Postage Stamp - Single Stamp	360	each	<input type="text"/>
OFFICE SUPPLIES				
107708	Dymo Costar Label (1" x 2 1/8") - 500 labels per roll	4	box	<input type="text"/>
AVE 08888	Magic Markers, Black	2	each	<input type="text"/>
UNV72220BX	Paper Clip, Jumbo	2	box	<input type="text"/>
84417-10053	SORTKWIK - Finger Moistener - 3/8 ounce container	2	bottle	<input type="text"/>
MMM 62003/4X1296	Tape, 3/4" Scotch	2	roll	<input type="text"/>
COMPONENT TOOLS				
53382	Compressed Air Spray Can - 10 Oz. 3M Dust Remover	1	bottle	<input type="text"/>
MEC SUPPLIES				
PC221670	Tissues - Facial	2	box	<input type="text"/>

* Expired Lot

Home Urine Collection

Non-Consumable

Part #	Description	Par	Unit	Count
OFFICE SUPPLIES				
MMM C38-BK	Tape Dispenser, Scotch	1	each	<input type="text"/>
COMPONENT TOOLS				
14783415	E-Z Fold Step Stool - 12" Black	1	each	<input type="text"/>
ML4042	Footstool - Metal	1	each	<input type="text"/>
08-952	Scissors - Lister bandage (5 1/2 " Stainless Steel)	1	each	<input type="text"/>
30.60	Temperature Strip Storage Box w/ Lid - 250mL	1	each	<input type="text"/>

* Expired Lot

Appendix E

Reminder Postcards English and Spanish

Reminder

NHANES Urine Collection Kit



You recently completed a health examination as part of the National Health and Nutrition Examination Survey ("NHANES"). You received urine collection kits before you left the exam center.

Please remember to mail the urine samples to the laboratory using the kits that were given to you.

Once you mail the urine samples, we will mail you a check for \$50.00, as promised.

We appreciate your participation and interest in NHANES.

If you have any questions, please call the toll-free number 1-888-458-4762.

NHANES
RP 5055F
1600 Research Blvd.
Rockville, MD 20850-9973
Project #8056.00.01.04



Recordatorio

Materiales de NHANES para recolectar muestras de orina



Recientemente usted completó un examen de salud como parte de la Encuesta Nacional de Examen de Salud y Nutrición ("NHANES"). Antes de irse del centro de examen, usted recibió los materiales para recolectar unas muestras de orina.

Por favor, recuerde enviar las muestras de orina al laboratorio usando los materiales que se le dieron.

Tan pronto como envíe las muestras de orina, le enviaremos un cheque por \$50 dólares, como le prometimos.

Agradecemos su participación e interés en NHANES.

Si desea hacer alguna pregunta, por favor llame gratis al número 1-888-458-4762.

NHANES
RP 5055F
1600 Research Blvd.
Rockville, MD 20850-9973
Project #8056.00.01.04

